RX NC SPRINTER®

Instructions For Use

RAPID EXCHANGE BALLOON DILATATION CATHETER

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Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway NE,
Minneapolis, MN55432, U.S.A.
www.medtronic.com

Manufactured At/EC Authorised Representative: or
Medtronic Ireland
Parkmore Business Park West,
Galway, Ireland.
Tel.: +353-91-708000. Fax.: +353-91-757524.

Made in: Medtronic Mexico
S. de R.L. de C.V.
Avenida Paseo del Cupacab 10510
Parque Industrial El Lago
Tijuana, B.C. 22570, Mexico
Tel: +52-664-903-1100. Fax.: +52-664-903-1175

0086

RBP

En.	Rated Burst Pressure
Ja.	最大拡張圧
Fr.	Pression de rupture théorique
De.	Garantierte Belastungsgrenze
It.	Pressione nominale di scoppio
Nl.	Nominale barstdruk
Pt.	Pressão de Ruptura Nominal
El.	Μέγιστη λειτουργική πίεση διάρρηξης
Es.	Presión máxima recomendada
Da.	Estimeret sprængningstryk
Sv.	Märkvärde för språngtryck
Fi.	Nimellinen murtumisaine
Cs.	Jmenovitý tlak protřzení
Et.	Nominaalne lõhkemisrõhk
Hu.	Előrelátható robbanási nyomás:
Lv.	Teorētiskais pārpīšanās spiediens.
Lt.	Nominalus plyšimo spaudimas
Pl.	Nominalne ciśnienie dopuszczalne
Sk.	Stanovený rozpnací tlak
Tr.	Nominal Patlama Basıncı
Zh.	额定爆裂压力

RBP

En.	Do not exceed rated burst pressure
Ja.	最大拡張圧を超えてはいけません
Fr.	Ne pas dépasser la pression de rupture théorique
De.	Garantierte Belastungsgrenze nicht überschreiten
It.	Non superare la pressione massima consigliata.
Nl.	De nominale barstdruk niet overschrijden
Pt.	Não exceda a pressão de ruptura nominal
El.	Μην υπερβείτε τη μέγιστη λειτουργική πίεση διάρρηξης
Es.	No superar la presión máxima recomendada
Da.	Det estimerede sprængningstryk må ikke overskrides
Sv.	Överskrid inte märkvärdet för språngtryck
Fi.	Älä ylittää nimellistä murtumisainetta
Cs.	Nepřekročte jmenovitý tlak protřzení
Et.	Mitte ületada nominaalset lõhkemisrõhku
Hu.	Ne haladjuk meg az előrelátható robbanási nyomást.
Lv.	Nepārsniedziet teorētisko pārpīšanās spiedienu.
Lt.	Neviršykite nominalaus plyšimo spaudimo.
Pl.	Nie przekraczać nominalnego ciśnienia dopuszczalnego
Sk.	Nepřekročte stanovený rozpnací tlak
Tr.	Nominal patlama basıncını aşmayın.
Zh.	不要超过额定的爆裂压力

En.	Maximum guidewire diameter
Ja.	最大適合ガイドワイヤー径
Fr.	Diamètre maximum du guide
De.	Maximaler Durchmesser des Führungsdrahtes
It.	Diametro massimo della guida
Nl.	Maximale diameter voerdraad
Pt.	Diâmetro máximo do fio-guia
El.	Μέγιστη διάμετρος οδηγού σύρματος
Es.	Diámetro máximo de la guía
Da.	Maksimal ledetrådsdiameter
Sv.	Maximal ledardiameter
Fi.	Ohjainvaijerin enimmäisläpimitta
Cs.	Maximální průměr zaváděcího drátu
Et.	Juhtetraadi maksimaalne läbimõõt
Hu.	A vezetődórt maximális átmérője:
Lv.	Katetra vadīklas maksimālais diametrs
Lt.	Maksimalus nukreipiančio stygos skersmuo
Pl.	Maksymalna średnica prowadnika
Sk.	Maximálny priemer vodiaceho drôtu
Tr.	Maksimum kılavuz kablo çapı
Zh.	最小导丝直径

! USA

En.	For U.S. Audiences only
Ja.	米国のクライアントのみ限定
Fr.	Uniquement pour les États-Unis
De.	Nur USA
It.	Solo per l'uso negli Stati Uniti
Nl.	Alleen voor klanten in de V.S.
Pt.	Apenas para utilização nos E.U.A.
El.	Μόνο για χρήση στις ΗΠΑ
Es.	Sólo para el mercado de EE.UU.
Da.	Gælder kun kunder i USA
Sv.	Endast för amerikansk publik
Fi.	Koskee vain Yhdysvaltoja
Cs.	Jen pro uživatele v USA
Et.	Ainult USA klientidele
Hu.	Csak az Egyesült Államokban alkalmazható
Lv.	Lietošanai tikai ASV
Lt.	Skirta tik JAV
Pl.	Tylko dla odbiorców w USA
Sk.	Len pre spotrebiteľov v USA
Tr.	Yalnızca A.B.D için
Zh.	仅限于美国听众

En.	For Post Dilatation of Balloon Expandable Stents
Ja.	バルーン拡張式ステントの拡張後に使用のこと
Fr.	Dilatation après mise en place de stents expansibles par ballonnet
De.	Für die nachträgliche Dilatation von ballonexpandierbaren Stents
It.	Per espansione post-rilascio di stent espansibili tramite palloncino
Nl.	Voor postdilatie van met een ballon expandeerbare stents
Pt.	Para Dilatação Posterior de Stents Expansíveis por Balão
El.	Για διαστολή, μετά την τοποθέτηση, stent που διαστέλλονται με μπαλόνι
Es.	Para la postdilatación de stents de balón expansibles
Da.	Til efterudvidelse af ballonexpanderbare stenter
Sv.	För dilatation av ballongexpanderande stentar efter insättningen
Fi.	Tarkoitettu pallolaajennettävien stenttien jälkilaajennukseen
Cs.	Pro následnou dilataci stentů roztažitelných balónek
Et.	Balloondilateeritavate stentide järelaiendamiseks
Hu.	A ballonnal tágítható stentek utólagos tágítására
Lv.	Ar balonu izplešamo stentu izplešanai pēc ievietošanas
Lt.	Skirta balioniniams išplečiamiems stentams išplėsti
Pl.	Do rozprężania uprzednio umieszczonych rozprężalnych protez naczyńnych
Sk.	Pre stenty expandovateľné dilatáciou balónika
Tr.	Balonla Genişletilebilen Stentlerin Dilatasyonu için
Zh.	用于可气囊扩张支架的后扩张

En.	Length
Ja.	長さ
Fr.	Longueur
De.	Länge
It.	Lunghezza
Nl.	Lengte
Pt.	Comprimento
El.	Μήκος
Es.	Longitud
Da.	Længde
Sv.	Längd
Fi.	Pituus
Cs.	Délka
Et.	Pikkus
Hu.	Hosszúság:
Lv.	Garums
Lt.	Ilgis
Pl.	Długość
Sk.	Dĺžka
Tr.	Uzunluk
Zh.	长度

En.	Do not use if package is damaged
Ja.	包装が破損している場合は、使用しないでください
Fr.	Ne pas utiliser si l'emballage est endommagé
De.	Bei Schäden an der Verpackung das Produkt nicht verwenden
It.	Non utilizzare il prodotto se la confezione è danneggiata.
Nl.	Niet gebruiken als de verpakking beschadigd is
Pt.	Não utilize se a embalagem estiver danificada
El.	Μη χρησιμοποιείτε εάν η συσκευασία είναι κατεστραμμένη
Es.	No utilizar si el envase está dañado
Da.	Må ikke bruges, hvis emballagen er beskadiget
Sv.	Använd inte om förpackningen är skadad
Fi.	Älä käytä, jos pakkaus on vioittunut
Cs.	Nepoužívejte, je-li obal poškozený
Et.	Mitte kasutada, kui pakend on kahjustatud
Hu.	Ne használja, ha a csomagolás sérült
Lv.	Nelietojiet, ja iepakojums ir bojāts
Lt.	Nenaudoti, jei pakuotė pažeista
Pl.	Nie używać, jeżeli opakowanie jest uszkodzone
Sk.	Nepoužívajte, ak je obal poškodený
Tr.	Ambalaj hasarlıysa kullanmayın
Zh.	如果包装破损, 请勿使用

En.	Consult instructions for use
Ja.	取扱説明書を参照してください
Fr.	Consulter le mode d'emploi
De.	Siehe Gebrauchsanweisung
It.	Consultare le istruzioni per l'uso
Nl.	Raadpleeg de gebruiksaanwijzing
Pt.	Consulte as instruções para utilização
El.	Συμβουλευτείτε τις οδηγίες χρήσης
Es.	Consultar las instrucciones de uso
Da.	Se brugsanvisningen
Sv.	Läs bruksanvisningen
Fi.	Lue käyttöohjeet
Cs.	Podívejte se do návodu k použití
Et.	Enne kasutamist lugege kasutusjuhiseid
Hu.	Olvassa el a használati utasítást
Lv.	Izmantojiet lietošanas instrukcijas
Lt.	Vadovaukites naudojimo instrukcijomis
Pl.	Przed użyciem zapoznać się z instrukcją obsługi
Sk.	Prečítajte si návod na použitie
Tr.	Kullanım talimatlarına bakın
Zh.	参阅使用说明

En.	Do Not Reuse
Ja.	再使用しないでください。
Fr.	À usage unique
De.	Nicht wieder verwenden
It.	Non riutilizzare
Nl.	Niet hergebruiken
Pt.	Não Reutilize
El.	Μην επαναχρησιμοποιείτε
Es.	No reutilizar
Da.	Kun til engangsbrug
Sv.	Får ej återanvändas
Fi.	Tuotetta ei saa käyttää uudelleen
Cs.	Jen pro jednorázové použití
Et.	Mitte kasutada korduvalt
Hu.	Ne használja újra!
Lv.	Atkārtoti nelietot
Lt.	Nenaudokite pakartotinai
Pl.	Nie używać ponownie
Sk.	Nepoužívajte na viackrát
Tr.	Tekrar Kullanmayın
Zh.	不可重复使用

En.	Use by
Ja.	使用期限
Fr.	Date de péremption
De.	Zu verwenden bis
It.	Data di scadenza
Nl.	Uiterste gebruiksdatum
Pt.	Use antes de
El.	Ημερομηνία λήξης
Es.	Utilizar antes de
Da.	Anvendes inden
Sv.	Använd före
Fi.	Viimeinen käyttöpäivämäärä
Cs.	Spotřebovat do
Et.	Kasutada kuni
Hu.	Felhasználható:
Lv.	Derīgs līdz
Lt.	Naudoti iki
Pl.	Do użytku przed
Sk.	Použite najneskôr do
Tr.	Son kullanma tarihi
Zh.	有效期

LOT

En.	Lot number
Ja.	ロット番号
Fr.	Numéro du lot
De.	Chargennummer
It.	Numero lotto
Nl.	Partijnummer
Pt.	Número do lote
El.	Αριθμός παρτίδας
Es.	Número de lote
Da.	Lotnr.
Sv.	Lottnummer
Fi.	Eränumero
Cs.	Číslo šarže
Et.	Partinumber
Hu.	Tételszám:
Lv.	Partijas numurs
Lt.	Serijos numeris
Pl.	Numer serii
Sk.	Šarža číslo
Tr.	Lot numarası
Zh.	产品批号

REF

En.	Product Code
Ja.	製品コード
Fr.	Référence
De.	Katalognummer
It.	Numero di catalogo
Nl.	Bestelnummer
Pt.	Código do Produto
El.	Κωδικός προϊόντος
Es.	Código del producto
Da.	Produktkode
Sv.	Produktkod
Fi.	Luettelonumero
Cs.	Kód výrobku
Et.	Tootekood
Hu.	Termékkód:
Lv.	Izstrādājuma kods
Lt.	Gaminio kodas
Pl.	Kod produktu
Sk.	Kód výrobku
Tr.	Ürün Kodu
Zh.	产品代码


En. Manufacturer
 Ja. 製造元
 Fr. Fabricant
 De. Hersteller
 It. Produttore
 Ni. Fabrikant
 Pt. Fabricante
 El. Κατασκευαστής
 Es. Fabricante
 Da. Producent
 Sv. Tillverkare
 Fi. Valmistaja
 Cs. Výrobce
 Et. Tootja
 Hu. Gyártó
 Lv. Ražotājs
 Lt. Gamintojas
 Pl. Producent
 Sk. Výrobca
 Tr. Üretici
 Zh. 制造商

EC|REP

En. EC Authorised Representative
 Ja. EC 認定代理店
 Fr. Représentant CE Autorisé
 De. Autorisierte EU-Vertretung
 It. Rappresentante autorizzato per l'Europa
 Ni. Erkend vertegenwoordiger in de EU
 Pt. Representante Autorizado na CE
 El. Εξουσιοδοτημένος Αντιπρόσωπος στην Ε.Ε.
 Es. Representante autorizado para la CE
 Da. Autoriseret EU-repræsentant
 Sv. Auktoriserad EU-representant
 Fi. Valtuutettu edustaja EY-alueella
 Cs. Autorizovaný zástupce pro ES
 Et. EÜ ametlik esindaja
 Hu. EU meghatalmazott képviselő
 Lv. EK autorizēts pārstāvis
 Lt. EC įgaliotasis atstovas
 Pl. Autoryzowany przedstawiciel we WE
 Sk. Autorizovaný zástupca pre EÚ
 Tr. AB Yetkili Temsilcisi
 Zh. EC 授权代表

STERILE|R

En. Sterilized using irradiation.
 Ja. 放射線滅菌済み
 Fr. Stérilisé par irradiation.
 De. Sterilisiert durch Bestrahlung.
 It. Il prodotto è stato sterilizzato per irradiazione.
 Ni. Gesteriliseerd door middel van bestraling.
 Pt. Esterilizado por irradiação.
 El. Αποστείρωση με ακτινοβολίαση.
 Es. Esterilizado mediante irradiación.
 Da. Steriliseret ved hjælp af bestråling.
 Sv. Steriliserat med stråling.
 Fi. Steriloitu vain säteilyttämällä.
 Cs. Sterilizovano zářením.
 Et. Steriliseeritud kiiritusega.
 Hu. Irradiáció útján sterilizált.
 Lv. Sterilizēts apstarojot.
 Lt. Sterilizuota naudojant iradiaciją.
 Pl. Sterylizacja w komorze radiacyjnej.
 Sk. Sterilizované ožiarením.
 Tr. Işinlama ile sterilize edilmiştir.
 Zh. 使用射线灭菌。

En. Rapid Exchange Balloon Dilatation Catheter
 Ja. ラピッド・エクスチェンジ・バルーン拡張カテーテル
 Fr. Cathéter de dilatation à ballonnet à échange rapide
 De. Rapid Exchange Ballondilatationskatheter
 It. Catetere dilatatore a palloncino a cambio rapido
 Ni. Ballondilatatiekatheter voor snelle uitwisseling
 Pt. Cateter de Dilatação por Balão de Troca Rápida
 El. Καθετήρας μπαλονιού διαστολής ταχείας αναλλαγής
 Es. Catéter balón de dilatación de intercambio rápido
 Da. Hurtigt udskifteligt ballondilatationskatheter
 Sv. Ballongdilatationskatheter för snabbt utbyte
 Fi. Pikavaihtojärjestelmällä varustettu pallolajennuskatetri
 Cs. Balónková dilatáční katétra umožňující rychlou výměnu
 Et. Kiirivahetus-ballondilatatsioonikateeter
 Hu. Rapid Exchange dilatációs ballonkatéter
 Lv. Ātri nomaināms balona paplašināšanas katetris
 Lt. Greitai pakeičiamas balioninis išplėtimo kateteris
 Pl. Cewnik z balonem typu Rapid Exchange
 Sk. Rýchlovýmenný dilatáčný katéter s balónikom
 Tr. Hızlı Değiştirilen Balonlu Dilatasyon Kateteri
 Zh. 快速交换气囊扩张导管

ENGLISH

Medtronic Balloon Catheter Description:

The NC Sprinter Rapid Exchange Balloon Dilatation Catheter is a Percutaneous Transluminal Coronary Angioplasty (PTCA) Rapid Exchange System. The balloon at the distal end of the catheter can be inflated to a defined diameter at a specific pressure (see labelling). The proximal end of the catheter has a female luer for attachment to an inflation device. The catheter provides a lumen which enables the use of a guide wire to position the catheter.

Radiopaque balloon marker(s) enable accurate placement. Shaft markers for brachial and femoral techniques are in place.

Contents:

- One (1) Balloon Dilatation Catheter
- One (1) Flushing Cannula*
- One (1) Re-Fold Tool*
- One (1) Looper™ Device*
- * Components contained in separate pouch

Indications:

- The balloon dilatation catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.
- The balloon dilatation catheter (balloon models 2.25 - 5.0mm) is also indicated for post-delivery expansion of balloon expandable stents.

Note: Bench testing was conducted with the NC Sprinter Rapid Exchange Balloon Dilatation Catheter and marketed balloon expandable stents. Consideration should be taken when this device is used with different manufacturers' stents due to difference in stent design.

Contra-Indications:

The catheter is contraindicated for use in:

- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of significant stenosis.

Warnings:

- For single patient, single procedure use only. Do NOT resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.**
- Do not use the catheter if its package has been opened or damaged.**
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.**
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.**
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.**
- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization.**
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.**
- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.**
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.**

Precautions:

- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the procedure for which it is to be used.
- If the device is kinked, it should not be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty.

- Appropriate anticoagulation, antiplatelet and vasodilator therapy should be administered to the patient.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn from the patient before removing any additional equipment.
- Care should be taken not to apply excessive force during preparation or use, as this may damage the device.

Adverse Effects:

Possible adverse effects include, but are not limited to the following:

- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel
- Hemorrhage or hematoma
- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast medium
- Hypo/hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Embolism

Materials To Be Used In Combination With A Balloon Catheter Include:

- Suitable guide wire, refer to label claim.
- 20 cc syringe for balloon preparation.
- Suitable guiding catheter, refer to label claim.
- 10 cc or smaller syringe for manual dye injections.
- Appropriate inflation medium (e.g. 50:50 sterile mixture of a contrast medium and saline.)
- Pressure-indicating inflation device.
- Hemostasis valve.

Directions For Use:

1. Preparation

- Select an appropriate balloon catheter for the target vessel.
- Remove the device from the sterile packaging.

Note: This device allows the guide wire lumen to be flushed in two(2) formats:

- With device contained within its hoop (Ref. In-Hoop Flushing section).
- With device removed from its hoop:
 - Remove the balloon protector (C) from around the balloon. Removing the balloon protector will also remove the balloon stylette.
 - Fill a 20cc syringe with 20cc of a saline solution.
 - Remove the flushing cannula (E) from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
 - Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
 - Insert the flushing cannula into the guidewire entry port (exchange joint).
 - Depress the syringe to flush the guide wire lumen until fluid emerges from the tip (distal end of the balloon catheter).
 - Remove the flushing cannula from the guide wire entry port (exchange joint).

Note: Do not discard the flushing cannula until the end of the angioplasty procedure as additional flushing may be required.

2. Balloon Purging

Purge air from the catheter using a 20 cc syringe filled with 2 to 3 cc of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a meniscus of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. A flow of bubbles will move from the balloon catheter into the inflation device. Balloon purging is complete when the flow of bubbles stops. Do not attempt Pre-Inflation technique to purge the balloon lumen.

3. Insertion Technique

- Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.
- Advance a guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit port (exchange joint).
- The hemostasis valve should be gradually tightened to control back flow. Excessive valve tightening may affect balloon inflation/deflation time as well as movement of the guide wire.
- Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion.



4. Balloon inflation

- Inflate the balloon to dilate the lesion using standard PTCA techniques.
- After each subsequent inflation, the distal blood flow should be assessed.
- If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do not exceed the rated burst pressure (see labelling).
- Confirm the results with fluoroscopy.

5. Removing the Catheter

- Apply negative pressure to the inflation device and confirm that the balloon is fully deflated.
- Withdraw the balloon catheter into the guiding catheter while preserving guide wire position and adjusting the hemostasis valve appropriately. Remove the balloon catheter from the system.
- After the deflated balloon catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
- Inspect the balloon catheter integrity.
- If reinserting the same balloon catheter, flush the guide wire lumen by using the flushing cannula to flush saline through the guide wire entry port (exchange joint) until fluid emerges from the tip (distal end of the balloon catheter). Prior to reinsertion, the balloon catheter should be wiped clean with gauze soaked with sterile normal saline. The balloon may be refolded using the Re-Fold Tool as described in the Re-Fold Tool Instructions for use.

Note: Do not discard the balloon stylette until the end of the angioplasty procedure as it may be required for balloon refolding.

In-Hoop Flushing

The in-hoop flushing cannula (E) enables the guidewire lumen to be flushed without removing the catheter from its packaging hoop

Instructions for use:

- With the balloon catheter contained within its hoop (A) manually detach the blue tubing (B).
- Remove the balloon protector (C) from around the balloon. Removing the balloon protector will also remove the balloon stylette.

Note: Do not discard the balloon stylette until the end of the angioplasty procedure as it may be required for balloon refolding.

- Fill a 20cc syringe with 20cc of a saline solution.
- Remove the flushing cannula from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
- Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
- Insert the flushing cannula into the tip (distal end of the balloon catheter).
- Depress the syringe to flush the guide wire lumen, until fluid emerges from the guide wire entry port (exchange joint). Droplets will be visible in the hoop.
- Remove the flushing cannula from the distal tip of the balloon.

Note: Do not discard the flushing cannula until the end of the angioplasty procedure as additional flushing may be required.

Re-Fold Tool

This Re-Fold Tool (F) is an accessory component that allows the balloon to be re-folded if required.

Instructions for use:


- Deflate the balloon by applying negative pressure to the inflation device and maintain under vacuum.
- Visually inspect the balloon to confirm that it is fully deflated.
- Remove the Re-Fold Tool from the accessories pouch (D).
- Ensure the balloon protector is not on the balloon stylette and load the non-flared end of the Re-Fold Tool onto the straight end of the stylette.
- Carefully load the stylette back through the distal tip of the catheter and past the proximal end of the balloon.
- While holding the catheter just proximal to the balloon, push the re-fold device over the balloon in a gentle twisting motion until the entire balloon is covered.
- Gently remove the Re-Fold Tool and stylette.
- Inspect the balloon for any potential damage. Discard the balloon catheter if there is any visible damage present on the balloon.

Looper™ Device

The LOOPER™ device (G) is an accessory component for use with Medtronic Rapid Exchange balloon catheters. The device allows the catheter to be fastened in a coiled configuration for ease of handling during use.

Instructions for use:

- Remove the LOOPER™ device from the accessories pouch (D).
- Form the catheter into a single or double loop when required.
- Press down the wings of the LOOPER™ device to open the clipping section for retention of the hypotube (H).
- Release the wings of the LOOPER™ device with the hypotube inside the clipping section to secure the catheter in a coiled configuration.

Disclaimer of Warranty

NOTE: ALTHOUGH THE CATHETER, THEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC INC. AND ITS AFFILIATES HAS NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC INC. AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC INC. AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC INC. AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

References

The physician should consult recent literature on current medical practice on balloon dilatation, such as published by American College of Cardiology / American Heart Association.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Patents:

Protected under one or more of the following patents:
U.S. Patents 7,001,358, 6,190,358 and 6,605,057. U.S. and foreign patents pending.

Compliance Chart:

The Compliance Chart is based on actual data from in vitro testing @ 37° C, rounded to two decimal places. See Page 29 for Compliance Chart.

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