# MITRA CLIP®

## INSTRUCTIONS FOR USE

### MITRA CLIP SYSTEM

Steerable Guide Catheter  
Clip Delivery System

Ref No. SGC01ST  
Ref No. CDS02ST

### MITRA CLIP SYSTEM ACCESSORIES

Stabilizer  
Lift  
Support Plate

Ref No. SZR01ST  
Ref No. LFT01ST  
Ref No. PLT01ST

**WARNING:** Read all instructions carefully. Failure to follow these instructions, warnings and precautions may lead to device damage or patient injury. Use of the MitraClip System should be restricted to those physicians trained to perform invasive endovascular and transseptal procedures and to those physicians trained in the proper use of the system.

### TABLE OF CONTENTS

INTENDED USE .....	2
CONTRAINDICATIONS .....	2
WARNINGS .....	2
PRECAUTIONS .....	2
POTENTIAL COMPLICATIONS AND ADVERSE EVENTS .....	3
PATIENT COUNSELING .....	4
HOW SUPPLIED .....	4
STORAGE .....	4
MITRA CLIP SYSTEM DIMENSIONS .....	4
GLOSSARY OF ACRONYMS .....	5
DEVICE DESCRIPTION .....	5
MITRA CLIP SYSTEM ACCESSORIES OVERVIEW .....	6
REQUIRED ACCESSORIES .....	8
ADDITIONAL REQUIRED EQUIPMENT NOT INCLUDED .....	8
MITRA CLIP PROCEDURE STEP-BY-STEP INSTRUCTIONS .....	8
1.0 DEFINITION OF TERMS .....	8
2.0 PATIENT PREPARATION .....	9
3.0 MITRA CLIP SYSTEM PREPARATION BEFORE USE .....	9
4.0 ACCESS TO THE MITRAL VALVE .....	12
5.0 STEERABLE GUIDE CATHETER INSERTION .....	12
6.0 CLIP DELIVERY SYSTEM INSERTION .....	13
7.0 INITIAL MITRA CLIP SYSTEM POSITIONING IN THE LEFT ATRIUM .....	13
8.0 FINAL MITRA CLIP SYSTEM POSITIONING .....	14
9.0 GRASPING THE LEAFLETS AND VERIFYING THE GRASP .....	14
10.0 CLOSING THE CLIP AND EVALUATING CLIP POSITION .....	15
11.0 MITRA CLIP DEVICE PRE-DEPLOYMENT CLIP ASSESSMENT .....	15
12.0 CLIP DEPLOYMENT .....	16
13.0 ADDITIONAL MITRA CLIP DEVICE PLACEMENT .....	17
14.0 MITRA CLIP SYSTEM REMOVAL .....	17
PATENTS .....	18
NON-HARMONIZED SYMBOLS USED FOR LABELING .....	168

RELEASED

**INTENDED USE**

The MitraClip System is intended for reconstruction of the insufficient mitral valve through tissue approximation.

**CONTRAINDICATIONS**

Patients with the following conditions should not be treated with the MitraClip System:

- Patients who cannot tolerate procedural anticoagulation or post procedural anti-platelet regimen
- Active endocarditis of the mitral valve
- Rheumatic mitral valve disease
- Evidence of intracardiac, inferior vena cava (IVC) or femoral venous thrombus

**WARNINGS**

The following anatomic considerations may result in suboptimal leaflet insertion and/or MitraClip implantation. Patients with these conditions who undergo the MitraClip procedure may have an increased risk of serious adverse events which may be minimized with proper device usage and echocardiographic assessment. These events include the potential for increased procedure time, tissue injury, and/or worsening MR which may require additional MitraClip intervention or surgical treatment.

- Severe leaflet calcification in the grasping area
- Severe calcification of the annulus and/or subvalvular apparatus, such as the chordae tendinae
- Severely restricted posterior leaflet
- Cleft or perforation in the grasping area
- Leaflet Flail gap  $\geq 10$ mm and/or Leaflet Flail width  $\geq 15$ mm
- Coaptation length  $< 2$ mm
- Intracardiac mass

The following anatomic considerations may result in serious adverse events, including tissue injury, worsening MR or in the case of the mitral valve area, mitral stenosis. If these events were to occur they can require additional percutaneous or surgical intervention for treatment. The occurrence of these events can be minimized with appropriate preoperative evaluation of the valve anatomy and regurgitant jet.

- Mitral valve orifice area  $< 4$  cm<sup>2</sup>
- Primary regurgitant jet outside of the A2-P2 area and/or presence of a clinically significant 2nd jet

Patients with the following considerations in whom the Steerable Guide Catheter is used may have an increased risk of having a serious adverse event which may be avoided with preoperative evaluation and proper device usage.

- Previous interatrial septal patch or prosthetic atrial septal defect (ASD) closure device which could result in significant difficulty in visualization or technical challenges during transseptal puncture and/or introducing the SGC into the left atrium.
- Known or suspected unstable angina or myocardial infarction within the last 12 weeks could increase the procedural morbidity and mortality, due to increased hemodynamic stress secondary to general anesthesia.
- Patients with active infection have an increased risk of developing an intraoperative and/or postoperative infection, such as sepsis or soft tissue abscess.
- Known or suspected left atrial myxoma could result in thromboembolism and tissue injury due to difficulty with device positioning.
- Recent cerebrovascular event (CVA) may increase the procedural morbidity associated with a transcatheter intervention, such as recurrent stroke.

The MitraClip Device should be implanted with sterile techniques using fluoroscopy and echocardiography (e.g. transesophageal [TEE] and transthoracic [TTE]) in a facility with on-site cardiac surgery and immediate access to a cardiac operating room.

Read all instructions carefully. Failure to follow these instructions, warnings and precautions may lead to device damage, user injury, or patient injury. Use universal precautions for biohazards and sharps while handling the MitraClip System to avoid user injury.

Use of the MitraClip System should be restricted to those physicians trained to perform invasive endovascular and transseptal procedures and those trained in the proper use of the system.

For the Clip Delivery System and Steerable Guide Catheter only: this device is designed for single use only. Cleaning, re-sterilization and/or re-use may result in infections, malfunction of the device and other serious injury or death.

**PRECAUTIONS**

NOTE the "Use by" date specified on the package.

Inspect all product prior to use. Do not use if package is opened or damaged.

The inside of the outer pouch is not a sterile barrier. The inner pouch within the outer pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

**POTENTIAL COMPLICATIONS AND ADVERSE EVENTS**

The following ANTICIPATED EVENTS have been identified as possible complications of the MitraClip procedure.

Allergic reaction (anesthetic, contrast, Heparin, nickel alloy, latex)	Hypotension/hypertension
Aneurysm or pseudo-aneurysm	Infection and pain at insertion site
Arrhythmias	Infection and pain at incision site
Atrial fibrillation	Lymphatic complications
Atrial septal defect requiring intervention	Mesenteric ischemia
Arterio-venous fistula	Mitral stenosis
Bleeding	Mitral valve injury
Cardiac arrest	MitraClip erosion, migration or malposition
Cardiac perforation	MitraClip Device thrombosis
Cardiac tamponade/Pericardial Effusion	MitraClip System component(s) embolization
Chordal entanglement/rupture	Multi-system organ failure
Coagulopathy	Myocardial infarction
Conversion to standard valve surgery	Nausea/vomiting
Death	Pain
Deep venous thrombus (DVT)	Peripheral ischemia
Dislodgement of previously implanted devices	Prolonged angina
Dizziness	Prolonged ventilation
Drug reaction to anti-platelet/anticoagulation agents/contrast media	Pulmonary congestion
Dyskinesia	Pulmonary thrombo-embolism
Dyspnea	Renal insufficiency or failure
Edema	Respiratory failure/atelectasis/pneumonia
Emboli (air, thrombus, MitraClip Device)	Septicemia
Emergency cardiac surgery	Shock, Anaphylactic or Cardiogenic
Endocarditis	Single leaflet device attachment (SLDA)
Esophageal irritation	Skin injury or tissue changes due to exposure to ionizing radiation
Esophageal perforation or stricture	Stroke or transient ischemic attack (TIA)
Failure to deliver MitraClip to the intended site	Urinary tract infection
Failure to retrieve MitraClip System components	Vascular trauma, dissection or occlusion
Fever or hyperthermia	Vessel spasm
Gastrointestinal bleeding or infarct	Vessel perforation or laceration
Hematoma	Worsening heart failure
Hemolysis	Worsening mitral regurgitation
Hemorrhage requiring transfusion	Wound dehiscence

**PATIENT COUNSELING**

Patients undergoing any procedures known to potentially be associated with bacteremia after implantation of the MitraClip Device should be prescribed prophylactic antibiotic therapy prior to such procedures.

Short-term anticoagulation therapy may be necessary after mitral valve repair with the MitraClip Device. Prescribe anticoagulation and other medical therapy per institutional guidelines.

After placement of a MitraClip Device, the Implant Identification Card should be filled out and the patient should be instructed to carry it at all times.

All patients should be advised to limit strenuous physical activity for at least the first month post-procedure or longer if warranted.

Physicians should consider the following in counseling patients about the MitraClip Device:

- Discuss the risks associated with MitraClip Device placement.
- Discuss the risk/benefit considerations for the patient.

**HOW SUPPLIED**

**Contents:**

CDS02ST: One (1) MitraClip Device, one (1) Clip Delivery System, one (1) MitraClip Device Implant Card, one (1) Instructions for Use.

SGC01ST: One (1) Steerable Guide Catheter, one (1) Dilator, one (1) Silicone Pad, three (3) Fasteners, and one (1) Instructions for Use.

**Sterile:**

For the Clip Delivery System and Steerable Guide Catheter only: these devices are provided sterile, in a thermoformed tray with lid, in sealed pouches. Parts of the devices that are in either direct or indirect contact with circulating blood are non-pyrogenic. Confirm the "Use by" date specified on the package. DO NOT use the MitraClip System if the "Use by" date has passed. These devices are intended for single use only. DO NOT re-sterilize. DO NOT use if the package is opened or damaged. Cleaning, re-sterilization and/or re-use may result in infections, malfunction of the device and other serious injury or death.

The white Guide tip shape retainer and transparent protective tubing are provided sterile and pre-installed on the distal tip of the Steerable Guide Catheter. The Fasteners and the Silicone Pad used with the Stabilizer are provided sterile with the Steerable Guide Catheter. The Dilator, Fasteners and the Silicone Pad are intended for single use only. DO NOT re-sterilize. Cleaning, re-sterilization and/or re-use may result in infections, malfunction of the device and other serious injury or death.

**Non-sterile:**

The Stabilizer, Support Plate and Lift are provided non-sterile. Follow the cleaning and sterilization instructions provided with the Stabilizer, Support Plate and Lift.

**STORAGE**


Handle with care. Store in original packaging. Keep dry. Keep away from sunlight.

**MITRACLIP SYSTEM DIMENSIONS**

Component	Dimension
<i>Delivery Catheter</i>	
Extended Length (from Sleeve curved at 90 degrees)	45 mm – 70 mm
Catheter Shaft Outer Diameter	3.4 mm (10 Fr)
<i>Steerable Sleeve</i>	
Working Length	1095 mm
Catheter Distal Shaft Outer Diameter	5.3 mm (16 Fr)
<i>MitraClip Device</i>	
Closed Clip Length (Figure 1A)	15 mm maximum
Grasping Width at 120 degrees (Figure 1B)	17 mm minimum
Clip Width at 180 degrees (Figure 1C)	20 mm maximum
Arm Width (Figure 1D)	5 mm maximum
Arm Length (Coaptation Length) (Figure 1E)	9 mm maximum
<i>Steerable Guide Catheter</i>	
Working Length	800 mm
Catheter Shaft Inner Diameter	5.5 mm (16 Fr)
Catheter Shaft Outer Diameter	8.1 mm (24 Fr)
Catheter Distal Tip Diameter	7.7 mm (23 Fr)
Catheter Septal Crossing Diameter	7.4 mm (22 Fr)
<i>Dilator</i>	
Working Length	1220 mm
Shaft Inner Diameter	1.0 mm (3 Fr)
Shaft Outer Diameter	5.4 mm (16 Fr)
Distal Tip Outer Diameter	1.5 mm (4 Fr)

RELEASED

Figure 1: MitraClip Device Dimensions


**GLOSSARY OF ACRONYMS**

Guide:	Steerable Guide Catheter	RA:	Right Atrium
CDS:	Clip Delivery System	LA:	Left Atrium
Sleeve:	Steerable Sleeve	LV:	Left Ventricle
DC:	Delivery Catheter	RO:	Radiopaque
Clip:	MitraClip Device (Implant)	MR:	Mitral Regurgitation

**DEVICE DESCRIPTION**

The MitraClip System consists of two parts: 1) the Clip Delivery System and 2) the Steerable Guide Catheter. The Clip Delivery System consists of three major components: 1) the Delivery Catheter 2) the Steerable Sleeve and 3) the MitraClip Device. The Clip Delivery System is introduced into the body through a Steerable Guide Catheter which includes a dilator. The Clip Delivery System and Steerable Guide Catheter constitute the MitraClip System.

The Clip Delivery System (Figures 2, 4 and 5) is used to advance and manipulate the implantable MitraClip Device for proper positioning and placement on the mitral valve leaflets. The Clip Delivery System is designed to deploy the implant in a way that requires multiple steps to ensure safe delivery of the device.

The outer surfaces of the Delivery Catheter and the Steerable Guide Catheter have a hydrophilic coating.

The MitraClip Device (Figure 6) is a percutaneously implanted mechanical Clip. The MitraClip Device grasps and coapts the mitral valve leaflets resulting in fixed approximation of the mitral leaflets throughout the cardiac cycle. The MitraClip Device is placed without the need for arresting the heart or cardiopulmonary bypass. The implantable MitraClip Device is manufactured with metal alloys and polyester fabric (Clip cover) that are commonly used in cardiovascular implants.

The MitraClip Device arms can be adjusted to any position from fully opened, fully inverted and fully closed. These positions are designed to allow the MitraClip Device to grasp and approximate the leaflets of the mitral valve using controls on the Delivery Catheter Handle. The MitraClip Device can be locked, unlocked and repeatedly opened and closed. The Grippers can be raised or lowered repeatedly.

The MitraClip Device is magnetic resonance conditional. It is considered safe for patients undergoing magnetic resonance imaging procedures with:

- Static magnetic field up to 3 Tesla;
- Maximum spatial gradient in static field of 2500 gauss/cm or less;
- Maximum whole-body averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning.

Magnetic resonance image quality may be compromised if the area of interest is in the exact same area, or relatively close to the MitraClip Device. A maximum image artifact of 60 x 70 mm was measured in testing conducted in a 3T magnetic resonance system. It may be necessary to optimize the magnetic resonance imaging parameters due to the presence of the implant.


The MitraClip Device can be removed using standard surgical techniques (and/or according to applicable study protocols) and can be disposed of according to institutional guidelines (and/or applicable study protocols).

The Steerable Guide Catheter (Figure 3a) is used to introduce the Clip Delivery System into the left side of the heart through the interatrial septum. The Steerable Guide Catheter is also used to position and orient the Clip Delivery System to the appropriate location above the mitral valve. The Dilator (Figure 3b) is used for the introduction of the Steerable Guide Catheter into the femoral vein and left atrium.

**MITRACLIP SYSTEM ACCESSORIES OVERVIEW**

Several accessories are used in conjunction with the MitraClip System including: 1) a Stabilizer, 2) a Lift 3) a Support Plate, 4) a Silicone Pad and 5) Fasteners. The Stabilizer is provided separately as a non-sterile reusable device and must be cleaned and sterilized prior to each use. The Stabilizer is used on the sterile field to support and position the Steerable Guide Catheter and Clip Delivery System during the procedure. The Lift and Support Plate are provided separately as non-sterile reusable devices and must be cleaned prior to each use. The Lift and Support Plate are used outside the sterile field to provide a stable platform for the Stabilizer and MitraClip System during the procedure. Follow the cleaning and sterilization instructions provided with the Stabilizer, Support Plate and Lift. The Silicone Pad and Fasteners are single use accessories and are provided sterile with the Steerable Guide Catheter packaging. The Silicone Pad is used on the sterile field under the Stabilizer to prevent incidental movement of the Stabilizer during the procedure. The Fasteners are used on the sterile field to secure the Steerable Guide Catheter and Clip Delivery System to the Stabilizer.

Legend of Figure Labels	
<p><b>Figure 2: Clip Delivery System (CDS)</b></p> <p>1 Delivery Catheter Handle 2 Delivery Catheter Fastener 3 A/P Knob 4 M/L Knob 5 Steerable Sleeve Handle 6 Clip Introducer 7 MitraClip Device</p> <p><b>Figure 3a: Steerable Guide Catheter</b></p> <p>8 Hemostasis Valve 9 Alignment Marker 10 Flush Port 11 +/- Knob 12 Proximal Shaft 13 Distal Shaft 14 Radiopaque Tip Ring</p> <p><b>Figure 3b: Dilator</b></p> <p>15 Rotating Hemostatic Valve 16 Flush Port 17 Echogenic Spiral Groove</p> <p><b>Figure 4: CDS Handles</b></p> <p>18 Actuator Knob 19 Release Pin 20 Arm Positioner 21 Lock Lever Cap 22 Gripper Lever Cap 23 Lock Lever 24 Gripper Lever</p>	<p>25 Delivery Catheter Top Flush Port (Bottom Flush Port Not Shown) 26 Delivery Catheter Handle 27 Delivery Catheter Fastener 28 Sleeve Flush Port 29 A/P Knob 30 M/L Knob 31 Steerable Sleeve Handle</p> <p><b>Figure 5: CDS Distal End</b></p> <p>32 Longitudinal Alignment Marker 33 Key 34 Steerable Sleeve Shaft 35 Radiopaque Alignment Markers 35a Proximal 35b Distal 36 Sleeve Radiopaque Tip Ring 37 Delivery Catheter Shaft 38 Delivery Catheter Radiopaque Ring 39 MitraClip Device</p> <p><b>Figure 6: MitraClip Device Positions</b></p> <p>A Clip fully closed (low profile) B Clip opened to 180 degrees C Clip closed to 120 degrees D Clip closed to 60 degrees E Clip closed to 20 degrees F Clip inverted G Clip fully inverted</p>

RELEASED

Figure 4: CDS Handles


Figure 5: CDS Distal End


Figure 6: MitraClip Device Positions

A: Clip fully closed

B: Clip opened to 180 degrees

C: Clip closed to 120 degrees

D: Clip closed to 60 degrees

E: Clip closed to 20 degrees

F: Clip inverted

G: Clip fully inverted

The Steerable Guide Catheter and Clip Delivery System (Steerable Sleeve, Delivery Catheter and Clip) are steered and actuated by the use of control knobs, levers and fasteners located on the handles.

Table 1: MitraClip System Handle Controls

Device	Control	Function
Steerable Guide Catheter	+/- Knob	Tip deflection
Steerable Sleeve	M/L Knob	Tip deflection
	A/P Knob	Tip deflection
Delivery Catheter	Lock Lever	Locks-Unlocks Clip via Lock Line
	Gripper Lever	Raises-Lowers Grippers via Gripper Line
	Arm Positioner	Opens-Closes Clip Arms
	DC Fastener	Locks-Unlocks DC translation/torque
	Lock Lever Cap	Controls removal of Lock Line
	Gripper Lever Cap	Controls removal of Gripper Line
	Release Pin Actuator Knob	Prevents Actuator Knob turning Turns actuator shaft for Clip deployment

RELEASED

**REQUIRED ACCESSORIES**

SZR01ST: One (1) Stabilizer  
 LFT01ST: One (1) Lift  
 PLT01ST: One (1) Support Plate  
 One (1) Silicone Pad, three (3) Fasteners (All are included sterile with the Steerable Guide Catheter)

**ADDITIONAL REQUIRED EQUIPMENT NOT INCLUDED**

Transseptal sheath and guidewire  
 Transseptal needle  
 Step-up dilators  
 260 cm of 0.9 mm (0.035") super stiff exchange length guidewire  
 High pressure three way stopcocks (5)  
 Arterial high pressure extension tubing (3)  
 50-60 cc syringes with luer fitting (2)  
 1000 ml pressure bags (2)  
 Sterile IV tubing with thumbwheel occluders (2)  
 Heparinized sterile saline solution (2, 1 liter bags)  
 Rolling IV Pole  
 Sterile Basin

**MITRACLIP PROCEDURE STEP-BY-STEP INSTRUCTIONS**

**1.0 DEFINITION OF TERMS**

*Defined Terms are in italics throughout document.*

TERM	DEFINITION AND RELATED TECHNIQUE
<i>Lock the Clip</i>	<ol style="list-style-type: none"> <li>1. Rotate the Lock Lever outward.</li> <li>2. Fully advance the Lock Lever.</li> <li>3. Rotate the Lock Lever inward to engage the lever.</li> </ol>
<i>Unlock the Clip and Open the Clip Arms</i>	<ol style="list-style-type: none"> <li>1. Rotate the Lock Lever outward and then retract the lever until the mark on the lever is fully exposed.</li> <li>2. Rotate the Lock Lever inward to engage the lever.</li> <li>3. Turn the Arm Positioner at least 1/2 turn in the "Close" (clockwise) direction.</li> <li>4. Turn the Arm Positioner in the "Open" (counter-clockwise) direction until the desired <i>Clip Arm Angle</i> is achieved.</li> <li>5. <i>Lock the Clip</i></li> </ol> <p><b>WARNING:</b> <b>Always Lock the Clip immediately after the desired Clip Arm Angle is achieved. Clip could be damaged and not unlock or open. Inability to open the Clip may result in valve injury, single leaflet device attachment (SLDA), worsening of mitral regurgitation, deployment of the Clip in an unintended location, or conversion to surgical intervention.</b></p> <p><b>NOTE 1:</b> If Clip does not open smoothly, retract the Lock Lever farther, then repeat steps 2 – 5.</p> <p><b>NOTE 2:</b> If the Clip Arms fail to open visibly (as observed under fluoroscopic guidance), use the following techniques in the order provided, as needed:</p> <ol style="list-style-type: none"> <li>A. Stop and return <i>Arm Positioner to Neutral</i>. Retract Lock Lever farther, then turn the Arm Positioner farther in the "Close" direction before turning in the "Open" direction.</li> <li>B. Turn the <i>Arm Positioner to Neutral</i>, then incrementally iterate the amount of Arm Positioner rotation in the "Close" direction followed by rotation in the "Open" direction. Iterate until Clip opens or until it is no longer possible to rotate the Arm Positioner in the "Close" direction.</li> <li>C. Turn the <i>Arm Positioner to Neutral</i>, iterate the amount of Lock Lever retraction past the mark in 5 mm increments, and rotate the Arm Positioner fully in the "Close" direction, before rotating in the "Open" direction, until Clip opens.</li> <li>D. Advance the Gripper Lever and repeat NOTE 2, Step C. Retract the Gripper Lever after Clip opens.</li> <li>E. If in the LA and free of tissue, release the DC Fastener, then release the Sleeve curves and repeat NOTE 2, Step C.</li> </ol> <p><b>WARNING:</b> <b>Failure to release the DC Fastener before releasing Sleeve curves may result in device damage and/or embolization.</b></p> <ol style="list-style-type: none"> <li>F. If the Clip does not open after performing all steps in NOTE 2, DO NOT use the device.</li> </ol>
<i>Arm Positioner to Neutral</i>	Turn the Arm Positioner in the "Close" or "Open" direction until no resistance to turning is noted.
<i>Invert the Clip Arms</i>	<ol style="list-style-type: none"> <li>1. <i>Unlock the Clip</i>.</li> <li>2. Turn the Arm Positioner at least 1/2 turn in the "Close" direction.</li> <li>3. Turn the Arm Positioner in the "Open" direction until the Clip Arms invert (see Figure 6F). DO NOT over-invert the Clip Arms; stop turning the Arm Positioner when resistance is first noted.</li> <li>4. <i>Lock the Clip</i>.</li> </ol>
<i>Raise the Grippers</i>	<ol style="list-style-type: none"> <li>1. Rotate the Gripper Lever outward.</li> <li>2. Slowly retract the Gripper Lever (under fluoroscopic observation) until the mark on the lever is just exposed.</li> </ol> <p><b>NOTE:</b> If pulling beyond the mark is required, advance the Gripper Lever back to the mark once the Grippers are fully raised.</p> <ol style="list-style-type: none"> <li>3. Rotate the Gripper Lever inward to engage the lever.</li> </ol>
<i>Lower the Grippers</i>	<ol style="list-style-type: none"> <li>1. Rotate the Gripper Lever outward.</li> <li>2. Fully advance the Gripper Lever.</li> <li>3. Rotate the Gripper Lever inward to engage the lever.</li> </ol>
<i>Clip Arm Angle</i>	<ul style="list-style-type: none"> <li>• Angle between the inner edges of both Clip Arms.</li> <li>• All <i>Clip Arm Angles</i> are measured using fluoroscopy with optimal view allowing clear observation of the tip of the Clip and both arms in the same plane so they appear as a "V" (see Figure 6).</li> </ul>

TERM	DEFINITION AND RELATED TECHNIQUE
<i>Grasping Arm Angle</i>	A <i>Clip Arm Angle</i> of approximately 120 degrees. <b>NOTE:</b> Establish <i>Grasping Arm Angle</i> after closing the Clip from a larger <i>Clip Arm Angle</i> .
<i>Fully Close the Clip Arms</i>	Turn the Arm Positioner in the "Close" direction until the Clip Arms contact the DC. <ul style="list-style-type: none"> <li>Under direct visualization, the Clip is fully closed when the Clip Covering contacts the DC.</li> <li>Under fluoroscopic observation, the Clip is fully closed when the inner edges of the Clip Arms are parallel.</li> </ul> <p><b>WARNING:</b> Never close the Clip while the Lock Lever is in an unlocked state. Clip could be damaged and not unlock or open. Inability to open the Clip may result in valve injury, single leaflet device attachment (SLDA), worsening of mitral regurgitation, deployment of the Clip in an unintended location, or conversion to surgical intervention.</p>
<i>Establish Final Arm Angle</i>	Pre-deployment <i>Clip Arm Angle</i> that reflects the <i>Clip Arm Angle</i> post-deployment. <ol style="list-style-type: none"> <li>With the Lock Lever fully advanced, turn the Arm Positioner in the "Open" direction until resistance is first noted. The Clip Arms may open slightly and then remain in a stable position.</li> <li>Confirm that the Clip is locked by observing slight Delivery Catheter shaft deflection using fluoroscopy.</li> </ol> <p><b>NOTE:</b> If continued opening of the Clip Arms is noted, reconfirm that the Lock Lever is completely advanced. Close the Clip Arms, and <i>Establish Final Arm Angle</i>.</p>

**2.0 PATIENT PREPARATION**

- Prepare the patient per institution's standard practice for transseptal catheterization.
- Place support plate under patient's leg in the region between the area of the upper leg and the knee and place the Lift over the ipsilateral lower extremity prior to draping the patient.
- Place the Lift on the Support Plate such that the front edge (i.e., the edge that corresponds with the shorter legs of the Lift) is approximately 80 cm from the patient's mid sternum.
- Adjust the height of the Lift so that the front edge of the Lift is close to the patient's leg, but is not impinging on it. Adjust the back legs to be 2 or 3 notches above the front legs (i.e., the back legs of the Lift are taller than the front legs).
- Ensure the Lift and Support Plate are covered completely by sterile drape during the procedure. Use towels as necessary to minimize direct contact between the patient and all surfaces of both the Lift and Support Plate.
- Prepare the patient for invasive hemodynamic monitoring.

**3.0 MITRACLIP SYSTEM PREPARATION BEFORE USE**

**WARNING:** DO NOT use the MitraClip System after the "Use By" date stated on the package label, and never reuse or re-sterilize the system. Use of expired, reused, or re-sterilized devices may result in infection, endocarditis, and/or sepsis.

**WARNING:** Always inspect the MitraClip System and its packaging to verify no damage has occurred as a result of shipping and handling and that the sterile barrier has not been compromised. DO NOT use the device if damage is detected. Use of product with a compromised sterile barrier may result in infection, endocarditis, and/or sepsis. Use of damaged product may result in patient injury.

- DO NOT remove the protective cover placed over the Clip.

**WARNING:** DO NOT handle the Clip directly; leave it in the protective cover to avoid potential contamination. Removal of the protective cover may result in infection, endocarditis, and/or sepsis. Removal of the protective cover may result in damaged product which may result in patient injury.

- The preparation is most easily accomplished with the aid of an assistant.

**3.1 STEERABLE GUIDE CATHETER PREPARATION**

**WARNING:** All lumens contain air when shipped. Use proper de-airing techniques before and during use to minimize the risk of air embolism.

- Carefully remove the white Guide tip shape retainer and transparent protective tubing from the Guide tip.
- Inspect Steerable Guide Catheter and Dilator to verify they are undamaged.

**WARNING:** DO NOT use if damage is detected. Use of damaged product may result in air embolism, vascular and/or cardiac injury.

- Remove the sterile package containing Fasteners and Silicone Pad from the Steerable Guide Catheter tray.
- Fill a basin with 1000 cc of heparinized saline.
- Flush and de-air the Guide and Dilator with heparinized saline:
  - Connect 3-way stopcocks to the Guide and Dilator flush ports.
  - De-air the Dilator, then close the stopcock and the Rotating Hemostatic Valve.
  - Hydrate 5-10 cm of the distal end of the Dilator with heparinized saline.
  - Insert the Dilator approximately 10 cm into Guide then remove.
  - Connect high pressure tubing and a 50-60 cc syringe filled with heparinized saline to the Guide flush port.

- 3.1.5.6 De-air the Guide.
  - 3.1.5.6.1 With the tip raised, displace all air from the Guide while tapping along the length of the catheter shaft.
  - 3.1.5.6.2 Cover the Guide tip with finger once heparinized saline exits the Guide.
  - 3.1.5.6.3 Close the Guide stopcock.
- 3.1.6 Submerge the Guide tip in the basin of heparinized saline.
- 3.1.7 While the Guide tip is submerged in the basin of heparinized saline, remove finger from Guide tip and check the Guide valve for leaks by raising the handle to a vertical position for a minimum of 30 seconds.
- 3.1.8 Hydrate 5-10 cm of the distal end of the Dilator with heparinized saline.
- 3.1.9 Cover the Guide tip with finger and insert the Dilator into the Guide while Guide tip remains submerged in the basin of heparinized saline.
  - 3.1.9.1 While advancing the Dilator, continually watch for air in the Guide Hemostasis Valve housing. If needed, remove finger from Guide tip and aspirate while assuring the Guide tip is submerged.
  - 3.1.9.2 Remove finger from tip of Guide when the Dilator tip approaches the Guide tip.
  - 3.1.9.3 Advance the Dilator until the curve is extended from the Guide tip.

### 3.2 STEERABLE GUIDE CATHETER FUNCTIONAL INSPECTION

- The functional inspection is most easily accomplished with the aid of an assistant.
- The Guide functional inspection should be performed with the Guide tip and Dilator tip submerged in a basin of heparinized saline to prevent air from entering the lumens. If the Guide tip and/or Dilator tip fails to remain submerged during inspection, flush the Guide and/or Dilator with heparinized saline to completely remove air.

**WARNING: Failure to completely remove air may result in air embolism.**

**WARNING: All catheter manipulations should be done with care. DO NOT continue to rotate or manipulate any of the handle controls if significant resistance is noted. Use of damaged product may result in air embolism, vascular and/or cardiac injury.**

#### Guide Inspection

- 3.2.1 Inspect all Guide parts to verify they are undamaged.

**WARNING: DO NOT use device if damage is detected. Use of damaged product may result in air embolism, vascular and/or cardiac injury.**

- 3.2.2 To confirm proper tip deflection with "+" knob rotation:
  - 3.2.2.1 Rotate the +/- Knob in the "+" direction until the Guide is curved to approximately 80 degrees.
  - 3.2.2.2 Remove hand from the +/- Knob and check that the knob does not slip.
  - 3.2.2.3 Return the +/- Knob to the neutral position.
  - 3.2.2.4 Repeat steps 3.2.2.1 through 3.2.2.3.
- 3.2.3 To confirm proper tip deflection with "-" knob rotation:
  - 3.2.3.1 Rotate the +/- Knob in the "-" direction until the Guide curve is substantially straightened.
  - 3.2.3.2 Remove hand from the +/- Knob and check that the knob does not slip.
  - 3.2.3.3 Return the +/- Knob to the neutral position.
  - 3.2.3.4 Repeat steps 3.2.3.1 through 3.2.3.3.
- 3.2.4 Retract the Dilator until the tip is 3-5 cm beyond the Guide tip. Position the Dilator to create a smooth transition.

### 3.3 STABILIZER PREPARATION

- 3.3.1 Assemble the sterilized Stabilizer by placing the two Fasteners in the Stabilizer. Ensure that the Fasteners can be fully threaded into the Stabilizer holes. Set the Stabilizer aside in a protected sterile environment for later use.

### 3.4 CLIP DELIVERY SYSTEM PREPARATION

- 3.4.1 Inspect the Clip, DC shaft, and Sleeve tip to verify they are undamaged.

**WARNING: DO NOT use the device if damage is detected. Use of damaged product may result in air embolism, device or device component embolization, vascular and/or cardiac injury.**

#### Sleeve Preparation

- 3.4.2 Connect 3-way stopcocks to the Sleeve flush port and bottom DC flush port.
- 3.4.3 Remove the cap from the Clip Introducer.
- 3.4.4 Place the cap on the top flush port of the DC Handle.
- 3.4.5 Connect a 3-way stopcock to the Clip Introducer flush port.
- 3.4.6 Connect one high pressure tube to each drip line from the pressurized bags with sterile heparinized saline; flush and de-air the lines.
- 3.4.7 Connect one high pressure tube to the 3-way stopcock on the bottom flush port of the DC Handle and one high pressure tube to the 3-way stopcock on the flush port of the Sleeve Handle.
- 3.4.8 Flush and de-air the Sleeve with heparinized saline.
  - 3.4.8.1 With the tip raised and the shaft held taut, displace all air from the Sleeve lumen while tapping along the length of the catheter shaft.
  - 3.4.8.2 While flushing, release the DC Fastener, retract and advance the DC Handle to remove residual air from the lumen.

**WARNING: Using excessive force when pulling the DC Radiopaque Ring against the Sleeve tip, while translating the DC shaft, may result in device damage including distal tip embolization.**

- 3.4.8.3 Secure the DC Fastener with DC Handle fully advanced.

#### Delivery Catheter Preparation

**WARNING: DO NOT handle the Clip directly; leave in the protective cover to avoid potential contamination. Removal of the protective cover may result in damaged product which may result in patient injury.**

- 3.4.9 Attach a 50-60 cc syringe filled with heparinized saline to the 3-way stopcock on the Clip Introducer.
- 3.4.10 De-air the Clip Introducer, then close the stopcock.
- 3.4.11 Temporarily remove the cap from top flush port of the DC Handle.
- 3.4.12 Flush and de-air DC Handle and all lumens of the DC with heparinized saline.
- 3.4.13 After de-airing the DC Handle chamber, replace the cap to close off top flush port of the DC Handle.
- 3.4.14 Retract and advance the Lock Lever several times to remove residual air from the lumens.
- 3.4.15 Loosen the Lock Lever and the Gripper Lever Caps to de-air. DO NOT turn lever caps more than 1/2 turn in the "Open" direction. After de-airing, tighten the lever caps.
- 3.4.16 With the tip raised and the shaft held taut, displace all air from the DC while tapping along the length of the catheter shaft.
- 3.4.17 Confirm continuous flow from the distal end of the DC.

### 3.5 CLIP DELIVERY SYSTEM FUNCTIONAL INSPECTION

- The functional inspection is most easily accomplished with the aid of an assistant.
- 3.5.1 Inspect all Clip Delivery System parts, including the Clip, to verify they are undamaged.

**WARNING: DO NOT use device if damage is detected. Use of damaged product may result in air embolism, device or device component embolization, vascular and/or cardiac injury.**

**WARNING: All catheter manipulations should be done with care. DO NOT continue to, rotate or manipulate any of the handle controls if significant resistance is noted. Use of damaged product may result in air embolism, vascular and/or cardiac injury.**

#### Sleeve Inspection

**WARNING: DO NOT deflect the Sleeve more than 90 degrees during the inspections below. Use of damaged product may result in cardiac injury.**

- 3.5.2 To confirm proper tip deflection with "A" knob rotation:
  - 3.5.2.1 With the DC handle fully advanced and the shaft held taut, rotate the A/P Knob approximately 3/4 turn in the "A" direction from neutral to confirm that the distal tip deflects
  - 3.5.2.2 Remove hand from the A/P Knob and check that the knob does not slip.
  - 3.5.2.3 Return the A/P Knob to the neutral position.
  - 3.5.2.4 Repeat steps 3.5.2.1 through 3.5.2.3.
- 3.5.3 To confirm proper tip deflection with "P" knob rotation:
  - 3.5.3.1 With the DC handle fully advanced and the shaft held taut, rotate the A/P Knob approximately 3/4 turn in the "P" direction from neutral to confirm that the distal tip deflects.
  - 3.5.3.2 Remove hand from the A/P Knob and check that the knob does not slip.
  - 3.5.3.3 Return the A/P Knob to the neutral position.
  - 3.5.3.4 Repeat steps 3.5.3.1 through 3.5.3.3.
- 3.5.4 To confirm proper tip deflection with "M" knob rotation:
  - 3.5.4.1 With the DC Handle fully advanced and the shaft held taut, rotate the M/L Knob in the "M" direction until the distal tip deflects to approximately 90 degrees to confirm distal tip deflection.
  - 3.5.4.2 Remove hand from the M/L Knob and check that the knob does not slip.
  - 3.5.4.3 Return the M/L Knob to the neutral position.
  - 3.5.4.4 Repeat steps 3.5.4.1 through 3.5.4.3.

#### Delivery Catheter and Clip Inspection

**WARNING: DO NOT handle the Clip directly, leave in the protective cover to avoid potential contamination. Removal of the protective cover may result in infection, endocarditis, and/or sepsis. Removal of the protective cover may result in damaged product which may result in patient injury.**

**NOTE:** If *Clip Arm Angle* is greater than *Grasping Arm Angle*, close the Clip to *Grasping Arm Angle*; if *Clip Arm Angle* is less than *Grasping Arm Angle*, *Unlock the Clip and Open the Clip Arms* to 180 degrees then close the Clip to *Grasping Arm Angle*.

- 3.5.5 Carefully inspect the Grippers to confirm the cover is intact and not damaged.

**WARNING: DO NOT use the device if damage is detected. Use of damaged product may result in cardiac injury and/or may lead to inability to reduce MR.**

- 3.5.6 *Raise the Grippers.*

**CAUTION: Raising the Grippers more often than needed, retracting the Gripper Lever forcefully, or retracting the Gripper Lever more than 1.5 cm beyond the mark may damage the Gripper cover and impair CDS performance.**

- 3.5.7 *Unlock the Clip.*

**WARNING: Retracting the Lock Lever forcefully may result in the inability to lock or *unlock the Clip*. Damage could occur causing the Clip to not unlock or open. Inability to open the Clip may result in valve injury or lead to deployment of the Clip in an unintended location.**

- 3.5.8 *Invert the Clip Arms.*

**WARNING: DO NOT continue turning the Arm Positioner if resistance is felt. Use of damaged product may result in cardiac injury.**

- 3.5.9 *Lock the Clip.*
- 3.5.10 Close the Clip to *Grasping Arm Angle*.
- 3.5.11 *Lower the Grippers* once to de-air the lumens.
- 3.5.12 Release the DC Fastener and torque the DC Handle clockwise and counterclockwise 1/4 turn while translating the shaft.

**WARNING:** Using excessive force when pulling the DC Radiopaque Ring against the Sleeve tip, while translating the DC shaft, may result in device damage including distal tip embolization.

- 3.5.13 Secure the DC Fastener with DC Handle fully advanced.
- 3.5.14 Close the Clip to a *Clip Arm Angle* of approximately 20 degrees.
- 3.5.15 *Establish Final Arm Angle.*
- 3.5.16 Return the *Arm Positioner to Neutral.*
- 3.5.17 *Unlock the Clip.*
- 3.5.18 Open the Clip to *Grasping Arm Angle.*
- 3.5.19 *Lock the Clip.*
- 3.5.20 Return the *Arm Positioner to Neutral.*
- 3.5.21 Release the DC Fastener and retract the DC fully against the Sleeve.
- 3.5.22 Secure the DC Fastener.
- 3.5.23 Temporarily discontinue heparinized saline flushes.

**The following steps should be performed just before use of the CDS:**

- 3.5.24 Re-start heparinized saline flushes.
- 3.5.25 *Raise the Grippers.*
- 3.5.26 *Fully Close the Clip Arms.*
- 3.5.27 *Lower the Grippers.*
- 3.5.28 Without removing the protective cover, carefully slide the Clip Introducer over the Clip.

**WARNING:** DO NOT compress the Clip Arms. Compressing the Clip Arms may result in inability to open the Clip. Inability to open the Clip may result in valve injury or lead to deployment of the Clip in an unintended location.

- 3.5.29 Stop when the tip of the Clip is just proximal to the tip of the Clip Introducer.
- 3.5.30 Turn the *Arm Positioner to Neutral.*

**WARNING:** Failure to Fully Close the Clip Arms, before insertion or retraction into the Clip Introducer may result in difficulty or inability to advance or retract the Clip, which may result in vascular and/or cardiac injury, air embolism, and/or the need for surgical intervention.

**WARNING:** Heparinized saline flush should be continuous throughout the procedure. Ensure flow is visible through the drip chamber, that the tubing is free from kinks and/or obstruction and appropriate pressure of 300 mm Hg is maintained. Discontinuing flush may result in air embolism and/or thrombus formation.

#### 4.0 ACCESS TO THE MITRAL VALVE

**NOTE:** This is a suggested sequence for the procedure. Variations may be used based upon patient anatomy.

- 4.1 Access the LA to accommodate the Guide tip using transvenous, transseptal techniques and equipment.
- 4.2 Heparinize the patient.

**WARNING:** Failure to administer heparin once transseptal access has been achieved may result in thrombus formation.

- 4.3 Carefully place a 260 cm super stiff 0.9 mm (0.035") exchange length guidewire in the left upper pulmonary vein or LA. Dilate the subcutaneous tissue and femoral vein to accommodate the Guide shaft using standard dilation technique.

#### 5.0 STEERABLE GUIDE CATHETER INSERTION

**WARNING:** Confirm a smooth transition between the Dilator and the tip of the Guide to minimize the risk of vascular and/or cardiovascular injury.

**CAUTION:** Always use pressure monitoring, echocardiography and fluoroscopy for guidance and observation during use of the MitraClip System.

**WARNING:** Always use a careful, deliberate, and iterative approach to positioning the MitraClip System. It is recommended to make multiple small adjustments rather than single large adjustments. Large adjustments may result in vascular and/or cardiac injury.

- 5.1 Rotate the +/- Knob in the "-" direction until the Guide curve is substantially straightened.
- 5.2 Wet the surface of the Guide shaft with sterile saline.
- 5.3 Insert the Guide-Dilator assembly over the stationary guidewire into the femoral vein.

**WARNING:** DO NOT use excessive force to advance or manipulate the Guide-Dilator assembly. If resistance is encountered, use echocardiography and/or fluoroscopy to assess before proceeding. Use of excessive force may result in arrhythmias, vascular and/or cardiac injury, including creation of a clinically significant atrial septal defect.

- 5.4 Advance the Guide-Dilator assembly to the RA. Rotate the +/- Knob to Neutral, then place tip of the Dilator partially across the atrial septum.
- 5.5 Slowly dilate the atrial septum by gradually advancing the tip of the Guide-Dilator assembly.

**WARNING:** DO NOT rapidly advance the Guide-Dilator assembly across the atrial septum. Rapid advancement may result in vascular and/or cardiac injury.

- 5.6 Advance the Guide-Dilator assembly until the tip of the Guide extends approximately 3 cm in the LA.
- 5.7 Adjust Guide deflection and torque to position the tip away from adjacent tissues.
- 5.8 Place the Silicone Pad on the sterile drape over the Lift. Place the Stabilizer onto the Silicone Pad.
- 5.9 Secure the Guide in the Stabilizer slot using the Fastener. Ensure the Fastener engages the metallic tube on the Guide shaft. The Guide handle should be immediately adjacent to the Stabilizer, such that they are in contact with each other.

- 5.10 Retract the Dilator approximately 5 cm into the Guide, leaving the guide wire in the left upper pulmonary vein or LA.

**CAUTION:** Always loosen the Fastener before torquing the Guide to prevent stripping the screw.

- 5.11 Retract the guidewire into the tip of the Dilator. Remove the Dilator and guidewire while gently aspirating the Guide (starting when the Dilator is approximately halfway retracted into the Guide, approximately 40 cm) using a 50-60 cc syringe. Cover Guide Hemostasis Valve with finger upon Dilator removal.

**NOTE:** Avoid contacting tissue or creating a vacuum in the Guide lumen. If necessary, position the Guide handle below the level of the LA to allow blood to fill the Guide lumen.

**WARNING:** DO NOT create a vacuum while removing the dilator from the Guide; air may enter the lumen of the Guide which may result in air embolism.

**WARNING:** Failure to fully retract guidewire into the Dilator may result in air embolization.

## 6.0 CLIP DELIVERY SYSTEM INSERTION

- 6.1 Confirm the Guide lumen is completely de-aired.

**WARNING:** To minimize the potential of air embolism, DO NOT introduce the CDS into the Guide until the Guide lumen has been completely de-aired.

- 6.2 Confirm there is a slow, continuous heparinized saline flush through both the Sleeve and the DC.

**CAUTION:** Failure to continuously flush the CDS with heparinized saline may reduce device performance.

**WARNING:** Heparinized saline flush should be continuous throughout the procedure. Ensure flow is visible through the drip chamber and that tubing is free from kinks and/or obstruction and pressure of 300 mm Hg is maintained. Discontinuing flush may result in air embolism and/or thrombus formation.

- 6.3 Confirm tip of the Clip is just proximal to the tip of the Clip Introducer.

- 6.4 Carefully remove the protective cover surrounding the Clip and the Clip Introducer.

- 6.5 Confirm that the stopcock on the Clip Introducer flush port is closed and that the Clip Introducer is de-aired.

- 6.6 While flushing heparinized saline on the Guide Hemostasis Valve, place the tip of the Clip Introducer against the Guide Hemostasis Valve and advance the Clip Introducer straight into the valve in a continuous motion while rotating the Clip Introducer in small clockwise and counterclockwise motions until the Clip can be observed distal to the valve.

**WARNING:** DO NOT continue to advance the Clip Introducer if resistance is felt; the Guide Hemostasis Valve, Clip Introducer or the Clip may be damaged. Damage to these components may result in air embolism, vascular or cardiac injury.

**WARNING:** To minimize the potential of air embolization, ensure proper de-airing when inserting the Clip Introducer into the Guide Hemostasis Valve.

- 6.7 Leave the Clip Introducer fully inserted in the Guide Hemostasis Valve throughout the procedure.

- 6.8 Align the Longitudinal Alignment Marker on the Sleeve shaft with the Alignment Marker on the Guide Hemostasis Valve.

- 6.9 Turn the +/- Knob to neutral then carefully advance the CDS through the Guide under fluoroscopic guidance. Stop when the tip of the Clip is even with the tip of the Guide.

**NOTE:** If resistance to CDS advancement is felt, reduce Guide deflection.

- 6.10 Under echocardiographic guidance, advance the CDS and retract the Guide iteratively as needed while maintaining the Guide in the LA. Stop when the Guide RO Tip Ring is between the RO Alignment Markers of the Sleeve, as confirmed under fluoroscopic guidance.

- 6.11 Position the Sleeve Handle in the Stabilizer slot.

- 6.12 Confirm that the Clip is free from the left atrial wall and valve tissue.

**WARNING:** Failure to confirm that the Clip is free from the left atrial wall and valve tissue may result in cardiac injury.

## 7.0 INITIAL MITRACLIP SYSTEM POSITIONING IN THE LEFT ATRIUM

**NOTE:** Positioning is achieved with iterative adjustments of the Guide and CDS using torque, translation and knob adjustments. The goals of positioning are:

- Positioning the Clip centrally over the valve with respect to anterior-posterior and medial-lateral directions.
- Aligning the Clip so the DC Shaft is perpendicular to the plane of the mitral valve.
- Positioning the distal tip of the Clip at least 1 cm above the leaflets.

**WARNING:** Excessive torque on the Guide and translation of the MitraClip System may inadvertently displace the tip of the Guide from the LA, which may result in arrhythmias or cardiac injury.

**WARNING:** DO NOT continue to rotate or manipulate any of the handle knobs if significant resistance is noted; device damage may occur and result in cardiac injury.

- 7.1 Adjust the Guide position as necessary to maintain that the Clip is free from adjacent tissue.

- 7.2 Adjust Sleeve deflection using the M/L Knob and/or the A/P Knob to deflect the Clip towards the apex. Retract the DC Radiopaque Ring against the Sleeve tip as necessary.

- 7.3 During Sleeve deflections confirm that the Guide RO Tip Ring is between the RO Alignment Markers of the Sleeve prior to making maximum Sleeve deflections.

**WARNING:** DO NOT deflect the Sleeve tip more than 90 degrees as device damage may occur. Use of damaged product may result in cardiac injury.

- 7.4 Secure the Sleeve handle in the Stabilizer using the Fastener.
- 7.5 To reposition the MitraClip System, move the Stabilizer and the system together until positioning is adequate.
- 7.6 Adjust the MitraClip System position to maintain adequate height above the mitral valve in the LA.

**WARNING:** Maintain the Clip above the leaflets until ready to grasp to minimize the risk of Clip entanglement in the chordal apparatus. Clip entanglement may result in cardiac injury, worsening mitral regurgitation, difficulty or inability to remove the Clip and conversion to surgical intervention.

**8.0 FINAL MITRACLIP SYSTEM POSITIONING**

- 8.1 *Raise the Grippers*

**CAUTION:** Raising the Grippers more often than needed, retracting the Gripper Lever forcefully, or retracting the Gripper Lever more than 1.5 cm beyond the mark may damage the Gripper cover and impair CDS performance.

- 8.2 *Unlock the Clip and Open the Clip Arms* to approximately 180 degrees.

**WARNING:** Retracting the Lock Lever forcefully may result in the inability to unlock Clip. Inability to open the Clip may result in valve injury or lead to deployment of the Clip in an unintended location.

- 8.3 *Lock the Clip.*

**CAUTION:** Failure to immediately advance the Lock Lever after Clip Arm opening may affect DC shaft straightness.

- 8.4 Adjust the MitraClip System to reposition the Clip as necessary. Confirm that the distal tip of the Clip is at least 1 cm above the leaflets.

- 8.5 Rotate the DC handle to align the Clip Arms perpendicular to the line of coaptation. DO NOT rotate the Clip more than 90 degrees in each direction.

- 8.6 Carefully translate the DC shaft multiple times to release stored torque. Fully retract the DC.

**WARNING:** Failure to fully release stored torque may result in unwanted Clip Arm orientation changes during grasping. Torque of the DC Handle more than 180 degrees may result in DC damage and cardiac injury.

- 8.7 Complete final MitraClip System positioning in the LA using multiple imaging planes. Re-secure the Guide and Sleeve Fasteners.

**9.0 GRASPING THE LEAFLETS AND VERIFYING THE GRASP**

- 9.1 Advance the DC distally to position the Clip approximately 2 cm below the valve. Ensure that the Clip Arms are oriented perpendicular to the line of coaptation.

**WARNING:** Failure to confirm that the Clip Arms are perpendicular to the line of coaptation may result in loss of leaflet capture and insertion. Loss of leaflet capture and insertion may result in inadequate reduction of mitral regurgitation and may result in a single leaflet device attachment (SLDA).

**WARNING:** DO NOT make substantial Clip Arm orientation adjustment in the LV. Clip entanglement in sub-valvular apparatus may result in cardiac injury and worsening mitral regurgitation; and may result in difficulty or inability to remove the Clip, and conversion to surgical intervention.

**WARNING:** Always ensure that either the Grippers are raised or that the Clip is closed while in the LV to avoid potential cardiac injury.

- 9.2 Close the Clip to the *Grasping Arm Angle*.

- 9.3 Without using excessive force, retract the DC to grasp both anterior and posterior leaflets.

**WARNING:** An improper grasp will allow one or both leaflets to move freely. Closing and deploying the Clip in this situation may result in loss of leaflet capture and insertion. Loss of leaflet capture and insertion may result in inadequate reduction of mitral regurgitation and may result in a single leaflet device attachment (SLDA).

- 9.4 If the grasp appears satisfactory, *Lower the Grippers* onto the leaflets.

**WARNING:** Failure to confirm that both Grippers have been lowered onto the leaflets prior to closing the Clip may result in loss of leaflet capture and insertion. Loss of leaflet capture and insertion may result in inadequate reduction of mitral regurgitation and may result in a single leaflet device attachment (SLDA).

**WARNING:** DO NOT adjust the position of the MitraClip System after grasping the leaflets, valve injury may occur.

- 9.5 Close the Clip until the *Clip Arm Angle* is approximately 60 degrees. Release tension on the DC and secure the DC Fastener.

- 9.6 Use echocardiographic imaging to verify insertion of both leaflets and satisfactory grasp by observation of:

- Leaflet immobilization
- Single or multiple valve orifice(s)
- Limited leaflet mobility relative to the tips of both Clip Arms
- Adequate MR reduction.

- 9.6.1 If grasping fails to hold both leaflets and the Clip retracts to the LA, reposition the MitraClip System.

9.6.1.1 *Unlock the Clip and Open the Clip Arms* to approximately 180 degrees and reorient the Clip Arms in the LA, as needed, then repeat grasping steps.

9.6.1.1.1 If significant repositioning is necessary, *Fully Close the Clip Arms and Lower the Grippers* then repeat positioning and grasping steps.

RELEASED

9.6.2 If the Sleeve limits DC travel during grasping, an inadequate grasp may require repositioning of the MitraClip System.

9.6.2.1 *Raise the Grippers, Unlock the Clip and Open the Clip Arms* to 180 degrees, and advance the DC handle. Repeat positioning and grasping steps as necessary.

## 10.0 CLOSING THE CLIP AND EVALUATING CLIP POSITION

10.1 Slowly close the Clip just until the leaflets are coapted and MR is sufficiently reduced. The Clip should maintain a distinct "V" shape.

**WARNING:** **DO NOT use excessive force to close the Clip further than is necessary to adequately reduce MR. Leaflet injury may occur.**

**WARNING:** **Closing the Clip too tightly may result in inability to deploy the Clip. Inability to deploy the Clip may result in worsening mitral regurgitation, cardiac injury, a single leaflet device attachment (SLDA), and/or conversion to surgical intervention.**

10.2 Use echocardiographic imaging to verify valve function, satisfactory coaptation, and insertion of both leaflets by observation of:

- Leaflet immobilization
- Single or multiple valve orifice(s)
- Limited leaflet mobility relative to the tips of both Clip Arms
- Adequate MR reduction.

If the Clip position is not satisfactory, *Raise the Grippers and Invert the Clip Arms*.

10.2.1 Retract the inverted Clip into the LA.

10.2.2 Confirm both leaflets move freely.

10.2.3 Repeat positioning steps, as necessary, then repeat grasping steps.

## 11.0 MITRACLIP DEVICE PRE-DEPLOYMENT CLIP ASSESSMENT

11.1 Confirm DC Handle is secure.

**WARNING:** **Failure to secure the DC Handle may result in leaflet injury or loss of leaflet insertion with resultant worsening mitral regurgitation, single leaflet device attachment (SLDA), and/or conversion to surgical intervention.**

11.2 *Establish Final Arm Angle.*

**WARNING:** **DO NOT turn the Arm Positioner more than 1/2 turn in the "Open" direction once initial resistance is felt to prevent device deployment. Inability to deploy the Clip may result in worsening mitral regurgitation, cardiac injury, a single leaflet device attachment (SLDA), and/or conversion to surgical intervention.**

11.3 Turn the Arm Positioner to the "closed" side of the neutral position.

11.3.1 Use echocardiographic imaging to verify valve function, satisfactory coaptation, and insertion of both leaflets by observation of:

- Leaflet immobilization
- Single or multiple valve orifice(s)
- Limited leaflet mobility relative to the tips of both Clip Arms
- Adequate MR reduction.

11.4 Perform mean pressure gradient assessment prior to proceeding to deployment.

11.5 Establish Gripper Line Removability

**WARNING:** **Failure to Establish Gripper Line Removability prior to deployment of the Clip may result in inability to remove the Gripper Line. Intervention may be required.**

11.5.1 Confirm the Gripper Lever is fully advanced.

11.5.2 Increase the flush rate to the DC and Sleeve. Remove the Gripper Lever Cap and "O" ring. Unwrap the two ends of the Gripper Line. Remove the plastic cover from the lines and separate the two ends so that no twists or knots are present.

11.5.3 With one free end of the Gripper Line in each hand, confirm that the Gripper Line is removable by pulling slowly on one end until the other end of the Gripper Line moves approximately 3-5 cm. If the Gripper Line is confirmed to be removable, continue to Clip Deployment.

**NOTE:** If excessive resistance is noted, stop, and pull on the other free end.

**WARNING:** **While Establishing Gripper Line Removability, ensure that both ends of the Gripper Line remain exposed. Failure to maintain exposure of both Gripper Line ends may result in an inability to remove the Gripper Line in its entirety and could lead to conversion to surgical intervention.**

**WARNING:** **Pulling the Gripper Line too quickly or with excessive force may raise the Grippers, break the Gripper Line and/or disturb leaflet capture and insertion. This may result in worsening mitral regurgitation, and could lead to a single leaflet device attachment (SLDA), and/or conversion to surgical intervention.**

11.5.3.1 If excessive resistance is noted at both ends of the Gripper Line (resulting in failure to Establish Gripper Line Removability), stop and remove the Clip Delivery System.

**NOTE:** The removal of the Clip Delivery System is most easily accomplished with the aid of an assistant.

11.5.3.1.1 Hold both free ends of the Gripper Lines together and apply tension to maintain the Grippers in a raised position through Step 11.5.3.1.4.

11.5.3.1.2 *Invert the Clip Arms and then Lock the Clip.*

11.5.3.1.3 Release the DC Fastener and retract the inverted Clip into the LA. Retract DC shaft until the DC Radiopaque Ring is fully against the tip of the Sleeve.

11.5.3.1.4 *Fully Close the Clip Arms.*

11.5.3.1.5 Continue to Section 14.2.1.5: MITRACLIP SYSTEM REMOVAL WITH CLIP ATTACHED to remove the Clip.

## 12.0 CLIP DEPLOYMENT

### 12.1 Deployment Step 1: Lock Line Removal

- 12.1.1 Remove the Lock Lever Cap and "O" ring. Unwrap the two ends of the Lock Line. Remove the plastic cover from the lines so that no twists or knots are present.
- 12.1.2 Grasp one of the free ends of the Lock Line, confirm the line moves freely, and slowly remove the Lock Lever. Pull the Lock Line coaxial to the Lock Lever. If resistance is noted, stop and pull on the other free end to remove the Lock Line.
- 12.1.3 *Establish Final Arm Angle.*

**NOTE:** The Clip Arms may open slightly before remaining in a stable position. If Arms open more than slightly, close the Clip to the desired Arm position and re-*Establish Final Arm Angle.*

- 12.1.4 Turn the *Arm Positioner to Neutral.*

### 12.2 Deployment Step 2: Delivery Catheter Shaft Detachment

- 12.2.1 Confirm that the Arm Positioner is Neutral and that the two ends of the Gripper Line have been unwrapped from under the cap and are not twisted or knotted. Remove the release pin from the DC Handle.
- 12.2.2 Turn the Actuator Knob of the DC approximately 8 turns counterclockwise. If it is difficult to turn the Actuator Knob, confirm that the Arm Positioner moves freely.

**WARNING:** **Failure to stop turning the Actuator Knob when resistance is felt or turning the Actuator Knob in the clockwise direction may result in inability to deploy the Clip. Inability to deploy the Clip may result in worsening mitral regurgitation, cardiac injury, a single leaflet device attachment (SLDA), and/or conversion to surgical intervention.**

- 12.2.3 Release the DC Fastener then retract the Actuator Knob after it is fully unthreaded.
- 12.2.4 Confirm that the Release Pin groove is fully exposed.
- 12.2.5 Retract the DC Handle such that the Clip has separated at least 1 cm from the DC tip.
- 12.2.6 Secure the DC Fastener.
- 12.2.7 Allow several minutes after catheter shaft detachment before proceeding to the final Clip deployment step. Use echocardiographic imaging to verify valve function, satisfactory coaptation, and insertion of both leaflets by observation of:
  - Leaflet immobilization
  - Single or multiple valve orifice(s)
  - Limited leaflet mobility relative to the tips of both Clip Arms
  - Adequate MR reduction.

**WARNING:** **If Clip placement and/or MR reduction is not satisfactory after Deployment Step 2: Delivery Catheter Shaft Detachment, DO NOT proceed to Deployment Step 3: Gripper Line Removal. Intervention may be required to remove the Clip.**

### 12.3 Deployment Step 3: Gripper Line Removal

- 12.3.1 Grasp one of the free ends of the Gripper Line, confirm the line moves freely and slowly remove the line. Pull the Gripper Line coaxial to the Gripper Lever. If resistance is noted, stop and pull on the other free end to remove the Gripper Line. Maintain at least 1 cm separation between the DC tip and the Clip while slowly removing the Gripper Line.

**WARNING:** **If less than 1 cm separation is present between the DC tip and the Clip before or during Gripper Line retraction, it may be difficult to remove the Gripper Line in its entirety.**

**WARNING:** **Pulling the Gripper Line too quickly or with excessive force may raise the Grippers, resulting in device damage and/or compromise leaflet capture and insertion. This may result in worsening mitral regurgitation, and could lead to a single leaflet device attachment (SLDA).**

- 12.3.1.1 If the Gripper Line does not move easily, release the DC Fastener and incrementally release Sleeve curves (M/L Knob and A/P Knob). Secure DC Fastener once Sleeve curves are released.

**WARNING:** **Failure to release the DC Fastener before releasing Sleeve curves may result in device damage and/or embolization.**

- 12.3.1.2 If the Gripper Line still does not move easily, partially release Guide curves.
- 12.3.1.3 If the Gripper Line still does not move easily, the CDS may also be partially retracted into the tip of the Guide, or completely removed by pulling only on the Sleeve Handle, to facilitate Gripper Line removal.

**WARNING:** **Retracting the CDS by pulling on the DC Handle may result in device damage and/or embolization.**

- 12.3.2 Release the DC Fastener and retract the DC Handle until the DC Radiopaque Ring is against the tip of the Sleeve.
- 12.3.3 Secure the DC Fastener.

#### 12.3.4 Confirm that the Clip position is stable.

- 12.3.5 Use echocardiographic imaging to verify valve function, satisfactory coaptation, and insertion of both leaflets by observation of:
  - Leaflet immobilization
  - Single or multiple valve orifice(s)
  - Limited leaflet mobility relative to the tips of both Clip Arms
  - Adequate MR reduction.

**12.4 If placing an additional Clip proceed to Section 13.0. If not placing an additional Clip proceed to section 14.0.**

**13.0 ADDITIONAL MITRACLIP DEVICE PLACEMENT**

**WARNING:** Use caution not to displace or dislodge an implanted Clip when placing an additional Clip; Clip detachment from leaflet(s) may occur which may result in a single leaflet device attachment (SLDA) or device embolization.

**13.1 When placing an additional Clip, the following are recommended:**

- 13.1.1 In the LA, ensure Clip Arms are oriented perpendicular to the line of coaptation and Grippers are raised.
- 13.1.2 Cross into the LV with a *Clip Arm Angle* of <90 degrees.
- 13.1.3 Use both fluoroscopy and echocardiography when crossing into the LV and during grasping.
- 13.1.4 *Unlock the Clip and Open the Clip Arms* to 180 degrees. Ensure that the Clip Arms are oriented perpendicular to the line of coaptation then Close the Clip to the *Grasping Arm Angle*.

**WARNING:** DO NOT use excessive force or retraction distance during grasping. This may compromise leaflet capture and insertion. Loss of leaflet capture and insertion may result in inadequate reduction of mitral regurgitation and may result in a single leaflet device attachment (SLDA).

**14.0 MITRACLIP SYSTEM REMOVAL**

**WARNING:** During MitraClip System removal always retract the CDS by pulling only on the Sleeve Handle. Retracting the CDS by pulling on the DC Handle may result in device damage and/or device or component embolization, and may result in vascular and/or cardiac injury.

**WARNING:** Failure to release the DC Fastener before releasing Sleeve curves may result in device damage and/or device or component embolism.

**WARNING:** Failure to utilize echocardiographic guidance while releasing Sleeve deflection may result in cardiac injury.

**14.1 MITRACLIP SYSTEM REMOVAL AFTER CLIP DEPLOYMENT**

- 14.1.1 Removal of the CDS While Leaving the Guide in Place.
  - 14.1.1.1 Release the DC Fastener.
  - 14.1.1.2 Slowly release Sleeve deflection by rotating the M/L Knob and the A/P Knob to neutral.
  - 14.1.1.3 Secure DC Fastener once Sleeve curves are released.
  - 14.1.1.4 Straighten the Guide with the +/- Knob when the Delivery Catheter tip is free from the left atrial wall and the mitral valve.
  - 14.1.1.5 Release the Sleeve Fastener and retract the CDS approximately 10 cm into the Guide by pulling only on the Sleeve Handle.
  - 14.1.1.6 Confirm that the Clip Introducer is still fully advanced in the Guide Hemostasis Valve.
  - 14.1.1.7 Retract the CDS by pulling only on the Sleeve Handle and position the Delivery Catheter tip inside the Clip Introducer. Begin gently aspirating the Guide (starting when the CDS is approximately halfway into the Guide, approximately 40 cm retracted) using a 50-60 cc syringe.
  - 14.1.1.8 Remove the CDS and the Clip Introducer simultaneously from the Guide by pulling on the Sleeve shaft and Clip Introducer. Ensure the Delivery Catheter tip is inside the Clip Introducer by visualizing the Proximal Sleeve alignment marker just outside the Clip Introducer. Aspirate the Guide during removal of the CDS and Clip Introducer. Cover Guide Hemostasis Valve with finger upon CDS removal. If necessary, position the Guide Handle below the level of the LA to allow blood to fill the Guide Lumen.

**WARNING:** DO NOT remove the tip of the CDS from the Guide without removing the Clip Introducer simultaneously. Failure to remove the Clip Introducer simultaneously may result in air embolism.

**WARNING:** DO NOT create a vacuum while removing the CDS from the Guide; air may enter the lumen of the Guide which may result in air embolism.

- 14.1.1.9 Aspirate using a 50-60 cc syringe to remove any remaining air from the Guide.
- 14.1.2 Removal of the CDS and Guide simultaneously.
  - 14.1.2.1 Release the DC Fastener.
  - 14.1.2.2 Slowly release Sleeve curves by rotating the M/L Knob and the A/P Knob to neutral.
  - 14.1.2.3 Secure the DC Fastener once Sleeve curves are released.
  - 14.1.2.4 Straighten the Guide with the +/- Knob when the Delivery Catheter tip is free from the left atrial wall and the mitral valve.
  - 14.1.2.5 Release the Sleeve Fastener and retract the CDS approximately 10 cm into the Guide by pulling only on the Sleeve Handle.
  - 14.1.2.6 Carefully retract the Guide tip into the RA. The Guide may be straightened further with the +/- Knob if desired.
  - 14.1.2.7 Remove the MitraClip System from the femoral vein, while providing hemostasis.

**14.2 MITRACLIP SYSTEM REMOVAL WITH CLIP ATTACHED**

- 14.2.1 Removal of the CDS while leaving the Guide in place.
  - 14.2.1.1 Confirm Clip is locked.
  - 14.2.1.2 Fully Close the Clip Arms.

**WARNING:** Failure to Fully Close the Clip Arms prior to retraction into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.

- 14.2.1.3 Lower the Grippers.
- 14.2.1.4 Release the DC Fastener and retract the DC Handle until the DC Radiopaque Ring is fully against the tip of the Sleeve.
- 14.2.1.5 Slowly release Sleeve deflection by rotating the M/L Knob and the A/P Knob to neutral.
- 14.2.1.6 Secure the DC Fastener once Sleeve curves are released.

- 14.2.1.7 Straighten the Guide with the +/- Knob when the tip of the MitraClip Device is free from the left atrial wall and the mitral valve.

**WARNING: Failure to straighten the Guide prior to retracting the Clip into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.**

- 14.2.1.8 Release the Sleeve Fastener and retract the CDS into the Guide by pulling only on the Sleeve Handle.

**NOTE:** If resistance is noted, advance and rotate the Clip by rotating the DC Handle then retract the CDS into the Guide. The Guide and/or Sleeve position may also be adjusted to facilitate Clip entry into the Guide. If necessary, retract the Sleeve or advance the Clip to create a 2-3cm separation to facilitate Clip entry into the Guide.

**WARNING: Failure to utilize fluoroscopic guidance while retracting the CDS into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.**

- 14.2.1.9 Confirm that the Clip Introducer is still fully advanced in the Guide Hemostasis Valve.

- 14.2.1.10 Retract the CDS by pulling only on the Sleeve Handle and position the Clip inside the Clip Introducer. Begin gently aspirating the Guide (starting when the CDS is approximately halfway into the Guide, approximately 40 cm retracted) using a 50-60 cc syringe.

- 14.2.1.11 Remove CDS and Clip Introducer simultaneously from the Guide by pulling on the Sleeve shaft and Clip Introducer. Ensure the Clip is inside the Clip Introducer by visualizing the Proximal Sleeve alignment marker just outside the Clip Introducer. Aspirate the Guide during removal of the CDS and Clip Introducer. If necessary, position the Guide Handle below the level of the LA to allow blood to fill the Guide lumen.

**WARNING: DO NOT remove the tip of the CDS from the Guide without removing the Clip Introducer simultaneously and with the Clip inside the Clip Introducer. Failure to remove the Clip Introducer simultaneously may result in air embolism.**

**WARNING: DO NOT create a vacuum while removing the CDS from the Guide; air may enter the lumen of the Guide which may result in air embolism.**

**WARNING: DO NOT re-use the CDS after removal. Replace the CDS with a new device. Reinserting the CDS after removal may result in inability to open the Clip. Inability to open the Clip may result in valve injury or lead to deployment of the Clip in an unintended location.**

- 14.2.1.12 Aspirate using a 50-60 cc syringe to remove any remaining air from the Guide.

- 14.2.2 Simultaneous removal of CDS and Guide.

- 14.2.2.1 Confirm Clip is locked.

- 14.2.2.2 Fully Close the Clip Arms.

**WARNING: Failure to Fully Close the Clip Arms prior to retraction into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.**

- 14.2.2.3 Lower the Grippers.

- 14.2.2.4 Release the DC Fastener and retract the DC Handle until the DC Radiopaque Ring is fully against the tip of the Sleeve.

- 14.2.2.5 Slowly release Sleeve deflection by rotating the M/L Knob and the A/P Knob to neutral.

- 14.2.2.6 Secure the DC Fastener once Sleeve curves are released.

- 14.2.2.7 Straighten the Guide with the +/- Knob when the tip of the MitraClip Device is free from the left atrial wall and the mitral valve.

**WARNING: Failure to straighten the Guide prior to retracting the Clip into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.**

- 14.2.2.8 Release the Sleeve Fastener and retract the CDS approximately 10 cm into the Guide by pulling only on the Sleeve Handle.

**NOTE:** If resistance is noted, advance and rotate the Clip by rotating the DC Handle then retract the CDS into the Guide. The Guide and/or Sleeve position may also be adjusted to facilitate Clip entry into the Guide. If necessary, retract the Sleeve or advance the Clip to create a 2-3cm separation to facilitate Clip entry into the Guide.

**WARNING: Failure to utilize fluoroscopic guidance while retracting the CDS into the Guide may result in device damage, inability to remove the CDS and/or vascular and cardiac injury.**

- 14.2.2.9 Carefully retract the Guide tip into the RA. The Guide may be straightened further with the +/- Knob if desired.

- 14.2.2.10 Remove the MitraClip System from the femoral vein, while providing hemostasis.

#### PATENTS

This device and its use are covered by one or more of the following US patents: 7,736,388; 7,682,369; 7,666,204; 7,655,015; 7,608,091; 7,604,646; 7,563,267; 7,464,712; 7,288,097; 7,226,467; 7,048,754; 6,770,083; 6,752,813; 6,629,534; 6,461,366; 6,269,819. Additional US and Foreign patents pending.

© 2015, Evalve, Inc. All rights reserved.

MitraClip is a trademark of the Abbott Group of Companies.

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