## Instructions For Use

---

- DA** Brugsanvisning
- DE** Gebrauchsanweisung
- EL** Οδηγίες χρήσης
- ES** Instrucciones de uso
- FR** Mode d'emploi
- IT** Istruzioni per l'uso
- KO** 사용 지침
- NL** Gebruiksaanwijzing
- PL** Instrukcja obsługi
- PT** Instruções de utilização
- SV** Bruksanvisning
- ZH** 使用说明

## EN | Retrieval Devices

**CAUTION/RX ONLY**

Federal Law (USA) restricts sale of this device by or on the order of a physician.

**DEVICE DESCRIPTION**

The device consists of a flexible wire cable and snare loop with net which can be extended and retracted from the flexible outer sheath using a three ring handle. The device length is 230 cm; device diameter is 2.5 mm.

**INTENDED USE**

This device is intended to be used to retrieve excised polyps, tissue samples foreign bodies and calculi during flexible and rigid endoscopy procedures.

**WARNINGS**

- This device is designed for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure which may in turn result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Do not use this device for any purpose other than the stated intended use.
- Visually inspect the package prior to use for any tears, irregularities or damage. If an abnormality is detected, the sterile condition of the device may have been compromised. Do not use the device.
- Before use, confirm that the net, catheter and handle are not damaged. If any abnormalities such as significant deformations or excessive bends are found, do not use the device.
- Confirm that there are no significant deformations, sharp parts, burrs or edges.
- Gently run your fingertips over the entire length of the catheter to check for any crushes, excessive bends, broken areas or other damage.
- Make sure that there are no cracks on the handle
- Do not store the package containing the device in places where it may become damaged, wet or improperly sealed. Otherwise, the sterility of the device may be compromised and pose an infection control risk or cause tissue irritation.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

**CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired site.

**PRECAUTIONS**

- Refer to the package label for minimum channel size required for this device.
- The endoscopic retrieval of foreign objects, food bolus, tissue samples or excised polyps should only be performed by persons having adequate training in clinical endoscopic technique.
- Avoid blindly passing the device past any foreign object, particularly if the entire lumen is blocked.
- Use care to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea.
- This device is not recommended for the retrieval of sharp foreign objects.
- After extraction, if the retrieved material does not exit the net upon deployment, the material should be removed from the net by rinsing the net in a basin of water; do not "pick" retrieved material from the net with a finger or tool as this may damage the net.

**INSTRUCTIONS FOR USE**

Read the instructions for use in its entirety prior to use.

1. Prior to advancing the device into the endoscope, uncoil the device and drape in a "U" shaped configuration. The device can be damaged if it is coiled while the handle is operated.
2. Holding the proximal end of the device in one hand and the distal sheath in the opposite hand, slowly retract the finger rings until the net is completely withdrawn into the sheath. Confirm that there are no significant deformations or abnormalities.
3. Once the foreign object, food bolus, tissue sample or excised polyp for retrieval has been identified and with the net fully retracted within the catheter, advance the catheter into the accessory channel using short, 1.0 – 1.5 cm strokes until the distal end of the sheath is endoscopically visualized.

4. Advance the distal end of the sheath slightly past the foreign object, food bolus, tissue sample or excised polyp, however, avoid blindly passing the device past any foreign object, food bolus, tissue sample or excised polyp if the entire lumen is blocked. Open the net by advancing the finger rings forward. Confirm that the net is fully open via endoscopic observation.
5. Endoscopically manipulate the net over the foreign object, food bolus, tissue sample or excised polyp.
6. Retrieve the foreign object, food bolus, tissue sample or excised polyp by slowly retracting the finger rings on the handle. Gentle traction should be applied on the handle to keep the net closed.
7. Retract the catheter retrograde toward the endoscope so that the net is endoscopically observed yet not obscuring the endoscopic view. Endoscopic observation is necessary during extubation so as not to lose sight of the foreign object, food bolus, tissue sample or excised polyp during removal.
8. Once the endoscope and the net have been extubated, the retrieved foreign object, food bolus, tissue sample or excised polyp can be removed from the net by advancing the finger rings forward to open the net. Specimens should be prepared according to standard technique for histologic evaluation.

**PRODUCT DISPOSAL**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

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**EndoChoice, Inc.**
**Atlantico Systems, Ltd.**