Solyx™

SIS System

Solyx Blue

SIS System

Directions for Use	2
Instrucciones de uso	10
Mode d'emploi	19
Gebrauchsanweisung	28
Istruzioni per l'uso	37
Gebruiksaanwijzing	46
Instruções de Utilização	55

50526171-01

2017-05

TABLE OF CONTENTS

WARNING 3

DEVICE DESCRIPTION..... 3

INTENDED USE/INDICATIONS FOR USE 3

CONTRAINDICATIONS..... 3

GENERAL WARNING..... 4

POST PROCEDURAL WARNING 4

PRECAUTIONS..... 4

ADVERSE EVENTS 5

HOW SUPPLIED..... 5

OPERATIONAL INSTRUCTIONS..... 5

 Prior to Use..... 5

 Figure 1: Solyx™ SIS System and Solyx Blue SIS System..... 6

WARNING 6

 Steps for Use..... 6

 Incision and Dissection..... 6

 Figure 2: Dissection Pathway..... 6

 Sling Placement 6

PRECAUTION 7

 Figure 3: Mesh Assembly Placement onto Delivery Device 7

 Figure 4: Mesh Orientation on Delivery Device..... 7

WARNING 7

 Figure 5..... 8

 Figure 6..... 8

 Figure 7..... 8

 Figure 8: Contralateral Delivery Device Placement..... 9

WARRANTY 9

Solyx™

SIS System

Solyx Blue

SIS System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician trained in use of surgical mesh for repair of stress urinary incontinence.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

WARNING

This product is intended for use only by clinicians with adequate training and experience in the surgical treatment of stress urinary incontinence (SUI). The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.

DEVICE DESCRIPTION

The Solyx and Solyx Blue SIS (Single Incision Sling) Systems are sterile single use systems each consisting of one (1) delivery device and one (1) mesh assembly. The mesh assembly is comprised of a polypropylene knitted mesh with polypropylene carriers at each end of the distal mesh. The carrier is designed to be placed on the tip of the delivery device. The disposable delivery device consists of a handle, a stainless steel shaft and a deployment mechanism. The delivery device is designed to facilitate the passage of the mesh assembly through bodily tissues for placement into the obturator internus muscle.

INTENDED USE/INDICATIONS FOR USE

The mesh implant is intended for use as a suburethral sling for the treatment of stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

CONTRAINDICATIONS

A mesh implant is contraindicated in the following patients:

- Pregnant patients, patients with the potential for future growth or patients who are considering future pregnancies.
- Any patients with soft tissue pathology into which the implant is to be placed.
- Patients with any pathology which would compromise implant placement.
- Patients with any pathology that would limit blood supply or infections that would compromise healing.

GENERAL WARNING

The risks and benefits of performing a suburethral sling procedure in the following patients should be carefully considered:

- Careful consideration should be given to performing this procedure for patients with untreated coagulopathies or who are being treated with either anticoagulants or antiplatelet agents.
- Patients with hypertonic bladders or vesico ureteral reflux
- Take special care in cases of bladder prolapse because of anatomical distortion. If the patient requires a cystocele repair, it should be done prior to the suburethral sling procedure
- Vaginal and urinary tract infection should be treated prior to a suburethral sling implantation procedure.
- User should be familiar with surgical procedures and techniques involving nonabsorbable meshes.
- Good surgical practices should be followed for management of contamination or infected wounds.
- Mesh is considered a permanent implant. Removal of mesh or correction of mesh related complication may involve multiple surgeries.
- Complete removal of mesh may not be possible and additional surgeries may not always fully correct the complications.

POST PROCEDURAL WARNING

- If subsequent infection occurs, follow appropriate medical intervention practices.
- The patient should be advised that future pregnancies may negate the effects of this procedure and the patient may again become incontinent.

PRECAUTIONS

- The use of polypropylene mesh in urogynecologic procedures such as the treatment of stress urinary incontinence, regardless of the route of delivery (transvaginal, suprapubic or transobturator), has been associated with cases of erosion. Erosion has been reported in bladder, vagina, urethra and ureter, and bowel. Treatment of the erosion may require surgical removal.
- As with all surgical procedures, certain risk factors are known to impact patient outcomes in the pelvic floor which include, but are not limited to, impaired vascularity (e.g. diabetes, smoking status, estrogen status, pelvic floor radiation exposure, etc.), age, pelvic floor myalgia, impaired wound healing (e.g. diabetes, steroid usage, etc.), or active infection in or near the surgical site. The above pathophysiologic conditions must be considered when determining whether the patient is an appropriate candidate for mesh implantation, either by transvaginal, suprapubic or transobturator route.
- Standard surgical practices should be followed for the suburethral sling procedure as well as for the management of contaminated or infected wounds.
- Bleeding can occur. Check carefully before releasing patient from the hospital.
- Ensure the mesh is placed tension free under the mid-urethra.
- Use of this device should be done with the understanding that subsequent infection may require removal of the mesh.
- Physician should determine when it is suitable for each patient to return to normal activities.
- Patients should be counseled when to resume vigorous activities (heavy lifting, exercise), and intercourse after the procedure.
- Should dysuria, bleeding or other problems occur, the patient should be instructed to immediately contact their physician.

- Do not use any mechanical means of contact with the mesh (such as clips, staples, etc.) within the urethral support region of the mesh as mechanical damage to the mesh may occur.
- Avoid excessive tension on the mesh during handling.

ADVERSE EVENTS

The following adverse events have been reported due to suburethral sling placement, but are not limited to:

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the mesh implant could include:
 - erosion/exposure/extrusion of the mesh through the vaginal or urethral mucosa, bladder wall or other surrounding tissue
 - scarring/scar contracture
 - device migration
 - fistula formation and inflammation

The occurrences of these events may require surgical intervention and possible removal of the entire mesh.

- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Allergic reaction has been reported.
- Known risks of surgical procedures for the treatment of incontinence include:
 - pain, ongoing pain (pelvic, vaginal, groin/thigh, dyspareunia)
 - infection
 - detrusor instability
 - complete failure of the procedure
 - voiding dysfunction (incontinence, mild to moderate incontinence due to incomplete urethral support or due to overactive bladder)
 - bruising, bleeding (vaginal, hematoma formation)
 - abscess
 - vaginal discharge
 - dehiscence of vaginal incision
 - edema and erythema at the wound site
 - Perforation or laceration of vessels, nerves, bladder or urethra may occur during placement.

The occurrence of these events may require surgical intervention. In some instances the response to these events may persist as a permanent condition after the intervention.

HOW SUPPLIED

The Solyx™ and Solyx Blue SIS Systems are sterile, single use systems each consisting of one (1) delivery device and one (1) mesh assembly. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Store in a cool, dry, dark place. Rotate inventory so that devices are used prior to the expiration date displayed on package label.

OPERATIONAL INSTRUCTIONS

Prior to Use

The Solyx and Solyx Blue SIS Systems are supplied sterile and are intended for single patient use only. Carefully examine the system to verify that neither the contents nor the sterilized package has been damaged in shipment. **DO NOT USE** if the sterile barrier or product is damaged. Immediately return damaged product to Boston Scientific. The design of both systems allows the operator a transvaginal route of delivery. **See Figure 1** for parts description.


Figure 1: Solyx™ SIS System and Solyx Blue SIS System

WARNING

Ensure that the bladder is empty prior to initiating the use of this product. Ensure that the bladder, urethra and other important landmarks are properly identified.

Steps for Use

Incision and Dissection

1. Prepare and drape the patient using standard surgical practice. Ensure the bladder is empty.
2. Make a 1.0 cm to 1.5 cm vertical midline incision on the anterior vaginal wall at the level of the mid-urethra. Dissect bilaterally to the interior portion of the inferior pubic ramus at the 45° angle off the midline creating a pathway for delivery device placement.
3. Dissect bilaterally to the interior portion of the inferior pubic ramus at a 45° angle off the midline creating a pathway for delivery device placement. (See Figure 2).


Figure 2: Dissection Pathway

Sling Placement

4. Place the mesh assembly onto the delivery device by placing the delivery device tip into the mesh carrier (See Figure 3). Note: The mesh carrier should be pushed onto the delivery device tip until it is flush with the end of the delivery device shoulder (See Figure 3). While placing the mesh carrier onto the delivery device, make sure that the mesh is oriented so that it lies on the outside curve of the delivery device (See Figure 4). The mesh is now ready for placement.

PRECAUTION

Hold each side of the mesh carrier (See Figure 3) to avoid possible glove puncture.


Figure 3: Mesh Assembly Placement onto Delivery Device


Figure 4: Mesh Orientation on Delivery Device

5. Insert the delivery device into the dissection pathway targeting placement of the carrier at a 45° angle off the midline. (See Figure 5) Advance the delivery device towards the obturator foramen just lateral to the inferior pubic ramus (See Figure 5) until the midline mark on the delivery device is approximately at the midline position under the urethra. Deposit the carrier by gripping the deployment mechanism with one hand and pulling the delivery device handle back with the other hand (See Figure 6). This action will deposit the carrier into surrounding muscle tissue releasing it from the delivery device.

NOTE: Once the carrier is deposited in tissue, it is not designed to be reconnected onto the shaft tip for additional tension/adjustment.

WARNING

If excessive resistance is encountered during advancement/withdrawal, stop and determine remedial action prior to proceeding.


Figure 5.


Figure 6.


Figure 7

6. Remove the delivery device from the incision. Bring the deployment mechanism back to start position, flush with the handle (See Figure 7).

NOTE: Once the carrier is deposited in tissue, it is not designed to be reconnected onto the shaft tip for additional tension/adjustment.

7. Place the second mesh carrier onto the delivery device following instructions from Step 4. Make sure that the mesh is not twisted so that it will lie flat under the urethra when placement is complete.

8. Insert the delivery device into the dissection pathway on the contralateral side targeting placement of the carrier at a 45° angle off the midline (See Figure 8). Advance the delivery device towards the obturator foramen just lateral to the inferior pubic ramus (See Figure 8) until appropriate mesh placement is achieved. Deposit the carrier by gripping the deployment mechanism with one hand and pulling the delivery device handle back with the other hand (See Figure 6). This action will deposit the carrier into surrounding muscle tissue, releasing it from the delivery device.


Figure 8: Contralateral Delivery Device Placement

9. Remove the delivery device from the incision.
10. Cystoscopy may be performed at this time, to be determined at the physician's discretion.
11. Close incisions according to usual methods.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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All rights reserved.
**Boston Scientific Limited
**Boston Scientific (Australia) Pty Ltd
© 2017 Boston Scientific Corporation or its affiliates.