DJO, LLC
1430 Decision Street
Vista
CA 92081-8553
USA

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CLASS I **(10-15 mmHg)**

Indications

- Diseases of the superficial venous system
- Functional disorders (tired, aching, heavy legs)
- Superficial varicosis without oedema
- Prevention (travel & long journeys, phlebitis and oedema)

CLASS II **(15-20 mmHg)**

Indications

- Diseases of the superficial & fascial venous system
- Functional disorders (tired, aching, heavy legs)
- Varicosis with oedema
- Superficial thrombophlebitis
- Post-sclerotherapy
- Prevention (travel & long journeys, phlebitis and oedema)
- Overweight
- Professions requiring long periods of standing

CLASS III **(20-36 mmHg)**

Indications

- Prevention of deep vein thrombosis in patients confined to bed rest
- Diseases of the deep venous system
- Post-operatively (1 month)
- Varicose veins complicated by atrophy
- Severe postphlebotic cases
- Postphlebotic syndrome

Contra-indications

- Advanced peripheral obstructive arterial disease
- Severe neuropathy (diabetic)
- Arterial ulcers
- Systemic oedema

Composition

Microtrans/Fantaisie/Transparent:

Elastane/Polyamide

Coton:

Elastane/Polyamide/Cotton


For Single Patient Use Only

Our knee high stockings, hold ups and tights are manufactured to a high quality, attaining an unmatched level of transparency. Nevertheless, it is important to take care when handling the stockings to ensure that they do not become damaged.

HOW DO YOU PUT ON YOUR STOCKINGS?

Prior to applying the stockings, it is advisable to elevate the legs above the heart for at least 5 minutes, to ensure blood return. Ideally, the stockings should be applied in the morning before getting up.

1. Slip the hand into the stocking.
2. Turn the stocking inside out, keeping hold of the heel with your hand.
3. Slip your foot in gently up to the top of the stocking. Ensure that the heel is positioned correctly avoiding any creases.
4. Gently slide the stocking up the leg and adjust accordingly.
5. Smooth out any wrinkles along the leg.

While putting on your stockings, take care not to catch the material with your nails or rings. To avoid tearing, be careful not to pull on the silicone band.

Care Instructions:

1. Delicate wash up to 40°C.
2. Use mild soap, do not use chlorine bleach.
3. Rinse well with tepid water.
4. Do not tumble dry.
5. Dry flat, away from the sun or direct heat.

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DJO Asia-Pacific Limited
DJO Global, Inc
U.S.A.