
ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: For immobilization and protected range of motion associated with posttraumatic treatment. ACL, PCL, LCL, and MCL surgeries, patella realignment and meniscal or articular repairs. This product was designed to complement the variety of medical treatments common to the above afflictions. The range of motion settings and frequency and duration of use should be determined by your prescribing Healthcare professional.

CONTRAINDICATIONS: None

WARNINGS AND PRECAUTIONS: If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

APPLICATION INFORMATION:

1. Unfasten strap ends and open brace flat.
2. Place leg on top of flat brace. Shorter upright is applied to the thigh.
3. Telescope paddles to desired length by depressing the slider button until the desired length is reached. Ensure that both sides are equal in length.
4. Align the upright with the attached foam liner to midline of leg. Wrap foam pads snugly around leg, leaving the opposite upright unattached, and attach anteriorly as shown. Position and press into place the remaining upright. Make sure medial and lateral hinge height are equal.
5. Loop the straps through the D-rings and fasten the two straps closest to the knee first, then fasten the remaining straps.
6. For easy, one piece removal, unfasten the straps and open the foam pads. The brace can then be removed and reapplied as a single unit.

INSTRUCTIONS FOR HINGE ADJUSTMENT:

7A) Adjust the hinge by depressing the Flexion and Extension buttons inward and rotate stop until the desired angle aligns with button center. Release push button and ensure it is fully retracted to the extend position.

7B) Utilize the quick-lock feature at 0, 15, or 30 degrees, by sliding the red quick lock button inward until the lock is engaged. Ensure uprights do not rotate. If another lock setting is required, see Step 5A and ensure both Flexion and Extension buttons are set at the same angle setting.

OTHER ADJUSTMENTS:

- Bend each bar by holding thigh/calf bar firmly against a solid surface and apply gentle and constant pressure to the hinge in the direction desired.
- Bend each side bar an equal amount above and below the hinge.
- For optimal ease of application post-operatively, pre-fit the brace prior to surgery if possible.

CLEANING INSTRUCTIONS:

Hand wash foam liners in water (30°C) with mild detergent. Rinse thoroughly.
AIR DRY only, do not heat dry.
Regular cleaning of the brace is recommended.
Always consult your physician or therapist before altering the brace.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

COMPOSITION: Polypropylene 30%, Aluminum 40%, Nylon 25%, Delrin 3%, Stainless Steel 2%.

INTENDED FOR USE ON A SINGLE PATIENT.

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