




ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INDICATIONS: Management of edema and fracture healing • Acute ankle sprains • Soft tissue injuries of the lower leg • Stress and stable fractures of the lower leg and ankle • Stable delayed union or non-union fractures of the distal tibia and fibula • Achilles tendon repairs

CONTRAINDICATIONS: This device is contraindicated for tibial fractures and unstable fractures of the lower leg, ankle and foot. Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity. The determination of when to apply the brace rests strictly at the discretion of the treating physician.

This product was designed to complement the variety of medical treatments common to the above afflictions. The frequency and duration of use would be determined by your prescribing Healthcare professional.

Inflate and deflate bladder once before use to ensure integrity of air chambers.

- 1) Select correct size of walker for the patient. The end of the patient's toes should be within 1 in (2.5 cm) from, without extending past, the front end of the walker. Loosen straps and remove liner from boot.
- 2) Place foot in liner and secure with contact closure. Make sure the heel fits snugly into the posterior portion of the liner. Secure the liner firmly around the foot, ankle, and lower leg.
- 3) Lift the end of the boot inner sole and position one end of the toe wrap on the contact closure strip.
- 4) Step into boot with heel seated firmly to the back of the shell. Fold the toe wrap over toes and secure on liner.
- 5) Position the front shell so the distal edges (toe region) fit inside the rear shell and the proximal edges (shin region) fit outside the rear shell. Fasten the front shell to the back shell by securing the straps, starting with the instep strap first, followed by all other straps. Top strap must be secured outside of shell.
- 6) To inflate walker liner press and release the air bulb several times until liner is providing comfortable cushion and support. DO NOT OVER INFLATE.
- 7) To release air from the walker liner, turn the inflation/deflation knob several times counter-clockwise. When re-inflating, turn knob clockwise until snug. Release air if liner feels too tight.

- Accessory pads may be used where necessary to aid in patient comfort
- Do Not over tighten straps. This may result in reduction of blood flow or sensation.
- Not for use on patients with compromised sensation in the foot and ankle region.
- Patients must use care when walking on slippery or wet surfaces while wearing the walker.

USE AND CARE:

- To clean liners, hand wash in cold water using mild detergent and air dry. Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

Always consult with your physician or therapist before making changes to the brace.

CAUTION: Patients with diabetes, vascular deficiencies and neuropathies should be advised to perform circulation checks.

WARNING: This product is to be used under the supervision of a medical professional. this device is not intended for unsupervised public use. If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

MATERIAL COMPONENTS:

Shell: POLYPROPYLENE
Outer Sole: NATURAL RUBBER
Inner Sole: EVA
Straps: NYLON
Liner: NYLON AND FOAM

NOT MADE WITH NATURAL RUBBER LATEX. FOR SINGLE PATIENT USE ONLY. RX PRODUCT

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT IN THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL. THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.

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