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INSTRUCTIONS FOR USE – SMR SHOULDER SYSTEM

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g. product literature, surgical technique).

1. PRODUCT INFORMATION

The SMR Shoulder System is intended for partial or total primary or revision shoulder joint replacement. The components are intended for use in cemented and uncemented applications as specified in the table of section 2.1. The main goal of the SMR Anatomic joint prosthesis is to reproduce the articular anatomy, in part or totally, while in the SMR Reverse shoulder the geometry of the natural gleno-humeral joint is inverted. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depend in part on the pre-operative situation, intra-operative options and the post-op rehabilitation. SMR shoulder System components include humeral stems, humeral bodies, adaptor tapers, humeral heads, and glenoid components. SMR shoulder components should not be used with components from other systems or other manufacturers.

NOTE: Some components of SMR System are allowed for the use with SMR Stemless System.

1.1. MATERIALS

The materials of the SMR Shoulder System are reported in the table of Section 2.1. Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon. Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone cement and applying bone cement for prosthesis fixation.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50°C/32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened, ensure that both the model and size of the implant correspond to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged.

Components removed from the package should not be used if they are dropped or suffer other accidental trauma. Devices should not be modified in any way.

The device's code and lot number should be recorded in the patient's case history by using the labels included in the component packaging. The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

The SMR Shoulder System is intended for partial or total, primary or revision shoulder joint replacement in skeletally mature patients. The SMR Anatomic Shoulder System is indicated for partial or total, primary or revision shoulder joint replacement in patients suffering from disability due to: • non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis; • inflammatory degenerative joint disease such as rheumatoid arthritis; • treatment of acute fractures of the humeral head that cannot be treated with other fracture fixation methods; • revision of a failed primary implant; • cuff tear arthropathy (CTA Heads only).

The Large Resection Stems are indicated for oncology applications.

The SMR Reverse Shoulder System is indicated for primary, fracture or revision total shoulder replacement in a grossly rotator cuff deficient joint with severe arthropathy (disabled shoulder). The patient's joint must be anatomically and structurally suited to receive the selected implants and a functional deltoid muscle is necessary to use the device. The Modular SMR Shoulder System allows the assembly of components in various humeral and glenoid constructs. The constructs are intended for cemented and uncemented use as specified in the following table. In the Anatomic shoulder the humeral construct consists of the humeral stem, the humeral body, the adaptor taper and the humeral head. In the Reverse shoulder the humeral construct consists of the humeral stem, the reverse humeral body and the reverse liner. On the humeral side the fixation of the humeral stem determines if the construct is cemented or uncemented. The Anatomic glenoid construct consists of an all polyethylene glenoid or a metal back assembled with a liner while the Reverse glenoid construct consists of the metal back, the connector and the glenosphere. On the glenoid side, the fixation of the all polyethylene glenoid or the metal back determines if the construct is cemented or uncemented. The intended use of the Glenoid Bone Graft Instruments (not approved in Canada) is to enable the preparation of the bone graft from the patient's humeral head or allograft from the bone bank, to restore the glenoid anatomy in case of glenoid deficiency (e.g. glenoid type B2 or C according to Walch's classification). The Glenoid Bone Graft Instruments are intended to be used only with Axioma TT Metal Back and L1 TT Metal Back.

System		Components	Material	Use		Available in Canada
A	R			C	NC	
•	•	SMR Stems (Cemented, Cemented Revision)	Ti6Al4V	X		•
•	•	SMR Stems (Cementless Finned, Cementless Revision)	Ti6Al4V		X	•
•	•	SMR Large Resection stems	Ti6Al4V	X		•
•	•	SMR Modular Spacers	Ti6Al4V	X		•
•		SMR Humeral Bodies (Trauma, Finned)	Ti6Al4V	X	X	•
•	•	SMR Reverse Humeral Body	Ti6Al4V	X	X	•
•	•	SMR Reverse HA Coated Humeral Body	Ti6Al4V+HA	X	X	•
	•	Humeral Extension	Ti6Al4V	X	X	•
•		SMR Humeral Heads (Standard, CTA)	CoCrMo	X	X	•
			Ti6Al4V	X	X	
•		SMR Adaptor Tapers (Neutral, Eccentric)	Ti6Al4V	X	X	•
•		SMR CTA Head Adaptor for Reverse Humeral Body	Ti6Al4V	X	X	•

	•	SMR Glenospheres	CoCrMo		X	•
			Ti6Al4V		X	•
			UHMWPE X-Lima +Ti6Al4V		X	
	•	SMR Connectors	Ti6Al4V		X	•
	•	Reverse Liners	UHMWPE	X	X	•
			UHMWPE X-Lima	X	X	
			CoCrMo	X	X	
			Alumina	X	X	
•		SMR Cemented Glenoids	UHMWPE	X		•
•		SMR 3 Pegs Cemented Glenoids	UHMWPE X-Lima	X		
			UHMWPE *	X		•
•	•	SMR Metal Back Glenoids	Ti6Al4V+PoroTi		X	•
			Ti6Al4V+PoroTi+HA		X	•
•	•	SMR TT Metal Back Baseplate	Ti6Al4V		X	•

•	•	SMR TT Metal Back Peg	Ti6Al4V		X	•
•		SMR Metal Back Liner	UHMWPE		X	•
			UHMWPE X-Lima		X	
	•	SMR Glenoid Plates	Ti		X	

Material Standards

Ti6Al4V (ISO 5832-3 – ASTM F1472) – **CoCrMo** (ISO 5832-12 – ASTM F1537) – **Ti** (ASTM F67) – **UHMWPE** (ISO 5834-2 – ASTM F648) – **Alumina** (ISO 6474) – **PoroTi Titanium Coating** (ASTM F1580) – **HA Hydroxyapatite Coating** (ISO 13779)

A = Anatomic; R = Reverse; C = Cemented; NC = Not Cemented

* Not available in EU.

2.2. CONTRAINDICATIONS

Absolute contraindications include: • local or systemic infection; • septicaemia; • persistent acute or chronic osteomyelitis; • confirmed nerve lesion compromising shoulder joint function; • deltoid muscle insufficiency.

Relative contraindications include: • vascular or nerve diseases affecting the concerned limb; • poor bone stock (for example due to osteoporosis or extended previous revision surgery) compromising the stability of the implant; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials.

In cases of bone tumors, use an appropriate system designed to treat cases requiring large bone resections (SMR Large Resections stems). The use of primary or revision implants not designed and intended for use in cases of bone resection may result in a poor outcome and / or failure of the implant or implant fixation.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • wrong size of components; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or

resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • use our combinations with products, prosthesis or instruments of another manufacturer; • errors of operative technique.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques. Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

In selecting patients for surgery, the following factors can be critical to the eventual success of the procedure:

- **Partial Shoulder replacement:** in cases with a deficient and unreconstructable rotator cuff, a CTA-head is indicated;
- **Total Shoulder replacement:** the rotator cuff must be intact or reconstructable. In cases with a deficient and unreconstructable rotator cuff, a hemiprosthesis with a CTA head or a Reverse Total Shoulder Arthroplasty is indicated;
- **Reverse Shoulder replacement:** the bone stock of the glenoid and humerus must be able to support the implant. In cases with significant bone loss and in which adequate fixation on the glenoid side cannot be obtained, a hemiarthroplasty with a CTA-head should be performed.
- **Bone Grafting (not approved in Canada):** once removed the graft should be inspected to ensure the bone quality is adequate for use with the glenoid bone grafting technique. The glenoid bone grafting technique should never be used with poor quality bone, as it may compromise bone healing.

NOTE: with CTA heads the use of Trauma Humeral Bodies is recommended to avoid possible impingement between the head and the body when using the Finned Humeral Body and the eccentric tapers in the lower position.

NOTE: the size large Metal Back Glenoid is not suitable for coupling with the 36 mm Glenospheres and the 40 mm Glenospheres.

NOTE: For glenoid bone grafting techniques (not approved in Canada) it is important to use a peg size Medium, Long or X-Long in order to enable a minimum section of the peg into the native bone, providing component stability. The following table identifies the allowed (✓) / not allowed (X) combinations between the bone graft thickness and the pegs dimensions:

		TT METAL-BACK - PEG SIZE							
		SMALL-R				SMALL / STD			
		Short	Medium	Long	X-Long	Short	Medium	Long	X-Long
Bone Graft Thickness	5 mm GRAFT	X	✓	✓	✓	X	✓	✓	✓
	10 mm GRAFT	X	X	✓	✓	X	✓	✓	✓
	15 mm GRAFT	X	X	✓	✓	X	X	✓	✓
	15° SLOPED GRAFT	X	✓	✓	✓	X	✓	✓	✓
	20° SLOPED GRAFT	X	✓	✓	✓	X	✓	✓	✓

NOTE: The Dia. 50, 52 and 54 mm Humeral Heads with + 3mm increased height cannot be coupled to the Long Adaptor Tapers (both concentric and eccentric).

The Dia. 52 and 54 mm Humeral Heads with + 2mm increased height cannot be coupled to the Long Adaptor Tapers (both concentric and eccentric).

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future.

The possible impact of the factors mentioned in section 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and survival rate of the implant. The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications. Use only instruments and trial components specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery. Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues.

The surgeon should make the patient aware of the limitations of limb function after shoulder arthroplasty and that the reconstructed joint must be protected from full load bearing for a period of time. Excessive physical activity or trauma to the replaced shoulder joint can lead to premature failure of the shoulder arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear. In particular the following precautions should be presented to the patient by the surgeon:

- avoid repeated high weight lifting;
- keep body weight under control;
- avoid sudden peak loads (consequences of activities such as contact sports, playing tennis) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation, such as for example excessive degrees of abduction and/or internal rotation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure. Care must be taken when handling Trabecular Titanium devices; they should not get in contact with fabrics that can release particles inside the trabecular structure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects can occur in shoulder arthroplasty include:

- loosening of the prosthetic components, for example due to osteolysis;
- prosthesis dislocation and instability;
- breakage or damage of the device;

- stem subsidence; • instability of the system because of inadequate soft tissue balancing; • dissociation due to incorrect coupling of the devices; • infection; • local hypersensitivity; • local pain; • periprosthetic fractures, including intraoperative fractures; • bone resorption; • temporary or permanent nerve damage; • fretting of modular junctions; • excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles; • metallosis; • heterotopic ossification; • scapular notching; • rotator cuff failure; • glenoid erosion (only in hemi procedures); • loss of range of motion; • additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, haematomas, systemic allergic reactions, systemic pain, blood loss, donor site morbidity, temporary or permanent paralysis.

4. STERILITY

a. Implants

All implantable components of the SMR Shoulder System are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by EtO or radiation (E-beam or Gamma rays) while UHMWPE components are sterilized by EtO. HA coated implants are sterilized only by radiation (E-beam or Gamma rays). Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the “Instrument Care, Cleaning, Disinfections and Sterilization” brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

The SMR System has not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near the implant side.

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Limacorporate S.p.A.