Delta Acetabular System

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INSTRUCTIONS FOR USE – DELTA ACETABULAR SYSTEM

Before using a LimaCorporate product, the surgeon must carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g.: product literature, surgical technique).

For each product of DELTA System refer to the specific surgical technique.

1. PRODUCT INFORMATION

DELTA System is intended to replace the natural acetabulum in total hip arthroplasty in cases of primary and revision surgery. The components are intended for use in cementless applications.

The main goal of joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give back articular mobility to the patient. The degree of pain reduction and mobility depend, in part, on the pre-operative situation, intra-operative options and the post-op rehabilitation.

DELTA acetabular system includes acetabular metal shells (TT, PF, Fins, ST-C, One TT, Revision, Revision TT and Multihole TT cups), liners and other accessories (spacers, hemispheric modules and bone screws). The components of DELTA acetabular system must not be used outside of the allowed combinations or with components belonging to other manufacturers.

Information on the compatibility of the DELTA acetabular system components with components from other LimaCorporate hip systems is available from LimaCorporate.

1.1. MATERIALS

Components	Material
DELTA TT acetabular cup	Ti6Al4V
DELTA PF acetabular cup	Ti6Al4V+PoroTi+HA
	Ti6Al4V+PoroTi
DELTA Fins acetabular cup	Ti6Al4V+PoroTi+HA
DELTA ST-C acetabular cup	Ti6Al4V+PoroTi
	Ti6Al4V+PoroTi+HA
DELTA One TT acetabular cup	Ti6Al4V

DELTA Multihole TT acetabular cup	Ti6Al4V
DELTA Revision TT acetabular cup	Ti
DELTA Revision acetabular cup	Ti+PoroTi
Hemispheric modules	Ti6Al4V
Ceramic liners	Bilox®Delta
Polyethylene liners	UHMWPE+Ti6Al4V
	UHMWPE X-Lima*+Ti6Al4V
	LimaVit™ (UHMWPE X-Lima* + Vitamin E) +Ti6Al4V
Dual mobility liners	CoCrMo
Mobile liner	UHMWPE
	LimaVit™ (UHMWPE X-Lima* + Vitamin E)
Spacers	Ti6Al4V
Bone screws	Ti6Al4V
Material Standards	
Ti6Al4V (ISO 5832-3 - ASTM F1472) – Ti (ISO 5832-2 – ASTM F67) – UHMWPE (ISO 5834-2- ASTM F648) - PoroTi Titanium Coating (ASTM F1580) – HA Hydroxyapatite Coating (ISO 13779) – CoCrMo (ISO 5832-12- ASTM F1537) - Bilox®Delta (ISO 6474-2)	

*crosslinked UHMWPE

Some patients may be hypersensitive or allergic to the implant materials; this must be appropriately considered and evaluated by the surgeon prior the surgery.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50° C / 32-122°F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature. Once the package is opened, ensure that both the model and size of the implant correspond to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged.

Components removed from the package must not be used if they are dropped or suffer other accidental trauma. Devices must not be modified in any way.

The device's code and lot number need to be recorded in the patient's case history by using the labels included in the component packaging.

The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g.: taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g.: HIV, hepatitis); • immune system response/rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

The DELTA System is indicated for use in total hip arthroplasty. DELTA PF, DELTA ST-C, DELTA Fins, DELTA TT, DELTA One TT and DELTA Multihole TT acetabular cups are intended to be used for reduction or relief of pain and/or improved hip function in skeletally mature patients with the following conditions:

- non-inflammatory degenerative joint disease such as osteoarthritis or avascular necrosis;
- hip dislocation using protruded liners, spacers (which allow better protection, multiple offset and soft tissue tensioning) and using coupling for dual mobility (which reduces jumping distance);
- rheumatoid arthritis;
- post-traumatic arthritis;
- correction of functional deformity in case of acetabulum verticalization, anteversion and retroversion;
- fractures of femoral neck.

DELTA TT, DELTA One TT, DELTA Multihole TT, DELTA Revision and DELTA Revision TT cups with spacers and hemispheric modules are indicated also in:

- revision of previous unsuccessful femoral head replacement, cup arthroplasty or other procedure;
- clinical management problem where arthrodesis or alternative reconstruction techniques are less likely to achieve satisfactory results;
- where bone stock is of poor quality or is inadequate for other reconstructive techniques as indicated by deficiencies of the acetabulum.

The DELTA system is intended for cementless use.

2.2. CONTRAINDICATIONS

Absolute contraindications include:

- local or systemic infection;
- septicemia;
- persistent acute or chronic osteomyelitis;
- confirmed nerve or muscle lesion compromising hip joint function.

Relative contraindications include:

- vascular or nerve diseases affecting the concerned limb;
- poor bone stock (for example due to osteoporosis) compromising the stability of the implant;
- metabolic disorders which may impair fixation and stability of the implant;
- any concomitant disease and dependence that might affect the implanted prosthesis;
- metal hypersensitivity to implant materials.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • medical disabilities which can lead to an unnatural gait and loading of the hip joint; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

COMBINATIONS ALLOWED/NOT ALLOWED

The allowed and not allowed combinations are also reported in the DELTA System surgical techniques.

- DELTA PF cups outer diameter 50 and 52 mm for liners size Large must not be coupled with BIOLOX® DELTA liners;
- All angled spacers are not allowed for coupling with DELTA TT, DELTA PF, DELTA Fins and DELTA ST-C cups;
- Extra Large spacers can be coupled with DELTA Multihole TT cup only, using liners size Large;
- BIOLOX® DELTA liners must be coupled only with BIOLOX® DELTA femoral heads;
- Bone screws are intended to increase primary stability in acetabular reconstruction surgery if needed. The following LimaCorporate acetabular components allow for additional bone screw fixation: DELTA TT, DELTA PF, DELTA Fins, DELTA One TT, DELTA Multihole TT, DELTA Revision, DELTA Revision TT;
- Any Liner size must match the size of any DELTA cup or internal size of the spacer;
- The Spacer size has to match the DELTA cup size used;
- The DELTA Dual Mobility with CoCrMo liner is allowed with DELTA TT, PF, Fins, ST-C, One TT, Multihole TT, Revision, Revision TT;
- The DELTA Dual Mobility with ceramic liner is allowed with DELTA PF, DELTA ST-C, DELTA Fins and DELTA TT; with DELTA Multihole TT in combination with spacer size Extra-Large only;
- DELTA Dual Mobility System can be coupled only with femoral heads dia. 28mm size S, M and L;
- DELTA One TT, DELTA Multihole TT, DELTA Revision, DELTA Revision TT use with BIOLOX® DELTA liner is allowed only in combination with spacers.

Use of Skirted Femoral Heads:

In primary surgery, due to reduced flexion - extension range of motion, skirted femoral heads (#XL, #XXL, #XXXL) should not be used with protruded liners.

In revision surgery where the hip femoral stem is left in place and the acetabulum is reconstructed, surgeons may consider the use of a skirted femoral head with the DELTA system Cup protruded liner to achieve appropriate joint stability. However, this may lead to impingement between the head and the cup liner that could potentially cause damage to the implants requiring further surgery. If the surgeon believes that the use of a skirted head with a protruded liner is necessary to achieve joint stability, the patient should be warned of the chance of impingement of the components as well as the possible consequences, including the risks and possible complications. In such circumstances, the patient should also be advised to limit flexion - extension joint movements to help minimize the potential risk for impingement and associated complications.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

The surgeon should carefully plan the surgery, it is essential for the safety of the device that the covering angle of the cup does not exceed 45° (an angle of 40° is better) and that the angle of ante-version is between 10° and 20°.

Complications or failures of the total hip replacement are more likely to occur in heavy and highly active patients. The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing hip replacement.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future.

The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and the survival rate of the implant.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications.

Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants must be protected from scratches or any other damage. All component coupling surfaces must be clean and dry before assembly. The stability of component couplings must be verified as described in the surgical technique.

Care must be taken when handling Trabecular Titanium devices; they must not get in contact with fabrics as this can release particles inside the Trabecular structure.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues. The surgeon should make the patient aware of the limitations of limb function after hip arthroplasty and that the reconstructed joint must be protected from full load bearing for a period of time. Excessive physical activity or trauma to the replaced hip joint can lead to premature failure of the hip arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon: • avoid repeated high weight lifting; • avoid excessive or repeated stair climbing; • keep body weight under control, overweight conditions may adversely affect the outcomes of joint replacement; • avoid sudden peak loads (consequences of activities such as running and skiing) or movements which can lead to sudden stops or twisting; • avoid positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that most commonly and frequently occur in hip arthroplasty include: • excessive wear of UHMWPE components (i.e. damaged articular surfaces or presence of particles); • loosening of the prosthetic components; • prosthesis dislocation and instability; • damage to the prosthetic implant; • instability of the system because of inadequate soft tissue balancing; • dissociation due to incorrect coupling of the devices; • infection; • local hypersensitivity; • local pain; • temporary or permanent nerve damage; • fractures of the devices; • limb length discrepancies; • additional surgery; • blood loss; • squeaking.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular disturbances, haematomas, systemic allergic reactions and systemic pain.

4. STERILITY

a. Implants

All implantable components of the DELTA Acetabular system are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by radiation or EtO and UHMWPE components by EtO. HA coated and Trabecular Titanium implants are sterilized by E-Beam and Gamma radiation. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to the component geometry and material, as well as the MR power, duration and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants. The DELTA System has not been evaluated for safety and compatibility in the MR environment. The DELTA System has not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy. These systems include components which are passive metallic devices and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities, including image distortion in MR and X-ray scanner

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Limacorporate S.p.A.