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Instructions for use

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PTI0001-006-0720
(Ausgabedatum: 23.07.2020)

Inhaltsverzeichnis / Content

Deutsch	4
English	36
Español	66
Português	102
Italiano	138
Français	172
USA	206

Inhalt

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English

Contents

1. Symbols	38
2. Description	39
2.1 Intended Purpose.....	39
2.2 Indications	39
2.3 Contraindications.....	40
2.4 Factors Interfering With Success	40
2.5 Possible Side Effects.....	41
2.6 System Compatibility	43
2.6.1 Primary Surgery.....	43
2.6.2 Revision Surgery	43
2.6.3 BioBall™ Adapter Neck Length Sizes	44
3. Use	46
3.1 General Instructions.....	46
3.2 Instructions on Patient Information	48
3.3 Features.....	49
3.4 Pre-Operative Planning	52
3.5 Implant Handling.....	53
3.5.1 Instructions for Handling	53
3.5.2 Sterile Implants.....	57
3.5.3 Resterilisation.....	58
3.5.4 Storage and Handling.....	58

3.6 Handling of Instruments	60
3.6.1 Cleaning, disinfection and sterilisation instructions for reusable surgical instruments	60
3.6.2 Warnings	60
3.6.3 Limitations on Reprocessing	61
3.6.4 Location of Use	61
3.6.5 Preparation for Cleaning/Pre-Cleaning.....	61
3.6.6 Automatic Cleaning and Disinfection.....	61
3.6.7 Manual Cleaning.....	62
3.6.8 Instrument Maintenance.....	62
3.6.9 Checking Functionality	62
3.6.10 Cleaning and Disinfecting Empty Trays.....	63
3.6.11 Packing	63
3.6.12 Sterilisation.....	63
3.6.13 Storage.....	63
3.6.14 Processing Instructions in accordance with DIN EN ISO 17664.....	64
3.6.15 Reusable Instruments.....	64
4. Explantation	64
5. Information Updates.....	64

1. Symbols

 0482	CE mark		Temperature limit
	Manufacturer		Do not re-use
	Date of manufacture		Consult instructions for use
	Use-by date		Caution
	Batch code		
	Catalogue number		Quantity
	Sterilized using ethylene oxide		
	Sterilized using irradiation		
	Do not re-sterilize		
	Non-sterile		
	Do not use if package is damaged		
	Keep dry		

2. Description

2.1 Intended Purpose

Use of implants contrary to intended purpose

- Risk of injury due to implant fracture!
- Implants must only be used in accordance with intended purpose.

BioBall™ Adapters 12/14, 14/16, MS 8/10, MS 10/12, MSZI, MSV4, MST1, MSBG, MSPC, MSSR and MSSY are for use as a replacement part in hip revision operations in combination with a BioBall™ head. They serve to preserve the existing anchored hip stem or total hip endoprosthesis (Hip TEP). The BioBall™ Adapter 12/14 can also be used during the primary operation for correcting positioning with only the approved stems of the Merete GmbH.

2.2 Indications

- Bearing couple revisions
- Intraoperative correction of offset, neck length, lateralisation and anteversion/retroversion with anchored prosthesis stem
- 12/14 adapter: intraoperative correction of offset, neck length, lateralisation and anteversion/retroversion during primary operation as well

2.3 Contraindications

- Acute or chronic infections in the hip joint or the immediate vicinity
- Patients with joint diseases that may be successfully treated with another, joint salvage treatment
- Any comorbidities that could pose a risk to the function or success of the implant, especially severe muscular, nervous or vascular disorders with specific effects on the limb to be operated upon
- Severely damaged in-situ stem tapers (visible changes in shape, or palpable defects, such as localised wear, abrasion/material loss, or scratches/ridges) or implants which cannot be clearly identified
- Allergies to any of the materials used

2.4 Factors Interfering With Success

- Obesity or pre-obesity
- Local bone tumours
- Osteoporosis or osteomalacia
- Malformations, congenital hip dislocation, severe axial displacement of the knee
- Systemic or metabolic disorders
- Smoking, Alcohol or drug abuse
- Physical activities associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads (e.g. hard physical labour, certain types of sport)
- Patients with pathological mental or neurological conditions, or patients who are not capable of following the necessary post-operative care instructions

2.5 Possible Side Effects

MRI Safety Information/Indications for Use

Non-clinical testing has demonstrated that the Merete Hip Implant System (consisting of cemented or non-cemented hip stem, taper adapter, metal or ceramic head ball, inlay and cup from the materials TiAl6V4 ELI (ISO 5832-3), Vivium® (ISO 5832-9), CoCrMo (ISO 5832-4/5832-12), BIOLOX® delta ceramic (ISO 6474-2), PE/XPE (ISO 5834-2)) is MR conditional.

A patient with the entire assembled Merete Hip Implant System can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla.
- Maximum spatial gradient field of 3,000 Gauss/cm (30 T/m).
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) at 1.5 Tesla or 3.0 Tesla of 1 W/kg for 15 minutes of scanning. Under the scan conditions defined above, the Merete Hip Implant System is expected to produce a maximum temperature rise of less than 6°C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Merete Hip Implant System extends at least 1 cm and up to approximately 8 cm from the device and exhibits geometric distortion in the image when imaged with a gradient echo pulse sequence or a fast-spin echo pulse sequence and a 1.5 Tesla MRI system or a 3.0 Tesla MRI system.

The side effects listed below are among the most common adverse effects of implantation procedures:

- Early or late onset infections
- Dislocation, subluxation, insufficient range of motion, undesired shortening or lengthening of the affected extremities as a result of suboptimal implant positioning
- Bone fractures due to unilateral overuse or weakened bone substance
- Reduced bone density due to stress shielding or bone resorption as a tissue response to abrasion particles
- Tissue reactions, osteolysis and loosening of the implant due to metal corrosion or accumulation of abrasive particles or loose cement
- Corrosion with local tissue reactions or pain
- Aseptic loosening
- Noise development ("squeaking")
- Corrosion and fretting
- Local tissue reactions and hypersensitivity
- Dissociation of modular components
- Temporary or permanent nerve damage due to pressure or haematoma
- Wound haematoma and delayed wound healing
- Vascular disorders, including vein thrombosis, pulmonary embolism and heart failure
- Heterotopic ossification
- Nerve damage as a result of surgical trauma

With ceramic components, a risk of fracture can never be ruled out entirely. The following factors can increase this risk:

- Obesity or pre-obesity
- Alcohol or drug abuse
- Physical activities associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads (e.g. hard physical labour, certain types of sport)

2.6 System Compatibility

WARNING

Combination of BioBall™ Adapters with hip stems with neck insertion system

- Risk of injury due to premature implant failure!
- BioBall™ Adapters must not be combined with hip stems that use a neck insertion system.

2.6.1 Primary Surgery

Insofar as the BioBall™ Adapter label does not indicate otherwise, BioBall™ Adapters may be used in combination with either metal or ceramic BioBall™ heads. In primary surgery, BioBall™ Adapter 12/14 is only intended for use together with the Merete-brand hip stems approved for that purpose.

2.6.2 Revision Surgery

Insofar as the BioBall™ Adapter label does not indicate otherwise, BioBall™ Adapters may be used in combination with either metal or ceramic BioBall™ heads (possible combination see 2.6.3)! Surgeons wishing to perform revisions using BioBall™ Adapters with hip stems from other manufacturers must check taper (adapter-stem) compatibility prior to the operation. If using a 12/14 taper in such cases, this must adhere to the applicable CeramTec BIOLOX® specifications. The taper may not display any kind of shape-altering damage, severe abrasion/material loss, or deep scratches/burrs or similar surface defects. Use the BioBall™ AdapterSelector™ to check the taper geometry. If desired, Merete can provide information regarding suitable tapers. No biomechanical testing information is available on the use of BioBall™ Adapters with hip stems from other manufacturers. Consequently, only manufacturer-approved extensions may be used.

NOTE

BioBall™ DELTA™ ceramic heads must not be used with BioBall™ special adapters (MS 8/10, MSZI, MSV4 Offset 2XL and 3XL, MST1, MSBG, MSPC, MSSR and MSSY).

2.6.3 BioBall™ Adapter Neck Length Sizes

BioBall™ Adapter	Standard	Offset	BioBall™ Head
12/14	S-5XL	M-5XL	Ceramic or Metal
14/16	M-5XL	2XL-5XL	
MS 10/12	S-3XL	M-3XL	
MSV4	M-3XL	M-XL	
MS 8/10	S-2XL	M-2XL	Metal only
MSZI	S-3XL	–	
MST1	M-3XL	M-3XL	
MSV4	–	2XL-3XL	
MSBG	M-2XL	–	
MSPC	M-L	–	
MSSR	M-XL	–	
MSSY	S-XL	–	

BioBall™ Adapters may be used with stems whose tapers are made of the following materials:

12/14, 14/16, MS 8/10, MS 10/12, MSZI, MSV4, MST1, MSBG, MSPC and MSSY:

- TiAl6V4 ELI
- TiAl6Nb7
- CoCr alloys
- Stainless steel

MSSR:

- TiAl6V4 ELI
- TiAl6Nb7

Based on their dimensions, sliding pairs may be formed using only the following combinations of materials:

- BioBall™ metal heads may only be combined with UHMWPE inlays or cups.
- BioBall DELTA™ heads may only be combined with BIOLOX® delta inlays, or with UHMWPE inlays or cups.

Implant Materials

Materials used in manufacturing implants are indicated on the label.

BioBall™ heads may be made of the following materials:

- Metal head: Vivium™ (High Nitrogen Stainless Steel)
- Ceramic head: BIOLOX® delta

BioBall™ Adapters are made of:

- Ti6Al4V ELI alloy

Acceptable combinations may be made with cups or inlays made of the following materials:

- BIOLOX® delta
- UHMWPE

Manufacturers' labels indicate the materials used for the production of a particular inlay or cup.

Additional information on the chemical composition and mechanical properties of the materials used is available from Merete on request.

3. Use

3.1 General Instructions

WARNING

Use of damaged or defective implants

- Risk of injury due to premature implant failure!
- Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

WARNING

Use of damaged or defective instruments

- Risk of injury due to premature implant failure!
- Instruments with identifiable damage may not be used.

WARNING

Use of implant/instrument contrary to intended use

- Damage to/destruction of instrument/implant and injury to patient!
- Ensure correct handling of implant/instrument. Do not misuse.

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure (e.g. implant loosening, fretting or corrosion)!
- BioBall™ Adapters may only be combined with stem tapers after taper specifications have been clearly identified and matched.

WARNING

Use of implants which have been previously used

- Risk of injury due to premature implant failure!
 - Risk of sepsis!
- Implants are only approved for single use, not repeated use.

WARNING

Foreign bodies (e.g. cement residues, tissue, bones) between implant components

- Risk of injury due to implant failure!
- Thoroughly clean any foreign bodies from implant components.

WARNING

Use of implants contrary to intended purpose

- Risk of injury due to implant failure!
- Implants must only be used in accordance with intended purpose.

Prior to using Merete products, surgeons and assisting staff must study in detail the safety information in this product information sheet as well as the product-specific guidelines listed in the surgical technique. The relevant documentation is available from Merete on request. Surgeons must also be aware of any remaining risks associated with the products he or she intends to use, and must inform patients of such risks in advance. Implant operations must only be performed by surgeons who are not only qualified to carry out such operations, but also have extensive verified knowledge and experience in this field. The surgeon bears all responsibility for adverse effects or complications arising from misdiagnosis, improper surgical technique, incorrect implant selection or handling, or failure to observe the safety instructions provided in this product information sheet. Neither the manufacturers nor authorised Merete product representatives may be held liable in such cases. Before operating, study the techniques outlined in the surgical technique carefully. Training in the correct handling of implants and instruments by an

authorised Merete representative is essential.

The implantation of a joint endoprosthesis should only be considered if, after careful evaluation, this represents the best solution compared with other possible nonsurgical or joint salvage alternatives. The substitution of a severely altered, diseased joint through prosthesis can reduce pain and allow a good mobility and weight bearing. An implanted artificial joint prosthesis can never reach the full operability of a natural, healthy joint. Various factors may, in the longer term, lead to the loosening or wear of the even stably implanted artificial joint. The revision surgery may therefore be required, however, it is associated with significant health risk. Severe morphological changes in the supporting bone structure can affect implant performance, as can wear and tear on the mobile components. In extremely rare cases, especially if there is trauma or overuse, such changes may cause a component to fail. Patients should always be made aware of these potential risks so that they can take appropriate precautions to reduce overuse and mechanical stress on the implant.

3.2 Instructions on Patient Information

Surgeons must make their patients aware of information mentioned in Section 2, i.e. any factors that could hinder the success of the operation as well as potential complications that may occur as a result of specific indications. There is a risk of the implant failing prematurely, for example due to implant luxation, dislocation, or loosening caused by changes in load bearing behavior, fatigue or fracture of the cement bed, and/or tissue reactions to the implant and its abrasion particles; see also 2.5. Further risks exist in connection with the general surgical procedures to be performed, as well as risks of an operation under anesthesia. With ceramic components, a risk of fracture can never be ruled out entirely. This risk must be explained to the patient. Patients also need to be informed of any measures they can take to help prevent potential consequences of these factors. Furthermore, the patient has to be informed about alternative methods of treatment and other systems with adequate intended uses. The operating surgeon should keep a written record of all information given to the patient.

3.3 Features

MRI Safety Information/Indications for Use

Non-clinical testing has demonstrated that the Merete Hip Implant System (consisting of cemented or non-cemented hip stem, taper adapter, metal or ceramic head ball, inlay and cup from the materials TiAl6V4 ELI (ISO 5832-3), Vivium® (ISO 5832-9), CoCrMo (ISO 5832-4/5832-12), BIOLOX® delta ceramic (ISO 6474-2), PE/XPE (ISO 5834-2)) is MR conditional.

A patient with the entire assembled Merete Hip Implant System can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla.
- Maximum spatial gradient field of 3,000 Gauss/cm (30 T/m).
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) at 1.5 Tesla or 3.0 Tesla of 1 W/kg for 15 minutes of scanning. Under the scan conditions defined above, the Merete Hip Implant System is expected to produce a maximum temperature rise of less than 6°C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Merete Hip Implant System extends at least 1 cm and up to approximately 8 cm from the device and exhibits geometric distortion in the image when imaged with a gradient echo pulse sequence or a fast-spin echo pulse sequence and a 1.5 Tesla MRI system or a 3.0 Tesla MRI system.

WARNING

Damage to taper connection

- Risk of implant failure!
- Ensure careful implantation.
- Do not use damaged implants.

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure (e.g. implant loosening, fretting or corrosion)!
- BioBall™ Adapters may only be combined with stem tapers after taper specifications have been clearly identified and matched.

WARNING

Foreign bodies (e.g. cement residues, tissue, bones) between implant components

- Risk of injury due to failure of implant!
- Thoroughly clean any foreign bodies from implant components.

WARNING

Breaking of ceramic components

- Risk of injury due to implant failure!
- When performing revision surgery following breakage of a ceramic component, do not use metal heads.
- Replacement component must also be a ceramic head.

The BioBall™ System is made up of a titanium adapter which may only be used with a BioBall™ head. The primary indication is the revision of existing prosthetic stems. The BioBall™ Adapter 12/14 can also be used during the primary operation for purposes of intraoperative correction with approved stems of the Merete GmbH. In general, it is possible to vary offset, neck length, lateralisation and anteversion. BioBall™ Adapters change hip stem offset, thereby changing the biomechanics of the joint. As a result of these biomechanical alterations, excessive strain on the implant may cause it to loosen, break, or wear down, thus shortening its lifespan. BioBall™ heads may only be combined with BioBall™ Adapters. Under no circumstances may BioBall™ implants be reused. When performing revision surgery following breakage of a ceramic component, do not use metal heads – replacement component must also be a ceramic head. The individual components of the implant system are assembled with secure taper connections. Only an exact match of the taper dimensions makes this possible. Taper sizes are indicated on product labels – and, where possible, on the implants themselves – in the form of two numbers and an angle (e.g., 12/14 5°42' or 14/16 – 6°0'). Heads may only be attached to new, unused BioBall™ Adapter tapers. BioBall™ Adapters may only be combined with stem tapers after taper specifications have been clearly identified and matched (see 2.6.3). It is particularly important to keep contact surfaces clean and dry during assembly (see 3.5)! Only combinations using Merete-approved bearing couples (see 2.6.3) are permitted. Using large heads may produce larger quantities of abrasion particles, depending on the materials and sizes of the mating components used. This can negatively impact implant durability. For best results, it is important to choose mating component materials and head sizes with care, tailoring the selection to the individual patient's profile.

Additional information (surgical technique) is available from Merete on request. There are no known interactions with any medications.

3.4 Pre-Operative Planning

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure (e.g. implant loosening, fretting or corrosion)
- BioBall™ Adapters may only be combined with stem tapers after taper specifications have been clearly identified and matched.

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Implantation of trial implants

- Risk of injury due to failure of trial implant!
- Only use trial implants in order to select a suitable permanent implant.
- Trial implants are not suitable for permanent implantation.

Operation planning should be done on the basis of in-depth evaluations of patient X-rays, which provide the information necessary for determining the appropriate prosthesis type, size, and possible combinations. If desired, X-ray templates for pre-operative planning are available from Merete. It is also important to carry out pre-operative tests on the patient to rule out allergic reaction to the implant materials. Use trial components during the operation in order to check correct implant size and positioning. The components planned for implantation must be on hand in all available sizes. It is absolutely essential to check the implant label (type, item number, material and size) against the information on the package before insertion. Use the additional enclosed labels for documentation in surgery reports and for re-ordering. Check that all implant surfaces are free of soiling.

Important: Always check the sterilisation expiration date!

When combining with other Merete components, observe the instructions in the relevant user instructions, instrumentation guides and surgery videos, if any. In case of doubt, contact Merete GmbH.

3.5 Implant Handling

3.5.1 Instructions for Handling

WARNING

Use of implants contrary to intended purpose

- Risk of injury due to implant failure!
- Implants must only be used in accordance with intended purpose.

WARNING

Use of implants which have been previously used

- Risk of injury due to premature failure of implant!
- Risk of sepsis!
- Implants are only permitted for single use, not for repeated use.

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Combination with over-long heads

- Risk of injury due to implant failure!
- Impaired component safety due to higher lever forces.

WARNING

Use of damaged or defective implants

- Risk of injury due to premature implant fracture!
- Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

WARNING

Damage to head

- Risk of implant failure!
- Never strike the head or the adapter directly with a hammer.
- It is advisable to secure the head in place with light hammer blows in an axial direction on the plastic impactor.

WARNING

Damage to taper connection

- Risk of implant failure!
- Ensure careful implantation.
- Do not use damaged implants.

WARNING

Foreign bodies in the taper connection

- Risk of implant failure!
- Thoroughly clean all foreign bodies from the taper connection.

WARNING

Breaking of ceramic components

- Risk of injury due to implant fracture!
- When performing revision surgery following breakage of a ceramic component, do not use metal heads.
- Replacement component must also be a ceramic head.

WARNING

Use of instruments with electrical energy

- Risk of injury due to implant failure!!
- Do not damage the surfaces of the implants under any circumstances.

Implants must be inserted with the appropriate Merete instruments. The use of Merete instruments for any other purpose is not permitted.

Do not re-use an implant that has been removed under any circumstances. Do not use any components that have been damaged during handling under any circumstances. When inserting and repositioning implants, the operating surgeon must ensure that implant surfaces have not been scratched or dented in any way. Even a tiny scratch can significantly reduce the life of an implant. Silicone covers for intraoperative taper protection are available from Merete upon request. All components must be checked intraoperatively for function. Do not use stems with damaged or deformed taper regions, and do not use stems/acetabular components with non-standard configuration or geometry. In revision surgeries, BioBall™ heads must only be used in combination with unused BioBall™ Adapters. Adhere to the instructions in 3.4 "Pre-Operative Planning", and 2.6 "System Compatibility", when working with previously used tapers. Use the BioBall™ AdapterSelector™ to check the prosthesis taper and the compatibility to the BioBall™ Adapter. The wound must be thoroughly cleaned before closure. In particular, cement residue and bone splinters can significantly impair the bearing performance and thus lead to premature wear of the bearing surfaces.

Please observe the following when implanting BioBall™ components:

- Rinse and dry the stem taper to ensure that all foreign bodies (including bone fragments, soft tissues, bone cement and other substances) are completely removed
- Before positioning BioBall™ components, check all components as well as the stem taper for damage, deformation, wear or contamination
- Apply the BioBall™ Adapter to the prepared stem taper using axial pressure while simultaneously turning it slightly wards the right. Employ axial pressure only at the selected angle (CCD or anteversion). Then place the head under axial pressure on the BioBall™ Adapter. After that, check both head and adapter to ensure that they are set correctly into place.

The intended purpose of the product may lead to limitations in technical range of motion (ROM) arising from the following factors:

- Distance from stem shoulder to centre of ball
- Exterior contour of the stem neck (with S through XL adapters)
- CCD angle of stem in situ (medialisation and lateralisation limit ROM)
- Cup unfavourably skewed in acetabulum and/or usage of luxation-reducing cup systems
- Cup geometry

Never use a ceramic head which has fallen onto a hard surface or otherwise been damaged.

When performing revision surgery following breakage of a ceramic component, do not use metal heads – replacement components must also be made of ceramic.

Implantation in children:

- Complete bone growth is a prerequisite for the application of the BioBall™ System.

- Use the BioBall™ System for children in exceptional cases. The decision is up to the user.
- Adjust the size and type of the implant as best as possible to the age, height, weight and bone development of the child.
- Due to growth, a subsequent operation for removal and/or revision is more likely in children.
- Perform follow-up checks at shorter intervals.

Training in the correct handling of implants by an authorised Merete representative is essential.

3.5.2 Sterile Implants

WARNING

Risk of infection due to non-sterile implants!

- Do not use implants whose packaging is damaged.
- Do not use implants whose expiry date has passed.

WARNING

Use of soiled implants

- Risk of sepsis!
- Use only implants without identifiable soiling.
- Handle implants only with sterile surgical gloves.

WARNING

Use of damaged or defective implants

- Risk of injury due to premature implant fracture!
- Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

NOTE

Observe symbol on packaging:
"Do not re-use".

Implants which are delivered sterile are clearly labelled "STERILE". Sterile implants have been sterilised with 25 to 42 kGy (2.5 to 4.2 Mrad) gamma rays. All sterile implants should be stored unopened in their original packaging in a cool – but frost-free – dark, dry place until they are to be used. Before using any implants, check the sterilisation expiration date on the product label, and check the protective packaging for damage. The red sterilization dot on the package serves as an indicator that the product is sterile. Do not use products in damaged packaging. Observe the rules of asepsis when removing products from their protective packaging.

3.5.3 Resterilisation

WARNING

Resterilisation of implants

- Risk of injury due to premature implant fracture caused by adverse material changes!
- Implants delivered sterile by Merete GmbH must not be resterilised and/or repacked.
- Products whose expiry date has passed may be returned to Merete GmbH.

Merete GmbH wishes to provide patients with cutting-edge technology products at all times. As such, we do not offer product resterilisation.

3.5.4 Storage and Handling

WARNING

Risk of infection due to non-sterile implants!

- Do not use implants whose packaging is damaged.
- Do not use implants whose expiry date has passed.

WARNING

Use of soiled implants

- Risk of sepsis!
- Use only implants without identifiable soiling.
- Handle implants only with sterile surgical gloves.

WARNING

Use of damaged or defective implants

- Risk of injury due to premature implant fracture!
- Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

Store implants in their unopened original packaging at room temperature or below (but frost-free) in a clean and dry environment. Check the sterilisation expiration date on the product label and ensure that the protective packaging is intact before using any implants. Do not use products in damaged packaging. Observe the rules of asepsis when removing products from their protective packaging. Implants must be handled with the greatest of care. Even the slightest surface defects increase wear and tear, which could lead to complications. Therefore, do not mark the prosthesis surfaces in any way and avoid touching them with metallic or hard objects (e.g. instruments). Any mechanical manipulation or alteration of implants negates their approval for use and is not permitted. Do not implant any components of the prosthesis that have been handled incorrectly or show signs of damage (e.g. scratches – see also 3.5.2 “Sterile Implants”).

3.6 Handling of Instruments

3.6.1 Cleaning, disinfection and sterilisation instructions for reusable surgical instruments

Procedure

- Cleaning
- Disinfection
- Sterilisation with hot steam (DIN EN ISO 17665-1)

3.6.2 Warnings

Instruments which are supplied non-sterile are clearly labelled "NON-STERILE". Non-sterile instruments must be cleaned, disinfected (if applicable), and sterilised prior to use. Only qualified personnel may prepare instruments for use. Use only approved cleaners and disinfectants (RKI, DGHM/VHA, FDA) (pH < 12 for metal instruments). Plastic instruments and instruments with plastic components must not be sterilised using dry heat under any circumstances.

NOTE	Sterilisation of instruments supplied non-sterile If Merete products are sterilised by the user, this must be noted in the surgical report. All relevant labels and user instructions must be retained. → Observe current RKI guidelines. → Observe the standard preparation instructions provided.
-------------	---

3.6.3 Limitations on Reprocessing

Other manufacturing information and recommendations must be followed.

3.6.4 Location of Use

It is advisable to process instruments for re-use as soon as possible after using them. Remove surface soiling with a disposable/paper towel. Immediately after use, the instrument may be laid in a disinfectant solution or in hot water (80°C), in order to reduce the risk of infection for the user and to facilitate cleaning.

3.6.5 Preparation for Cleaning/Pre-Cleaning

Dismantle all instruments as far as possible. Take care not to lose any small components or screws. Pre-cleaning of instruments: completely immerse the instrument in an enzymatic or alkaline cleaning solution (pH ≤ 12) and soak for 10 minutes. Clean the instrument with a soft plastic brush. Then rinse the instrument with running water for at least 1 minute.

3.6.6 Automatic Cleaning and Disinfection

Mechanised cleaning is preferable to manual cleaning, if this is an option. Manual pre-cleaning improves automatic cleaning performance. The machine should offer a suitable thermal disinfection programme: A0 value > 3000, or at least 10 minutes at 93°C in older machines. Alternatively, if using a chemical disinfection method, bear in mind the risk of residue being left on the instruments. When choosing a cleanser, make sure that it is compatible with instrument materials. Follow manufacturer instructions when loading cleaning machines. Place instruments in such a way so as to allow complete, thorough rinsing of all ducts and cavities. Use deionised water for the final rinse. Be sure to allow sufficient drying time. Immediately after the programme has finished, remove instruments from the machine and, if necessary, dry them with a soft, absorbent, lint-free cloth.

3.6.7 Manual Cleaning

Pre-clean instruments manually to remove heavy surface contamination. Never use abrasive cleaning agents or metal brushes. Lay the instrument into the cleaning solution, following manufacturer's instructions as regards concentration, soaking time, and compatibility with the instrument materials. It must be ensured that the instrument is fully covered in cleaning solution and that all cavities, lumens and openings can be vented. Clean lumens and drill holes using appropriate brushes. After cleansing, rinse using deionised water, and dry thoroughly. For subsequent ultrasonic cleaning, make sure that the ultrasonic bath is preheated according to device manufacturers or cleaning agent manufacturer's instructions. When loading the bath, make sure that the cleaning solution completely covers the instruments and that all cavities, lumens and openings are fully vented. Cleaning of the instruments proceeds at 35-40 kHz for 5 minutes. After UV cleansing has finished, rinse instruments thoroughly with clear, running water, making sure to flush out cavities, lumens and openings wherever applicable.

3.6.8 Instrument Maintenance

Allow instruments to cool down to room temperature. Oil moving parts lightly with sterilisable, steam-penetrable surgical lubricating oil.

3.6.9 Checking Functionality

WARNING

Use of damaged or defective instruments

- Risk of injury due to premature implant fracture!
- Instruments with identifiable damage may not be used.

WARNING

Use of soiled instruments

- Risk of sepsis!
- Use only instruments without identifiable soiling.
- Handle instruments only with sterile surgical gloves.

After each cleaning/disinfection, the instruments must be inspected for cleanliness, functionality and damage (e.g., bent, broken, worn or missing parts). Damaged instruments must be discarded and may not be used again. The instrument trays must be checked for completeness on the basis of the attached tray inserts. Training in the correct handling of implants and instruments by an authorised Merete representative is essential.

3.6.10 Cleaning and Disinfecting Empty Trays

Clean and disinfect empty trays using the same procedure, and under the same conditions, as for instruments. Be sure that the tray is completely dry before inserting instruments.

3.6.11 Packing

Prior to steam sterilisation, cleaned, disinfected instruments should be packed in suitable containers or sterilisation packs (DIN EN ISO 11607-1).

3.6.12 Sterilisation

Instrument sterilisation is to be performed using an EN ISO 17665-validated procedure. We recommend a fractionated vacuum method when using a steam steriliser. The maximum allowable sterilisation temperature is 134°C (273°F) plus tolerances as per DIN EN ISO 17665-1. Sterilisation time (exposure time at the sterilisation temperature) should be at 5 minutes at 134°C (273°F). The recommended pressure level is 3 bar.

Valid in Switzerland only: We recommend sterilising reusable, invasive instruments for 18 minutes at 134°C in saturated, confined steam!

3.6.13 Storage

After sterilisation, the instruments must be kept in their sterilisation packaging and stored in a dry, dust-free place.

3.6.14 Processing Instructions in accordance with DIN EN ISO 17664

The processor is responsible for ensuring that the processing proce-

cedure actually employed (using the materials, equipment, and personnel available in the processing facilities) achieves the desired results. Normally, this means that the procedure must be validated and subject to routine monitoring. Likewise, any deviation from the instructions provided should be carefully evaluated by the processor to determine its effectiveness and possible negative consequences. As regards the processing procedure, we would also like to point out the vital importance of following the Robert Koch Institute guidelines (published in the 2012, 55 issue of the German Federal Health Gazette) as well as any relevant national regulations.

3.6.15 Reusable Instruments

Handle instruments with care in order to ensure their long-term functionality. Insofar as Merete has not specified otherwise, instruments may be reused indefinitely, provided that they are fully functional (see all points in 3.6). Instruments with a defined limited period of usage must be labelled as follows: e.g. Q1/2022. On the expiry of the defined quarter, the instrument must be inspected and approved by Merete GmbH to verify its continued functionality. Operators are to use instruments in accordance with RKI guidelines and MPBetreibV regulations.

4. Explanation

For explanation purposes, always request instrument sets from Merete GmbH.

There are no special requirements for explanation of components.

5. Information Updates

If more than two years have elapsed between the issue/revision date and the consultation date, please contact your authorised Merete representative for updated product information, or contact Merete GmbH at www.merete.de or +49 (0)30 77 99 80-0.

Subject to technical change without notice.

Further information available from Merete GmbH or your authorised Merete representative.

New clinical experience may result in changes to certain procedures.
Merete GmbH recommends regular attendance at refresher courses.

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USA English

IMPORTANT PRODUCT INFORMATION FOR THE USA IntraBlock™ BioBall™ Hip System (IBS)

Symbols

0482

CE mark

Keep dry

Manufacturer

Temperature limit

Date of
manufacture

Do not reuse

Use-by date

Consult instructions
for use

Batch code

Caution, consult
accompanying
documents

Catalogue
number

Sterilized using
ethylene oxide

Quantity

Sterilized using
irradiation

Caution: Federal
law (USA) restricts this
device to sale by, or
on order of a physician

Do not resterilize

Non-sterile

Do not use if
package is
damaged

Important Product Information (Package Insert)

Caution

Federal Law restricts this device to sale by or on the order of a physician.

Description

The IntraBlock™ BioBall™ Hip System (IBS) is composed of femoral stems, heads, taper adapter sleeves and acetabulum cups. The System is designed for uncemented use for either primary or revision hip arthroplasty.

The IBS stem is straight, collarless and flat with a tapered design. The stem is manufactured from titanium alloy (Ti-6Al-4V ELI) conforming to ASTM F-136 with a proximal titanium plasma spray (TPS) coating. A lateralized version of the device is available. Both versions have a neck shaft angle of 135° and a 12/14 - 5°42' morse taper trunion (tolerance 5°37'-5°42'). The Ceramic Head is made of BioloX® delta¹ (Alumina Matrix Composite ISO 6474-2). The taper adapter sleeves are made from titanium alloy (Ti-6Al-4V ELI) with a 12/14 - 5° 42' morse taper trunion tolerance 5°37'-5°42' in the four various neck lengths. The cementless MultiCup Locking PressFit is manufactured from titanium alloy (Ti-6Al-4V ELI) conforming to ASTM F-136 with a titanium plasma spray (TPS) coating on the outer surface. It is used with a UHMWPE Inlay. Holes with threaded plugs give the option for an additional fixation with locking screws. The MultiCup is available in various sizes.

Intended Use

The IntraBlock™ BioBall™ Hip System (IBS) is intended for use in total hip arthroplasty. The IBS is indicated for:

- Patients suffering from severe hip pain and disability due to rheumatoid arthritis, osteoarthritis, traumatic arthritis
- Avascular necrosis of the femoral head
- Correction of functional deformity
- Treatment of non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement, unmanageable by other techniques
- Revision procedures where other treatment or devices have failed

The IBS is indicated for cementless use only.

Contraindications

- Acute or chronic infections, either local or systemic,
- Severe muscular, nervous or vascular disease which endanger the specific extremities,
- Defective bone structures which would impede adequate anchoring of the implant,
- Patients who are younger than 60 years whose joint disease is such, that good results are indicated by using another reconstructive procedure (e.g. osteotomy),
- All accompanying diseases which could endanger the functioning and the success of the implant,
- Stem taper situated in situ shows evidence of severe damage (noticeable damage as erosion, abrasion/loss of material, scratches/ridges or similar injuries or visible deformation) or cannot be clearly identified,
- Patients with mental or neurological disease conditions or patients who are not capable of following the necessary postoperative treatment instructions,

- unfavorable femur morphology,
- Revision with extensive metaphyseal bone defect, i.e. insufficient proximal rotational stability of the implant,

Warnings

- All implants are for single use only.
- Heads must not be sterilized on the hip stem.
- Autoclaved heads should not be cooled rapidly.
- **Important notice:** In revision surgery before using the BioBall™ system you must use the AdapterSelector™ to check the state of the prosthesis cone.
- BioBall™ heads must only be used in combination with unused BioBall™ Adapters.
- The stem cone and head bore should be dry and free of contamination.
- A ceramic head should not be implanted if the head or the cone of the stem are possibly damaged (e.g., if the head is dropped on the floor or if the stem cone is scratched by an instrument or if the head and stem cone are attached then detached).

Cautions

The head should be placed on the stem gently while keeping the head and cone in alignment, then firmly attached by sharply hitting the head with a soft plastic hammer.

Factors which could impair the success of the operation:

- Allergies to the implanted material, especially metal (e.g. cobalt, chrome, nickel etc.) as well as polyethylene and bone cement,
- Adiposity or overweight of the patient,
- Local bone tumors,
- Osteoporosis or osteomalacia,
- Malformations, congenital hip luxations, significant axial malposition in the knee area,
- Systemic diseases or metabolic disturbances,
- Alcohol and drug abuse,
- Physical activity associated with strong shocks, in which the implant is exposed to blows and/or excessive strain (e.g. hard physical work, certain kinds of sports, etc.).

Adverse Effects

Possible adverse effects:

The below enlisted adverse effects are counted among the most common consequences after an implantation:

- Loosening of the implant resulting from changed conditions of the load transfer or wearing and fracture of the cement bed and/or tissue reaction of the implant and the associated abrasion products,
- Failure or subsidence of the implant in case of over-loading,
- Early and late infection,
- Dislocation, subluxation, insufficient range of motion, undesired shortening or lengthening of the extremity as a result of sub-optimal positioning of the implant,
- Bone fracture resulting from unilateral over-bearing or weakened bone substance,
- Reduction of bone density as a result of stress shielding or bone resorption as a tissue reaction to abrasion particles,

- Aseptic loosening,
- Noise emission ("squeaking"),
- Temporary or chronic neural damage resulting from pressure or hematoma,
- Wound hematoma and delayed wound healing,
- Vascular disease including venal thrombosis, pulmonary embolism and cardiac arrest,
- Heterotopic ossification,
- Nerve damage due to surgical trauma.

MRI Safety Information

The BioBall System have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the BioBall System in the MR environment is unknown. Scanning a patient who has this screws may result in patient injury.

Further information for use

Cleaning, disinfection, and sterilization instructions for reusable surgical instruments and non-sterile implants.

Process

- Cleaning
- Disinfection
- Sterilization with hot steam (DIN EN ISO 17665-1)

Warnings

Instruments that are supplied non-sterile are clearly marked with “NON STERILE” and must be cleaned, disinfected and sterilized prior to use. The Instruments may only be processed by qualified personnel. Only approved cleaning agents and disinfectants (FDA) are to be used ($\text{pH} \leq 12$ pH for Instruments made of metal). Instruments made from synthetic materials or ones that containing components made from synthetic materials are not to be sterilized using dry heat.

Restrictions regarding reprocessing

The presented cleaning processes have been validated. Other methods of cleaning may be suitable but need to be validated by the user of the device. Differing instructions and recommendations of the manufacturer are to be observed.

Point-of-use-processing

It is advisable to prepare Instruments for reuse as soon as possible after having used them. Macroscopic surface contamination can be difficult to remove by automated cleaning procedures. Prior to cleaning remove macroscopic contamination with a disposable towel/paper towel.

Keep Instruments moist after use to prevent contamination from drying. Instruments may be placed in a disinfectant solution or hot water 176°F (80°C) immediately after use in order to facilitate cleaning and to reduce risk of infection.

Preparation for cleaning

Dismantle all Instruments as far as possible. Ensure to keep all small components and screws. Pre-cleaning of instruments:

- Completely immerse the Instruments in an enzymatic or alkaline cleaning solution ($\text{pH} \leq 12$) and soak for 10 min
- Clean the Instruments with a surgical scrub brush
- Then rinse the instrument/implant for at least 1 min with deionized water

Automatic cleaning

Automatic cleaning is preferable to manual cleaning, if this is an option. The machine should offer a suitable thermal disinfection program. Minimum cycle steps and parameters are:

- Rinse 1 min with cold water < 109,4°F (< 43°C)
- Cleaning 5 min with cleaning agent (131 °F (55°C) or follow the manufacturer's instructions)
- 1 min neutralization with warm water
- 1 min rinse (note: final rinsing is to be carried out with deionized water)

Automatic disinfection

Choose a program for an A_0 value > 3000, or at least 10 minutes at 200°F (93°C) in older machines. Alternatively, if using a chemical disinfection method, bear in mind the risk of residue being left on the Instruments.

Manual cleaning

Begin by removing major surface contamination from the Instruments using a soft nylon brush or a soft, lint-free cloth, along with either clear running water or a cleaning solution. Never use abrasive cleaning agents or metal brushes. Place the Instruments in the cleaning solution, following manufacturer's instructions as regards concentration, soaking time, and compatibility with the Instruments materials. Ensure that the Instruments is completely submerged in the cleaning solution. Be sure to vent all cavities, lumens and openings. Clean lumens and drill holes using appropriate brushes. After cleansing, rinse using deionized water, and dry thoroughly. Subsequent ultrasonic cleaning. Be sure that the ultrasonic bath is preheated according to device manufacturer or cleaning-agent manufacturer instructions. When loading the bath, make sure that the cleaning solution completely covers the Instruments, and that all cavities, lumens and openings are fully vented. Clean Instruments at 35-40 kHz for five minutes. After ultrasonic cleansing has finished, rinse Instruments thoroughly with deionized water, making sure to flush out cavities, lumens and openings wherever applicable.

Instruments maintenance

Allow instruments to cool down to room temperature. Lubricate moving parts lightly with sterilizable, steam-penetrable surgical lubricating oil.

Checking functionality

After each cleaning / disinfection, inspect the Instruments for cleanliness, functionality, and damage (e.g., bent, broken, worn or missing parts). Never use damaged instruments or implants.

Cleaning and disinfecting empty trays

Clean and disinfect empty trays using the same procedure and under the same conditions as for Instruments. Be sure that the tray is completely dry prior to inserting Instruments.

Packing

Prior to steam sterilization, cleaned, disinfected Instruments should be inserted into suitable containers or sterilization packages (DIN EN ISO 11607-1).

Sterilization

BioBall™ implants are delivered in sterile packing, sterilized by Gamma Radiation:	
Method	Gamma Radiation
Radiation Dose	2.5 Mrad
Sterility Assurance Level	10 ⁻⁶
Sterility Validation Method	ISO 11137-2 / VD _{max} 25
Packaging	Double Peel PE Pouch
Pyrogenicity	Not labeled pyrogen free

Instruments are delivered non-sterile. Following sterilization method is recommended	
Method	Steam Sterilization
Cycle	Pre-Vacuum
Temperature	270°F [132°C]
Expose Time	4 minutes
Drying Time	20 – 30 minutes
Sterility Assurance Level (SAL)	10 ⁻⁶
Sterility Validation Method	Overkill method per ISO 17665-1:2006
Wrapping	FDA cleared wrap

The trays containing implants and instruments are only provided for packaging, not for sterilization.

Storage

After sterilization, the Instruments must be kept in their sterilization packaging and stored in a dry, dust-free place.

Preparation instructions in accordance with DIN EN ISO 17664

The preparer is responsible for ensuring that the preparation procedure actually employed (using the materials, equipment, and personnel available in the preparation facilities) achieves the desired results. Normally, this means that the procedure must be validated and subject to routine monitoring. Likewise, any deviation from the instructions provided should be carefully evaluated by the preparer to determine its effectiveness and possible negative consequences.

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Merete Technologies Inc.