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114800	Discovery Ulna Bearing Revision Kit (includes ulna bearing and two locking pins)	-
--------	--	---

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01-50-0991
Date: 05/07

BIOMET® ELBOW AND SHOULDER JOINT REPLACEMENT PROSTHESES

ATTENTION OPERATING SURGEON

DESCRIPTION

Biomet manufactures a variety of elbow and shoulder joint replacement prostheses intended for primary and revision joint arthroplasty for use in cemented applications.

Elbow joint replacement components include: humeral and ulnar components, and a hinge component. Components are available in a variety of surface finishes including bond coat (a thin layer of titanium plasma spray) and Interlok® finish. Small diameter cement plugs are available as specialty components.

Shoulder joint replacement components include humeral stems, humeral heads, and glenoid components. Components are available in a variety of designs and size ranges for both primary and revision applications. Specialty components include glenoid screws, centering sleeves, bipolar heads, and intercalary segments.

MATERIALS

Elbow

Humeral Stem	CoCrMo alloy or titanium alloy
Ulnar Stem	CoCrMo alloy or titanium alloy
Bearing Components	Ultra-High Molecular Weight Polyethylene (UHMWPE)
Axles	CoCrMo alloy
Surface Coating	Titanium alloy
Screws	Titanium alloy
Small Diameter Cement Plugs	UHMWPE

Shoulder

Humeral Stems	CoCrMo alloy or titanium alloy
Attachment Sleeves	CoCrMo alloy
Humeral Head	CoCrMo alloy
Glenoid Components	Ultra-High Molecular Weight Polyethylene (UHMWPE) / titanium alloy / 316 LVM stainless steel
Glenoid Screws	Titanium alloy
Centering Sleeves	Polymethylmethacrylate (PMMA)
Bipolar Heads	CoCrMo alloy / UHMWPE / titanium alloy
Intercalary Segments	CoCrMo alloy or titanium alloy
Porous Coating	Titanium alloy

INDICATIONS

1. Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis.
2. Rheumatoid arthritis.
3. Revision where other devices or treatments have failed.
4. Correction of functional deformity.
5. Treatment of acute or chronic fractures with humeral epicondyle (elbow) involvement or humeral head (shoulder), which are unmanageable using other treatment methods.
6. Oncology applications.

Patient selection factors to be considered include: 1) need to obtain pain relief and improve function, 2) ability and willingness of the patient to follow instructions, including control of weight and activity levels, 3) a good nutritional state of the patient, and 4) the patient must have reached full skeletal maturity.

CONTRAINDICATIONS

Absolute contraindications include: infection, sepsis, and osteomyelitis.

Relative contraindications include: 1) uncooperative patient or patient with neurologic disorders who are incapable of following directions, 2) osteoporosis, 3) metabolic disorders which may impair bone formation, 4) osteomalacia, 5) distant foci of infections which may spread to the implant site, and/or 6) rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram.

WARNINGS

Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components. Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure. Inadequate preclosure cleaning (removal of surgical debris) can lead to excessive wear. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Thoroughly clean and dry connecting segments, including taper, prior to attachment of components to avoid crevice corrosion and improper seating. Use clean gloves when handling implants. Laboratory testing indicates that implants subjected to body fluids, surgical debris or fatty tissue has lower adhesion strength to cement than implants handled with clean gloves. Do not modify implants. The surgeon is to be thoroughly familiar with the implants and instruments, prior to performing surgery.

1. Properly align and completely seat connecting components including tapers. Failure to properly align and completely seat the components together can lead to disassociation. Thoroughly clean and dry all connectors, including tapers prior to attachment of modular components to avoid crevice corrosion and improper seating.
2. Disassociation of the humeral head component from the humeral stem component has been reported. Failure to properly align and completely seat the components together can lead to disassociation. Thoroughly clean and dry tapers prior to attachment of modular head component to avoid crevice corrosion and improper seating.
3. Dislocation of the bipolar shoulder component has been reported. Closed reduction should be attempted with caution to prevent disassociation of the bipolar component. Do not use excessive force during closed reduction. The bipolar component may be impinged against the glenoid.

4. Care is to be taken to assure complete support of all parts of the device embedded in bone cement to prevent stress concentrations that may lead to failure of the procedure. Complete preclosure cleaning and removal of bone cement debris, metallic debris and other surgical debris at the implant site is critical to minimize wear of the implant articular surfaces. Implant fracture due to cement failure has been reported.

Biomet® joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals they cannot be expected to withstand the activity levels and loads of normal healthy bone and joint tissue.

Accepted practices in postoperative care are important. Failure of the patient to follow postoperative care instructions involving rehabilitation can compromise the success of the procedure. The patient is to be advised of the limitations of the reconstruction and the need for protection of the implants from full load bearing until adequate fixation and healing have occurred. Excessive activity, trauma and excessive weight have been implicated with premature failure of the implant by loosening, fracture, and/or wear. Loosening of the implants can result in increased production of wear particles, as well as accelerate damage to bone, making successful revision surgery more difficult. The patient is to be made aware and warned of general surgical risks, possible adverse effects as listed, and to follow the instructions of the treating physician including follow-up visits.

PRECAUTIONS

Specialized instruments are designed for Biomet® joint replacement systems to aid in the proper implantation of the prosthetic components. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, excessive wear and device failure. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments that have experienced extensive use or excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Biomet recommends that all instruments be regularly inspected for wear and disfigurement.

Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient.

- Patient must avoid placing excessive loads on the implant.
- Patient must avoid lifting more than 5lbs with the operated arm after surgery.
- Patient must avoid putting full body weight on the operated arm when rising from a seated position.
- Patient must avoid sudden or strenuous pulling activities after surgery, as these can produce excessive stress on the operated arm.

POSSIBLE ADVERSE EFFECTS

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process. Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis or osteolysis may result in loosening of the implant.
2. Early or late postoperative infection and allergic reaction.
3. Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
4. Infection is a rather common problem in elbow procedures
5. Impairment due to injury of the ulnar nerve is a major concern in elbow procedures.
6. Loosening, migration, and/or fracture of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, and/or excessive activity.
7. Periaricular calcification or ossification, with or without impediment of joint mobility.
8. Inadequate range of motion due to improper selection or positioning of components.
9. Undesirable shortening or lengthening of limb.
10. Dislocation and subluxation due to inadequate fixation, improper positioning, trauma, excessive range of motion, and/or excessive activity. Muscle and fibrous tissue laxity can also contribute to these conditions.
11. Fatigue fracture of component can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, or excessive weight.
12. Fretting and crevice corrosion can occur at interfaces between components.
13. Wear and/or deformation of articulating surfaces.
14. Accelerated wear of glenoid articular cartilage.
15. Intraoperative or postoperative bone fracture and/or postoperative pain.
16. Axle or bearing components may disassociate causing the elbow to disarticulate.
17. Revision and post-traumatic patients are susceptible to higher wear rates if varus/valgus constraints are compromised.

STERILITY

Prosthetic components are sterilized by exposure to a minimum dose of 25 kGy of gamma radiation. Do not resterilize. Do not use any component from an opened or damaged package. Do not use implants after expiration date has passed.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Comments regarding this device can be directed to Attn: Regulatory Dept., Biomet Inc., P.O. Box 587, Warsaw, IN 46581 USA, Fax: 574-372-1683.

Biomet and Interlok are registered trademarks of Biomet Manufacturing Corp. in the United States.

Authorized Representative: Biomet U.K., Ltd.
Waterton Industrial Estate
Bridgend, South Wales
CF31 3XA, U.K.

DISCOVERY™ SYSTÈME DE COUDE

P.O. Box 587
56 East Bell Drive
Warsaw, Indiana 46581 USA

Date: 08/07

Biomet manufactures a variety of elbow joint replacement prostheses intended for primary and revision joint arthroplasty for use in cemented applications. Elbow joint replacement components include humeral and ulnar components, and in some instances, hinge components. Components are available in a variety of surface finishes including bond coat (a thin layer of titanium plasma spray), porous titanium plasma spray and Interlok™ finish.

Humeral stem	CoCrMo alloy or titanium alloy
Ulnar stem	CoCrMo alloy or titanium alloy
Bearing components	Ultra-High Molecular Weight Polyethylene (UHMWPE)
Axles	CoCrMo alloy
Connectors	CoCrMo alloy
Surface coating	Titanium alloy
Locking clips/screws	Titanium alloy

1. Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis.
2. Rheumatoid arthritis.
3. Revision where other devices or treatments have failed.
4. Correction of functional deformity.
5. Treatment of acute or chronic fractures with humeral epicondyle involvement, which are unmanageable using other treatment methods.

In addition to the above indications, BiAxial™ elbow components are also indicated for joints with both intact and limited soft tissue structure about the elbow.

Absolute contraindications include: infection, sepsis, and osteomyelitis.

Relative contraindications include: 1) uncooperative patient or patient with neurologic disorders who is incapable of following directions, 2) osteoporosis, 3) metabolic disorders which may impair bone formation, 4) osteomalacia, 5) distant foci of infections which may spread to the implant site, and/or 6) rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram.

Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components. Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure. Inadequate preclosure cleaning (removal of surgical debris) can lead to excessive wear. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Thoroughly clean and dry connecting segments, including taper, prior to attachment of components to avoid crevice corrosion and improper seating. Use clean gloves when handling implants. Laboratory testing indicates that implants subjected to body fluids, surgical debris or fatty tissue have lower adhesion strength to cement than implants handled with clean gloves. Do not modify implants. The surgeon is to be thoroughly familiar with the implants and instruments prior to performing surgery.

Patient selection factors to be considered include: 1) need to obtain pain relief and improve function, 2) ability and willingness of the patient to follow instructions, including control of weight and activity levels, 3) a good nutritional state of the patient, and 4) the patient must have reached full skeletal maturity.

Elbow joint replacement prostheses have not received FDA clearance for non-cemented application (USA).

1. Properly align and completely seat connecting components including tapers. Failure to properly align and completely seat the components together can lead to disassociation. Thoroughly clean and dry all connectors, including tapers prior to attachment of modular components to avoid crevice corrosion and improper seating.
2. Care is to be taken to assure complete support of all parts of the device embedded in bone cement to prevent stress concentrations that may lead to failure of the procedure. Complete preclosure cleaning and removal of bone cement debris, metallic debris and other surgical debris at the implant site is critical to minimize wear of the implant articular surfaces. Implant fracture due to cement failure has been reported.
3. BiAxial Elbow: Insertion of the axle clip must be performed properly. Complete seating of the clip using a new, unused clip is necessary to prevent disassociation. If an axle clip is removed for any reason, do not reuse.

Biomet® joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals, they cannot be expected to withstand the activity levels and loads of normal healthy bone and joint tissue.

Accepted practices in postoperative care are important. Failure of the patient to follow postoperative care instructions involving rehabilitation can compromise the success of the procedure. The patient is to be advised of the limitations of the reconstruction and the need for protection of the implants from full load bearing until adequate fixation and healing have occurred. Excessive activity, trauma and excessive weight have been implicated with premature failure of the implant by loosening, fracture, and/or wear. Loosening of the implants can result in increased production of wear particles, as well as accelerate damage to bone, making successful revision surgery more difficult. The patient is to be made aware and warned of general surgical risks, possible adverse effects as listed, and to follow the instructions of the treating physician including follow-up visits.

Specialized instruments are designed for Biomet® joint replacement systems to aid in the proper implantation of the prosthetic components. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, excessive wear and device failure. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments

are subject to wear with normal usage. Instruments that have experienced extensive use or excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Biomet recommends that all instruments be regularly inspected for wear and disfigurement.

Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient.

- Patient must avoid placing excessive loads on the implant.
- Patient must avoid lifting more than 5lbs with the operated arm after surgery.
- Patient must avoid putting full body weight on the operated arm when rising from a seated position.
- Patient must avoid sudden or strenuous pulling activities after surgery, as these can produce excessive stress on the operated arm.

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process. Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis or osteolysis may be a result of loosening of the implant.
2. Early or late postoperative infection and allergic reaction.
3. Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
4. Infection is a rather common problem in elbow procedures.
5. Impairment due to injury of the ulnar nerve is a major concern in elbow procedures.
6. Loosening, migration, and/or fracture of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, and/or excessive activity.
7. Periarticular calcification or ossification, with or without impediment of joint mobility.
8. Inadequate range of motion due to improper selection or positioning of components.
9. Undesirable shortening or lengthening of limb.
10. Dislocation and subluxation due to inadequate fixation, improper positioning, trauma, excessive range of motion, and/or excessive activity. Muscle and fibrous tissue laxity can also contribute to these conditions.
11. Fatigue fracture of component can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, or excessive weight.
12. Fretting and crevice corrosion can occur at interfaces between components.
13. Wear and/or deformation of articulating surfaces.
14. Intraoperative or postoperative bone fracture and/or postoperative pain.
15. Axle or bearing components may disassociate causing the elbow to disarticulate.
16. Revision and post-traumatic patients are susceptible to higher wear rates if varus/valgus constraints are compromised.

Prosthetic components are sterilized by exposure to a minimum dose of 25 kGy of gamma radiation. Do not resterilize. Do not use any component from an opened or damaged package. Do not use implants after expiration date has passed.

Caution: Federal law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

Comments regarding this device can be directed to Attn: Regulatory Dept., Biomet, Inc., P.O. Box 587, Warsaw, IN 46581 USA, Fax: 574-372-3968.

Biomet® and all other trademarks herein are the property of Biomet, Inc. or its subsidiaries.

Authorized Representative: Biomet U.K., Ltd.
Waterton Industrial Estate
Bridgend, South Wales
CF31 3XA, U.K.

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Lima Orthopaedics South Africa
14 6Th Avenue
Maraisburg - Johannesburg - Gauteng
South Africa 1709

Lima SM spa
Strada Borrana 38
47899 Serravalle, Republic of San Marino
Tel.: +378 0549 961911
Fax: +378 0549 961912
E-mail: info@lima-sm.com

Discovery™ Trademark is licensed to Limacorporate spa for use only in the European Economic Area (E.E.A.) and Switzerland.

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, see the package insert or Biomet's website.

limacorporate.com

This publication is not intended for distribution in the U.S.
Ce support est destiné aux professionnels de santé. Lire attentivement les instructions et notices

The DISCOVERY™ Elbow System is manufactured by Biomet.
The DISCOVERY™ Elbow System is distributed in the European Economic Area (E.E.A.) and Switzerland by Limacorporate spa

B.1401.2F.000.1

011600