# Discovery Elbow System

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## INSTRUCTIONS FOR USE – DISCOVERY ELBOW SYSTEM

### Intended for cemented fixation

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g.: product literature, surgical technique).

### 1. PRODUCT INFORMATION

The main goal of an elbow joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and on the post-op rehabilitation. Discovery elbow system is a semi-constrained elbow replacement system intended for primary and revision joint arthroplasty for use in cemented applications. Discovery elbow system includes humeral stem, ulnar stems, humeral condyle components and ulnar bearing kit.

#### 1.1. MATERIALS

| Components                                                                                          | Material                                         |
|-----------------------------------------------------------------------------------------------------|--------------------------------------------------|
| <b>Humeral stem</b>                                                                                 | Ti6Al4V<br>Ti6Al4V + Bond Coat                   |
| <b>Ulnar stem</b>                                                                                   | UHMWPE + Ti6Al4V<br>UHMWPE + Ti6Al4V + Bond Coat |
| <b>Humeral Condyle Set</b>                                                                          | CoCrMo + Ti6Al4V                                 |
| <b>Condyle Lock Screw</b>                                                                           | Ti6Al4V                                          |
| <b>Ulnar Bearing Kit</b>                                                                            | UHMWPE + Ti6Al4V                                 |
| <b>Material standards</b>                                                                           |                                                  |
| <b>CoCrMo</b> (ISO 5832-12, ASTM F1537) – <b>UHMWPE</b> (ASTM F648-98) – <b>Ti6Al4V</b> (ASTM F136) |                                                  |

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon.

Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone and applying the bone cement for prosthesis fixation.

## **1.2. HANDLING AND STORAGE**

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50° C / 32-122° F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. Components removed from the package should not be used if they are dropped or suffer other accidental impacts. Devices should not be modified in any way. The device's code and lot number should be recorded in the patient's case history by using the labels included in the components packaging.

The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

## **2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS**

### **2.1. INDICATIONS**

Discovery Elbow system is indicated for patients with: • Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis, • Rheumatoid arthritis, • Revision where other devices or treatments have failed, • Correction of functional deformity, • Treatment of acute or chronic fractures with humeral epicondyle involvement, which are unmanageable using other treatment methods.

All Discovery Elbow components are intended for use with bone cement.

### **2.2. CONTRAINDICATIONS**

Absolute contraindications include: • infection, • sepsis, • osteomyelitis.

Relative contraindications include: • uncooperative patient or patient with neurologic disorders who is incapable of following directions, • osteoporosis, • metabolic disorders which may impair bone formation, • osteomalacia, • distant foci of infections which may spread to the implant site, • rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram.

## 2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • obesity<sup>1</sup>; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • muscle deficiencies; • refusal to modify postoperative physical activities; • patient history of infections; • systemic diseases and metabolic disorders; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • patient's resistance to disease generally weakened (HIV, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • use our combinations with products, prosthesis or instruments of another manufacturer; • patient smoking.

## 3. WARNINGS

### 3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques. Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of devices required based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

### COMBINATIONS ALLOWED/NOT ALLOWED

|                                                                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| X-small humeral components should only be used with x-small ulnar components.<br>X-small humeral condyles should only be used with x-small humeral components.<br>X-small bearing kits should only be used with x-small ulnar components. |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future. The possible impact of the factors mentioned in section 2 and 3.4. should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors. Patient selection factors to be considered include: 1) need to obtain pain relief and improve function, 2) ability and willingness of the patient to follow instructions, including control of weight and activity levels, 3) a good

---

1. According to the definition of the World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 30

nutritional state of the patient, and 4) the patient must have reached full skeletal maturity.

**Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.** After use, the device may be a potential biohazard. Reuse of devices labeled for single-use may result in product contamination, patient infection and/or failure of the device to perform as intended. Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

### 3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants be available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may negatively affect system performance and survival rate of the implant. Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure. Inadequate preclosure cleaning (removal of surgical debris) can lead to excessive wear. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Care is to be taken to assure complete support of all parts of the device embedded in bone cement to reduce the risk of stress concentrations that may lead to failure of the procedure. Complete preclosure cleaning and removal of bone cement debris, metallic debris and other surgical debris at the implant site is critical to minimize wear of the implant articular surfaces. Implant fracture due to cement failure has been reported.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labeled indications.

Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken to protect the surfaces involved in the coupling between components (tapers); articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique. Use clean gloves when handling implants. Laboratory testing indicates that implants subjected to body fluids, surgical debris or fatty tissues have lower adhesion strength to cement than implants handled with clean gloves.

### 3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect any changes in position of the implant or surrounding tissues.

The surgeon should make the patient aware of the limitations of limb function after elbow arthroplasty and that the reconstructed joint must be protected from full load bearing until adequate fixation and healing have occurred. Excessive physical activity or trauma to the replaced elbow joint can lead to premature failure of the elbow arthroplasty through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon: • avoid placing excessive loads on the implant; • avoid lifting more than 2,25kg (5 lbs) with the operated arm after surgery; • avoid repeated high weight lifting; • avoid putting full body weight on the operated arm when rising from seated position; • keep body weight under control; • avoid sudden peak loads (consequences of activities such as contact sports, playing tennis) or movements which can lead to sudden stops or twisting; • avoid positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

Patient smoking may result in delayed healing, non-healing and/or compromised stability in or around the placement site.

### 3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that can occur in elbow arthroplasty include: • Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process. Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis or osteolysis may be a result of loosening of the implant.

• Early or late postoperative infection and/or allergic reaction. • Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device. • Infection is a rather common problem in elbow procedures. • Impairment due to injury of the ulnar nerve is a major concern in elbow procedures. • Loosening, migration, and/or fracture of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, and/or excessive activity. • Periarticular calcification or ossification, with or without impediment of joint mobility. • Inadequate range of motion due to improper selection or positioning of components. • Undesirable shortening or lengthening of limb. • Dislocation and subluxation due to inadequate fixation, improper positioning, trauma, excessive range of motion, and/or excessive activity. Muscle and fibrous tissue laxity can also contribute to these conditions. • Fatigue fracture of component can occur as a result of loss

of fixation, strenuous activity, malalignment, trauma, non-union, or excessive weight. • Fretting and crevice corrosion can occur at interfaces between components. • Wear and/or deformation of articulating surfaces. • Intraoperative or postoperative bone fracture and/or postoperative pain. • Bearing and/or condyle components may disassociate causing the elbow to disarticulate. • Revision and post-traumatic patients are susceptible to higher wear rates if varus/valgus constraints are compromised.

#### **4. STERILITY**

##### **a. Implants**

All implantable components of the Discovery Elbow system are provided sterile with a Sterility Assurance Level (SAL) of  $10^{-6}$ . All components are sterilized by radiation. Do not use any component from a package that has been previously opened or appears to be damaged. Do not use implants after the expiration date printed on the label.

##### **b. Instruments**

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from [www.limacorporate.com](http://www.limacorporate.com) in the *Products section*). Users should validate their specific cleaning, disinfection and sterilization processes and equipments.

#### **5. MAGNETIC RESONANCE IMAGING (MRI)**

Discovery Elbow System components have not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. The risks associated with a passive implant in an MR environment have been evaluated for other prosthetic systems and are known to include heating, migration and image artifacts at or near the implant site.

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**Limacorporate S.p.A.**