Multigen Plus System

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INSTRUCTIONS FOR USE – MULTIGEN PLUS SYSTEM

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

Multigen Plus components are knee prosthesis components designed for long-term implantation inside the human body and must be exclusively utilized by orthopaedic surgeons having a good knowledge of the specific operative technique. The main goal of a knee joint prosthesis is to reproduce the articular anatomy reducing pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and on the post-op rehabilitation. Multigen Plus system is a modular system used to replace knee joint during surgery. Both primary design (CR and PS versions) and revision design (CCK and H versions) are available for Multigen Plus system.

- CR version is used when the posterior cruciate ligament (PCL) is intact or in which the ligament is repaired at the same time as the joint is replaced or, in absence of the posterior cruciate ligament, only with the proper ultra-congruent liner.

- PS version is used when the PCL is not preserved;

- CCK is a semi-constrained system used when the PCL is not preserved; in this case the post is higher and thicker and it is intended to be coupled with a box in the femoral component: this coupling ensures a greater stability to the joint compared to the PS version.

- H versions is a constrained (or hinged) system used only in cases of revisions in which severe instability, laxity of the knee and great bone loss are encountered.

Multigen Plus components should not be used alone as single components, or with implants from other manufacturers.

1.1. MATERIALS

The materials of Multigen Plus System are reported in the table of Section 2.1.

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon.

Multigen Plus components can be inserted cementless or implanted using cement (PMMA). Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone and applying the bone cement for prosthesis fixation.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50° C / 32-122° F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package should not be used if they are dropped or suffer other accidental impacts. Devices should not be modified in any way.**

The device's code and lot number should be recorded in the patient's case history by using the labels included in the components packaging.

The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

The components forming original LimaCorporate systems must be rigorously coupled following the given indications. The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

Multigen Plus knee system includes femoral components, tibial components, tibial liners and patellar components: it is a modular system indicated for use in total knee replacement, both primary and revision.

Multigen Plus total knee system is indicated for use in knee arthroplasty in skeletally mature patients with the following conditions: • Advanced articular destruction generated by primary degenerative or post-traumatic arthrosis or rheumatoid arthritis; • traumatic events with articular pain; • avascular necrosis; • congenital or acquired deformity; • failures of previous operations, articular reconstruction, arthrodesis or total arthroplasty.

Multigen Plus components are intended for both cemented and cementless use

Components	Material	Use	
		C	NC
Femoral component	CoCrMo	X	
	CoCrMo + PoroTi		X
	CoCrMo + CFR Peek + Ti6Al4V (for H design)	X	
Tibial plate	Ti6Al4V + UHMWPE	X	
	Ti6Al4V + PoroTi		X
	CoCrMo	X	
	CoCrMo + PoroTi		X
	CoCrMo + CFR Peek (for H design)	X	
Tibial liner	UHMWPE (for CR, PS and UC design)	X	X
	UHMWPE + Ti6Al4V (for CCK design)	X	
	UHMWPE + CoCrMo + Ti6Al4V (for H design)	X	
Tibial stem and modules	Ti6Al4V		X
Femoral stems and modules	Ti6Al4V		X
Femoral and tibial augments	Ti6Al4V	X	
Patella	UHMWPE	X	
Material Standards			
Ti6Al4V (ISO 5832-3 - ASTM F1472) – CoCrMo (ISO 5832-4 – ASTM F75) – UHMWPE (ISO 5834-2 – ASTM F648) – Poroti Titanium Coating (ASTM F1580) – CFR Peek (ASTM F2026)			

C = Cemented; NC = Not Cemented

2.2. CONTRAINDICATIONS

Contraindications include: • acute or chronic infections, local or systemic infections; • serious muscular, neurological or vascular diseases affecting the concerned limb; • bone destruction or poor bone quality, which could compromise implant's stability; • any concomitant disease and dependence that might affect the implanted prosthesis; • allergy to material.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight; • strenuous physical activities (active sports, heavy physical work); • fretting of modular junctions; • incorrect implant positioning (e.g. varus positioning); • insufficient bone to support the femoral and/or tibial components; • medical disabilities which can lead to an unnatural gait and loading of the knee joint; • Improper fixation of the device; • Bone resorption (e. g. due to osteolysis); • errors of operative technique; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumor, infections); • severe deformity leading to impaired anchorage or improper positioning of implants.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques. The responsibility for negative outcomes and complications resulting from a failure to follow these instructions for use, from incorrect material handling and from incorrect clinical indications resides with the surgeon. In such circumstances the manufacturer accepts no liability. The components of Multigen Plus system should not be used outside the allowed combinations or with components belonging to other manufacturer.

The permitted combinations of femur and tibia sizes for Multigen Plus systems are shown in the following tables:

For CR/PS systems with mobile plates (the size of the mobile liner must match the size of the femoral component)

		Femoral Component = Tibial Liner					
		#0	#1	#2	#3	#4	#5
Tibial Plate	#0	✓	✓				
	#1	✓	✓	✓			
	#2	✓	✓	✓	✓		
	#3	✓	✓	✓	✓	✓	
	#4	✓	✓	✓	✓	✓	✓
	#5	✓	✓	✓	✓	✓	✓

For CR/PS systems with fixed coupling (the size of the fixed liner is the same of the size of the fixed tibial plate)

		Femoral Component					
		#0	#1	#2	#3	#4	#5
Tibial Plate = Tibial Liner	#0	✓					
	#1		✓	✓	✓	✓	✓
	#2		✓	✓	✓	✓	✓
	#3		✓	✓	✓	✓	✓
	#4		✓	✓	✓	✓	✓
	#5		✓	✓	✓	✓	✓

For ultracongruent (UC) liners with fixed coupling (the size of the fixed liner is the same of the size of the fixed tibial plate)

		Femoral Component					
		#0	#1	#2	#3	#4	#5
Tibial Plate = Tibial Liner	#0	✓					
	#1		✓				
	#2		✓	✓			
	#3		✓	✓	✓		
	#4		✓	✓	✓	✓	
	#5		✓	✓	✓	✓	✓

For CCK system (the size of the tibial liner is the same of the size of the tibial plate)

		Femoral Component				
		#1	#2	#3	#4	#5
Tibial Plate = Tibial Liner	#1	✓	✓			
	#2	✓	✓	✓		
	#3		✓	✓	✓	
	#4			✓	✓	✓
	#5				✓	✓

For H system (the size of the tibial liner must match the size of the femoral component)

		Femoral Component = Tibial Liner				
		#1	#2	#3	#4	#5
Tibial Plate	#1	✓	✓			
	#2	✓	✓	✓		
	#3	✓	✓	✓	✓	
	#4	✓	✓	✓	✓	✓
	#5	✓	✓	✓	✓	✓

Note: the information about allowed model combinations are reported also in the surgical techniques. Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices based on the anatomy and specific conditions of each patient. Inadequate preoperative planning can lead to improper selection of the implants and/or incorrect implant positioning. Complications or failures of the total knee replacement are more likely to occur in heavy and highly active patients. The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing knee replacement. LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can fail or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time, in fact, can be influenced by the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection, pain, stiffness, reduced range of motion). The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors. Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person. Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery. It is extremely important to inform the patient before the surgery on the factors that may prejudice the success of the operation. The information to the patient has to be documented before the surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants are available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may affect system performance and survival of the implant. The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labelled indications. Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery. Care must be taken to protect the surfaces involved in the coupling between components (e.g. taper coupling between the tibial plate and the tibial stem and coupling between the tibial plate and the tibial liner); the articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique. Particles of bone cement or bone can cause an increased rate of wear of the polyethylene insert (so-called third body wear). To avoid this problem the joint and the surrounding bone-implant interfaces should be carefully washed and cleaned and all foreign material must be removed.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect possible changes in position or condition of the implant or surrounding tissues. The surgeon should make the patient aware of the limitations of limb function after knee arthroplasty and that caution is needed, especially in the post-operative period. Excessive physical activity or trauma to the involved knee can lead to premature failure through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon:

- avoid repeated high weight lifting;
- keep body weight under control, overweight conditions may adversely affect the outcomes of joint replacement;
- serious morphological alterations of the supporting bone structures or articular damage of the sliding coupling might modify the performance expected by the prosthetic implant;
- avoid sudden peak loads (consequences of activities, such as running and skiing), movements which can lead to sudden stops or twisting;
- avoid movements associated with excessive range of motion and positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that can occur in knee arthroplasty include: • loosening of the prosthetic components; • prosthesis dislocation and instability; • damage to the prosthetic implant; • instability of the system because of inadequate soft tissue balancing; • dissociation due to incorrect coupling of the devices; • infection; • tissue reaction to the implant material or abrasion; • local hypersensitivity; • bone fracture; • local pain; • periprosthetic fractures; • temporary or permanent nerve damage; • flexion contracture; • limited range of motion; • lengthening or shortening of the leg; • anesthetic and intraoperative events; • patellofemoral disorders; • arterial injuries; • wound problems; • fracture of the devices; • excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles; • additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, hematomas, systemic allergic reactions and systemic pain.

4. STERILITY

a. Implants

All implantable components of the Multigen Plus knee system are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . Metal components are sterilized by radiation or EtO and UHMWPE components by EtO. Do not use any component from a package that has been previously opened or appears to be damaged. **Do not use implants after the expiration date printed on the label.**

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the *Products* section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment's.

5. MAGNETIC RESONANCE IMAGING (MRI)

The Multigen Plus System has not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near the implant site.

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Limacorporate S.p.A.