G1/G3

IT

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Main characteristics

Intended use

G1 and **G3** are radiopaque bone cements for surgical use to perform arthroplasty procedures, such as hip replacement, knee replacement, ankle replacement, shoulder replacement and other joint replacements, specifically formulated to allow fixation of prosthetic devices on the living bone.

G1 and **G3** acrylic-based bone cements demonstrate a short setting time and fast hardening, extremely high mechanical performance, high intrusion, low heat generation rate (versus ISO 5833:2002 – “Implant for surgery-acrylic resin cements”).

G1 standard viscosity allows manual spatula application and comes in a 20 g and 40 g powder package (**G1 20** and **G1 40**).

G3 low viscosity allows syringe application and comes in a 20, 34.4, 60 g powder package (**G3 20**, **G3 40** and **G3 60**).

G1 and **G3** bone cements must be prepared immediately before application, using the two sterile components.

Preparation, handling and application of **G1** or **G3** bone cements must be performed only by qualified healthcare professionals, specifically trained in the procedure and under the direct supervision of the physician responsible for the procedure.

The manufacturer declines any liability in case of the use of **G1** or **G3** bone cements not strictly corresponding to the specified intended use, or the use of **G1** or **G3** bone cements by personnel not adequately qualified and trained.

Composition

Powder components
PMMA (Polymethyl Methacrylate)
BPO (Benzoyl Peroxide)
Barium Sulfate
Liquid components
MMA (Methyl Methacrylate)
DMPT (N,N dimethyl-p-toluidine)
HQ (Hydroquinone)

Content

	G1		G3		
	20	40	20	40	60
Powder weight (g)	20	40	20	34.4	60
Liquid volume (ml)	9.1	18.3	11.6	20.0	35.0 ¹

Packaging and sterility

The manufacturing and packaging process of **G1** and **G3** is performed under strict quality procedures in a controlled environment, which conforms to the most stringent international standards.

The liquid component is sterilized by filtration, the powder component is sterilized by ethylene oxide. The sterile liquid is contained in an amber glass vial packed in a blister pack in sterile conditions and subsequently sterilized by ethylene

oxide. The powder is packed in two bags in sterile conditions. The internal bag in medical-grade paper and polyethylene containing the powder component is inserted in another bag in Tyvek and polyethylene and both are sterilized by ethylene oxide. The two bags are packed in a non-sterile protective aluminium wrapping.

Before using **G1** or **G3**, carefully check the integrity of the packaging. If the package is incomplete, damaged or unsealed, the cement cannot be used and must be discarded. After opening the package, it is mandatory, and the responsibility of the operator, to use an aseptic handling technique. Any error in handling and during transfer into the sterile field might affect bone cement sterility, the sterility of the surgical intervention and imply a risk of severe complications for the patient, such as infections and sepsis.

Clinical indications

Bone cements **G1** and **G3** are intended to be used during surgical interventions caused by:

- Arthroplasty revision
- Joint replacement after a traumatic event
- Osteoarthritis
- Osteoporosis hangover

Contraindications

Do not use **G1** and **G3** when a patient is known to be hypersensitive to the constituents of bone cement or to the contrast medium (barium sulfate).

Absolute contraindications are: local or systemic infections not completely resolved; application where there is no possibility of bone regeneration.

Side effects

- It is commonly known that all bone cements may potentially and rarely directly or indirectly cause adverse patient reactions affecting the cardiovascular system. Hypotensive reactions have occurred between 10 and 165 seconds following the application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- Rare severe adverse events, associated with the use of bone cements include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism, thrombophlebitis, haemorrhage, haematoma, infection of surgical wounds.
- Use of bone cement during pregnancy and breastfeeding: no contraindication known, but use of bone cement should be avoided.
- The patient must be made aware and accept potential risks and complications connected with the use of bone cement. The patient needs to be informed about the measures that could be taken to reduce possible negative effects.

¹ In two 17.5 ml vials to be used at the same time.

Familiarization with G1 and G3 bone cements

Before using **G1** or **G3**, the surgeon should become familiar with their properties, handling, application, and with the handling of other devices which come into contact with the bone cement, or are going to be used. The surgeon shall perform simulations of the procedure following correct medical practice and the instructions of the manufacturers of involved devices. It is also recommended for surgeons to practice mixing, handling and application prior to use. The characteristics of the bone cement and the timing of the different phases depend on the ambient temperature, humidity and mixing technique. The surgeon shall simulate the entire process of mixing, handling and application in similar environmental conditions to the operating room.

Warnings when using G1 or G3 bone cements

- Monitor the patient carefully for any change in blood pressure during and immediately following the application of the bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements. Hypotensive reactions have occurred between 10 and 165 seconds following the application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- Rare severe adverse events, associated with the use of bone cements include myocardial infarction, cardiac arrest, cerebrovascular accident and pulmonary embolism.
- In the case of cardiovascular or pulmonary complications it is necessary to control and re-establish an appropriate level of volemia.
- Be careful during the bone cement mixing phase and strictly follow the instructions.
- Avoid direct skin or eye contact with the liquid component of dough and reduce as far as possible exposure to monomer vapors, which may cause irritation of the airways, eyes and, in rare cases, affect the liver.
- Ventilate the room to eliminate monomer vapors. The liquid component is volatile and flammable. In the presence of monomer vapors do not use electrocautery instruments or other high temperature sources.
- Do not use latex gloves or other latex devices. The liquid component is a lipid solvent which can cause glove perforation and may damage exposed tissues. PVP gloves (three layers: polyethylene, vinyl copolymer, polyethylene) or Viton®-butyl gloves give adequate protection for a long period. If surgical synthetic rubber gloves are going to be used, it is advisable to wear a second pair of gloves on top, suitable for bone cement handling.
- Operators wearing contact lenses shall not mix

- bone cement or be exposed to monomer vapors.
- The polymerization of bone cement is an exothermic reaction which finishes only when the bone cement becomes hard inside the body. Respect the waiting time before application (fig. 1) and consider that generated heat may damage bone tissue or other tissues in contact with bone cement.
- **Do not add any other material to G1 or G3 bone cements.** The addition of non-approved ingredients (powder, water solutions,...) severely compromises physical and chemical characteristics of bone cement both in the preparation phase and after the implant.
- Use always the entire content of one powder bag mixed with the entire content of one single vial (two vials in the case of **G3 60**). It is not allowed to mix more than one powder dose at a time.
- The operator should be instructed on bone cement application techniques, and should respect them strictly. The implant bone cavity should be carefully cleaned, irrigated and dried before bone cement application to avoid contamination of the bone cement through blood or bone marrow. The graft of a foreign body in contact with biological tissues increases the risk of infection in the post-surgical phase.
- Single use device. This device cannot be re-used.
- The safety and effectiveness of G1 and G3 bone cement in pregnant women or children has not yet been established. Bone cement must not be used during the first trimester of pregnancy. During all other stages of pregnancy, the product must only be used in the case of life-threatening circumstances.
- When treating children, G1 and G3 bone cement can be used to save limbs only in the event that other treatment procedures will not have a good probability success.

Application instructions

When using **G1** or **G3**, respect normal application rules for bone cements to limit unpredicted effects. The bone cement implant needs to be performed by a specialist physician, specifically trained.

Preliminary work process

To the prepare bone cement, the following is necessary:

- a sterile working surface;
- a sterile bowl made of ceramic, stainless steel, polypropylene or other material specifically approved for contact with bone cement dough;
- a sterile spoon or spatula made of ceramic, stainless steel, polypropylene or other material specifically approved for contact with bone cement dough.

Only after checking the integrity of the entire packaging can the external non-sterile aluminium laminated envelope be opened by an assistant who will then open only the contained outer pouch, which is the sterility barrier, and then the internal pouch, containing the powder, can be laid on the sterile working surface using an aseptic technique. Similarly the assistant will open the blister containing the vial and then place the glass vial in the sterile field using an aseptic technique.

Any error in handling and during transfer into the sterile field might affect bone cement sterility, the sterility of the surgical intervention and imply a risk of severe complications for the patient, such as infections and sepsis.

A dose is prepared by pouring the entire content of the liquid component vial (two vials in case of **G3 60**) into a bowl containing the entire quantity of powder included in one bag.

Preparation and application

The preparation and application of **G1** or **G3** bone cements take place in four subsequent phases:

- I. mixing;
- II. waiting (in this phase the bone cement, when touched, sticks to the gloves);
- III. application (in this phase the bone cement has the aspect of a fluid dough, it does not stick and can be handled, injected by a syringe and a canula);
- IV. setting (in this phase the bone cement should have already been implanted and hardening completed. Non-implanted bone cement has in this phase becomes too viscous, it cannot be used anymore and must be discarded).

The duration of phases II to IV depends on the ambient temperature and humidity. A higher temperature accelerates hardening, while a lower temperature slows it down (see fig. 1).

Prepare the mix choosing the appropriate type of bone cement and the dose which is necessary for the kind of surgical intervention and the volume of the cavity to be filled. Always open the phial at a safe distance from the container (risk of shards of glass). All the liquid must be poured into the container first.

Carefully open the bag and pour the entire contents into the container.

The powder must always be added to the liquid.

Mix the entire quantity of powder and liquid in the container well, until you obtain an even mixture that is free of air bubbles.

Mix by hand inside the sterile bowl using a sterile spoon or spatula conforming to materials' specifications.

If you use a mixing or application device for surgical cements, follow the manufacturer's instructions.

Use a regular stirring action, not too fast, and continue for one minute. Do not exceed the mixing time.

Viscosity increases progressively as a consequence of the polymerization reaction, during phases II to IV.

Keep waiting until the completion of phase II (waiting phase) and then proceed with the application (phase III, the cement has become dough).

Dose

Choose the dose considering the intervention to be performed. As a preventive measure, it is advisable to keep at hand at least a second package of **G1** or **G3** bone cement to be able to prepare an extra quantity of bone cement, if required during the procedure.

The maximum recommended dose of cement is as follows:

Pack	No. of packs
G1 20	6
G1 40	3
G3 20	6
G3 40	3
G3 60	2

Other uses

G1 and **G3** are not permitted for any use, other than the specified intended use.

Storage

G1 and **G3** bone cements must be stored in the sealed original package, inside a dry, clean storage room, at a storage temperature not below 5°C and not above 25°C.

Expiry date, sterility, disposal

The expiry date is written on the external box and on the secondary outer blister labels.

G1 and **G3** cannot be used after the expiry date. **G1** and **G3** cannot be used if the package is unsealed or the secondary blister is not present or damaged.

Opened or damaged **G1** and **G3** packages must be disposed of with all their content. **G1** and **G3** are sterilized by ethylene oxide, aseptic filtration and radiation and cannot be re-sterilized. A yellowish colour of the powder or of the liquid is not normal, in this case the cement cannot be used and must be discarded. Dispose of the remaining **G1** or **G3** bone cements and the content of partially utilized, expired, un-usable or damaged packages, following rules and procedures applicable to this kind of hospital waste.

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Limacorporate S.p.A.