Medtronic

SYNCHROMED® II

8637

Programmable pumps

SYNCHROMED® II

Pompes programmables

SYNCHROMED® II

Bombas programables

Implant manual • Manuel d'implantation • Manual de
implantación

USA Rx only

2003

Filename Date Time
UC200xxxxxx EN
4.6 x 6.0 inches (140 mm x 197 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version: 03-03-2010

M221311A028 Rev A

2009-10
Printing instructions: *Refer to the "Implant Manual" category Table 1 in doc# A00002 for Neuro Core European Printing Instructions.*

Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Do not reuse

Sterilized using ethylene oxide

Caution, consult accompanying documents


Use by

Date of manufacture

Manufacturer

Temperature limitation

Keep away from magnets

Serial number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.

For USA audiences only

Authorized representative in the European community

Medtronic® and SynchroMed® are registered trademarks of Medtronic, Inc.

Table of contents

- Description 5**
- Package contents 6**
- Patient identification card 6**
- Device specifications 7**
- Device longevity 10**
- Flow rate accuracy 11**
 - Measurement error 11
 - Fluid volume 11
 - Environmental conditions 12
- Declaration of Conformity 14**
- Instructions for use 15**
 - Preparing for pump implant 15
 - Sterile procedure 16
 - Emptying the pump 16
 - Preparing to fill the pump 16
 - Filling the pump 17
 - Replacing an implanted pump 17
 - Preparing the pump pocket 18
 - Implanting the pump 19
 - Programming the pump 20
 - Refilling the pump 21
 - Accessing the catheter port 21

Refer to the Indications, Drug Stability, and Emergency Procedures reference manual for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, and component disposal.

Refer to the appropriate drug labeling for indications, contraindications, warnings, precautions, dosage and administration information, and screening procedures.

Appendix A—Refilling the pump 22

- Contraindications 22
- Preliminary procedures 22
- Emptying the pump 23
- Performing a reservoir rinse 24
- Refilling the pump 25
- Programming the pump 27
- After the refill procedure 27

Appendix B—Accessing the catheter access port 28

- Contraindications 28
- Preliminary procedures 29
- Preparing the access site 29
- Intraspinal applications only 29
- Vascular applications only 31
- Flushing the catheter access port 32
- Clearing a catheter occlusion 32

Description

The implantable Medtronic Model 8637 SynchroMed II programmable pump is part of an infusion system that stores and delivers a prescribed drug to a specific site. The implanted infusion system consists of a Model 8637 SynchroMed II pump and a catheter.

The catheter connects to the pump catheter port. The pump is anchored in the pump pocket using the suture loops on the outside of the pump (Figure 1).


Figure 1. Pump exterior view.

The drug is stored in the pump reservoir (Figure 2). Per a programmed prescription, the drug moves from the pump reservoir, through the pump tubing, catheter port, and catheter, to the infusion site. The catheter access port (CAP) allows injection of drug directly into the implanted catheter for drug administration and diagnostic purposes. Drug injected into the CAP bypasses the pump mechanism and goes directly through the catheter port into the implanted catheter to the infusion site. The CAP allows entry of a 24-gauge noncoring needle to prevent accidental injection during refill procedures (which use the 22-gauge noncoring needle supplied in the refill kit).

The manufacturer and model code recorded on a radiopaque identifier are visible using standard x-ray procedures.


Figure 2. Pump interior view.

Package contents

- Pump
- Needle, 22-gauge (black sheath)
- Needle, 24-gauge (purple sheath)
- Product literature
- Registration form
- Patient identification card
- **[! USA]** Warranty card

Patient identification card

A patient identification card is packaged with this device. Advise the patient to carry the identification card at all times.

[! USA] The patient identification card packaged with the device is temporary; a permanent card is mailed to the patient when Medtronic receives the registration form.

[! USA] The implant registration form registers the device warranties and creates a record of the device in the Medtronic implant data system.

Device specifications

Table 1. Shipping and operating values for the Model 8637 SynchroMed II pump

	Shipping	Operating
Fluid in reservoir	Sterile water	—
Shipping flow rate	0.006 mL/day	—
Infusion modes	Simple continuous	Single bolus Priming bolus Bridge bolus Simple continuous Flex Minimum rate Stopped pump
Alarms		
Critical alarm	Disabled	Enabled with an interval programmed
Non-critical alarm	Disabled	Enabled with an interval programmed

Table 2. Device specifications for the Model 8637 SynchroMed II pump^a

	8637-20	8637-40
Pump		
Thickness (including septum)	19.5 mm	26.0 mm
Weight (empty/full)	165/185 g	175/215 g
Displacement volume	91 mL	121 mL
Diameter (including CAP)	87.5 mm	87.5 mm
Pump reservoir		
Volume	20.0 mL	40.0 mL
Residual volume	1.4 mL	1.4 mL
Fill volume at shipping	17.5 mL	37.5 mL
Pump tubing		
Volume ^b	0.25 mL	0.25 mL
Reservoir fill port		
Septum puncture life	500 punctures	500 punctures
Catheter access port		
Prime volume	0.14 mL	0.14 mL
Septum puncture life	500 punctures	500 punctures
Flow rate		
Maximum programmable ^c	24 mL/day	24 mL/day
Minimum programmable ^c	0.048 mL/day	0.048 mL/day
Stopped pump maximum leakage	0.030 mL/day	0.030 mL/day
Bacterial retentive filter		
Pore size	0.22 µm (micron)	0.22 µm (micron)
Power source		
Battery	Lithium hybrid cathode	Lithium hybrid cathode
Longevity	Rate dependent (Figure 3)	Rate dependent (Figure 3)

Table 2. Device specifications for the Model 8637 SynchroMed II pump^a (continued)

	8637-20	8637-40
Radiopaque identifier	NGP	NGV
Reservoir pressure	20.68 kPa to 34.75 kPa	20.68 kPa to 34.75 kPa

^a All measurements are approximate.

^b **If the pump has been replaced and the catheter has not been replaced and has not been aspirated**, use a priming bolus of 0.300 mL to fill the pump tubing with drug **before** connecting the catheter and implanting the pump.

^c Actual limits depend on pump calibration constant and selected infusion mode.

Table 3. Material of components in the Model 8637 SynchroMed II sterile package

Component	Material	Material contacts human tissue	Material contacts drug
Pump			
Exterior	Titanium	Yes	No
Reservoir	Titanium	No	Yes
Reservoir valve	Titanium	No	Yes
Tubing	Silicone rubber	No	Yes
Reservoir fill port septum	Silicone rubber	Yes	Yes
Catheter access port septum	Silicone rubber	Yes	Yes
Catheter port	Titanium	Yes	Yes
Bacterial retentive filter	Polyvinylidene fluoride	No	Yes
Suture loops	Titanium	Yes	Yes
Propellant	Inert gas	No	No
Needles	Stainless steel	Yes	Yes

Device longevity

Device longevity is a function of flow rate. Flow rates affect the battery voltage and motor revolutions (Figure 3).


Figure 3. Typical device longevity based on flow rate.

Device longevity is the calculated number of service months remaining based on actual usage rates.¹ An elective replacement indicator (ERI) message displays on the programmer when the pump nears the end of its service life (EOS). At ERI, the pump continues to operate within specifications. The ERI thresholds allow the pump to operate for a minimum of 90 days, at rates up to 1.5 mL/day, between ERI activation and EOS (Figure 4). When activated, ERI is date stamped and displayed by the programmer after interrogating the pump. The EOS activation indicates the pump has reached the end of its service life. At EOS, the pump stops, but telemetry is available until the pump battery is depleted.

¹ Device longevity sources include battery life (voltage), device life (years), and motor life (revolutions).


Figure 4. ERI and EOS.

Flow rate accuracy

The flow rate accuracy of the pump is within $\pm 14.5\%$ of the programmed flow rate at 0.048-24 mL/day, 37 °C, 50% reservoir volume, and 300 meters above sea level. Measurement error, fluid volume, and changes in environmental conditions (eg, body temperature and atmospheric pressure) all affect the flow rate. The effects of these changes on flow rate are cumulative if the conditions exist simultaneously.

Measurement error

The apparent flow rate based on clinical measurements can vary due to measurement error (eg, syringe measurement accuracy, human error, and the volume of fluid in the extension tubing and filter).

Fluid volume

The flow rate of the pump varies slightly with the volume of fluid in the pump reservoir. The pump flow rate decreases as the reservoir volume approaches 1 mL. The pump flow rate decreases rapidly and then stops as the reservoir volume decreases from 1 mL to 0 mL. Therefore, the pump should be refilled prior to reaching 1 mL or less. Typically, the flow rate decreases by about 4% as the volume is reduced from the half-full volume to a volume of 1 mL. The usable volume is the reservoir volume minus 1 mL (Figure 5).


Figure 5. Flow rate accuracy as a function of fluid volume in reservoir.

Environmental conditions

Body temperature

The flow rate of the pump varies with body temperature. The flow rate increases as the temperature increases above 37°C and decreases as the temperature decreases below 37°C (Figure 6).


Figure 6. Flow rate accuracy as a function of temperature (typical effect).

Atmospheric pressure

Patients living or traveling (eg, airline flights, mountain climbing) at altitudes above sea level are exposed to lower atmospheric pressures. Within days of exposure to the lower pressures, the flow rate of the pump can increase and then stabilize at the higher flow rate. In circumstances where a potential increase in flow rate may pose a risk to a patient, reprogramming the infusion prescription offsets this higher flow rate (Figure 7).

In rare instances, exposure to the lower atmospheric pressure can cause the pump to deliver more than 14.5% of the programmed flow rate while the patient is exposed to the lower pressure. Consider changes in drug concentrations or changes to pump programming for patients exposed to lower pressures.


Figure 7. Flow rate accuracy as a function of altitude (typical effect).

Declaration of Conformity

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment and Directive 90/385/EEC on Active Implantable Medical Devices.

For additional information, contact the appropriate Medtronic representative listed on the inside back cover of this manual.

Instructions for use

Implanting physicians should be experienced in pump and catheter implant procedures and should be thoroughly familiar with all product labeling.

Cautions:

- Do not implant a pump that was dropped onto a hard surface or shows signs of damage. Implanting a pump that has been dropped or damaged can result in lack of intended therapy, and require additional surgery to replace the pump.
- Do not implant the pump unless pump operation has been confirmed. Failure to confirm pump operation before implant can result in additional surgery to replace the pump.
- Do not prematurely activate the pump reservoir valve. Activation of the pump reservoir valve seals the pump reservoir valve closed. Unusual resistance or the inability to inject the entire fill volume may indicate activation of the pump reservoir valve. If the valve closes, a portion of the reservoir contents must be delivered or removed before filling can be completed. Procedural delays can occur. To prevent activation of the pump reservoir valve during emptying and filling procedures:
 - completely aspirate all contents of the pump reservoir before filling.
 - do not allow air into the pump reservoir through an open needle in the septum or an unclamped extension.
 - do not exceed the maximum reservoir volume indicated in the pump labeling.
- To use the nonsterile clinician programmer in a sterile field, place a sterile barrier between the patient and the programming head to prevent infection. Do not sterilize any part of the clinician programmer. Sterilization may damage the programmer.

Preparing for pump implant

1. Assemble equipment and supplies.

Sterile items

- The pump package containing the pump, 22-gauge noncoring needle (for filling the pump), and 24-gauge noncoring needle (for flushing the catheter access port)
- Empty 20-mL syringes (for emptying the pump)
- 0.22- μ m (micron) filter
- Syringe containing prescribed fluid (volume not to exceed the reservoir volume of the pump)
- 10-mL syringe with 1–2 mL of sterile, preservative-free saline (for flushing the catheter access port)

Nonsterile items

- Medtronic clinician programmer
- ### 2. Before opening the shelf package, use the clinician programmer to interrogate the pump and verify pump battery status and current settings.

- a. Confirm that there are no active alarm events.

Note: If the pump is still in Shelf Mode, audible alarms are disabled. The pump must be interrogated to determine if an alarm has been activated.

- b. Confirm that the pump calibration constant displayed on the screen matches the calibration constant printed on the shelf package.

Warning: The calibration constant displayed on the programmer screen after reading the pump status must match the calibration constant printed on the shelf package. If calibration constants differ, contact the appropriate Medtronic representative listed on the inside back cover of this manual. Using an incorrect calibration constant can result in a clinically significant or fatal drug underdose or overdose.

3. Attach a "FOR YOUR RECORDS" label (enclosed in the shelf package) to the patient's record. This label displays the pump model number, reservoir size, calibration constant, and serial number.

Note: Updating the pump with the new parameters can be performed at this time or after the implant procedure. Refer to "Programming the pump" on page 20 for instructions.

Sterile procedure

1. Open the sterile pump package and remove the pump.
2. Remove the protective cap from the catheter port (a small amount of water might be present in the protective cap).

Emptying the pump

1. Assemble the 22-gauge noncoring needle and the empty syringe.
2. Insert the needle into the reservoir fill port until the needle touches the metal needle stop.
3. Withdraw the sterile water from the pump into the empty syringe (the pump is shipped nearly full).

Note: If the volume of fluid in the pump reservoir exceeds the volume of the syringe used for emptying, remove the filled syringe and needle. Attach an empty syringe and needle, and repeat until the pump reservoir is empty.

4. Empty the pump reservoir until air bubbles no longer appear in the syringe, ensuring all water and air is removed from the pump reservoir.
5. Remove the syringe and needle from the reservoir fill port.

Preparing to fill the pump

1. **If using preservative-free ziconotide sterile solution**, refer to the preservative-free ziconotide sterile solution drug labeling for instructions for use of this drug with the pump.
2. **For all indicated drugs except preservative-free ziconotide sterile solution**, refer to Table 4 to determine the fill method.

Notes:

- A change in concentration is not recommended at the time of replacement.

- The pump reservoir capacity is 20 mL or 40 mL. Because some sterile water remains in the pump reservoir, the final concentration of drug varies based on the fill method.

Table 4. *Expected concentration of drug in pump reservoir based on fill method*

Pump reservoir capacity	Filling without rinsing	Rinsing with 3 mL of drug	Rinsing with 10 mL of drug
8637-20	93%	98%	99%
8637-40	97%	99%	100%

3. If you are rinsing the pump reservoir before filling, rinse and discard the appropriate volume based on the fill method shown in Table 4.

Filling the pump

1. Attach the filter to the syringe containing the prescribed fluid.
2. Attach the needle to the syringe containing the prescribed fluid and filter, and purge the air from the fluid pathway.
3. Read the actual fill volume in the syringe.
4. Insert the needle into the reservoir fill port, and inject the prescribed fluid slowly into the pump reservoir.
5. If the reservoir valve is activated before the pump is filled completely, discontinue injection, remove the needle from the reservoir fill port, and return to "Emptying the pump" on page 16, step 4.
6. When filling is complete, remove the needle from the reservoir fill port.
7. Flush the catheter access port using a 24-gauge noncoring needle and a syringe filled with 1 to 2 mL of saline (or a heparinized solution for vascular applications, if not contraindicated).
 - a. Gently insert the needle into the catheter access port until the needle touches the metal needle stop.
 - b. Inject fluid into the catheter access port until fluid is observed at the catheter port.
 - c. Remove the needle from the catheter access port.
8. If implanting a new pump, go to "Preparing the pump pocket" on page 18.

Replacing an implanted pump


1. Remove the suture at the catheter connector and disconnect the implanted pump from the implanted catheter. To avoid damage to the pump connector, leave the connector attached to the catheter.

Note: If you are replacing a SynchroMed EL Model 8626, 8626L, 8627, or 8627L Implantable Pump, interrogate the replaced pump for catheter volume information that may be stored in the pump. If you are replacing a SynchroMed II Model 8637 Implantable

Pump, interrogate the replaced pump for catheter volume information. Enter the catheter volume information into the clinician programmer.

2. If not replacing the catheter, slowly aspirate 1 to 2 mL of fluid from the catheter using a 1-mL tuberculin syringe. Leave the syringe in place to avoid CSF loss. Aspirating directly from the catheter clears the catheter of drug and confirms catheter patency.

Note: Conditions might exist under which the catheter is not patent or is not aspirated. If the catheter is not patent it must be replaced. Refer to the technical manual packaged with the catheter for catheter replacement instructions.

 **Warning:** During vascular applications, do not aspirate blood through the catheter access port or catheter. Blood sampling or aspiration through the catheter access port is contraindicated in vascular applications. Residual blood from aspiration or blood sampling can occlude the catheter or pump and inhibit drug delivery, resulting in a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose, and require surgical revision or replacement.

3. Continue the implant with one of the two following procedures:
 - **If the catheter has been replaced or aspirated**, proceed to "Implanting the pump" on page 19.
 - **If the catheter has not been replaced and has not been aspirated**, use a priming bolus of 0.300 mL to fill the pump tubing with drug **before** connecting the catheter and implanting the pump. Refer to the programming guide packaged with the programmer software for information on how to calculate and program this bolus. The pump internal tubing prime must be complete before attaching the catheter to the pump. If not, drug present in the catheter is bolused into the intrathecal space. Proceed to "Implanting the pump" on page 19.

 **Warnings:**

- If this is a pump replacement and the catheter has not been replaced and has not been aspirated, the pump tubing should be primed **before** connecting the catheter and implanting the pump. **Do not** program a postoperative priming bolus after the catheter has been connected to the pump. Programming a postoperative priming bolus in this situation can result in a clinically significant or fatal overdose.
- Use the catheter length recorded at implant or catheter revision when calculating catheter volume. The actual implanted catheter length and catheter model number are required to accurately calculate catheter volume. **A universal value does not exist that can be used as a substitute for this knowledge.** An inaccurate calculation of the catheter volume can result in a clinically significant or fatal drug underdose or overdose.

Preparing the pump pocket

Prepare the subcutaneous pocket using an incision in the lower abdomen.

 **Caution:** Select a location in the lower abdomen that is:

- away from bony structures (eg, 3 to 4 cm) to minimize discomfort at the pump site.
- away from areas of restriction or pressure to minimize the potential for skin erosion and patient discomfort.
- a minimum of 20 cm away from another programmable device to minimize telemetry interference and incorrect or incomplete programming.
- in an area accessible to the patient for proper operation of a patient control device (if applicable).

In the pediatric population, care must be taken to select an appropriate location by taking into consideration:

- available body mass.
- presence of ostomies.
- growth and development.

Implanting the pump

1. Connect the implanted catheter to the pump by pushing the pump connector onto the catheter port and suturing.
2. Place the filled pump into the prepared pocket.

 **Cautions:**

- Implant the pump no more than 2.5 cm from the surface of the skin in order to maintain access to the reservoir and catheter access ports. Implantation of the pump is contraindicated if the pump can not be implanted 2.5 cm or less from the surface of the skin.
- Place the pump in the prepared pocket so:
 - the reservoir fill port is anteriorly oriented and the reservoir fill port and catheter access port will be easy to access after implant.
 - no sutures to the skin will be directly over the reservoir fill port or the catheter access port.
 - the catheter is not kinked or twisted and is secured well away from the pump ports.

Improper component placement can result in inaccessible pump ports, inadequate drug delivery, component damage, or procedural delays, and require surgical revision or replacement.

3. Suture the pump in the subcutaneous pocket using the following steps:
 - a. Suture first to the fascia in the bottom of the subcutaneous pocket.
 - b. Use these two sutures and the lower suture loops on the pump to draw the pump into the pocket.
 - c. Tie the sutures.

- d. Suture the remaining two loops at the top of the pump pocket.
 - e. Tie the sutures, securing the pump into the pocket.
4. Irrigate the pump pocket.
 5. Close the incisions per normal procedure and apply dressing.

Programming the pump

1. Enter the following into the clinician programmer: patient information, catheter model number, implanted catheter length (in centimeters), drug name and concentration, and the volume of prescribed fluid placed in the pump reservoir at implant.

Note: If you are replacing a SynchroMed EL Model 8626, 8626L, 8627, or 8627L Implantable Pump, interrogate the replaced pump for catheter volume information that may be stored in the pump. If you are replacing a SynchroMed II Model 8637 Implantable Pump, interrogate the replaced pump for catheter volume information. Enter the catheter volume information into the clinician programmer.

Warning: Use the catheter length recorded at implant or catheter revision when calculating catheter volume. The actual implanted catheter length and catheter model number are required to accurately calculate catheter volume. **A universal value does not exist that can be used as a substitute for this knowledge.** An inaccurate calculation of the catheter volume can result in a clinically significant or fatal drug underdose or overdose.

2. If the catheter is new or has been aspirated, program a postoperative priming bolus to advance the drug from the reservoir to the catheter tip.
3. If the catheter has not been replaced or aspirated, program the pump to deliver the prescribed infusion.

Warning: If this is a pump replacement and the catheter has not been replaced and has not been aspirated, the pump tubing should be primed **before** connecting the catheter and implanting the pump. **Do not** program a postoperative priming bolus after the catheter has been connected to the pump. Programming a postoperative priming bolus in this situation can result in a clinically significant or fatal overdose.

4. Set the Low Reservoir Alarm (to at least 1 mL).
5. Program the pump with new parameters.
Note: Refer to the programming guide supplied with the programmer software for instructions on programming the pump.
6. Print out the patient's prescription and pump settings (pump status).
7. Place the prescription and pump settings (pump status) in the patient's records.
8. Determine the refill date from the printout.
9. Schedule a refill appointment.

Refilling the pump

Medtronic recommends that refill procedures be performed using the appropriate Medtronic refill kit. Refer to "Appendix A—Refilling the pump" on page 22 for pump refill instructions.

Accessing the catheter port

Medtronic recommends that catheter access port procedures be performed using the appropriate Medtronic CAP kit. Refer to "Appendix B—Accessing the catheter access port" on page 28 for catheter access port instructions.

Appendix A—Refilling the pump

Contraindications

Medtronic refill kits are contraindicated for all catheter access port procedures.

Be certain you are accessing the correct port when injecting fluids into the reservoir fill port or accessing the catheter access port of an implanted pump. ALWAYS:

- identify the pump model and reservoir volume.
- identify the location of the reservoir fill port and the catheter access port.
- use the instructions, noncoring needles, and other accessories provided in the appropriate kit.
- verify the location of the correct port during needle insertion, using other medical procedures as appropriate.
- refer to the appropriate drug labeling for indications, contraindications, warnings, precautions, adverse events, and dosage and administration information.

Improper injection into the pump pocket or catheter access port allowing drug leakage into surrounding tissue can result in significant tissue damage or a loss of or change in symptom control, drug withdrawal symptoms, or a clinically significant or fatal drug overdose. Improper injection into the catheter access port can result in a clinically significant or fatal drug overdose. Observe the patient after the pump refill procedure for any signs or symptoms that could indicate a pocket fill or any other drug-related adverse event due to the refill procedure.

Cautions:

- Use strict aseptic technique when accessing the reservoir fill port or the catheter access port. Failure to use aseptic technique can contaminate fluids or tissues and result in local or systemic infection.
- Use extreme caution when accessing the reservoir fill port or catheter access port of the implanted pump if local or systemic infection is suspected. Avoid contaminating the system or further spreading the infection. Local or systemic infection may require surgical revision or replacement of infusion system components.
- Do not spill or leak vesicant or cytotoxic drug into adjacent tissue during pump procedures. Spillage or leakage of vesicant drug into adjacent tissue can result in significant local tissue damage.

Preliminary procedures

1. Gather the following sterile equipment:

From the appropriate refill kit:

- Extension tubing set with clamp
- 0.22- μ m (micron) filter
- 22-gauge noncoring needle

- 20-mL empty syringe
- Fenestrated drape
- Template

Locally supplied:

- Syringe containing prescribed fluid (volume of fluid not to exceed the reservoir volume of the pump)
 - Cleansing agent¹
 - Sterile gloves³
 - Alcohol pads or swabs³
2. Prepare the programmer for use. Refer to the Programming Guide packaged with the pump software for instructions.

Emptying the pump

1. Identify the pump model, reservoir volume, and location of the reservoir fill port and catheter access port.
2. Confirm that the refill volume of the prescribed fluid does not exceed the reservoir volume of the pump.
3. Prepare the injection site by cleansing the area.
4. Open the kit and glove packages. Put on the sterile gloves.
5. Place the drape, exposing the pump site.
6. Using sterile procedures, assemble the needle, extension tubing, and empty syringe (Figure 8).
7. Locate the reservoir fill port septum by placing the template over the pump and aligning the edges of the template with the edges of the pump.
8. Close the clamp (Figure 8).
9. Gently insert the 22-gauge needle through the template center hole and into the center of the reservoir fill port septum until the needle touches the needle stop (Figure 8).

Caution: Do not use excessive force when accessing the reservoir fill port. Excessive force can damage the pump, requiring surgical replacement.

¹ Included in the Model 856X Refill Kit


Figure 8. Close clamp and insert needle into reservoir fill port.

10. Open the clamp and slowly withdraw the fluid from the reservoir into the empty syringe.
11. If the syringe maximum capacity is reached before the reservoir is completely empty, more than one syringe will be needed to empty the pump.
 - a. Close the clamp.
 - b. Remove the full syringe.
 - c. Attach an empty syringe.
 - d. Repeat steps 8–11 until the reservoir is empty.
12. Empty the reservoir completely, until air bubbles are present in the extension tubing and syringe. The amount withdrawn from the pump should approximately equal the reservoir volume from the current pump status or the calculated expected volume.
13. Close the clamp and remove the syringe.

Note: Keep the needle in the reservoir fill port septum and the clamp closed for the pump refill procedure that follows.
14. If decreasing concentrations or changing the prescribed fluid, go to "Performing a reservoir rinse" on page 24. If a reservoir rinse is not required, go to "Refilling the pump" on page 25.

Performing a reservoir rinse

To prevent drug overdose when decreasing concentrations or changing the prescribed fluid in the pump reservoir, always rinse the reservoir twice between solutions to remove the drug that remains in the reservoir after emptying the pump.

Warning: Rinse the pump reservoir between solutions when changing drug or decreasing drug concentrations in the pump reservoir. A significant amount of drug may be present in the pump reservoir after emptying the pump. This residual volume can not be removed by emptying the pump. Rinsing the pump reservoir between solutions minimizes the amount of drug in this residual volume but does not eliminate it. Failure to account for residual drug in the pump reservoir can result in a concentration that is different than intended and a clinically significant or fatal drug underdose or overdose.

For programmable pumps, program a bridge bolus after rinsing the pump reservoir. The bridge bolus advances the old drug (the drug remaining in the pump tubing, catheter access port, and catheter after emptying and refilling the pump) to the catheter tip at the specified flow rate. Refer to the programming guide for specific procedures.

1. Empty the pump completely.
2. Fill the pump with sterile preservative-free Sodium Chloride Injection, USP.
3. Empty the pump completely.
4. Repeat steps 2 and 3.

Refilling the pump

1. Confirm that the refill volume of the prescribed fluid does not exceed the reservoir volume of the pump.
2. Attach the filter to the syringe containing the prescribed fluid and purge the air from the fluid pathway.
3. Attach the syringe with the prescribed fluid and filter to the extension tubing set (Figure 9).
4. Open the clamp and slowly (maximum 1 mL per 3 seconds) depress the plunger on the syringe to inject the prescribed fluid into the pump reservoir.


Figure 9. Open clamp and inject into pump reservoir.

- Caution:** Do not prematurely activate the pump reservoir valve. Activation of the pump reservoir valve seals the pump reservoir valve closed. Unusual resistance or the inability to inject the entire fill volume may indicate activation of the pump reservoir valve. If the valve closes, a portion of the reservoir contents must be delivered or removed before filling can be completed. Procedural delays can occur. To prevent activation of the pump reservoir valve during emptying and filling procedures:
- completely aspirate all contents of the pump reservoir before filling.
 - do not allow air into the pump reservoir through an open needle in the septum or an unclamped extension.
 - do not exceed the maximum reservoir volume indicated in the pump labeling.
5. If you have activated the reservoir valve, complete steps a through f below. Otherwise, proceed to step 6.
- a. Discontinue injection, close the clamp, and remove the needle from the septum.
 - b. Remove the syringe with prescribed fluid and attached filter.
 - c. Attach an empty 20-mL syringe to the extension tubing set (needle attached).

- d. Insert the needle into the reservoir fill port, open the clamp, and aspirate until all fluid/air is removed.
- e. Close the clamp and remove the syringe containing the aspirate from the extension tubing set and discard the syringe.
- f. Repeat steps 1 through 4.

Note: For pumps with a reservoir valve, the amount of time before the valve will release is dependent on the duration and the amount of pressure applied after the valve is first activated. The more pressure exerted, the longer it may take to release the valve.

6. When filling is complete, close the clamp and carefully remove the needle from the reservoir fill port septum.
7. Remove the cleansing agent from the patient's skin using an alcohol pad.
8. Apply an adhesive bandage, if desired.
9. Discard all components of the kit.

Programming the pump

1. If the drug concentration or drug has been changed, program a bridge bolus.
2. If any prescription information has changed, enter the changed information into the clinician programmer: for example drug name, drug concentration, infusion rate, or volume of prescribed fluid in the pump reservoir.
3. Update the pump.

After the refill procedure

1. Print out the patient's prescription and pump settings (pump status).
2. Place the prescription and pump settings (pump status) information in the patient's records.
3. Determine the refill date from the printout.
4. Schedule a refill appointment.

Appendix B—Accessing the catheter access port

Contraindications

Medtronic catheter access port kits are contraindicated for all refill procedures.

Be certain you are accessing the correct port when injecting fluids into the reservoir fill port or accessing the catheter access port of an implanted pump. ALWAYS:

- identify the pump model and reservoir volume.
- identify the location of the reservoir fill port and the catheter access port.
- use the instructions, noncoring needles, and other accessories provided in the appropriate kit.
- verify the location of the correct port during needle insertion, using other medical procedures as appropriate.
- refer to the appropriate drug labeling for indications, contraindications, warnings, precautions, adverse events, and dosage and administration information.

Improper injection into the pump pocket or catheter access port allowing drug leakage into surrounding tissue can result in significant tissue damage or a loss of or change in symptom control, drug withdrawal symptoms, or a clinically significant or fatal drug underdose. Improper injection into the catheter access port can result in a clinically significant or fatal drug overdose. Observe the patient after the pump refill procedure for any signs or symptoms that could indicate a pocket fill or any other drug-related adverse event due to the refill procedure.

Cautions:

- Use strict aseptic technique when accessing the reservoir fill port or the catheter access port. Failure to use aseptic technique can contaminate fluids or tissues and result in local or systemic infection.
- Use extreme caution when accessing the reservoir fill port or catheter access port of the implanted pump if local or systemic infection is suspected. Avoid contaminating the system or further spreading the infection. Local or systemic infection may require surgical revision or replacement of infusion system components.
- Do not use excessive force, do not use syringes smaller than 10 mL, and do not exceed an injection rate of more than 5 mL per minute when accessing the catheter access port. Failure to follow these instructions can create excessive pressure and cause catheter damage or disconnection allowing leakage into surrounding tissue resulting in tissue damage or a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose. In addition, do not rock the needle sideways while injecting. Rocking can break the needle inside the pump, requiring surgical revision or replacement.
- Always use a bacterial-retentive filter when injecting through the catheter access port. The catheter access port does not contain a filter. Injecting into the catheter

access port without using a filter can contaminate fluids or tissues and result in local or systemic infection.

Preliminary procedures

Gather the following sterile equipment:

From the appropriate catheter access port kit:

- Extension tubing set with clamp
- 0.22- μ m (micron) filter
- 10-mL empty syringe (for aspirating in intraspinal applications)
- Fenestrated drape
- Appropriate template and noncoring needle
- The template and needle are pump-specific. Match the name on the template with the name of the pump. Match the color and gauge of the needle to the color and gauge indicated on the template.

Locally supplied:

- Syringe containing prescribed fluid (10 mL minimum)
- Syringe with 5 mL saline (for flushing if desired)
- Cleansing agent
- Sterile gloves
- Alcohol pads or swabs
- Adhesive bandage, optional

Preparing the access site

1. Identify the pump model, reservoir volume, and location of the reservoir fill port and catheter access port.
2. Prepare the injection site by cleansing the area.
3. Open the kit and glove packages. Put on the sterile gloves.
4. Place the drape, exposing the pump site.
5. Place the template over the pump, aligning the edges of the template with the edges of the pump. Locate the catheter access port septum. The catheter access port is located near the edge of the pump, not in the center of the pump.

Intraspinal applications only

Warning: During intraspinal applications, before injecting fluids through the catheter access port, aspirate approximately 1- 2 mL from the catheter. A significant amount of drug may be present in the catheter access port and catheter, and failure to remove the drug during catheter access port injections can result in a clinically significant or fatal drug overdose.

1. Using sterile procedures, assemble the needle, extension tubing, and empty syringe (Figure 10).
2. Close the clamp and gently insert the needle into the catheter access port septum until the needle touches the needle stop (Figure 10).


Figure 10. Close clamp and insert needle into catheter access port.

3. Open the clamp and aspirate approximately 1–2 mL to ensure removal of drug from the catheter access port and catheter. Refer to the specifications in the appropriate catheter technical manual for the catheter volume.
4. Close the clamp and remove the syringe.
Note: Keep the needle in the catheter access port septum and the clamp closed for the procedure that follows.
5. Attach the filter to the syringe containing the prescribed fluid and purge the air from the fluid pathway.
6. Attach the syringe with the prescribed fluid and filter to the extension tubing set (Figure 11).
7. Open the clamp and inject directly at an infusion rate not greater than 5 mL per minute (Figure 11).


Figure 11. *Open clamp and inject into catheter access port.*

8. During injection, check the needle puncture site for swelling.
9. When the drug has been injected, close the clamp, and carefully remove the syringe.
10. Proceed to "Flushing the catheter access port" on page 32.

Vascular applications only

1. Using sterile procedures, assemble the needle, extension tubing, filter, and syringe containing the prescribed fluid.
2. Purge the air from the fluid pathway.
3. Close the clamp and gently insert the needle into the catheter access port septum until the needle touches the needle stop.
4. Open the clamp and inject directly at an infusion rate not greater than 5 mL per minute (Figure 12).

⚠ Warning: During vascular applications, do not aspirate blood through the catheter access port or catheter. Blood sampling or aspiration through the catheter access port is contraindicated in vascular applications. Residual blood from aspiration or blood sampling can occlude the catheter or pump and inhibit drug delivery, resulting in a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose, and require surgical revision or replacement.


Figure 12. Open clamp and inject into catheter access port.

5. During injection, check the needle puncture site for swelling.
6. When the drug has been injected, close the clamp, and carefully remove the syringe, leaving the filter attached to the extension tubing.
7. Proceed to "Flushing the catheter access port".

Flushing the catheter access port

1. For vascular applications, flushing the catheter access port is required. For intraspinal applications, flushing the catheter access port is optional.
To flush the catheter access port, attach the syringe containing 5 mL of sterile saline to the extension tubing.
Note: For vascular applications, use a heparinized solution (if heparin is not contraindicated). For vascular applications where heparin is contraindicated and intraspinal applications, use a preservative-free solution.
2. Open the clamp and inject directly at an infusion rate not greater than 5 mL per minute.
3. During injection, check the needle puncture site for swelling.
4. Flushing the catheter access port also flushes the drug from the catheter. Drug therapy from the pump will be lost until the catheter is refilled by the drug flow from the pump. If you want to minimize the amount of time without therapy, program a purge bolus.
5. Upon completion of the flush, close the clamp, and carefully remove the needle from the catheter access port septum.
6. Remove the cleansing agent from the patient's skin using an alcohol pad.
7. Apply an adhesive bandage, if desired.
8. Discard all components of the kit.

Clearing a catheter occlusion

For vascular applications only:

1. Assemble a 10-mL syringe, stopcock, and the correct catheter access port needle supplied in the catheter access port kit.

2. Fill the syringe with 5 mL of a heparinized solution (if not contraindicated).
3. Enter the catheter access port as described in "Vascular applications only" on page 31.
4. Gently alternate irrigation and aspiration to open the catheter.
5. If resistance continues, aspirate as much fluid as possible from the catheter access port. Use a 1-mL tuberculin syringe, stopcock, and catheter access port needle and inject 0.5 mL of heparinized solution (1000 U/mL) to lightly pressurize the system. Close the stopcock, and remove the needle.

Caution: When clearing a catheter occlusion, do not inject more than 0.5 mL of fluid into the catheter access port. Injection of more than 0.5 mL may cause catheter damage or disconnection, allowing leakage into surrounding tissue, resulting in tissue damage or a loss of or change in therapy.

6. Check for catheter patency in 10 minutes using a 10-mL syringe filled with 5 mL of heparinized solution (1000 U/mL). Lightly pressurize the syringe and wait 10–15 seconds to observe for loss of resistance. Do not force infusion.
7. If the catheter remains occluded, aspirate through the catheter access port to remove the original 0.5 mL of heparinized solution.
8. Repeat step 5 or steps 4 and 5 using 0.5 mL of heparinized solution or streptokinase (refer to drug labeling for complete instructions) and wait 10 to 30 minutes before checking for catheter patency. Attempt to dislodge the clot before flushing.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.
© Medtronic, Inc. 2009