Medtronic

INTRATHECAL CATHETER

8731SC

with sutureless connector

CATHÉTER INTRATHÉCAL

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INTRATHEKAL- KATHETER

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INTRATHECALE KATHETER

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Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Do not reuse

Sterilization: ethylene-oxide gas


Caution, consult accompanying documents

Use by

Manufacturing date

Storage temperature

Lot number

Volume

Conformité Européenne (European Conformity).
This symbol means that the device fully complies with
European Directive AIMD 90/385/EEC.

For USA audiences only

Medtronic®, IsoMed®, SynchroMed®, and MiniMed® are registered trademarks of Medtronic, Inc.

Table of contents

Description 5

Package contents 5

Registration form 5

Device specifications 5

Instructions for use 7

Positioning the patient 7

Placing the catheter 7

Tunneling the catheter 11

Connecting and anchoring the catheter 11

Attaching the catheter to the pump 13

Calculating catheter length and volume 19

Refer to the Indications, Drug Stability, and Emergency Procedures reference manual for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, and component disposal.

Refer to the appropriate drug labeling for indications, contraindications, warnings, precautions, dosage and administration information, and screening procedures.

Description

The implantable Medtronic Model 8731SC Intrathecal Catheter (Figure 1) is part of an infusion system that stores and delivers parenteral drugs to the intrathecal space. The implanted infusion system components consist of a Medtronic pump and a Model 8731SC catheter. The catheter connects to the pump at the catheter port.

Note: The Medtronic Model 8731SC Intrathecal Catheter is designed for use with Medtronic implantable infusion pumps with the exception of Medtronic MiniMed infusion pumps.

Package contents

- Spinal segment with inserted guide wire
- Pump segment with attached sutureless pump connector, connector pin, and strain-relief sleeve
- Separate connector pin
- 15 T-gauge introducer needles (2) (11.4 cm and 9.3 cm)
- Transparent spinal segment strain-relief sleeves (2)
- Opaque pump segment strain-relief sleeves (2)
- V-wing anchors (2)
- Ruler (nonsterile)
- Product literature
-  Warranty card
-  Registration form

Registration form

 The implant registration form registers the device warranties and creates a record of the device in the Medtronic implant data system.

Device specifications

The two-piece catheter has a thin-walled spinal segment and a thicker-walled pump segment. These segments connect to each other with a connector pin and strain-relief sleeves. The pump segment connects to the pump with a sutureless pump connector.


Figure 1. Model 8731SC Intrathecal Catheter.

Table 1. Device specifications for the Model 8731SC^a catheter

Description	Value
Length	
Total catheter	104.1 cm
Spinal segment	38.1 cm
Pump segment	66.0 cm
Spinal segment	
Outer diameter	1.4 mm (4.2 French)
Inner diameter	0.53 mm
Marker spacing	1 cm starting at 14 cm
Catheter tip	Closed with 6 side holes
Pump segment	
Outer diameter	2.2 mm (6.5 French)
Inner diameter	0.53 mm
Catheter volume	0.0022 mL/cm
Connector pin inner diameter	0.41 mm
Transparent strain-relief sleeve outer diameter	4.1 mm
Opaque strain-relief sleeve outer diameter	4.0 mm
Guide wire outer diameter	0.46 mm
Introducer needle (long)	15 T-gauge, 11.4 cm
Introducer needle (short)	15 T-gauge, 9.3 cm
Trimnable segments	Spinal and pump segments
Pump segment to spinal segment separation force	> 5.0 N
Pump connector to pump separation force	> 10.0 N

^a All measurements are approximate.

Table 2. Material of components in the Model 8731SC sterile package

Component	Material	Material contacts human tissue	Material contacts drug
Catheter	Silicone	Yes	Yes
Catheter tip	Titanium	Yes	Yes
Pump connector	Silicone	Yes	Yes
	Silicone (external component)	Yes	No
	Titanium	Yes	No
	Nylon	Yes	No
Connector pin	Titanium	Yes	Yes
Strain-relief sleeve	Silicone	Yes	No
V-wing anchor	Silicone	Yes	No
Guide wire	Stainless steel	Yes	No
	Guide wire handle	Acetal resin (plastic)	Yes

Table 2. Material of components in the Model 8731SC sterile package (continued)

Component	Material	Material contacts human tissue	Material contacts drug
Introducer needle	Stainless steel	Yes	No
Needle hub	Polymethyl pentene	Yes	No
Stylet	Stainless steel	Yes	No
Stylet cap	ABS	Yes	No

Instructions for use

Implanting physicians should be experienced in intrathecal-access procedures and should be thoroughly familiar with all product labeling.

Cautions:

To prevent damage to the catheter, which may require additional surgery to repair or replace the catheter:

- do not twist or stretch the catheter.
- use only rubber-tipped forceps on the catheter and exercise care when using sharp surgical instruments near the catheter.
- exercise care when using electrocautery near the catheter.
- keep the introducer needle in place while the incision is made in preparation for securing the catheter after catheter placement.
- do not kink or trim the guide wire and do not trim the catheter while the guide wire is in place.

To prevent damage to the catheter or injury to the patient caused by the tip of the guide wire protruding through a fluid outlet hole:

- do not trim the catheter before placement.
- do not trim the catheter before the catheter is positioned and secured to tissue.
- keep the guide wire fully inserted into the catheter during placement.

To prevent damage to surrounding tissue and injury to the patient, exercise care when using electrocautery near the introducer needle.

Positioning the patient

1. Place the patient in a lateral recumbent position with the lumbar region slightly flexed. Drape to allow fluoroscopic visualization of the spine in regions where the catheter will be inserted and where the catheter tip will be placed. Position the patient's hips, legs, and arms to avoid pressure points and interference with fluoroscopic imaging of the spine.
2. Administer the appropriate anesthesia.

Placing the catheter

1. Select the introducer needle of the desired length for the patient's needs.
2. Orient the bevel of the 15 T-gauge introducer needle cephalad or parallel to the dural fibers and insert while using fluoroscopy. Use a paramedian oblique approach to insert the needle into the intrathecal space at a shallow angle (approximately equal to 30°) (Figure 2 on page 8). The entry point of the needle into the skin (or fascia, if needle insertion is performed through an open incision) should be parallel to the vertebral pedicle, approximately 1 to 1.5 vertebral levels below the interlaminar space selected for dural puncture and 1–2 cm lateral to the midline on the side of the intended pump pocket (Figure 2 and Figure 3 on page 8).

Note: In retrospective analysis, low catheter complication rates were associated with the paramedian oblique technique.¹


Figure 2. Shallow paramedian oblique insertion side view.


Figure 3. Shallow paramedian oblique insertion top view.

3. Advance the needle until the dura is penetrated.
4. Remove the needle stylet and confirm needle location by observing cerebrospinal fluid (CSF) backflow.
5. Replace the needle stylet to stop CSF flow. Make certain the stylet handle key is properly seated in the slot in the needle hub.
6. Confirm the catheter guide wire is seated completely with the hub against the end of the catheter.
7. Verify that the needle has a cephalad orientation by making sure that the stylet handle key faces up.
8. Remove the stylet.
9. Thread the spinal segment tip (the end with fluid-outlet holes) through the needle and advance the catheter to the desired location.

Note: When the catheter tip reaches the curved point of the needle, a slight increase in advancement pressure will be noted. If using the 11.4 cm needle, the first centimeter marking from the catheter tip will be located at the needle hub. If using the 9.3 cm needle, the first centimeter marking

¹ Follett, K. *et al.* Prevention of intrathecal drug delivery catheter-related complications. *Neuromodulation*, vol. 6(1); 2003. pp. 32-41.

from the catheter tip will be approximately 2.3 cm outside of the needle hub.

△ **Caution:** If the catheter must be retracted during positioning, do not withdraw the catheter through the introducer needle. The needle tip can damage the catheter, requiring additional surgery to repair or replace the catheter. Withdraw the introducer needle and catheter together. Then, carefully retract the catheter from the introducer needle. Begin the procedure again.

10. Verify the catheter position using fluoroscopy.
11. Prepare the anchor site by making an incision at the needle site to expose an area of the fascia that is large enough to place the v-wing anchor and the catheter connection. The incision should be long enough to avoid introducing kinks or sharp bends when anchoring the catheter.

Notes:

- To protect the catheter, keep the introducer needle in place while the incision is made. Alternatively, this incision can be made prior to catheter placement.
 - Undermine the edges of the incision to develop a smooth fascial plane for the catheter anchor.
12. Carefully remove the needle from the fascia and grasp the catheter near the exit site to prevent catheter dislodgement.
 13. Carefully slide the needle to the end of the catheter to prevent the needle tip from damaging the catheter during guide wire removal.
 14. Hold the catheter securely at the spinal exit site and remove the needle and guide wire from the catheter simultaneously (Figure 4).

△ **Caution:** To prevent catheter damage or dislodgement during introducer needle and guide wire removal:

- hold the catheter straight.
- hold the catheter securely at the exit site (do not squeeze the catheter because excessive compression will increase the difficulty of guide wire removal).
- use minimal traction to avoid quick or sudden removal.


Figure 4. Removing the needle and guide wire.

Note: If the catheter becomes twisted or if the guide wire seizes up, **stop. Let the catheter relax** and return to its original shape. Then slowly start again.

△ **Caution:** If the catheter must be repositioned after the guide wire has been removed, remove the catheter from the patient's body, lay the catheter straight, and then carefully reinsert the guide wire into the catheter avoiding excessive force. Do not reinsert the guide wire into the catheter while the catheter is in the patient's body, which

can cause injury to the patient or damage the catheter, requiring additional surgery to repair or replace the catheter.

15. Confirm CSF backflow.
16. Tie a purse-string suture to close the ligament around the catheter.
17. Slide the v-wing anchor onto the spinal segment (Figure 5 on page 10) and continue sliding the anchor down as close as possible to the fascia entry point.


Figure 5. Attaching the v-wing anchor.

Cautions:

- Always use an anchor to secure the catheter to the surrounding tissue to prevent catheter dislodgement or kinking.
 - Do not tie sutures around the catheter, which may damage or occlude the catheter.
 - Do not overtighten sutures on the anchor, which may damage the component.
 - Do not use metallic sutures for ligation, which may damage the component.
18. Using rubber-tipped forceps, clamp the end of the spinal segment of the catheter to prevent CSF loss during tunneling and pocket formation.
 19. Using heavy nonabsorbable sutures, suture the anchor wings together at the notch with the wings flat against each other to properly engage and grasp the catheter (Figure 6).
 20. Using heavy nonabsorbable sutures, suture the anchor wings together through the suture holes and attach to the surrounding fascia (Figure 6).


Figure 6. Placing the v-wing anchor.

21. Unclamp the spinal segment and confirm CSF backflow.

22. Using rubber-tipped forceps, re-clamp the spinal segment.

Tunneling the catheter

1. Refer to the pump implant manual for instructions on preparing the subcutaneous pump pocket.

Warning: When placing the catheter, avoid making sharp bends or kinks, which may occlude the catheter. Occlusions can result in a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose and require surgical revision or replacement.

2. Tunnel subcutaneously from the spinal incision site toward the pump implant site with the Model 8583 or 8586 Catheter Passer. Refer to the catheter length to determine the appropriate passer tool to use.
3. Pass the pump segment from the pump pocket site to the spinal incision site. Leave enough catheter near the pump pocket to allow for one or two complete catheter coils behind the pump.

Connecting and anchoring the catheter

Warning: Firmly secure all connections. Failure to secure connections can allow drug or cerebrospinal fluid (CSF) leakage into surrounding tissue, resulting in tissue damage or a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose.

1. At the spinal incision site, trim the spinal segment of the catheter as necessary. Leave approximately 5 cm outside the fascia entry point to allow for patient movement and to prevent kinking (Figure 7).

Caution: When trimming, allow enough slack in the catheter to accommodate patient movement. Trimming the catheter too short or pulling the catheter taut may result in catheter damage or dislodgement, requiring additional surgery to repair or replace the catheter.


Figure 7. Trimming the spinal segment.

2. Save the trimmed catheter for catheter length measurements and volume calculations.
3. Slide the transparent strain-relief sleeve, small end first, onto the spinal segment (Figure 8 on page 12).


Figure 8. Sliding the transparent strain-relief sleeve onto the spinal segment.

4. Unless it is desirable to trim the pump segment at the spinal incision site, proceed to step 5. If it is desirable to trim the pump segment, perform the following steps:
 - a. Leaving a minimum of 10 cm of catheter (to allow for strain relief), trim the pump segment.
 - b. Save the trimmed catheter for catheter length measurements and volume calculations.
 - c. Slide an opaque strain-relief sleeve, small end first, onto the pump segment.
 - d. Insert the connector pin into the pump segment until the pump segment is against the closest large ring.
 - e. Slide the opaque strain-relief sleeve towards the connector pin until the sleeve snaps into place.
5. On the pump segment, grasp the strain-relief sleeve near the connector pin and insert the connector pin into the spinal segment until the spinal segment is against the closest large ring. Be careful not to disrupt the catheter placement in the spine (Figure 9).


Figure 9. Inserting the connector pin.

6. Unclamp the spinal segment.
7. Confirm catheter patency by verifying CSF flow through the pump connector.
8. Using rubber-tipped forceps, clamp the pump segment to prevent further CSF loss.

Note: If the catheter must be removed from the connector pin, grasp the spinal segment of the catheter next to the connector pin midsection. Then grasp the midsection of the connector pin with a rubber-tipped forceps and pull or twist the connector pin off of the spinal segment. **Do NOT pull on the catheter.** Trim the spinal segment of the catheter before reconnecting. Save the trimmed catheter for catheter length measurements and volume calculations.
9. Slide the transparent strain-relief sleeve of the spinal segment towards the connector pin until the sleeve snaps into place (Figure 10).

Note: Do not cover the center section of the connector pin and do not place sutures over the strain-relief sleeve.


Figure 10. Sliding the strain-relief sleeve over the connector pin.

10. Place a heavy nonabsorbable suture in each of the two grooves in the center of the connector pin and attach to the fascia (Figure 11 on page 13).


Figure 11. Anchoring the catheter.

Attaching the catheter to the pump

Proper alignment and full engagement of the sutureless pump connector to the pump catheter port is critical in ensuring the catheter is properly and completely connected to the pump.

Warnings:

- Firmly secure all connections. Failure to secure connections can allow drug or cerebrospinal fluid (CSF) leakage into surrounding tissue, resulting in tissue damage or a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose.
- When connecting the sutureless pump connector, keep the sutureless pump connector in line with the catheter port and do not angle the sutureless pump connector. If connected at an angle, the sutureless pump connector could detach after surgery or an occlusion could occur at the connection site.
 - Detachment of the sutureless pump connector could allow drug or cerebrospinal fluid (CSF) leakage into surrounding tissue resulting in tissue damage or a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose, and require surgical revision or replacement.
 - An occlusion will inhibit drug delivery resulting in a loss of or change in therapy, which may lead to a return of underlying symptoms, drug withdrawal symptoms, or a clinically significant or fatal drug underdose, and require surgical revision or replacement.


Figure 12. Identification of components.

1. With a thumb and forefinger, grasp the tapered portion of the sutureless pump connector (Figure 13).


Figure 13. Grasping the tapered portion.

2. At the pump pocket site, position the catheter port of the pump in line with the opening of the sutureless pump connector (see Figures 14 and 15 for SynchroMed II pump, Figures 16 and 17 for SynchroMed EL pump, and Figures 18 and 19 for IsoMed pump).


Figure 14. Attaching the sutureless pump connector to the SynchroMed II pump (Lateral View).


Figure 15. Attaching the sutureless pump connector to the SynchroMed II pump (Top View).


Figure 16. Attaching the sutureless pump connector to the SynchroMed EL pump (Lateral View).


Figure 17. Attaching the sutureless pump connector to the SynchroMed EL pump (Top View).


Figure 18. Attaching the sutureless pump connector to the IsoMed pump (Lateral View).


Figure 19. Attaching the sutureless pump connector to the IsoMed pump (Top View).

3. Attach the sutureless pump connector using one of the following methods (see Figures 14 and 15 for SynchroMed II pump, Figures 16 and 17 for SynchroMed EL pump, and Figures 18 and 19 for IsoMed pump):

- Method 1

Firmly press the pump connector onto the catheter port until the connector fully covers the catheter port. The connector snaps into place.

- Method 2

- a. Firmly squeeze precisely on the oval marks of the pump connector.
- b. While squeezing, carefully press the pump connector onto the catheter port until the connector fully covers the catheter port. The connector snaps into place.

Cautions:

- Exercise care when connecting the catheter or pump connector to metal connectors or fittings. These can cut or puncture the catheter or pump connector, requiring additional surgery to repair or replace the damaged component.
- When connecting the catheter to the pump, do not introduce air bubbles into the catheter, which can lead to spinal headache or affect infusion accuracy.

4. Release the thumb and forefinger.
5. Check that the sutureless pump connector is properly attached by using the following method:
 - a. Grasp the tapered portion of the connector.
 - b. Tug as if to remove the connector from the pump. The connector should feel firmly attached (Figure 20).
 - c. Rotate the sutureless pump connector at least 90° to the right and to the left (Figure 21).

Note: The catheter should not appear to 'wobble' when rotated. The catheter port should not be visible from any position.


Figure 20. Tug on connector to check for proper attachment to pump.



Figure 21. Rotate connector to check for proper attachment to pump.

Notes:

- If the sutureless pump connector does not feel firmly attached:
 - a. Place a thumb and forefinger precisely on the oval marks located on the two flat portions of the connector (Figure 14 on page 14).

- b. Firmly squeeze the connector to apply pressure precisely on the oval marks.
 - c. While squeezing on the ovals, pull the connector off the catheter port.
 - d. Repeat steps 1-5 of "Attaching the catheter to the pump" beginning on page 13.
 - If the sutureless pump connector is difficult to remove, try one or both of the following methods:
 - While squeezing on the ovals, gently rock the connector from side to side and slowly pull the connector off the catheter port.
 - If the finger pressure is inadequate, use a rubber-tipped forceps to carefully apply pressure precisely on the oval marks. While squeezing, pull the connector off the catheter port. Be careful not to damage the connector.
 - Do not use additional methods to secure this sutureless pump connector.
6. Unclamp the pump segment.
 7. Coil the excess catheter behind the pump and place the pump in the pocket.
 8. Refer to the pump implant manual for additional pump implant instructions.

Caution: Do not coil excess catheter in front of the pump. Coil excess catheter behind the pump to:

- minimize potential damage during pump replacement surgery.
- minimize potential kinking of the catheter.
- prevent damage to the catheter during refill and catheter access procedures.

Calculating catheter length and volume

Warning: Record the exact implanted catheter length, the catheter model number, and the calculated catheter volume in the patient's medical record. In the event of a catheter revision, record the changes made to the catheter length and the calculated catheter volume in the patient's medical records. Catheter measurement information will be required in later calculations. **A universal value does not exist that can be used as a substitute for this knowledge.** Attempting to estimate catheter volume in later calculations can result in a clinically significant or fatal drug underdose or overdose.

Note: If this is a partial catheter replacement, take into consideration the remaining implanted catheter volume when calculating total catheter volume.

1. Measure the trimmed lengths (cm) of catheters that were saved in step 2 and possibly step 4 and step 8 of "Connecting and anchoring the catheter" beginning on page 11.
2. Record the combined trimmed lengths in the patient's record.
3. Subtract the trimmed length from the original total catheter length. (Refer to Table 1 on page 6.)
Original total catheter length
– trimmed (discarded) catheter length
= total implanted catheter length
4. Multiply the implanted catheter length by the catheter volume per cm.
5. Record this information in the patient's record.

Example

Original total catheter length:	104.1 cm
Trimmed (discarded) catheter length:	5.1 cm
Model 8731SC catheter volume per cm:	0.0022 mL/cm
(104.1 cm – 5.1 cm) × 0.0022 mL/cm	= 0.218 mL

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.
© Medtronic, Inc. 2008