# Medtronic

**SPECIFY™ 2X8**

**39286**

Lead kit

**SPECIFY™ 2X8**

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**SPECIFY™ 2X8**

Elektrodenset

**SPECIFY™ 2X8**

Kit de electrodo

**SPECIFY™ 2X8**

Geleidingsdraadset

**SPECIFY™ 2X8**

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**SPECIFY™ 2X8**

Ledningssett

**SPECIFY™ 2X8**

Elektrodset

**SPECIFY™ 2X8**

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Rx only

2009

## Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

LEAD

Lead length

Do not reuse

STERILE

EO

Sterilized using ethylene oxide

Caution, consult accompanying documents

Date of manufacture

Use by

Temperature limitation

LOT

Lot number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

! USA

For USA audiences only

Medtronic® is a registered trademark of Medtronic, Inc.  
Specify™ is a trademark of Medtronic, Inc.

## Table of contents

### **Description 5**

### **Package contents 5**

### **Device specifications 5**

### **Instructions for use: spinal cord stimulation 9**

Preparing for surgery 9

Placing a surgical lead 9

Testing stimulation intraoperatively 10

Implant procedure: connecting the lead to two extensions 12

Tunneling the extensions 12

Creating a pocket for the lead-extension connections 12

Connecting the extensions 12

Creating strain-relief loops and suturing the connector boots 13

Implant procedure: connecting the lead directly to a neurostimulator 14

Anchoring the lead 14

Tunneling the lead for direct connection to a neurostimulator 15

Creating strain-relief loops when connecting the lead directly to a neurostimulator 16

**Refer to the indications sheet for indications and related information.**

**Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.**

**Refer to System Eligibility, Battery Longevity, Specifications reference manual packaged with the software application card for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.**

**[USA] Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.**

## Description


The Medtronic Specify 2x8 Model 39286 Surgical Lead is part of a neurostimulation system for pain therapy.

## Package contents

- Lead
- Connector boots, 1x8 low-profile (2 white radiopaque, 2 clear)
- Anchors (65 cm lead package):
  - two-wing (2)
  - three-wing (2)
  - EZ (2)
- Passing elevator
- Lead blank
- Tunneling tools: tunneling rod, tunneling tip (narrow), tube (passing straw)
- Wrench, torque
- Product literature
- [\[USA\]](#) Warranty card
- [\[USA\]](#) Registration form

## Device specifications

The lead has 16 electrodes on the distal end; the two proximal ends fit into eight-conductor connectors (ie, two Medtronic in-line 1x8 extensions or a Medtronic neurostimulator). One lead body has a white marker band to indicate electrodes 0-7.

**Figure 1.** Specify 2x8 Model 39286 Lead (electrode side shown).

**Table 1.** Device specifications<sup>a</sup> for lead Model 39286

<b>Description</b>	<b>Model 39286</b>
Connector	Octapolar, in-line
Shape	Contoured
Conductor resistance <sup>b,c</sup>	<36 $\Omega$ (30 cm) <77 $\Omega$ (65 cm)
Length (includes paddle)	30 cm, 65 cm
Diameter	1.3 mm
Distal end	
Number of electrodes	16
Electrode shape	Rectangular
Electrode size (width x length)	1.5 mm x 4.0 mm
Electrode stimulating area	6.0 mm <sup>2</sup>
Electrode spacing (edge to edge)	
in-line spacing	1.0 mm
row spacing	1.0 mm
Electrode paddle length	56.4 mm
Electrode paddle width	7.6 mm
Electrode paddle thickness	1.9 mm
Proximal end	
Lead contact length	1.5 mm
Lead contact spacing (center to center)	2.8 mm
Lead contact distance	22.5 mm

<sup>a</sup> All measurements are approximate.

<sup>b</sup> Electrical resistance of this device only.

<sup>c</sup> Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

**Table 2. Material of components in the Model 39286 package**

<b>Component</b>	<b>Material</b>	<b>Material contacts human tissue</b>
<b>Lead</b>		
Conductor wire	MP35N	No
Conductor wire insulation	Fluoropolymer	No
Electrodes	Platinum-iridium	Yes
Electrode paddle	Silicone rubber	Yes
Paddle mesh	Polyester	No
Insulation	Polyurethane	Yes
Marker band	Polyurethane, titanium dioxide	Yes
Proximal connector	MP35N	Yes
Adhesive	Epoxy, silicone rubber	Yes
<b>Connector boots</b>		
Clear	Silicone rubber	Yes
White radiopaque	Silicone rubber, barium sulfate	Yes
<b>Anchors</b>		
Two-wing, three-wing, EZ	Silicone rubber	Yes
Passing elevator	Acetal resin, silicone rubber, barium sulfate	Yes
Lead blank	Silicone rubber, barium sulfate	Yes
<b>Tunneling tools</b>		
Tunneling rod	Stainless steel	Yes
Tunneling tip	Stainless steel	Yes
Tube (passing straw)	Fluoropolymer	Yes
<b>Wrench, torque</b>		
Handle	Polymer	Yes
Shaft	Stainless steel	Yes

## Instructions for use: spinal cord stimulation

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

**Note:** Implant this lead after the patient completes a test stimulation period to determine the effectiveness of stimulation in managing pain.

**Warning:** Evaluate the geometry and suitability of the spinal canal (by using an appropriate imaging method) to ensure that the spinal canal can safely accommodate insertion and permanent placement of the surgical lead paddle and body. The following factors should be considered during this evaluation:

- Adequate space within the spinal canal
- Pre-existing conditions or surgical procedures that may have altered the patient's epidural anatomy
- Potential for post-implant tissue changes around the lead
- Changes in epidural space associated with typical spinal movement

Inadequate space to accommodate the lead may result in localized compression of the spinal cord, increasing the risk of spinal cord injury.

**Warning:** Safety and effectiveness of this lead has not been established for use with more than one neurostimulator. If this lead is connected to more than one neurostimulator, damage to nerve tissue may result.

**Caution:** Implanting the surgical lead paddle in anatomical locations that have a high degree of mobility (such as the cervical spine) is not recommended. The stresses associated with these locations may reduce the functional survival time of the lead, requiring reprogramming and/or early replacement of the lead to restore effective therapy.

**Cautions:**

- Do not bend, kink, or stretch the lead or extension, which may damage the component.
- Use only rubber-tipped forceps on the lead. Do not use sharp-edged instruments (eg, hemostat), which may nick or cut the insulation.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

## Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead length, and connector type.

## Placing a surgical lead

Consider the following suggestions when performing the laminotomy:

- Perform a laminotomy 1-2 vertebral spaces inferior to the desired final location of the surgical lead paddle.
- Use electrocautery to dissect the site.
- Avoid subperiosteal stripping to maintain a clear operative field.
- Inject additional local anesthesia just prior to stripping the paraspinal muscles from the spine.
- Additional anesthesia should be liberally administered on patient demand.
- Make the laminotomy width consistent with the lead paddle's width.

**Caution:** Ensure that the laminotomy is wide enough to accommodate the surgical lead paddle. Failure to do so may result in forcing the lead into the epidural space, which may damage the lead.

- If the interlaminar space is narrow, use instruments to reflect the ligamentum flavum and remove enough bone for the lead blank insertion. (A rongeur may be required during this exposure.)

- Use bipolar forceps to grasp and coagulate the epidural fat so that the fat can be stripped away, exposing the dura mater.

**Note:** A hemilaminotomy is preferred over a laminotomy because surgical leads are designed to fit in the epidural space under the vertebral body.

**Warning:** Introduce the passing elevator, the lead blank, and the lead at a shallow angle. Using an angle that is too deep could cause contusion to the spinal cord.

1. Carefully introduce the epidural passing elevator.
  - a. Slowly and with minimal force introduce the passing elevator into the epidural space along the midline, guiding the entry so that the approach angle is as shallow as possible.
  - b. Remove the passing elevator.
2. Carefully introduce the radiopaque lead blank:
  - a. Introduce the lead blank into the epidural space to ensure proper size and location of the lead site.
  - b. Remove and discard the lead blank.

**Caution:** Use bone wax or instrumentation to ensure that the bony edge of the laminotomy is smooth. Sharp edges may damage the lead, resulting in intermittent or loss of stimulation.

3. Position the surgical lead:
  - a. Using a rubber-tipped forceps to handle the lead paddle, carefully position the lead paddle in the epidural space.

**Note:** Make sure the stimulating electrodes face the dura mater.

  - b. Advance the lead paddle cephalad until the entire paddle is in the epidural space.
  - c. For bilateral pain, place the lead close to the midline.
4. Verify the lead position under fluoroscopy (AP and lateral views).


### Testing stimulation intraoperatively

**Note:** This procedure requires a snap-lid connector cable, which is available as an accessory kit.

**Caution:** To prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation):

- program parameter changes in small increments above the perception threshold (the amplitude value(s) at which the patient first perceives paresthesia).
  - decrease the amplitude(s) to 0.0 V before:
    - changing electrode polarities.
    - connecting or disconnecting the cable to the external neurostimulator.
    - replacing the external neurostimulator batteries.
    - turning ON the neurostimulator.
1. Use the lid latch to open the snap-lid connector cable's snap-lid connector (Figure 2a).

**Figure 2.** Snap-lid connector cable.

- Align the proximal end of one of the lead bodies against the inside end of an 8-contact lead groove of the connector (Figure 2a).
- Check that the lead-body contacts line up with the electrical contacts inside the lead groove.

**Note:** The lead bodies fit only one way into the snap-lid connector.

- Press the lead body gently into the lead groove (Figure 2b). Repeat steps 2-4 for the second lead body.

**Note:** Connect the lead body with the marker band to the 0-7 lead groove. The other lead body connects to the 8-15 lead groove.

- Close the lid until the latch snaps shut in order to seat the lead bodies in the lead grooves.

**Note:** Do not force the lid closed; it should close easily. If it does not, disassemble the components and repeat steps 1-5 in this section.

- After connecting the cable to the patient's implanted lead, pass the plug end out of the sterile field.
- After verifying that the external neurostimulator is turned off, connect the cable plug to the receptacle of the external neurostimulator.

**Warning:** Maintain adequate slack in the cable. If there is not enough slack and the cable is pulled, the implanted components may be dislodged.

- Identify optimal stimulation parameters, beginning at a pulse width of 210-240  $\mu$ s and a rate of 30 Hz.

**Note:** Ensure that the patient can provide immediate feedback.

- Increase the amplitude beginning from 0.0 V while asking the patient closed-ended questions to identify the perception threshold (the amplitude at which the patient first perceives paresthesia), the discomfort threshold (the amplitude at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

**Note:** If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.

- In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).
- Disconnect the cable from the lead by opening the snap-lid connector and gently lifting the lead bodies from the electrical contacts in the lead grooves. Leave the cable in the sterile field for additional parameter testing before closing.

**Caution:** Do not pull on the cable. Pulling on the cable may break a wire or dislodge the lead. A broken wire or dislodged lead may result in loss of stimulation and may require surgery to replace the lead.

## Implant procedure: connecting the lead to two extensions

Following intraoperative test stimulation, use this procedure to connect the surgical lead to two Medtronic in-line 1x8 extensions. Extensions are required for shorter lead lengths (ie, 30 cm) to reach the neurostimulator implant site, and the connector boots are used as the anchor points to secure the lead. Extensions are optional for longer leads (ie, 65 cm). Prior to connecting the extensions to a 65 cm lead, tunnel the 65 cm lead bodies (refer to "Tunneling the lead" within "Tunneling the lead for direct connection to a neurostimulator" on page 15 for instructions).

If connecting the surgical lead directly to a Medtronic neurostimulator, refer to "Implant procedure: connecting the lead directly to a neurostimulator" on page 14.

This procedure includes the following sections:

- Tunneling the extensions
- Creating a pocket for the lead-extension connections
- Connecting the extensions
- Creating strain-relief loops and suturing the connector boots

### Tunneling the extensions

For instructions on tunneling the extensions, refer to the appropriate extension manual.

**Note:** This procedure requires two Medtronic in-line 1x8 extensions, which are available in a separate kit.

### Creating a pocket for the lead-extension connections

For instructions on creating a pocket for the lead-extension connections, refer to the appropriate extension manual.

### Connecting the extensions

#### Cautions:

- Do not use saline or other ionic fluids at connections, which could result in a short circuit.
  - Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.
1. Slip a connector boot onto each lead body with the open end of the connector boots facing out, exposing the lead contacts.

#### Notes:

- Irrigate the connector boots with a nonionic antibiotic solution. Dry the interior of the connector boots.
  - A clear and a white radiopaque boot are provided to distinguish between the lead bodies.
2. Wipe the lead bodies and the extension setscrew connector junctions with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution. Dry all connections.
  3. Insert each lead connector fully into an extension setscrew connector junction.

#### Notes:

- Each lead contact must be aligned under each extension connector contact.
  - During insertion, some resistance is typical.
4. Use the torque wrench to tighten the setscrews to complete the electrical circuit.

#### Cautions:

- To avoid overtightening, do not use a hex wrench to tighten an extension setscrew. Overtightening an extension setscrew may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.

- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.
5. Using minimal force, pull on the connections to ensure they are secure.
  6. Confirm the integrity of the connections by rechecking the stimulation pattern.
  7. Slide the connector boots over the lead-extension connections and tie a ligature (using nonabsorbable silk) between the pairs of rings at both ends of each connector boot.

**Cautions:**


- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
- Do not tie ligatures around the lead or extension body, which may damage the insulation.
- Do not overtighten ligatures on the connector boot. Ligatures that are too tight may damage the component.
- Ensure that the end of the connector boot is secure to prevent fluid in the lead-extension connection, which may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

**Creating strain-relief loops and suturing the connector boots**

**Caution:** Use strain-relief loops as directed when implanting surgical leads. Implanting surgical leads without strain-relief loops may reduce the functional survival time of the lead, requiring reprogramming and/or early replacement of the lead to restore effective therapy.

1. While maintaining lead position, create a strain-relief loop with each of the lead bodies, by coiling the lead bodies near the connector boots (Figure 3).


**Figure 3.** Strain-relief loops made with the lead bodies.

**Cautions:**

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 4). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.

**Figure 4.** Coiling excess lead.

2. Insert the strain-relief loops into the lead-extension connection pocket, under the connection, leaving as much slack as possible in the lead bodies between the lead paddle and the lead-extension connections.
3. Suture each connector boot (using nonabsorbable silk) to the supraspinous ligament or to the deep fascia.

**Notes:**

- Do not suture the narrow end of the connector boots or anywhere on the lead bodies.
  - When using a longer lead, suturing the connector boots is at the clinician's discretion.
4. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.
  5. Refer to the neurostimulator implant manual for instructions on making the connection to the neurostimulator.

### Implant procedure: connecting the lead directly to a neurostimulator

Following intraoperative test stimulation, use this procedure to connect the surgical lead directly to a Medtronic neurostimulator. If using extensions to implant the lead, refer to "Implant procedure: connecting the lead to two extensions" on page 12.

This procedure includes the following sections:

- Anchoring the lead
- Tunneling the lead for direct connection to a neurostimulator
- Creating strain-relief loops when connecting the lead directly to a neurostimulator

#### Anchoring the lead

**Note:** The use of anchors is at the clinician's discretion. Proper use of anchors requires the use of strain-relief loops in order to reduce the stresses placed on the lead.

**Caution:** Do not tie ligatures around the lead or extension body, which may damage the insulation.

1. While maintaining lead position, slide the anchor onto the proximal end of each lead body and continue sliding the anchor down as close as possible to where the lead emerges from the vertebral column. (Leave slack in the lead bodies between the lead paddle and the anchors to allow for strain-relief loops.) Use care to maintain lead position.

**Notes:**

- Each lead body should be anchored separately.
  - The anchor can be lubricated with sterile water.
  - The EZ anchor can be trimmed to a shorter length.
2. Use 2-0 nonabsorbable suture (such as silk or some types of braided polyester mesh) to secure the anchor.

△ **Cautions:**

- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
- Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.

**For the EZ anchor:** Make a ligature, pass the suture around the anchor and make another ligature. Use a minimum of 4 sutures. Suture the anchor to the supraspinous ligament or deep fascia.

**For the winged anchor:** Tie ligatures around the grooves in the middle of the anchor to secure the anchor to the lead. Next, suture the anchor to the supraspinous ligament or deep fascia.

**Note:** Laboratory testing has shown that injecting silicone medical adhesive between the anchor and lead after tying the ligatures increases the anchor-to-lead holding force.

3. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

**Tunneling the lead for direct connection to a neurostimulator**

**Creating a subcutaneous pocket for the neurostimulator**

1. Make an incision equal to the length of the neurostimulator.
2. Use blunt dissection to create a subcutaneous pocket.

**Note:** Refer to the neurostimulator implant manual for the proper subcutaneous pocket depth.


**Tunneling the lead**

**Note:** If a longer tunneling tool is desired, a Medtronic tunneling tool accessory kit may be used.

1. After simulating the tunneling route, mark the patient's skin at the neurostimulator pocket. (The neurostimulator should be placed on the opposite side of the body from another active implanted device and should be placed preferably on the right side of the body to allow for future placement of cardiac devices on the patient's left side.)

△ **Cautions:**

- When routing the lead, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When routing two lead bodies, assure that the area between them is minimized (Figure 5). If the lead bodies are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).

**Figure 5.** Routing for two lead bodies.

2. After assembling the tunneling tools, tunnel subcutaneously between the lead incision site and the neurostimulator pocket.

△ **Caution:** Proceed slowly when the tunneling tool approaches the neurostimulator pocket. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.


**Notes:**

- Deep tunneling is not desirable.
  - Avoid the lower thoracic ribs.
  - If the tunneling tool does not extend between the lead incision site and the neurostimulator pocket, make an intermediate incision.
3. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.
  4. Gently insert the proximal end of both lead bodies through the passing straw.
  5. After carefully removing the passing straw, verify that the location of the lead has not changed.

**Creating strain-relief loops when connecting the lead directly to a neurostimulator**

△ **Caution:** Use strain-relief loops as directed when implanting surgical leads. Implanting surgical leads without strain-relief loops may reduce the functional survival time of the lead, requiring reprogramming and/or early replacement of the lead to restore effective therapy.


1. While maintaining lead position, create a strain-relief loop with each of the lead bodies, by coiling the lead bodies near the lead paddle (Figure 6).

**Figure 6.** Strain-relief loops made with the lead bodies (no connector boots).

△ **Cautions:**

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 7). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.

**Figure 7.** Coiling excess lead.

**Note:** Do not suture anywhere on the lead bodies.

2. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.
3. Refer to the neurostimulator implant manual for instructions on making the connection to the neurostimulator.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.
© Medtronic, Inc. 2008