Vessix[™] Catheter

Renal Denervation Catheter

Directions for Use	2
Instrucciones de uso	6
Mode d'emploi	10
Gebrauchsanweisung	14
Istruzioni per l'uso	18
Gebruiksaanwijzing	22
Instruções de Utilização	26

Vessix™ Catheter

Renal Denervation Catheter

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

This device should only be used by physicians trained in angiography and percutaneous interventional techniques.

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED IN THESE DIRECTIONS AND IN THE VESSIX™ GENERATOR OPERATOR'S MANUAL. FAILURE TO DO SO MAY RESULT IN SERIOUS COMPLICATIONS.

Refer to the instructions supplied with the interventional devices used in conjunction with the Vessix Catheter. (Figure 1 below)

WARNING

The Vessix Catheter is designed and intended to be used ONLY with the Vessix Generator

- Contents supplied STERILE using a Radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- No modification of this equipment is allowed.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- The system must be used under fluoroscopic observation with radiographic equipment that provides high quality images.
- The device must be used by trained physicians skilled in percutaneous, interventional techniques in a fully equipped catheterization laboratory.
- DO NOT retract the Vessix Catheter unless the balloon is fully deflated under vacuum. If significant resistance is met during manipulation, determine the cause of the resistance before proceeding. If excessive resistance is met during withdrawal of the Vessix Catheter, use fluoroscopy to verify the resistance, then remove the entire system simultaneously (catheter, guide wire, guiding sheath) using fluoroscopic guidance.
- To avoid the risk of patient injury or catheter failure, do not deliver more than five (5) treatment cycles combined per catheter. If more treatments are needed, use additional catheters.
- Use only the recommended balloon inflation medium (normal heparinized saline and contrast at ratios recommended in this instruction). Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the "Use By" (expiration) date specified on the package.
- Balloon pressure should not exceed the maximum inflation pressure of 5.0 atm.

DEVICE DESCRIPTION

The Vessix™ System consists of the Vessix Catheter, a non-compliant balloon catheter with electrical contacts mounted on the exterior of the balloon, the Vessix Generator, and connecting cables.

The Vessix Catheter is available in balloon diameters of **4.0 mm, 5.0 mm, 6.0 mm and 7.0 mm** with a balloon length of 25 mm and treatment length of 21 mm. The catheter has a total working length of approximately 90 cm. It is recommended that the Vessix Catheter be connected to the Vessix Generator by means of the non-sterile patient extension cable, which is approximately 3 meters in length. Optionally, the Vessix Catheter may be connected directly to the Vessix Generator.

The Vessix Catheter is inflated at a very low pressure until electrodes are properly apposed against the renal artery wall. Controlled low-powered radiofrequency (RF) energy is then delivered briefly to the electrodes by the Vessix Generator. This delivery of energy causes thermal heat to perfuse through the artery wall into the adventitia layer of the artery and results in denervation of the target renal nerves. During the treatment the RF generator delivers low average power at each electrode. The RF generator provides visible prompts to guide the operator through setup and treatment steps.


Figure 1

Contents

Qty Material

1 Vessix™ Catheter

INTENDED USE/INDICATIONS FOR USE

The Vessix Catheter is intended to be used to treat patients with uncontrolled hypertension.

CONTRAINDICATIONS

Contraindications to include those known but not limited to those known for conventional percutaneous interventional procedures.

- The Vessix Catheter is NOT intended for use in any artery other than the renal artery.
- The Vessix Catheter is NOT intended for use in a renal artery that has an implant, nor is it intended to be used in calcified renal arteries.

PRECAUTIONS

- Store in a cool, dry, dark place.
- Protect from direct sunlight and temperatures greater than 104 °F (40 °C) during storage.
- The device should only be used with 0.014 inch or 0.018 inch guide wires.
- Follow instruction for use supplied with the guiding sheath. A thorough understanding of the principles, clinical applications and risks associated with standard percutaneous procedures is necessary before using this product.
- Minimize torqueing of catheter shaft to avoid potential damage to the catheter.
- Do not reinsert the catheter once it has been removed from the sheath.
- If resistance is encountered during catheter insertion into sheath, do not force passage. Resistance may result in damage to the balloon circuitry if it is forced. If significant resistance is encountered, remove the catheter.

ADVERSE EVENTS

The following adverse events could possibly occur as a result of the use of the catheter, but are not limited to:

- Renal artery dissection or perforation
- Renal artery spasm
- Hypotension
- Acute closure of renal artery
- Thrombus in renal artery
- Stenosis in the renal artery
- Hematoma at access site
- Embolization
- Pyrogenic reaction
- Kidney malfunction
- Access site infection
- Hemorrhage
- Pain

HOW SUPPLIED

Sterile: The Vessix Catheter is supplied sterile. The device is sterilized by electron-beam radiation.

HANDLING AND STORAGE

Temperature:

- Operating temperature: 10 °C (50 °F) - 40 °C (104 °F)
- Storage temperature: 10 °C (50 °F) - 40 °C (104 °F)
- Transport temperature: -10 °C (14 °F) - 50 °C (122 °F)

Humidity:

- 30% - 75% RH, non-condensing

OPERATIONAL INSTRUCTIONS

A. Additional Materials and Equipment

- Guiding sheath in the appropriate size and configuration to access the renal artery
- Tuohy-Borst "y" adapter or hemostatic valve(s)
- Normal heparinized saline
- Standard contrast media
- Guide wires

B. Compatibility

- **Only use the Vessix™ Generator. Do not use any other RF generator.**
- Catheter Rated Voltage: 100 Vrms.
- Vessix Catheters must be used with either a 0.014 inch or 0.018 inch guide wire.
- All Vessix Catheter sizes, 4.0 mm, 5.0 mm, 6.0 mm, and 7.0 mm are recommended to be used with the following guiding sheath:
 - Cook Medical Flexor® Ansel Check-Flo® 8F Guiding Sheath with inner diameter of 0.113 inch (2.87 mm), Model # G49100 (KCFW-8.0-18/38-45-RB-ANL1-HC).

- The catheter lumen must be flushed prior to treatment with heparinized saline.

C. Renal Artery Access

1. Use a femoral puncture technique to insert an appropriately sized (5 or 6F) introducer sheath. Heparinize intravenously according to the standard practices of the institution.
2. Perform a diagnostic aortogram to obtain visual confirmation of all renal arteries. A selective arteriogram should be performed on each renal artery using a diagnostic catheter in order to determine the appropriate Vessix™ Catheter to be used for the treatment, per the table below:

Recommended Catheter Sizing Table				
Artery Size ->	3.0 – 4.0 mm	3.8 – 5.0 mm	4.7 – 6.0 mm	5.6 – 7.0 mm
4.0 mm balloon	Yes	x	x	x
5.0 mm balloon	x	Yes	x	x
6.0 mm balloon	x	x	Yes	x
7.0 mm balloon	x	x	x	Yes

3. The short sheath should be exchanged for an 8F guiding sheath.
4. The 8F guiding sheath should be placed near the ostium of the renal artery using a no-touch technique. Do not cannulate the renal ostium with the 8F guiding sheath.
5. Introduce a 0.014 inch or 0.018 inch guide wire into the renal vasculature via the guiding sheath. The wire should be placed beyond the bifurcation. Once access is secured, the Vessix Catheter should be prepared for treatment.

D. Preparation for Use

1. Using sterile techniques, open the sterile package, unfasten the electrical connector, free the lead, and then gently remove the catheter from the catheter hoop.
2. Slide the protective sheath off the balloon and inspect the device to ensure it has not been damaged. Remove the stainless steel shipping mandrel from the guide wire lumen.
3. Connect the catheter to the generator, which resides outside of the sterile field. The catheter-to-patient extension cable connection should remain outside of the sterile field.
4. Upon catheter to generator joining, the generator screen will confirm the catheter has been connected. (Figure 2 below)
5. Do not attempt to pull or reshape the electrical contacts mounted on the exterior of the balloon. Flush the guide wire lumen with normal heparinized saline in a routine manner. Care should be taken not to disrupt the circuitry on the balloon.
6. Attach a 10-20 cc syringe containing 2 parts normal heparinized saline to 1 part contrast to the inflation port of the catheter hub. Apply negative pressure for 20-30 seconds.
7. Release pressure slowly, allowing negative pressure to draw saline into the balloon lumen.
8. It is recommended to apply negative pressure again for 20-30 seconds prior to detaching the syringe.
9. Detach the syringe, leaving a meniscus of saline in the inflation port of the hub.
10. Prepare the inflation device with 2 parts normal heparinized saline to 1 part contrast, and purge to remove all air from the inflation device.
11. Attach the inflation device to the inflation port directly, ensuring no air bubbles remain at the connection.
12. Leave the inflation device on ambient pressure (neutral position). Do not pull negative pressure (vacuum) on the inflation device after balloon preparation or prior to delivering the catheter.

Figure 2

E. Placement of Vessix Catheter

1. Maintain neutral pressure on the inflation device. Insert the catheter over the guide wire through the sheath. If resistance is encountered, do not force passage. Resistance may indicate a problem and may result in damage to the balloon circuitry if it is forced. If significant resistance is encountered, remove and replace the catheter.
2. Advance the catheter over the guide wire until the balloon is fully inside the renal artery, and the distal electrode is located 3 - 5 mm proximally to the first bifurcation, using direct fluoroscopic visualization. Do not torque the catheter during insertion or placement.
3. Utilize the proximal and distal electrodes on the balloon as a reference to optimally position the catheter.

F. Vessix™ Catheter Activation and RF Treatment

1. Once connection to the generator is confirmed, the Vessix Catheter can be inflated moving in 0.5 atm increments up to a **pressure of 3.0 atm**. When proper apposition is confirmed by the generator upon balloon inflation, fluoroscopy should be used to demonstrate that the catheter is at the optimal position and that blood flow is occluded. Discontinue use of the catheter if excessive pressure is required to inflate the balloon.
2. Perform the denervation treatment. Follow the Vessix™ Generator Operator's Manual regarding the treatment procedure.
3. Upon completion of the treatment, deflate slowly and retract the balloon proximally towards the ostium of the renal artery.
4. If a second (full or partial) treatment in the same artery is possible and desired, use the distal electrode as a marker to place the balloon 5 mm proximal to the first treatment, so that there is no overlapping of treatment areas. (Note: To allow for a partial treatment it is acceptable for the balloon to be located partially outside of the renal artery in the aorta). Repeat steps F1 and F2 in order to perform a second treatment in the same artery.
5. When treatment is complete, deflate slowly and carefully retract the balloon back into the guiding sheath. If resistance is encountered, advance the balloon 5 to 10 mm back into the renal artery, and carefully retract the balloon again into the sheath.
6. Perform post procedural imaging on each renal artery treated.
7. Treat the contralateral renal artery at this time if the catheter is appropriately-sized. If treatment of the contralateral renal artery requires a different balloon size, remove the catheter (per section G), and insert an appropriately sized catheter (per sections D through F).

G. Vessix Catheter Deflation and Removal

1. Slowly deflate the balloon by 0.5 atm decrements and finally pull to negative pressure.
2. Carefully retract the catheter through the guiding sheath and out of the body. The guide wire should remain in the renal artery for a post-procedure selective renal arteriogram. If excessive resistance is met, using fluoroscopy carefully remove the entire system simultaneously (catheter, guide wire, and guiding sheath).
3. Upon completion, either manual compression or commercialized closure devices can be used to achieve hemostasis at the puncture site. Standard-of-care post-intervention monitoring procedures should be followed.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Cook Medical Flexor® Ansel Check-Flo® is a trademark of Cook Medical Incorporated

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or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd
Vessix Vascular, Inc.
© 2013 Boston Scientific Corporation