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EN Instructions for Use of Fistula needles and Safety Fistula Needles

Please read these instructions and warnings carefully before using the product.

Indications

Fistula needles and fistula needles with safety function are indicated for haemodialysis, haemofiltration and haemodiafiltration or other treatments requiring an extracorporeal blood circuit or vascular access when larger volumes of blood should be drawn from to the patient's circulation.

Contraindications

The use of a fistula needle and puncturing a vessel or a graft is in itself not a therapeutic procedure. There are no contraindications apart from risks generally associated with punctures by large bore needles in conjunction with extracorporeal blood circuits or due to known or suspected allergies against the materials used in the product. However, contraindications may be inherent in the procedures requiring the access to the blood. Therefore the use of this device is restricted to an order by a physician and should be used following the physician's instructions.

Warnings and precautions

1. Do not use if package is damaged or contaminated, protective caps are fallen off or missed or the bevel is unusually shaped or deformed. The blister packaging only assures sterility of the blood pathway.
2. Single use product. Do not reuse. Reuse may lead to infection or pyrogenic reactions.
3. Do not use adverse chemicals such as alcohol, acetone, ionic or non-ionic surfactants (e. g. quaternary ammonium compounds) or Freon as propellant. Povidone iodine is recommended as antiseptic solution.
4. Do not touch the cannula or the end surface of the Luer connector and use aseptic technique when preparing for venipuncture and during treatment. Adequate cleansing and disinfection of the vascular access area must be performed to avoid contamination during the needle insertion process.
5. Use an aseptic technique when connecting the needle to the bloodline. Lock the Luer connectors firmly and make sure all connections are visually inspected in intervals to detect leaks or air embolism. Connections should be visible all the time during treatment and checked carefully from time to time during the operation.
6. Luer connectors comply with ISO 594-1/2 standards. Safe and tight connection can only be assured when the matching Luer fitting is also in compliance with these standards. In order to avoid leakage or separation of the lines, taping assures connection. Using non-compliant adapters or excessive force may lead to stress cracks or leaks which are associated with risks of air embolism or blood loss.
7. The use of an air detector is mandatory to avoid air embolism. Do not use air for blood return after termination of the dialysis session.
8. Avoid kinking, twisting or occluding of the tubing as excessive pressure may damage the extracorporeal blood circuit.
9. Avoid needle sticks and direct contact with patient's blood. As exposure to blood or blood by-products carries the risk of transfer of infectious diseases, such as hepatitis or HIV, personal protective equipment such as gloves, gowns, masks and eye shields should be used routinely. Needles with safety function are recommended. Protect others against infective diseases by following adequate disposal procedures.
10. When using fistula needles with safety function always keep fingers behind front/back notch. Do not attempt to override the locking mechanism. Once retracted, do not push hub or tubing in the direction of the needle point. After use, hold the needle with point upward to prevent blood dripping and immediately discard into the designated container.
11. Do not unlock the device prior to cannulation because it is difficult to lock back if it is not centred properly. Visual inspection of locking mechanism should be done prior to cannulation. Tug slightly on the tubing to ensure that the lock is engage.
12. Buttonhole Needles feature a blunt tip. Do not use Buttonhole Needles for conventional needling procedures and make sure buttonhole formation is completed before use.
13. This product contains di-(2-ethylhexyl)phthalate (DEHP). Use for children, pregnant or nursing women only after careful risk evaluation to physician's instructions.

Mode of use

1. After thorough washing of the limb with soap and water the puncture area is carefully disinfected.
2. When the product is removed from the blister pack, it should be visually inspected for gross damage or particles in the lumen. After detaching the needle protector the bevel should be examined for visible deformity. The Luer cap should be tightened if loose. The needle may be flushed with saline or heparin solution. However, for puncture the needle should be air filled.
3. At a patient's request local anaesthesia may be used. With a tourniquet applied proximally a suitable puncture site is selected. In patients with limited area for punctures where two needles must be placed in close proximity, one or both wings can be easily removed by tearing along the perforation near the hub.
4. The markers on the hub reveal the bevel position. A black dot indicates the bevel up position for conventional puncture technique. A red dot indicates the reverse position for bevel down technique.
5. For grip during puncture both wings are folded over the hub. By its indentations they will form a compact handle to be grasped between thumb and index finger. In folded position the wings automatically block the hub and prevent rotation during puncture.
6. It is recommended to place the arterial needle first and then the venous needle. Each needle is inserted with clamp open but Luer cap tightened and at an angle of 25° to 45° to the skin. To prevent recirculation, the

bevel of the venous needle must always be 15 mm downstream from the arterial needle.

7. In inferior fistula conditions, rotation of the bevel or gradual withdrawal of the needle may lead to a better position of the bevel. The greatest amplitude of pulsation indicates the best position in the vessel. In fistulae with minimal pulsations in the tubing, correct bevel localisation should be checked by aspiration of blood by a 30 ml syringe. Aspiration and return of blood should be accomplished without resistance.
8. Place gauze underneath a safety needle if a steep angle is observed after cannulation. This is to reduce skin tensions and to prevent possible infiltration.
9. After being placed in the best position, the needles are secured to the skin. In the most commonly used method, tape is placed under tubing and hub, then over the hub, wing and skin on the opposite side at a 45° angle to the direction of the needle. The other end of the tape is crossed and secured in a similar fashion. A second tape is placed across the wings in a right angle to the needle.
In another method, tape is placed transversely under the tubing and over the rear edges of the wings. Two other stripes of tape are then placed longitudinally over the wings. This technique provides greater security for those who may be restless during treatment. Using any taping technique the needle entry site should remain visible.
10. Avoid taping on patient's clothing and covering of access site by pillow or blankets that may result in entanglement with the device, leading to dislodged needle or loose connection.
11. All air in the product must be eliminated before connecting it to the bloodline. This is accomplished by gradually loosening the Luer cap until the air is released in a controlled manner. The purged tubing is then occluded by the integrated tubing clamp. The Luer cap is removed and the bloodline is connected as prescribed. Connecting the arterial bloodline with the venous needle and vice versa should be avoided.
12. In rare case, cutting of small vessels in the subcutaneous tissue during puncture or reduced elasticity of the fistula vessel wall may cause bleeding from the puncture site with the needle in place. In such cases one can occlude the skin entry site by the cone-shaped hub of the needle. The needle is advanced until the hub compresses the skin exit site and is fixed by tape in this position.
13. When using a safety needle, secure the wing sheath at the access site prior to opening the external lock. Release the external lock first and then retract the needle completely by sliding the hub within the wing sheath, encapsulating the entire needle length while applying haemostatic pressure at the puncture site. Do not push the sliding hub towards the access site again.
14. An audible clicking sound is heard or locking force is felt when the needle is fully retracted into the final lock position.

Storage and shelf life

The product should be stored at 5°C to 30°C and 20 % r. H. to 80 % r. H. The product can be used until the expiry date if stored with triple protection in original cartons. Do not use if expiry date has lapsed. First-in / First-out principle should be observed.

Product disposal

Be careful after returning the blood to the patient and the end of the dialysis session when removing the needle. Promptly dispose the needles to a sharps container or when using a safety needle, activate the safety function before disposing the needle. In order to do so, clamp the tubing close to the hub before removal. Finally dispose the product according to internal or national regulations.

REACH information

This product contains DEHP in concentrations > 0.1 mass % according to article 33 and 59 (1, 10) of regulation (EG) Nr. 1907/2006 ("REACH")

Limitation of liability

The manufacturer shall not be liable for any misuse, improper handling, noncompliance with warnings and instructions, damage arising from results after the manufacturers release of the product, or failure or omission to inspect the product before use in order to insure that they are in proper condition, or any warranty given by independent distributors or dealers.

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Bionic Medizintechnik GmbH