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Instructions for use Renodrain Nephrostomy Catheter Sets with and without PC-coating

GANL Nr 06 Renodrain EN_12/30.04.2012

Please note: this medical product may be obtained and used only by trained medical personnel.

1. Description
RENODRAIN Puncture and Exchange-Sets with balloon-catheter or pigtail catheter for percutaneous nephrostomy. Catheter and puncture needles are sold separately, too.
Pigtail nephrostomy catheters are sold with Phosphorylcholine (PC) polymer coating on both inner and outer surface (Renodrain Yellow) or without coating (Renodrain).

2. Contents and packaging
- Nephrostomy catheter (biocompatible material PUR or 100% silicone).
- 2 dilators
- 1 dilator with split sheath (not in sets with pigtail FR/Ch 6 to 10)
- Guide wire, PTFE-coated
- Universal adapter
- Catheter plug
- 2-part or 3-part puncture needle (only in Puncture-Set)
- Double packaging and instructions for use
Sterilized with ethylene oxide.

3. Intended use
Note for user: Before the indication demands using the percutaneous nephrostomy technique, it should be examined if it is possible to remove the urinary retention by using a ureteral stent.

Diagnostic indications:
- Antegrade pyelography and ureteroscopy
- Withdrawal material for examination (histology, cytology, microbiology)
- Urodynamic measuring in the superior urinary tract
- Determination of renal function
Therapeutic indications:
- Decompression at urinary tract obstruction
- Urinary drainage at ureteral fistula
- Drug application (chemolitholysis, cytostatic agents)
- Methodes for antegrade manipulation in the superior part of the ureter (ureterembolisation, PNL, laser-radiation, ureterotomy)

4. Contraindications
- Coagulation defect
- Inflammation or dermatological disease in the puncture area
- Malformatted kidney and heavy spinal curvature
Contraindications arise from the constrictions of the corresponding basis application.

5. Side effects
Bleeding, vascular lesion, loss-parenchyma, haematuria, infection and pain in the puncture area, pain be connected with lack of urine production, pararenal extravasat formation, bacteraemia.

6. Instructions for use
General hygienic rules hast to be considered during the application of this device

Preparation, bedding and anaesthesia of the patient:
- Controlled stationary conditions, sterile conditions including washing off and covering, treated like an open operation or PNL
- Patient has to lie in prone position on the table.
- Antegrade PNS only with local infiltration anaesthesia, retrograde PNS or percutaneous operation depends from the patient's maximum resilience and the surgeons experience only with regional anaesthesia or intubation anaesthesia.

Technical performance:
- Localization of the renal pelvis under ultrasonic and radiological control.
- Afterwards precise initial puncture at the posterior axillary-line below the costal arch in direction of the inferior lateral kidney-circumference.
- Insertion of the guide wire through the needle into the renal pelvis. Then remove the needle.
- Dilation of the puncture channel over the wire with dilators in ascending size.

- Positioning the drainage catheter over the wire.
- Examination of the catheters position with contrast medium.
- **Catheter with pigtail:** Remove the guide wire in the case of an exact placement of the catheter. The distal catheter end forms itself into originally pigtail-design (not balloon catheter).

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- **Catheter with balloon:** Fill the balloon with 1.5-5 ml of sterile demineralised water. Assure correct placement of the catheter by gentle pulling. Remove guide wire not till then! After insertion of the catheter the funnel shall be connected with the counterpart of a urine bag and fixed.
- **Catheter without balloon:** Secure catheter by using sutures.
Exchange of the catheter:
- Insert guide wire with flexible tip first into established catheter.
- Deflate balloon (balloon-catheter) or remove fixation of catheter (catheter without balloon).
- Remove catheter via guide wire from the renal pelvis.
- Dilation of the puncture channel over the wire.
- Place new catheter over the guide wire into the renal pelvis. Assure correct placement of the catheter by gentle pulling. Remove guide wire not till then!
- During a long-term drainage a routine evaluation of the patient and catheter for correct catheter placement, infection and encrustation is required!

Magnusson Locking Pigtail
- After placing the catheter remove wire slowly and pull the movable part of the Magnusson Pigtail System simultaneously towards the fixed part of the connector. Lock both parts with a turn. Thus the pigtail is safely secured in the kidney.
- Attention: To secure the thread close to the catheter it is recommended to turn the movable part 360° before locking.
- To release the pigtail formation unlock the two parts due to a turn to the right and pull the movable part of the connector towards the body.

Removal of the catheter with balloon
Before removal the balloon shall be deflated completely by connecting an empty syringe at the valve and sucking until a vacuum is built.

7. Possible complications and/or risks
- Puncture of vessels in the renal pelvis wall may cause haematuria.
- Mispuncture.
- Incrustations, blood clots and kinks of the catheter could lead to partial or total obstruction of the catheter with the danger of an infected urinary congestion renal pelvis or an urosepsis.
- Dislocation of the catheter
- After removing of the catheter urine can outflow from the nephrostomy channel.
- Do not treat the products with disinfectant agents, since no information is available regarding possible material changes that this may cause.

8. Warnings / Precautions
- This product is only sterile if the packaging is undamaged or unopened. For single use only!
- All components must be carefully checked before use. Products that may be damaged must not be used.
- This product should be used only by physicians experienced in the use of such medical devices.
- Do not overflate balloon. Please inflate only to the maximum volume as indicated on catheter base or label.
- Before refilling balloon, empty completely and inflate to the maximum volume as indicated. Inflate balloon with little pressure. To reduce diffusion with silicone catheters, we recommend a solution of distilled water + 10% glycerin for inflation.
- Lubricate catheter only with water-soluble liquids. Do not use lubrication containing petroleum oil.
- In case of the use of a suture to fix the pigtail, the catheter must be removed latest after 30 days.
- Balloon catheters and pigtail catheter without PC-coating must be changed latest after 30 days.
- In exceptional cases the balloon will not unlock. In this case, the balloon can be deflated with several measures like cutting the valve, overflating the balloon or puncturing the balloon via the supraubric route. Which measure will be taken is under the responsibility of the applicant under consideration of the individual situation.

- Regular inspections of free passage and position of the catheter is required. Balloon-Catheter may change the balloon filling quickly. Filling should be checked every 1-2 weeks.

9. Interaction with other drugs
It is possible that the product may discolour due to the concomitant use of medications or even as a result of contact with urine. UROTECH has not carried out any tests in this regard; however our many years of experience have documented no indications that this impairs the product's performance in any way.

10. Disposal
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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Medizinische Technologie GmbH