RELEASED

PPL2104102 (2/25/15)

0086

# **XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL**

**Everolimus Eluting Coronary Stent System**

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RELEASED

# XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL

## Everolimus Eluting Coronary Stent System

### INFORMATION FOR PRESCRIBERS

#### Table of Contents

- 1.0 DEVICE DESCRIPTION
  - Table 1-1: Product Name and Sizes
  - Table 1-2: Drug Content in the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stents
  - Table 1-3: *In vitro* Device Specifications
- 2.0 HOW SUPPLIED
- 3.0 INDICATIONS
- 4.0 CONTRAINDICATIONS
- 5.0 WARNINGS
- 6.0 PRECAUTIONS
  - 6.1 Stent Handling – Precautions
  - 6.2 Stent Placement – Precautions
  - 6.3 Use in Conjunction with Other Procedures
  - 6.4 Stent / System Removal – Precautions
  - 6.5 Post Implantation – Precautions
  - 6.6 MRI Statement
  - 6.7 Drug Interactions
  - 6.8 Pregnancy

### 1.0 DEVICE DESCRIPTION

The XIENCE Xpedition, XIENCE Xpedition SV (Small Vessel) and XIENCE Xpedition LL (Long Lesion) Everolimus Eluting Coronary Stent System includes:

A pre-mounted L-605 cobalt chromium (CoCr) alloy XIENCE Xpedition stent with a coating that consists of a blend of the anti-proliferative drug everolimus and polymers. The product family consists of:

Table 1-1: Product Name and Sizes

Product Name	Stent Diameter (mm)	Stent Length (mm)
XIENCE Xpedition SV	2.0, 2.25	8, 12, 15, 18, 23, 28
XIENCE Xpedition	2.5, 2.75, 3.0, 3.25, 3.5, 4.0	8, 12, 15, 18, 23, 28
XIENCE Xpedition LL	2.5, 2.75, 3.0, 3.25, 3.5, 4.0	33, 38

The available dose of everolimus on the stent varies by size as follows:

Table 1-2: Drug Content in the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stents

Stent Diameter (mm)	Stent Length (mm)	Drug Dose (µg)
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	8	40
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	12	60
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	15	74
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	18	88
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	23	109
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	28	137
2.5, 2.75, 3.0, 3.25	33	157
2.5, 2.75, 3.0, 3.25	38	185

Stent Diameter (mm)	Stent Length (mm)	Drug Dose (µg)
3.5, 4.0	8	50
3.5, 4.0	12	75
3.5, 4.0	15	91
3.5, 4.0	18	116
3.5, 4.0	23	141
3.5, 4.0	28	174
3.5, 4.0	33	199
3.5, 4.0	38	232

- Two radiopaque markers, located underneath the balloon, which fluoroscopically mark the working length of the balloon and the expanded stent length.
- Two proximal delivery system shaft markers (95 cm and 105 cm proximal to the distal tip) indicate the relative position of the delivery system to the end of the brachial or femoral guiding catheter. Working catheter length is 145 cm.
- A shaft color change denotes the guide wire exit notch.

Table 1-3: *In vitro* Device Specifications

Stent Diameter (mm)	Stent Length (mm)	* Minimum Guiding Catheter Compatibility (ID)	** <i>in vitro</i> Stent Nominal Pressure		Rated Burst Pressure (RBP)		Stent Free Area (%)
			(atm)	kPa	(atm)	kPa	
2.0	8, 12, 15, 18, 23, 28	5F (0.056" / 1.42 mm)	10	1013	18	1824	81
2.25	8, 12, 15, 18, 23, 28	5F (0.056" / 1.42 mm)	10	1013	18	1824	83
2.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	84
2.75	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	86
3.0	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	87
3.25	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	88
3.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	10	1013	18	1824	85
4.0	8, 12, 15, 18, 23, 28, 33	5F (0.056" / 1.42 mm)	10	1013	18	1824	87
4.0	38	6F (0.066" / 1.68 mm)	10	1013	18	1824	87

\* See individual manufacturer specifications for (F) equivalent.

\*\* Assure full deployment of the stent (see section 9.5 *Deployment Procedure*). Deployment pressures should be based on lesion characteristics.

### 2.0 HOW SUPPLIED

**Sterile** – This device is sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if the package is open or damaged.

This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

**Contents** – One (1) XIENCE Xpedition, XIENCE Xpedition SV, or XIENCE Xpedition LL Everolimus Eluting Coronary Stent System; one (1) flushing tool.

**Storage** – Store at 15 – 30°C (59 – 86°F).

### 3.0 INDICATIONS

The XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent System are indicated for improving coronary luminal diameter in the following:

- Patients with symptomatic ischemic heart disease due to discrete *de novo* native coronary artery lesions.
- For restoring coronary flow in patients experiencing acute myocardial infarction who present within 12 hours of symptom onset.
- For the treatment of patients with concomitant diabetes, acute coronary syndrome, dual vessel lesions (two lesions in two different epicardial vessels), lesions residing within small coronary vessels; lesions where treatment results in the jailing of side branches (lesions with a side branch < 2 mm in diameter or an ostial stenosis < 50%); for the treatment of elderly patients (age ≥ 65), and for treatment of both men and women.
- For the treatment of patients presenting with in-stent restenosis in coronary artery lesions; chronic total occluded coronary artery lesions (defined as coronary artery lesions with TIMI flow 0 and lasting longer than 3 months); and coronary artery bifurcation lesions.

In all cases, the treated lesion length should be less than the nominal stent length (8 mm, 12 mm, 15 mm, 18 mm, 23 mm, 28 mm, 33 mm, or 38 mm) with a reference vessel diameter of ≥ 2.00 mm and ≤ 4.25 mm.

### 4.0 CONTRAINDICATIONS

The XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent System are contraindicated for use in:

- Patients in whom antiplatelet and / or anticoagulant therapy is contraindicated
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon
- Patients with a known hypersensitivity or contraindication to everolimus, cobalt, chromium, nickel, tungsten, acrylic, and fluoropolymers

### 5.0 WARNINGS

- Long-term outcome for this permanent implant containing polymers and everolimus is unknown at present.
- Judicious selection of patients is necessary, since the use of this device carries the associated risk of stent thrombosis, vascular complications, and / or bleeding events.
- Oral administration of everolimus in combination with cyclosporine has been associated with increased serum cholesterol and triglycerides.
- Persons allergic to L-605 cobalt chromium alloy, acrylic or fluoropolymers, or everolimus, may suffer an allergic reaction to this implant.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.

### 6.0 PRECAUTIONS

#### 6.1 Stent Handling – Precautions

- **For single use only.** Do not resterilize or reuse. Note the product "Use by" date.
- **The foil pouch is not a sterile barrier.** The inner header bag (pouch) within the foil pouch is the sterile barrier. **Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.**
- **Do not remove the stent from the delivery system,** as removal may damage the stent and / or lead to stent embolization. The stent system is intended to perform together as a system.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- **Do not manipulate, touch, or handle the stent with your fingers,** as this may cause coating damage, contamination, or dislodgement of the stent from the delivery balloon.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the stent.
- Implantation of the stent should be performed only by physicians who have received appropriate training.
- Coronary stent placement should only be performed at hospitals where emergency coronary artery bypass graft (CABG) surgery can be readily performed.
- Subsequent restenosis may require repeat dilation of the arterial segment containing the stent. The long-term outcome following repeat dilation of endothelialized stents is unknown at present.

#### 6.2 Stent Placement – Precautions

- **Do not prepare or pre-inflate the delivery system prior to stent deployment** other than as directed. Use balloon purging technique described under *Delivery System Preparation*.
- The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. Direct stenting in less complex coronary lesions with a predilate device for the XIENCE Xpedition EECSS (the XIENCE V EECSS) has been shown to be as effective and safe as stenting with pre-dilation for device lengths up to 28 mm in real-world settings. **If pre-dilation is performed,** limit the length of pre-dilation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the XIENCE Xpedition stent.
- When introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.
- Do not torque the catheter more than one (1) full turn.
- When treating multiple lesions within the same vessel, stent the distal lesion prior to stenting the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent, and reduces the chance of dislodging the proximal stent.
- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent, and may cause acute closure of the vessel requiring additional intervention (CABG, surgery, further dilation, placement of additional stents, or other).
- Do not expand the stent if it is not properly positioned in the vessel. (See section 6.4 *Stent / System Removal – Precautions*.)
- **Placement of a stent has the potential to compromise side branch patency.**
- The inflated balloon diameter of the system used to deploy the stent should approximate the diameter of the vessel. Oversizing of the stent can result in a ruptured vessel. To ensure full expansion of the stent, the balloon should be inflated to a minimum of nominal pressure.
- **Do not exceed RBP as indicated on product label.** Balloon pressures should be monitored during inflation. Use of pressures higher than specified on product label may result in a ruptured balloon with possible intimal damage and dissection.

- Stent retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- When multiple drug-eluting stents are required, only XIENCE Everolimus Eluting Coronary Stents with the identical cobalt-chromium stent substrate and identical drug-eluting polymer coating should be used (e.g., XIENCE V, XIENCE PRIME, XIENCE PRIME SV, XIENCE PRIME LL, XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL stents). Potential interaction with other drug-eluting stents or coated stents has not been evaluated, and should be avoided.
- If multiple stents are required, resulting in stent-to-stent contact, stent materials should be of similar composition. Placing multiple stents of different metals in contact with each other may increase the potential for corrosion *in vivo*, although *in vitro* corrosion tests using an L-605 CoCr alloy stent in combination with a 316L stainless steel alloy stent did not appear to increase corrosion.
- The extent of the patient's exposure to drug and polymer is directly related to the number of stents implanted. A patient can receive up to four XIENCE Xpedition, XIENCE Xpedition SV, or XIENCE Xpedition LL Everolimus Eluting Coronary Stents or other everolimus eluting coronary stents from the XIENCE family (i.e., XIENCE V, XIENCE PRIME) depending on the number of vessels treated and the lesion length. Those patients receiving bailout stenting will receive additional XIENCE family stents. The use of multiple XIENCE family stents will result in the patient receiving larger amounts of drug and polymer.
- The safety and effectiveness of the XIENCE Xpedition, XIENCE Xpedition SV, or XIENCE Xpedition LL Everolimus Eluting Coronary Stent in patients with prior brachytherapy of the target lesion or the use of brachytherapy to treat in-stent restenosis in a XIENCE Xpedition, XIENCE Xpedition SV, or XIENCE Xpedition LL Everolimus Eluting Coronary Stent have not been established. Both vascular brachytherapy and the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stents alter arterial remodeling. The synergy between these two treatments has not been established.
- Placement of a stent in the left main coronary artery has the potential to compromise blood flow to the distal anatomy.
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the stent may be damaged when retracting the unexpanded stent back into the guiding catheter.
- Should resistance be felt at any time during coronary stent system withdrawal, please follow the steps provided in section 6.4 Stent / System Removal – Precautions.

### 6.3 Use in Conjunction with Other Procedures

The safety and effectiveness have not been established with use of the following devices: mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent implantation.

### 6.4 Stent / System Removal – Precautions

#### Stent delivery system removal prior to stent deployment:

If removal of a stent system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guide catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guide catheter, the stent delivery system and the guide catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

#### Withdrawal of the stent delivery catheter from the deployed stent:

1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 – 15 seconds longer.
2. Position the inflation device to "negative" or "neutral" pressure.
3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

#### Note: If, during withdrawal of the catheter, resistance is encountered, use the following steps to improve balloon rewrap:

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

Failure to follow these steps and / or applying excessive force to the delivery system can potentially result in loss or damage to the stent and / or delivery system components.

If it is necessary to retain guide wire position for subsequent artery / lesion access, leave the guide wire in place and remove all other system components.

Stent retrieval methods (i.e., additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include, but are not limited to, bleeding, hematoma, or pseudoaneurysm.

### 6.5 Post Implantation – Precautions

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, a balloon catheter, or delivery system, to avoid disrupting the stent placement, apposition, and / or geometry.
- Antiplatelet therapy should be administered post-procedure (see section 8.1 Individualization of Treatment). Patients who require early discontinuation of antiplatelet therapy (e.g., secondary to active bleeding) should be monitored carefully for cardiac events. At the discretion of the patient's treating physician, the antiplatelet therapy should be restarted as soon as possible.

If the patient requires imaging, see section 6.6 – MRI Statement.

### 6.6 MRI Statement

Non-clinical testing has demonstrated that the XIENCE Xpedition stent, in single and in overlapped configurations up to 71 mm in length, is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 or 3 Tesla
- Spatial gradient field of 2500 Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for up to 15 minutes of scanning for each sequence

The XIENCE Xpedition stent should not migrate in this MRI environment. Non-clinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of the XIENCE Xpedition stent.

Stent heating was derived by using the measured non-clinical, *in vitro* temperature rises in a GE Excite 3 Tesla scanner and in a GE 1.5 Tesla coil in combination with the local specific absorption rates (SARs) in a digitized human heart model. The maximum whole body average SAR was determined by validated calculation. At overlapped lengths of up to 71 mm, the XIENCE Xpedition stent produced a non-clinical maximum local temperature rise of 3.3°C at maximum whole body average SAR of 2.0 W/kg (normal operating mode) for 15 minutes. These calculations do not take into consideration the cooling effects of blood flow.

The effects of MRI on overlapped stents greater than 71 mm in length or stents with fractured struts are unknown.

As demonstrated in non-clinical testing, an image artifact can be present when scanning the XIENCE Xpedition stent. MR image quality may be compromised if the area of interest is in the exact same area as, or relatively close to, the position of the XIENCE Xpedition stent. Therefore, it may be necessary to optimize the MR imaging parameters for the presence of XIENCE Xpedition stents.

### 6.7 Drug Interactions

Everolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the gut wall and liver, and is a substrate for the countertransporter P-glycoprotein. Therefore, absorption and subsequent elimination of everolimus may be influenced by drugs that affect these pathways. Everolimus has also been shown to reduce the clearance of some prescription medications when it was administered orally along with cyclosporine (CsA). Formal drug interaction studies have not been performed with the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Eluting Coronary Stent System. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall, when deciding to place the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent in a subject taking a drug with known interaction with everolimus.

Everolimus, when prescribed as an oral medication may interact with the following drugs or foods<sup>1,2</sup>:

- CYP3A4 / P-glycoprotein isozyme inhibitors
  - Antifungal agents (e.g., fluconazole, ketoconazole, itraconazole, posaconazole, voriconazole)
  - Macrolide antibiotics (e.g., erythromycin, clarithromycin, telithromycin)
  - Calcium channel blockers (e.g., verapamil, nifedipine, diltiazem)
  - Protease inhibitors (e.g., ritonavir, atazanavir, saquinavir, darunavir, indinavir, nelfinavir, amprenavir, fosamprenavir)
  - Other (e.g., cyclosporine, nefazodone, cisapride, metoclopramide, bromocriptine, cimetidine, danazol, sildenafil, terfenadine, astemizole, grapefruit / grapefruit juice, digoxin)
- CYP3A4 / P-glycoprotein isozyme inducers
  - Antibiotics (e.g., rifampin, rifabutin, ciprofloxacin, ofloxacin)
  - Anticonvulsants (e.g., carbamazepine, phenobarbital, phenytoin)
  - Non-nucleoside reverse transcriptase inhibitors (e.g., efavirenz, nevirapine)
  - Glucocorticoids (e.g., dexamethasone, prednisone, prednisolone)
  - HMGCoA reductase inhibitors (simvastatin, lovastatin)
  - Other (e.g., St. John's Wort)

### 6.8 Pregnancy

This product has not been tested in pregnant women or in men intending to father children. Effects on the developing foetus have not been studied. While there is no contraindication, the risks and reproductive effects are unknown at this time.

### 7.0 POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a stent in native coronary or peripheral arteries include:

- Abrupt closure
  - Allergic reaction to contrast
  - Aneurysm
  - Arterial perforation
  - Arterial rupture
  - Arteriovenous fistula
  - Arrhythmias, including atrial and ventricular
  - Bleeding complications, which may require transfusion
  - Coronary artery spasm
  - Coronary or stent embolism
  - Coronary or stent thrombosis
  - Death
  - Dissection of the coronary artery
  - Distal emboli (air, tissue, or thrombotic)
  - Drug reactions to antiplatelet agents / contrast medium
  - Embolization (stent or other)
  - Emergent or non-emergent coronary or peripheral artery bypass graft surgery
  - Fever
  - Hypersensitivity reactions
  - Hypertension / hypotension
  - Infection and pain at insertion site
  - Injury to the coronary or peripheral artery
  - Intervention due to:
    - Stent migration
    - Unintentional placement of stent
    - Partial stent deployment
    - Damaged stents
  - Ischemia, myocardial
  - Myocardial infarction
  - Nausea and vomiting
  - Palpitations
  - Peripheral vascular or nerve injury
  - Pseudoaneurysm
  - Renal insufficiency / failure
  - Restenosis of stented segment
  - Stroke / cerebrovascular accident
  - Thrombosis (stent or other)
  - Total occlusion of coronary artery
  - Unstable or stable angina pectoris
  - Vascular complications, including entry site, which may require vessel repair
  - Ventricular arrhythmias, including ventricular fibrillation and ventricular tachycardia
  - Vessel dissection
  - Vessel occlusion
  - Vessel spasm
- Adverse events associated with daily oral administration of everolimus (1.5 mg/day to 3.0 mg/day for at least 1 year, or 10 mg/day for at least 4 months):
- Abdominal pain
  - Acne / acneiform dermatitis
  - Anemia
  - Anorexia
  - Chest pain
  - Coagulopathy
  - Congestive heart failure
  - Conjunctivitis
  - Cough
  - Dehydration
  - Diarrhea
  - Dry mouth
  - Dry skin
  - Dysgeusia (abnormal taste)
  - Dyspepsia
  - Dysphagia
  - Dyspnea
  - Edema (including generalized edema, angioneurotic edema, and peripheral edema)
  - Epistaxis
  - Erectile dysfunction
  - Erythema
  - Fatigue
  - Fever
  - Graft thrombosis
  - Headache
  - Hemolysis
  - Hemoptysis
  - Hemorrhage
  - Hepatitis or hepatic dysfunction
  - Hypercholesterolemia
  - Hyperglycemia (possible new onset diabetes mellitus)
  - Hyperlipidemia
  - Hypersensitivity
  - Hypertension
  - Hypertiglyceridemia
  - Hypogonadism male
  - Hypophosphatemia
  - Impaired healing
  - Insomnia
  - Jaundice
  - Leukocytoclastic vasculitis
  - Leukopenia
  - Liver function test abnormalities
  - Lymphocele
  - Myalgia
  - Nail disorder
  - Nausea
  - Pain
  - Palmar-plantar erythrodysesthesia syndrome
  - Pancreatitis
  - Pericardial effusion
  - Pleural effusion
  - Pneumonia
  - Pneumonitis (including interstitial lung disease)
  - Proteinuria
  - Pruritis
  - Pulmonary alveolar proteinosis
  - Pulmonary embolism
  - Pyelonephritis
  - Rash
  - Renal failure (including acute renal failure)
  - Renal tubular necrosis
  - Sepsis
  - Skin exfoliation
  - Stomatitis / mouth ulceration and / or mucosal inflammation
  - Surgical wound complication
  - Thrombocytopenia
  - Thrombotic thrombocytopenic purpura / hemolytic uremic syndrome (TTP / HUS)
  - Urinary tract infection
  - Venous thromboembolism
  - Viral, bacterial, and fungal infections
  - Vomiting
  - Weakness
  - Wound infection

### 8.0 PATIENT SELECTION AND TREATMENT

#### 8.1 Individualization of Treatment

The risks and benefits described above should be considered for each patient before using the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent System. Patient selection factors to be assessed should include judgment regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Antiplatelet drugs should be used in combination with the XIENCE Xpedition, XIENCE Xpedition SV, XIENCE Xpedition LL Everolimus Eluting Coronary Stent. Physicians should use the information from the SPIRIT Clinical Trials, coupled with the current literature on drug-eluting stents and the specific needs of individual patients, to determine the specific antiplatelet / anticoagulation regimen to be used for their patients in general practice.

XIENCE V demonstrated low stent thrombosis rates in patients who either discontinued or interrupted dual antiplatelet therapy (DAPT) after 3 months post stent implantation. It is therefore recommended that patients treated with XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL stents remain on DAPT for at least 3 months after stent implantation.

It is very important that the patient comply with the post-procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could

<sup>1</sup> Cerican™ drug prescribing sheet, date of most recent amendment March 30, 2011. Novartis.

<sup>2</sup> Afinitor™ summary of product characteristics, date of authorization March 8, 2009. Novartis.

result in a higher risk of thrombosis, myocardial infarction, or death. Prior to percutaneous coronary intervention (PCI), if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventionalist and patient should carefully consider whether a drug-eluting stent and its associated recommended antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding, should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

## 9.0 CLINICIAN USE INFORMATION

### 9.1 Inspection Prior to Use

Prior to using the XIENCE Xpedition, XIENCE Xpedition SV, and XIENCE Xpedition LL Everolimus Eluting Coronary Stent System, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent does not extend beyond the radiopaque balloon markers. Do not use if any defects are noted.

### 9.2 Materials Required

- Appropriate guiding catheter(s)
- 2 – 3 syringes (10 – 20 cc)
- 1,000 u/500 cc heparinized normal saline (HepNS)
- 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
- Rotating hemostatic valve with 0.096 inch (2.44 mm) minimum inner diameter
- 60% contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

### 9.3 Preparation

#### 9.3.1 Packaging Removal

**Note:** The foil pouch is not a sterile barrier. The inner header bag (pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

1. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. When using a rapid exchange (RX) system, do not bend or kink the hypotube during removal.
2. Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product and replace with another. Follow product returns procedure for the unused device.

#### 9.3.2 Guide Wire Lumen Flush

1. Flush the guide wire lumen with HepNS using the flushing tool supplied with the product. Insert the flushing tool into the tip of the catheter and flush until fluid exits the guide wire exit notch.

**Note:** Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent on the balloon.

#### 9.3.3 Delivery System Preparation

1. Prepare an inflation device / syringe with diluted contrast medium.
2. Attach an inflation device / syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

**Note:** While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

**Note:** If air is seen in the shaft, repeat section 9.3.3 *Delivery System Preparation*, steps 3 through 5, to prevent uneven stent expansion.

### 9.4 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
2. The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. Direct stenting in less complex coronary lesions with a predicate device for the XIENCE Xpedition EECSS (the XIENCE V EECSS) has been shown to be as effective and safe as stenting with pre-dilation for device lengths up to 28 mm in real-world settings. **If pre-dilation is performed**, limit the length of pre-dilation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the XIENCE Xpedition stent.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel.

**Note:** If choosing between two stent diameters for tight lesions choose the smaller diameter stent and inflate. See product label for compliance information.

4. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
5. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion.
6. Carefully advance the delivery system into the guiding catheter and over the guide wire to the target lesion. Be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent system into the coronary artery.

**Note:** If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and the use of excessive force may result in stent damage or dislodgement. Maintain guide wire placement across the lesion and remove the delivery system and guiding catheter as a single unit.

7. Advance the delivery system over the guide wire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm stent position. If the position of the stent is not optimal, it should be carefully repositioned or removed (see section 9.6 *Removal Procedure*). The balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion.

**Note:** If removal of a stent system is required prior to deployment, ensure that the guide catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guide catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guide catheter, the stent delivery system and the guide catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

8. Tighten the rotating hemostatic valve. The stent is now ready to be deployed.

### 9.5 Deployment Procedure

**CAUTION: Refer to product label for *in vitro* stent inner diameter, nominal pressure, and RBP.**

1. Prior to deployment, reconfirm the correct position of the stent relative to the target lesion using the radiopaque balloon markers.
2. Deploy the stent slowly by pressurizing the delivery system in 2 atm increments, every 5 seconds, until stent is completely expanded. Fully expand the stent by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter ratio of about 1.1 times the reference vessel diameter (refer to product label for *in vitro* stent inner diameter, nominal pressure, and RBP).
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel and expand stent to nominal pressure at minimum. Maintain pressure for 30 seconds. If necessary, the delivery system can be pressurized or further pressurized to ensure complete apposition of the stent to the artery wall.
4. Maintain pressure for 30 seconds for full expansion of the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal stent expansion and proper apposition require that the stent be in full contact with the arterial wall.

**Note:** See section 9.6 *Removal Procedure* for instruction on withdrawal of stent delivery system.

5. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the stent to the artery wall.

**CAUTION: Do not exceed the labeled rated burst pressure (RBP) of 18 atm (1824 kPa).**

6. Fully cover the entire lesion and balloon treated area (including dissections) with the XIENCE Xpedition stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
7. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to guiding catheter position.

**Note:** See section 9.6 *Removal Procedure* for instruction on withdrawal of stent delivery system.

8. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
9. If the deployed stent size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent. If an initial angiographic appearance is suboptimal, the stent may be further expanded using a low profile, high pressure, noncompliant balloon dilation catheter. If this is required, the stented segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the stent geometry. Deployed stents should not be left underdilated.

**CAUTION: Do not dilate the stent beyond the following limits.**

Nominal Stent Diameter	Dilation Limit
2.0 – 2.5 mm	3.25 mm
2.75 – 3.25 mm	3.75 mm
3.5 – 4.0 mm	4.50 mm

10. If more than one XIENCE Xpedition stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE Xpedition stent should be positioned inside the deployed stent prior to expansion.
11. Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

### 9.6 Removal Procedure

**Withdrawal of the stent delivery catheter from the deployed stent:**

1. Deflate the balloon by pulling negative pressure on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 to 15 seconds longer.
2. Position inflation device on “negative” or “neutral” pressure.
3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

**If, during withdrawal of the stent delivery catheter, resistance is encountered, use the following steps to improve balloon rewrap:**

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

**Post stent delivery system withdrawal; Stent deployment confirmation**

1. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
2. If more than one XIENCE Xpedition stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped.
3. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE Xpedition stent should be positioned inside the deployed stent prior to expansion.
4. Reconfirm stent position and angiographic results to assess stented area. Repeat inflations until optimal stent deployment is achieved. If post-dilation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Ensure that the stent wall is in contact with the artery wall.

### 9.7 Post-deployment Dilation of Stent Segments

All efforts should be taken to ensure that the stent is not under-dilated. If the deployed stent size is still inadequate with respect to the vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and noncompliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guide wire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

**CAUTION: Do not dilate the stent beyond the following limits.**

Nominal Stent Diameter	Dilation Limit
2.0 – 2.5 mm	3.25 mm
2.75 – 3.25 mm	3.75 mm
3.5 – 4.0 mm	4.50 mm

## 10.0 SPIRIT AND XIENCE FAMILY OF CLINICAL TRIALS

The XIENCE Xpedition EECSS is based on the predicate devices XIENCE V EECSS and XIENCE PRIME EECSS.

The XIENCE Xpedition EECSS uses the identical stent platform, identical drug coating formulation, identical drug primer, identical nominal total drug content, and the identical stent contacting balloon materials as the XIENCE PRIME EECSS.

The XIENCE Xpedition EECSS differs from the XIENCE PRIME EECSS only in the stent delivery system. The XIENCE Xpedition stent delivery system utilizes the same principle of operation and materials as other Abbott Vascular RX stent systems and coronary dilatation catheters. Compared to the XIENCE V EECSS, the XIENCE Xpedition EECSS has the same stent security specification, similar stent placement on the balloon between the balloon markers, a similar tip entry profile, and a similar taper length for Xpedition stent sizes up to 28 mm in length. The crossing profile for XIENCE Xpedition is as good as or better than XIENCE V due to the smaller distal balloon seal area on the XIENCE Xpedition delivery system.

Based on the identical nature of the XIENCE Xpedition stent to the XIENCE PRIME stent, performance of the XIENCE Xpedition EECSS can be predicted to be similar to the performance of XIENCE V and XIENCE PRIME. Therefore, clinical trial data for XIENCE V and XIENCE PRIME are summarized in this section.

10.1 Pre-market Clinical Trials

Principal safety and effectiveness of the XIENCE V stent has been established from a series of pre-market clinical trials. SPIRIT III RCT was the pivotal randomized clinical trial (RCT) that established the non-inferiority of the XIENCE V stent to the TAXUS<sup>®</sup> Express<sup>®</sup> stent (TAXUS stent). The SPIRIT IV trial is a prospective, randomized, active-controlled, single-blinded, multi-center evaluation of the XIENCE V stent compared to the TAXUS Express stent (TAXUS stent) in the treatment of up to three *de novo* lesions  $\leq$  28 mm in length in native coronary arteries with RVD  $\geq$  2.5 mm to  $\leq$  4.25 mm. The SPIRIT IV trial established the superiority of the XIENCE V stent compared to the TAXUS Express stent. The SPIRIT Small Vessel (SV) Registry is a prospective, single-arm, open-label, US multi-center registry study that established the safety and effectiveness of the XIENCE V 2.25 mm diameter stent. SPIRIT PRIME is a prospective, open-label, multi-center nonrandomized clinical trial with two study arms using the core size XIENCE PRIME and XIENCE PRIME LL stent system that established the safety and effectiveness of the XIENCE PRIME and XIENCE PRIME LL stents. Tables 10.1-1 through 10.1-4 present the trial designs, angiographic results (for studies that required angiographic follow-up), and principal clinical outcomes at 1-year and from latest follow-up, respectively.

Table 10.1-1: SPIRIT Family of XIENCE V Clinical Trial Designs (Pre-market)

	SPIRIT III RCT	SPIRIT IV	SPIRIT Small Vessel Registry	SPIRIT PRIME Clinical Trial	
				Core Size Registry	Long Lesion Registry
<b>Study Type / Design</b>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Randomized</li> <li>Single-blinded</li> <li>Active-control</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Randomized</li> <li>Single-blinded</li> <li>Active-control</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Open-label</li> <li>Single-arm</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Open-label</li> <li>Single-arm</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Open-label</li> <li>Single-arm</li> </ul>
<b>Number of Subjects Enrolled</b>	Total: 1,002 XIENCE V: 668 TAXUS Express Control: 334	Total: 3,690 XIENCE V: 2,460 TAXUS Express Control: 1,230 <sup>1</sup>	Total: 150 2.25 mm XIENCE V	Total: 400 XIENCE PRIME	Total: 100 XIENCE PRIME
<b>Treatment</b>	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to three <i>de novo</i> lesions, maximum of two lesions per epicardial vessel	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels
<b>Lesion Size</b>	RVD: $\geq$ 2.5 $\leq$ 3.75 mm Length: $\leq$ 28 mm	RVD: $\geq$ 2.5 $\leq$ 4.25 mm <sup>2</sup> Length: $\leq$ 28 mm	RVD: $\geq$ 2.25 $<$ 2.50 mm Length: $\leq$ 28 mm	RVD: $\geq$ 2.25 $\leq$ 4.25 mm Length: $\leq$ 22 mm	XIENCE PRIME CS: RVD: $\geq$ 2.25 $\leq$ 4.25 mm Length: $\leq$ 22 mm XIENCE PRIME LL: RVD: $\geq$ 2.5 $\leq$ 4.25 mm Length: $>$ 22 mm and $\leq$ 32 mm
<b>Primary Endpoint</b>	In-segment late loss at 240 days	Ischemia-driven target lesion failure at 1 year (composite of cardiac death, target vessel MI or ischemia driven TLR)	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year
<b>Co-Primary Endpoint</b>	TVF at 270 days	None	None	None	None
<b>Clinical Follow-up</b>	30, 180, 240, 270 days, 1 to 5 years	30, 180, 270 days, 1 to 3 years	30 days, 240 days, 1 to 3 years	30, 180 days, 1 to 3 years	30, 180 days, 1 to 3 years
<b>Angiographic Follow-up</b>	240 days (N = 564)	None	240 days (N = 69)	None	None

<sup>1</sup> In the TAXUS stent arm, there was 1 subject who received 1 TAXUS<sup>®</sup> Liberté<sup>®</sup> stent

<sup>2</sup> RVD  $\geq$  2.5 mm to  $\leq$  3.75 mm and stent sizes up to 3.5 mm until 4.0 mm TAXUS is commercially available.

Table 10.1-2: SPIRIT Family of Clinical Trials Angiographic Results (Pre-market)

Angiographic Results	SPIRIT III RCT 240 Days		SPIRIT Small Vessel 240 Days
	XIENCE V (N = 376) (M = 427)	TAXUS (N = 188) (M = 220)	2.25 mm XIENCE V (N = 69) (M = 69)
In-Stent Late Loss (mm)	0.16 $\pm$ 0.41 (342)	0.30 $\pm$ 0.53 (158)	0.20 $\pm$ 0.40 (52)
In-Segment Late Loss (mm)	0.14 $\pm$ 0.39 (343)	0.26 $\pm$ 0.46 (158)	0.16 $\pm$ 0.41 (52)
In-Stent Binary Restenosis	2.3% (8/343)	5.7% (9/158)	3.8% (2/52)
In-Segment Binary Restenosis	4.7% (16/344)	8.9% (14/158)	9.6% (5/52)

Notes:

- Data are mean (mm)  $\pm$  SD or % (n/N).
- N is total number of patients. M is total number of lesions.
- SPIRIT III and SV 240-day includes follow-up window (240 + 28 days).

Table 10.1-3: SPIRIT Family of Clinical Trials Principal 1-Year Clinical Outcomes (Pre-market)

	SPIRIT IV		SPIRIT III RCT		SPIRIT Small Vessel	SPIRIT PRIME Clinical Trial	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	4.0% (97/2416)	6.8% (81/1195)	5.3% (35/655)	9.7% (31/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
TVF	5.5% (134/2416)	7.7% (92/1195)	8.5% (56/655)	11.6% (37/319)	11.0% (15/136)	N/A	N/A
MACE	4.1% (98/2416)	6.9% (82/1195)	6.0% (39/655)	10.3% (33/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
All Death	1.0% (25/2416)	1.3% (15/1195)	1.2% (8/657)	1.3% (4/320)	1.5% (2/136)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	0.4% (10/2416)	0.4% (5/1195)	0.8% (5/657)	0.9% (3/320)	1.5% (2/136)	0.3% (1/399)	0.0% (0/104)
MI	1.9% (45/2416)	3.1% (37/1195)	2.7% (18/655)	4.1% (13/319)	1.5% (2/136)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	2.2% (54/2416)	3.3% (39/1195)	3.4% (22/655)	4.7% (15/319)	2.9% (4/136)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	2.3% (56/2416)	4.6% (55/1195)	3.4% (22/655)	5.6% (18/319)	5.1% (7/136)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	2.2% (54/2416)	2.4% (29/1195)	3.2% (21/655)	4.7% (15/319)	5.9% (8/136)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.29% (7/2391)	1.10% (13/1181)	0.9% (6/650)	0.6% (2/316)	1.5% (2/136)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.3% (6/2385)	0.8% (10/1183)	0.8% (5/650)	0.3% (1/317)	0.7% (1/136)	0.5% (2/399)	0.0% (0/104)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- 1-year includes the follow-up window (365 + 28 days) for all trials.
- TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically-indicated TLR definition rather than ischemia-driven TLR.
- TVF includes cardiac death, MI, ischemia-driven TLR and TVR, non-target lesion. SPIRIT SV and PRIME used clinically-indicated TLR and TVR definition rather than ischemia-driven TLR and TVR definition, which was used for SPIRIT II, SPIRIT III, and SPIRIT IV.
- MACE includes cardiac death, MI, and ischemia-driven TLR.

Table 10.1-4: SPIRIT Family of Clinical Trials Principal Clinical Outcomes from Latest Follow-up (Pre-market)

	SPIRIT IV 3 Years		SPIRIT III RCT 5 Years		SPIRIT Small Vessel 2 Years	SPIRIT PRIME Clinical Trial 1 Year	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	9.5% (223/2348)	11.9% (138/1158)	13.4% (81/605)	20.6% (59/286)	8.3% (11/133)	4.5% (18/399)	7.7% (6/104)
TVF	13.3% (312/2348)	14.5% (168/1158)	20.3% (123/605)	26.6% (76/286)	12.0% (16/133)	N/A	N/A
MACE	9.8% (231/2348)	12.3% (142/1158)	14.4% (87/605)	22.0% (63/286)	8.3% (11/133)	4.5% (18/399)	7.7% (6/104)
All Death	3.4% (81/2348)	5.2% (60/1158)	6.0% (37/621)	10.3% (31/300)	1.5% (2/133)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	1.4% (34/2348)	1.9% (22/1158)	2.7% (17/621)	4.3% (13/300)	1.5% (2/133)	0.3% (1/399)	0.0% (0/104)
MI	3.1% (73/2348)	4.7% (55/1158)	4.6% (28/605)	7.0% (20/286)	1.5% (2/133)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	4.5% (105/2348)	6.0% (70/1158)	7.1% (43/605)	11.2% (32/286)	3.0% (4/133)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	6.3% (148/2348)	7.9% (92/1158)	8.9% (54/605)	12.9% (37/286)	5.3% (7/133)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	5.6% (132/2348)	5.4% (63/1158)	8.8% (53/605)	11.9% (34/286)	6.8% (9/133)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.62% (14/2263)	1.73% (19/1098)	1.5% (9/582)	1.9% (5/268)	1.5% (2/132)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.49% (11/2263)	1.28% (14/1098)	1.2% (7/582)	0.7% (2/268)	0.8% (1/132)	0.5% (2/399)	0.0% (0/104)

Notes:

- All counts presented in this table are SUBJECT counts. Subjects are counted only once for each event for each time period.
- Data includes the follow-up window of + 28 days for all trials.
- TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically-indicated TLR definition rather than ischemia-driven TLR.
- TVF includes cardiac death, MI, ischemia-driven TLR and TVR, non-target lesion.
- MACE includes cardiac death, MI and ischemia-driven TLR.

10.2 Post-Market Clinical Trials

The XIENCE V USA study is a prospective, multi-center, FDA mandated post-market study to evaluate the continued safety and effectiveness of the XIENCE V EECSS in real-world settings after it was commercialized in U.S., and also to support the FDA DAPT initiative. The objective of the SPIRIT V Single Arm Study (SAS) was to continue the assessment of the performance of XIENCE V EECSS in the treatment of patients with *de novo* coronary artery lesions. XIENCE V India is a prospective, open-label, multi-center, observational, single-arm registry to evaluate XIENCE V EECSS continued safety and effectiveness during commercial use in real-world settings in India. Tables 10.2-1 through 10.2-2 present clinical trial designs, principal clinical outcomes at 1 year and from latest follow-up, respectively.

The results from these post-marketing clinical trials demonstrate safety and effectiveness of XIENCE V in real-world settings. In addition, XIENCE V improved patient-reported outcomes (including better quality of life, reduced angina frequency, improved angina stability, and reduced physical limitation) at 6 months and the improvements were sustained for 1 year in coronary artery disease patients.

Table 10.2-1: SPIRIT and XIENCE V Family Clinical Trial Designs (Post-market)

	XIENCE V USA Phase I Cohort	XIENCE V USA Long-Term Follow-up Cohort	XIENCE V USA AV-DAPT Cohort	SPIRIT V (SAS)	XIENCE V India
Study Type / Design	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Prospective</li> <li>Single-arm</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Prospective</li> <li>Single-arm</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Randomized</li> <li>Double-blinded</li> <li>Placebo control</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Prospective</li> <li>Single-arm</li> </ul>	<ul style="list-style-type: none"> <li>Multi-center</li> <li>Prospective</li> <li>Single-arm</li> </ul>
Number of Subjects Enrolled	8040	4663	868	2663	977
Treatment	Per site standard care	Per site standard care	Patients were randomized to receive either thienopyridine or placebo treatment for additional 18 months along with aspirin	Maximum of one <i>de novo</i> , native target lesion per major epicardial vessel or side branch (no prior stent implant, no prior brachytherapy) maximum of 4 planned EES	Per site standard care
Lesion Size	No angiographic restrictions			RVD = $\geq 2.25 \leq 4.0$ mm Length $\leq 28$ mm by visual estimation.	No angiographic restrictions
Primary Endpoint	ARC definite and probable stent thrombosis up to 1 year	ARC definite and probable stent thrombosis from year 1 to 5	MACE (composite of all death, MI and stroke) 12-33 months	Composite rate of all death, MI, TVR at 30-day	ARC ST (Definite / probable) 1 year and yearly after through 3 years
Co-Primary Endpoint	Cardiac death or any MI at 1 year	Cardiac death or any MI from year 1 to 5	ARC definite and probable ST 12-33 months	None	Cardiac death and any MI at 1 year
Clinical Follow-up	14, 30, 180 days, and 1 year	2, 3, and 4 years	15, 24, 30, and 33 months	30 days and 1 and 2 years	14, 30, 180 days and 1, 2, and 3 years
Angiographic Follow-up	None	None	None	None	None

Table 10.2-2: SPIRIT and XIENCE V Family of Trials Principal Clinical Outcomes (Post-market)

	XIENCE V USA Phase I 1 Year	SPIRIT V (SAS)		XIENCE V India	
		1 Year	2 Years	1 Year	2 Years
	XIENCE V N = 8040	XIENCE V N = 2663	XIENCE V N = 2663	XIENCE V N = 990	XIENCE V N = 990
TLF (ARC)	9.4% (707/7522)	5.25% (138/2627)	7.49% (192/2562)	2.4% (24/986)	3.4% (32/942)
TLF	6.8% (513/7505)	N/A	N/A	N/A	N/A
All Death, MI (ARC) and TVR	N/A	7.04% (185/2627)	10.34% (265/2562)	N/A	N/A
Cardiac Death or MI (ARC)	7.2% (545/7522)	4.23% (111/2627)	5.74% (147/2562)	1.9% (19/986)	2.8% (26/942)
Cardiac Death or MI	3.3% (249/7505)	N/A	N/A	N/A	N/A
TLR	4.6% (349/7522)	1.90% (50/2627)	3.04% (78/2562)	1.2% (12/986)	1.5% (14/942)
TVR, non TLR	2.3% (176/7522)	1.45% (38/2627)	2.26% (58/2562)	0.1% (1/986)	0.1% (1/942)
All Death	2.6% (194/7522)	1.71% (45/2627)	2.97% (76/2562)	0.9% (9/986)	1.7% (16/942)
Cardiac Death	1.4% (108/7522)	1.10% (29/2627)	1.87% (48/2562)	0.9% (9/986)	1.7% (16/942)
All MI (ARC)	6.3% (475/7522)	3.54% (93/2627)	4.45% (114/2562)	1.3% (13/986)	1.6% (15/942)
All MI	2.2% (162/7505)	N/A	N/A	N/A	N/A
Stent Thrombosis					
ARC (Definite / Probable)	0.81% (60/7380)	0.65% (17/2607)	0.79% (20/2523)	0.51% (5/986)	0.53% (5/939)
ARC (Definite)	0.54% (40/7380)	N/A	N/A	0.41% (4/986)	0.43% (4/939)

- Notes:**
- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
  - XIENCE V USA 1-year data include 42-day window or through randomization date if occurred earlier than 407 days for the second enrollment phase. XIENCE V India 1-year data include 47-day window. SPIRIT V 1-year data include 28-day window.
  - SPIRIT V 2-year data include 28-day window. XIENCE V India 2-year data include 47-day window.
  - TLF (ARC) includes cardiac death, MI attributed to target vessel (per ARC definition), clinically-indicated TLR. TLF includes cardiac death, MI attributed to target vessel (per protocol definition), clinically-indicated TLR.

**11.0 PATENTS AND TRADEMARKS**

This product and / or its use may be covered by one or more of the following United States Patents: 5,514,154; 5,569,295; 5,636,641; 5,649,952; 5,665,772; 5,759,192; 5,780,807; 5,868,706; 6,131,266; 6,179,810; 6,369,355; 6,384,046; 6,419,693; 6,440,990; 6,482,166; 6,629,994; 6,656,220; 6,736,843; 6,746,423; 6,827,734; 6,887,219; 6,887,510; 6,890,318; 6,929,657; 6,939,373; 6,957,152. Other U.S. patents pending. Foreign patents issued and pending.

XIENCE V, XIENCE PRIME, and XIENCE Xpedition are trademarks of the Abbott Group of Companies.

Certican and Afinitor are registered trademarks of Novartis AG Aktiengesellschaft SWITZERLAND.

TAXUS Express and TAXUS Liberté are registered trademarks of Boston Scientific Corporation and its affiliates.

# XIENCE Xpedition, XIENCE Xpedition SV und XIENCE Xpedition LL

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