GRAFTMASTER RX

Coronary Stent Graft System

Koronarstentprothesensystem

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Coronary Stent Graft System

1.0 DEVICE DESCRIPTION

The GRAFTMASTER RX Coronary Stent Graft is constructed using a sandwich technique, whereby an ultra-thin layer of expandable polytetrafluoroethylene (PTFE) is placed between two GRAFTMASTER stents which are then pre-mounted on a balloon catheter delivery system. The stents are fabricated from medical grade 316L stainless steel.

The GRAFTMASTER delivery system is a rapid exchange co-axial design with the balloon and stent graft at the distal end of the catheter. The proximal lumen provides for inflation of the balloon with contrast medium. The central distal lumen permits a guide wire to be inserted through the lumen. The annular space between the distal outer member and the central distal lumen provides a fluid passage path from the proximal lumen to the balloon. The shaft of the catheter, the tip, and the tapers of the balloon are coated with HYDROCOAT hydrophilic coating.

Two radiopaque markers located on the distal end of the inner member are positioned to mark the working length of the balloon. The radiopaque markers fluoroscopically aid in positioning the stent pre-deployment and the delivery system for post-deployment dilatation. The balloon is designed to deliver an expandable stent of known diameter and length at specified pressures. Markers located on the proximal outer shaft help the physician gauge the delivery catheter position relative to the guiding catheter tip.

An adaptation arm on the proximal end of the catheter provides access to the inflation lumen. It is designed with a luer-lock fitting to facilitate connection to an inflation device.

Note: During stent graft deployment with the stent delivery system from crimped state to 4.80 mm, the stent graft may shorten up to 20%. Maximum post dilatation that can be achieved with a noncompliant post dilatation balloon is a maximum of 5.5 mm. With expansion to this diameter the system may shorten up to 25%. When choosing a GRAFTMASTER RX system for expansion in larger vessels a longer stent graft length is recommended to ensure treatment area is covered by the stent graft.

2.0 HOW SUPPLIED

Sterile – This device is sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if package is open or damaged.

This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning and / or re-sterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents – One (1) GRAFTMASTER RX Coronary Stent Graft System; one (1) protective sheath; one (1) Flexi-Clip; one (1) flushing tool.

Storage – Store in a dry, dark, cool place.

3.0 INDICATIONS

The GRAFTMASTER RX is a balloon-expandable pre-mounted coronary stent graft for intraluminal chronic placement in coronary arteries or aorto-coronary bypass grafts for the treatment of:

- Coronary artery aneurysm
- Coronary bypass-vein graft aneurysm
- Acute coronary artery perforation
- Acute coronary artery rupture

Patients being considered for stent graft implantation should be acceptable for coronary artery bypass graft surgery, and / or for coronary balloon angioplasty with ischemic heart disease due to discrete *de novo* or restenotic native coronary lesions.

4.0 CONTRAINDICATIONS

The GRAFTMASTER RX is contraindicated for use in:

- Patients in whom antiplatelet and / or anticoagulation therapy is contraindicated
- Patients who are judged to have a treatment area that prevents complete inflation of an angioplasty balloon or proper placement of the stent graft

5.0 WARNINGS

This device is intended for single-use only. Do not resterilize or reuse. Note the product "Use by" date specified on the package.

Judicious selection of patients is necessary since the use of this device carries the associated risk of subacute thrombosis, vascular complications, and / or bleeding events.

Persons allergic to 316L stainless steel (including the major elements iron, chromium, nickel, molybdenum), or PTFE may suffer an allergic reaction to this implant.

Implantation of the stent graft should be performed only by physicians who have received appropriate training.

Stent graft placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.

Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent graft. The long-term outcome following repeat dilatation of endothelialized stent grafts is unknown at present.

When multiple stents are required, stent materials should be of similar composition. Placing multiple stents of different metals in contact with each other may increase the potential for corrosion. The risk of *in vivo* corrosion does not appear to increase based on *in vitro* corrosion tests using an L-605 CoCr alloy stent (MULTI-LINK VISION Coronary Stent) in combination with a 316L stainless steel stent (MULTI-LINK TETRA Coronary Stent).

6.0 PRECAUTIONS

6.1 General Precautions

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

Employ aseptic techniques during removal from the package and during use.

Inspect all products prior to use. Do not use if the package is open or damaged, or if product is damaged.

The use of mechanical atherectomy devices or laser catheters is not recommended within the region that a stent graft has been placed.

6.2 Stent Graft Handling – Precautions

This device is intended for single-use only; do not reuse. Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing. Note the "Use by" (expiration) date on the product label.

Do not remove the stent graft from its delivery system. Removing the stent graft from the delivery system may damage the stent graft and / or lead to stent graft embolization.

The delivery system should not be used in conjunction with other stents.

Special care must be taken not to handle or in any way disrupt the stent graft position on the delivery system. This is most important during placement over the guide wire and the advancement through the hemostasis valve adapter and guiding catheter hub.

Do not manipulate, touch or handle the stent graft with your fingers, as this may cause contamination, or dislodgement of the stent graft from the delivery balloon.

Use only appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as it may cause uneven expansion and difficulty in deployment of the stent graft.

Use only diluted contrast medium.

6.3 Stent Graft Placement – Precautions

6.3.1 Stent Graft Preparation

- **Do not prepare or pre-inflate prior to stent graft deployment other than as directed.** Use the balloon purging technique described in section 8.2.3 *Delivery System Preparation*.
- While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent graft from the balloon.
- Use guiding catheters which have lumen sizes that are suitable to accommodate the stent graft delivery system (see package label).
- A large-bore Tuohy-Borst valve with a minimum internal diameter of ≥ 0.074 " (1.88 mm) is recommended.

6.3.2 Stent Graft Implantation

- The vessel should be pre-dilated with an appropriately sized balloon. Failure to do so may increase the difficulty of stent graft placement and cause procedural complications.
- Do not expand the stent graft if it is not properly positioned in the vessel (see section 6.4 *Stent Graft / System Removal – Precautions*).
- Implanting a stent graft may lead to dissection of the vessel distal and / or proximal to the stent graft and may cause closure of the vessel, requiring additional intervention (e.g., coronary artery bypass surgery, further dilatation, placement of additional stents, etc.).

- Placement of a stent graft has the potential to compromise side-branch patency.
- If more than one stent graft is required, the distal stent graft should be placed initially, followed by placement of the proximal stent graft. Stent graft placement in this order obviates the need to cross the proximal stent graft when placing the distal stent graft, and reduces the chances for dislodging the proximal stent graft.
- Do not expand the stent graft if it is not properly positioned in the vessel (see section 6.4 *Stent Graft / System Removal – Precautions*).
- **Do not exceed the rated burst pressure (RBP) as indicated on product label.** Monitor balloon pressures during inflation. Use of pressures higher than specified on the product label may result in a ruptured balloon with possible intimal damage and dissection.
- An unexpanded stent graft may be retracted into the guiding catheter one time only. An unexpanded stent graft should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the stent graft may be damaged when retracting the unexpanded stent graft back into the guiding catheter.
- Should resistance be felt at any time during coronary stent graft system withdrawal, please follow the steps provided in section 6.4 *Stent Graft / System Removal – Precautions*.
- Stent graft retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- Care should be taken to control the guiding catheter tip during stent graft delivery, deployment, and balloon withdrawal. Before withdrawing the stent graft delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.

6.4 Stent Graft / System Removal – Precautions

6.4.1 Stent Graft Removal Prior to Deployment

If removal of a stent graft system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent graft delivery system, and cautiously withdraw the stent graft delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent graft towards the guiding catheter, the stent graft delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

6.4.2 Withdrawal of the Delivery Catheter from the Deployed Stent Graft

1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 – 15 seconds longer.
2. Position the inflation device to “negative” or “neutral” pressure.
3. Stabilize guide catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stented segment.
4. Gently remove the stent graft delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

Note: If during withdrawal of the catheter resistance is encountered, use the following steps to improve balloon re-wrap.

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

Failure to follow these steps and / or applying excessive force to the delivery system can potentially result in loss or damage to the stent graft and / or delivery system components.

If it is necessary to retain guide wire position for subsequent artery / target site access, leave the guide wire in place and remove all other system components.

Stent graft retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.

6.5 Post-Implantation – Precautions

- Care must be exercised when crossing a newly deployed stent graft with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, a balloon catheter, or delivery system, to avoid disrupting the stent graft placement, apposition, and / or geometry.
- Antiplatelet therapy should be administered post-procedure (see section 7.0 *Concomitant Medication*). Patients who require early discontinuation of antiplatelet therapy (e.g., secondary to active bleeding) should be monitored carefully for cardiac events. At the discretion of the patient’s treating physician, the antiplatelet therapy should be restarted as soon as possible.

If the patient requires imaging, see section 6.6 *MRI Statement*.

6.6 MRI Statement

Nonclinical testing has demonstrated that the GRAFTMASTER RX Coronary Stent Graft, in single and in overlapped configurations up to 44 mm in length, is MR conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 or 3 Tesla
- Spatial gradient field of 2500 Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for up to 15 minutes of scanning for each duration of a sequence

The GRAFTMASTER stent graft should not migrate in this MRI environment. Nonclinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of the GRAFTMASTER stent graft.

Stent heating was derived by using the measured nonclinical, *in vitro* temperature rises in a GE Excite 3 Tesla scanner and in a GE 1.5 Tesla coil in combination with the local specific absorption rates (SARs) in a digitized human heart model. The maximum whole-body-averaged SAR was determined by validated calculation. At overlapped lengths of up to 44 mm, the GRAFTMASTER stent graft produced a nonclinical maximum local temperature rise of 1.8°C at a maximum whole-body-averaged SAR of 2.0 W/kg (normal operating mode) for 15 minutes. These calculations do not take into consideration the cooling effects of blood flow.

The effects of MRI on overlapped stents greater than 44 mm in length or stents with fractured struts are unknown.

As demonstrated in nonclinical testing, an image artifact can be present when scanning the GRAFTMASTER stent graft. MR image quality may be compromised if the area of interest is in the exact same area, or relatively close to, the position of the GRAFTMASTER stent graft. Therefore, it may be necessary to optimize the MR imaging parameters for the presence of the GRAFTMASTER stent graft.

It is suggested that patients register the conditions under which the implant can be safely scanned with the MedicAlert Foundation (www.medicalert.org) or an equivalent organization.

7.0 CONCOMITANT MEDICATION

This information is provided as a guide and is not intended to dictate medical practice. The conditions under which the device has been used to date include:

- Aspirin prior to the procedure
- Heparin and nitroglycerine during the procedure
- Anticoagulants / antiplatelets or antiplatelets alone post-treatment up to one month

Most experience to date with the coronary stents has been made with anticoagulant / antiplatelet therapy as compared to the use of antiplatelets alone. Administration of antiplatelets and anticoagulants will vary with patient conditions and medical histories. Prior anticoagulant / antiplatelet management is an essential adjunct to stent implantation. Proper heparin management before and after sheath removal is also critical to the overall success of the procedure.

8.0 OPERATOR'S MANUAL

8.1 Materials Required

- Appropriate guiding catheter(s)
- 2 – 3 syringes (10 – 20 cc)
- 1,000 u/500 cc heparinized normal saline (HepNS)
- 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
- Rotating hemostatic valve with 0.096 inch (2.44 mm) minimum inner diameter
- Contrast material diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device
- Guide wire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

8.2 System Preparation

Note: During stent graft deployment, with the stent delivery system from crimped state to 4.80 mm, the graft may shorten up to 20%.

8.2.1 Packaging Removal

Carefully remove the delivery system from its protective tubing for preparation of the delivery system. Do not bend or kink the hypotube during removal.

Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product and replace with another. Follow product returns procedure for the unused device.

Visually inspect the stent graft for uniformity, protruding coils, centering on the balloon and verify that the stent does not extend beyond the radiopaque balloon markers.

8.2.2 Guide Wire Lumen Flush

Flush the guide wire lumen with HepNS using the flushing tool supplied with the product. Insert the flushing tool into the tip of the catheter and flush until fluid exits the guide wire exit notch.

Note: Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent graft on the balloon.

8.2.3 Delivery System Preparation

1. Prepare an inflation device / syringe with diluted contrast medium.
2. Attach an inflation device / syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.
10. Moisten the stent graft with heparinized saline by submerging the stent graft into a sterile bowl containing the solution.

Note: Do not use gauze sponges to wipe down the stent graft, as fibers may disrupt the stent graft.

Note: While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent graft from the balloon.

Note: If air is seen in the shaft, repeat steps 3 through 5 to prevent uneven stent graft expansion.

8.3 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
2. **Pre-dilate the treatment area with a percutaneous transluminal coronary angioplasty (PTCA) catheter of appropriate length and diameter for the vessel / site to be treated.**
3. Predilatations of the vessel must take into account proximal atherosclerotic plaque beyond the treatment area which may prevent advancement of the device to the treatment area.
4. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
5. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the treatment area.
6. Carefully advance the delivery system into the guiding catheter and over the guide wire to the treatment area. Be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent graft system into the coronary artery.
7. Advance the delivery system over the guide wire to the treatment area under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent graft across the treatment area. Perform angiography to confirm stent graft position. If the position of the stent graft is not optimal, it should be carefully repositioned or removed. Expansion of the stent graft should not be undertaken if the stent graft is not properly positioned in the treatment area.

CAUTION: If resistance is encountered, do not force passage. Resistance may indicate damage to the device or movement of the stent graft on the balloon.

Note: If removal of a stent graft system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent graft delivery system and cautiously withdraw the stent graft delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent graft towards the guiding catheter, the stent graft delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

8. Tighten the rotating hemostatic valve. The stent graft is now ready to be deployed.

CAUTION: Avoid over-tightening the Tuohy-Borst valve, as this may restrict the flow of contrast medium in and out of the balloon, thereby slowing inflation / deflation.

8.4 Deployment Procedure

CAUTION: Refer to product label for *in vitro* stent graft outer diameter, nominal pressure, and RBP.

1. Prior to deployment, reconfirm the correct position of the stent graft relative to the treatment area using the radiopaque balloon markers.
2. Deploy the stent graft slowly by pressurizing the delivery system in 2-atm increments every 5 seconds until stent graft is completely expanded. Fully expand the stent graft by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a stent graft inner diameter ratio of about 1.1 times the reference vessel diameter (refer to product label for *in vitro* stent graft inner diameter, nominal pressure, and RBP).
3. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the stent graft to the artery wall. Maintain pressure for 30 seconds for full expansion of the stent graft. Fluoroscopic visualization during stent graft expansion should be used in order to properly judge the optimum stent graft diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal stent graft expansion and proper apposition require that the stent graft be in full contact with the arterial wall.

Note: See section 8.5 Removal Procedure for instruction on withdrawal of stent graft delivery system.

4. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the stent graft to the artery wall.

CAUTION: Do not exceed the labeled rated burst pressure (RBP) of 16 atm (1621 kPa).

5. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent graft delivery system withdrawal, pay particular attention to guiding catheter position.
6. Confirm stent graft position and deployment using standard angiographic techniques. Fluoroscopic visualization during stent graft expansion should be used in order to properly judge the optimum expanded stent graft diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent graft be in full contact with the artery wall. Stent graft wall contact should be verified through routine angiography.
7. If the deployed stent graft size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent graft. If the initial angiographic appearance is suboptimal, the stent graft may be further expanded using a low profile, high pressure, noncompliant balloon dilation catheter. If this is required, the stented segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the stent graft geometry. Deployed stent graft should not be left underdilated.

Note: Maximum post dilatation that can be achieved with a noncompliant post dilatation balloon is a maximum of 5.5 mm. With expansion to this diameter, the system may shorten by approximately 25%. When choosing a GRAFTMASTER RX system for expansion in larger vessels, a longer stent graft length is recommended to ensure treatment area is covered by the stent graft.

CAUTION: Do not dilate the stent graft beyond 5.5 mm.

8.5 Removal Procedure**8.5.1 Withdrawal of the Delivery Catheter from the Deployed Stent Graft:**

1. Deflate the balloon by pulling negative pressure on the inflation device. Confirm balloon deflation under fluoroscopy and wait 10 – 15 seconds longer.
2. Position inflation device on "negative" or "neutral" pressure.
3. Stabilize guiding catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stented segment.
4. Gently remove the stent graft delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

If during withdrawal of the stent graft delivery catheter resistance is encountered, use the following steps to improve balloon rewrap:

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

9.0 POTENTIAL ADVERSE EVENTS

Possible adverse events include, but are not limited to the following:

- Acute myocardial infarction
- Acute or subacute closure of the coronary artery
- Acute or subacute closure of the stent graft
- Allergic reaction to implanted stent material (316L) and / or PTFE
- Arteriovenous fistula
- Cardiac arrhythmia, including ventricular fibrillation
- Cardiac tamponade
- Cardiogenic shock
- Coronary artery dissection, perforation, rupture or injury
- Coronary artery spasm
- Death
- Deformation of stent graft symmetry on deployment or on recrossing with another catheter (IVUS, balloon, etc.)
- Distal emboli (air, particle, thrombogenic emboli)
- Disturbance of cardiac conductivity, rhythm
- Drug reactions, allergic reaction to contrast medium
- Emergency coronary artery bypass graft surgery
- Fever
- Hemorrhage or hematoma
- Hypo / hypertension
- Infection
- Inflammation
- Intimal tear
- Ischemia, myocardial
- Pseudoaneurysm
- Renal failure
- Restenosis of the treated artery
- Stent graft embolization
- Stroke / cerebrovascular accidents
- Subscribed tissue reaction
- Thrombosis of the coronary artery
- Total occlusion of the coronary artery
- Unstable angina

10.0 REFERENCES

The physician should consult recent literature on current medical practice on stent implantation, such as that published by the American College of Cardiology / American Heart Association (ACC / AHA).

11.0 PRODUCT INFORMATION DISCLOSURE

Abbott Vascular has exercised reasonable care in the manufacture of this device. Abbott Vascular excludes all warranties, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device as well as factors relating to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond Abbott Vascular's control directly affect this device and the results obtained from its use. Abbott Vascular shall not be liable for an incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

12.0 PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 5,916,264; 6,068,656. Other U.S. patents pending. Foreign patents issued and pending.

GRAFTMASTER, HYDROCOAT, and MULTI-LINK VISION are trademarks of the Abbott Group of Companies.

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