OMNILINK ELITE

Peripheral Stent System

CE
0086

EN	English	IT	Italiano	HU	Magyar
DE	Deutsch	FI	Suomi	BG	Български
SV	Svenska	FR	Français	RO	Română
DA	Dansk	ES	Español	SK	Slovensky
PL	Polski	RU	Русский	ET	Eesti keel
CS	Česky	PT	Português	LV	Latviešu
TR	Türkçe	NL	Nederlands	LT	Lietuvių
NO	Norsk	EL	Ελληνικά	SL	Slovenščina

Graphical Symbols for Medical Device Labeling; Графични символи за етикетиране на медицински уреди; Grafické symboly k označení zdravotnického prostředku; Grafiske symboler beregnet til mærkning af medicinske produkter; Grafische symbolen voor de etikettering van medische hulpmiddelen; Meditsiiniseadme märgistuse graafilised sümbolid; Lääketieteellisten laitteiden tuotetarroissa esiintyvät symbolit; Symboles graphiques pour l'étiquetage des dispositifs médicaux ; Grafische Symbole zur Beschriftung von Medizinprodukten; Γραφικά σύμβολα για τη σήμανση ιατρικών συσκευών; Grafikus szimbólumok orvosi eszközök címkézéséhez; Simboli grafici per l'etichettatura di dispositivi medicali; Grafiskie simboli medicīnisko ierīču apzīmēšanai; Medicinos prietaisų etikečių grafiniai simboliai; Symbole graficzne do oznaczania sprzętu medycznego; Símbolos Gráficos para Rotulagem de Dispositivos Médicos; Simboluri grafice pentru etichetarea dispozitivelor medicale; Графические символы на этикетках медицинских изделий; Grafické symboly na označovanie zdravotníckej pomôcky; Grafčni simboli za označevanje medicinskega pripomočka; Símbolos gráficos utilizados en el etiquetado de productos sanitarios; Grafiska symboler för märkning av medicinsk utrustning; Tibbi Cihaz Etiketleri için Grafik Semboller; Grafiske symboler for merking av medisinsk utstyr

INSTRUCTIONS FOR USE

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
Peripheral Stent System

1.0 DEVICE DESCRIPTION

The Omnilink Elite Peripheral Stent System (Omnilink Elite Stent System) is a flexible, balloon-expandable L605 cobalt-chromium stent pre-mounted on the balloon of an over-the-wire (OTW) stent delivery system. The OTW stent delivery system is compatible with a 0.035" (0.89 mm) guide wire and comes in lengths of 80 cm and 135 cm. The stent is mounted on the balloon between the two radiopaque markers. The delivery system can be utilized to optimize the stent wall apposition post stent deployment. See *Table 1* for stent dimensions.

Table 1: *In vitro Device Specifications**

Expanded Stent Diameter** (mm)	Stent Lengths (mm)	<i>In vitro</i> * Stent Deployment Pressure (bar)	Rated Burst Pressure RBP (bar)	Recommended Minimum Sheath / Introducer*** (F)
4.0	12, 16, 19	11	16	6
5.0	12, 16, 19, 29, 39, 59	11	16	6
6.0	12, 16, 19, 29, 39, 59	11	14	6
7.0	12, 16, 19, 29, 39, 59	11	14	6
8.0	19, 29, 39, 59	11	14	6
9.0	19, 29, 39, 59	11	14	6
10.0	19, 29, 39	11	14	6

* All data provided are based on *in vitro* testing. Assure full deployment of the stent. Deployment pressures should be based on lesion characteristics.

** The inflated balloon diameter of the system used to deploy the stent should approximate the diameter of the vessel.

*** See individual manufacturer specifications for (F) equivalent on box label and pouch label.

2.0 HOW SUPPLIED

Contents. One each: Omnilink Elite Peripheral Stent System.

Storage. Keep dry, keep away from sunlight, temperature limitation 15°C - 30°C (59°F - 86°F).

Non Pyrogenic.

Sterile. Sterilized with electron beam radiation. Do not use if package is open or damaged.

This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

3.0 INDICATIONS

The Omnilink Elite Peripheral Stent System is indicated for the treatment of *de novo* or restenotic atherosclerotic lesions in protected peripheral arteries and palliation of malignant strictures in the biliary tree.

4.0 CONTRAINDICATIONS

The Omnilink Elite Peripheral Stent System is contraindicated for use in:

Vessel Use:

- Patients with a contraindication for anti-platelet / anti-coagulant therapy
- Aneurysms immediately adjacent to the stent implantation site
- Presence of fresh, unlysed, unorganized thrombus or embolic material
- Target lesions that are resistant to complete balloon inflation
- Patients with a known hypersensitivity to cobalt or chrome
- Patients with a known hypersensitivity to silicone

Biliary use:

- Patients with bleeding disorders
- Use in perforated duct where the leakage can be enhanced by the prosthesis
- Severe ascites
- Patients with a known hypersensitivity to cobalt or chrome
- Patients with a known hypersensitivity to silicone

5.0 WARNINGS

- Do not use if package is open or damaged.

Should **unusual resistance** be felt **at any time** during lesion access or stent delivery system removal, the introducer sheath and stent delivery system **should be removed as a single unit**. If possible, retain the guide wire position for subsequent vessel access. Applying excessive force to the stent delivery system can potentially result in loss or damage to the stent and stent delivery system components. (See *Stent / System Removal – Precautions*)

- Since the use of this device carries the associated risk of subacute thrombosis, vascular complications and / or bleeding events, judicious selection of patients is necessary.
- Once fully deployed, the stent cannot be repositioned.
- This device should be used only by physicians trained in angiography and percutaneous transluminal angioplasty and / or biliary stent placement and cholangiography.

- The Omnilink Elite Peripheral Stent System is intended to perform as a system. The stent should not be removed for use in conjunction with other dilatation catheters, nor should the Omnilink Elite Peripheral Stent System be used in conjunction with other stents.
- When multiple stents are required, stent materials should be of similar composition.

6.0 PRECAUTIONS

Subsequent restenosis may require repeat dilatation of the vessel segment containing the stent. The long-term outcome following repeat dilatation of endothelialized Omnilink Elite peripheral stents is unknown at present.

Estimated vessel / duct reference size should not be less than 4.0 mm in diameter or bigger than 11.0 mm in diameter.

6.1 Stent Delivery System Handling – Precautions

- This device is designed for single use only. Do not resterilize or reuse.
- Use the stent system prior to the "Use By" date specified on the package.
- **Do not remove the stent from its delivery balloon**, as removal may damage the stent and / or lead to stent embolization.
- Carefully inspect the Omnilink Elite Peripheral Stent System prior to use to verify that the stent has not been damaged in shipment and that the device dimensions are suitable for the specific procedure. Take care to avoid unnecessary handling.
- Do not "roll" the mounted stent with your fingers, as this action may loosen the stent from the delivery balloon.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the stent.
- Refer to the instructions for use supplied with any interventional devices to be used in conjunction with Omnilink Elite Peripheral Stent System, for their intended uses, contraindications, and potential complications.
- Never advance delivery balloon catheter without the guide wire extending from the tip.

6.2 Stent Placement – Precautions

- Do not prepare or pre-inflate the balloon prior to stent deployment other than as directed. Use the balloon-purging technique described in the *Clinician Use Manual* in section 8.0.
- The expanded stent outer diameter should approximate the diameter of the vessel / duct. Oversizing of the stent can result in a ruptured vessel / duct. To ensure full expansion of the stent, the balloon should be inflated to a minimum of nominal pressure.
- Implanting a stent may lead to dissection or perforation of the vessel / duct distal and / or proximal to the stent and may cause acute closure of the vessel / duct requiring additional intervention (surgical intervention, further dilatation, placement of additional stents, or other).
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the

proximal stent in placement of the distal stent and reduces the chances for dislodging the stent.

- Do not expand the stent if it is not properly positioned in the vessel / duct. (See *Stent / System Removal – Precautions*)
- Fretting corrosion on overlapped Omniliink Elite stents in bile ducts has not been tested.
- Stenting across a major bifurcation may hinder or prevent future side branch access.
- Balloon pressures should be monitored during inflation. **Do not exceed the Rated Burst Pressure (RBP) as indicated on product label.** Use of pressures higher than specified on the product label may result in a ruptured balloon with possible vessel / duct damage or perforation.
- **Do not attempt to pull an unexpanded stent back through the introducer sheath; dislodgement of the stent from the balloon may occur.**
- Stent retrieval methods (use of additional wires, snares and / or forceps) may result in additional trauma to the duct / vasculature and / or the duct / vascular access site. Complications may include bleeding, hematoma or pseudoaneurysm.
- The Omniliink Elite Peripheral Stent System is intended for deployment and post-deployment dilatation of the stent only and should not be used to dilate other locations.

6.3 Stent / System Removal – Precautions

Should **unusual resistance** be felt at **any time** during either vessel / duct access or during removal of an undeployed stent, the stent delivery system and introducer sheath should be **removed as a single unit**. If possible, retain the guide wire position for subsequent vessel / duct access.

When removing the stent delivery system as a single unit:

- DO NOT retract the stent delivery system into the introducer sheath.
- Position the proximal balloon marker just distal to the tip of the introducer sheath.
- Advance the guide wire in the anatomy as far distally as safely possible.
- Secure the stent delivery system to the introducer sheath; then remove the introducer sheath and stent delivery system as a **single unit**.

Failure to follow these steps and / or applying excessive force to the stent delivery system can potentially result in loss or damage to the stent and / or stent delivery system components.

6.4 Post Implant – Precautions

Great care must be exercised when crossing a **newly deployed stent** with a guide wire or balloon catheter to avoid disrupting the stent geometry.

Magnetic Resonance Imaging (MRI) Information:

Non-clinical testing has demonstrated that the Omniliink Elite stent, in single and in overlapped configurations up to 100 mm in length, is MR Conditional as defined in ASTM F2503. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Spatial gradient field of 2500 Gauss/cm or less

- Maximum whole-body-average specific absorption rate (WB-SAR) of 2 W/kg for 15 minutes of scanning per sequence for patient landmarks above umbilicus (Total duration of all scans may exceed 15 minutes)
- Maximum WB-SAR of 1 W/kg for 15 minutes of scanning for patient landmarks below umbilicus
- Transmit RF body coil should be used in normal operating mode, as defined in IEC 60601-2-33

MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of the Omniliink Elite stent.

The Omniliink Elite stent should not migrate in this MRI environment. Magnetic force on the Omniliink Elite stent was tested according to ASTM F2052-06e. The deflection angle was measured to be 6° in a GE Signa 3 T HDx MR system. Stent heating was derived by using the measured non-clinical, *in vitro* temperature rise according to ASTM F2182-09 in a GE Signa HDx 3 Tesla scanner and in a GE 1.5 Tesla coil in combination with the local specific absorption rates (SARs) in a digitized human model. For the SAR conditions above, the maximum *in vivo* temperature rise was calculated to be 6°C at 64 MHz (1.5 T) and 128 MHz (3 T) for stent lengths 100 mm and less. These calculations do not take into consideration the cooling effects of blood flow.

The effects of MRI on overlapped stents greater than 100 mm in length or stents with fractured struts are unknown.

Image artifact may be present when scanning the Omniliink Elite stent as demonstrated in nonclinical testing performed according to ASTM F2119-07 in a GE Signa HDx 3 Tesla scanner.

The image artifact (both inside and outside the device lumen) extends approximately 5 mm from the device, using the spin echo sequence, and 10 mm from the device, using the gradient echo sequence. MR image quality may be compromised if the area of interest is in the exact same area, or relatively close to, the position of the Omniliink Elite stent. Therefore, it may be necessary to optimize the MR imaging parameters in the presence of Omniliink Elite stents.

7.0 POTENTIAL ADVERSE EFFECTS

Possible complications include, but are not limited to, the following:

For Vessel Use

- Abrupt closure
- Allergic reaction (contrast medium; drug; stent material)
- Aneurysm or pseudoaneurysm in vessel or at vascular access site
- Angina or coronary ischemia
- Arrhythmia (including premature beats, bradycardia, atrial or ventricular tachycardia, arterial or ventricular fibrillation)
- Arteriovenous fistula
- Bleeding complications requiring transfusion or surgical intervention
- Death
- Detachment of a system component or implantation in an unintended site
- Embolization (air, tissue, plaque, thrombotic material, stent)
- Emergent surgery

- Fever
- Hematoma or hemorrhagic event, with or without surgical repair
- Hyperperfusion syndrome
- Hypotension / hypertension
- Infection
- Ischemia or infarction not covered under other adverse effects
- Myocardial infarction
- Pain (leg, foot, and / or insertion site)
- Partial stent deployment
- Peripheral nerve injury
- Pulmonary embolism
- Renal failure or insufficiency
- Restenosis
- Shock
- Stent malapposition or migration, which may require emergency surgery to remove stent
- Stent strut fracture
- Stroke
- Target limb loss (amputation of toe, foot, and / or leg)
- Thrombosis or occlusion
- Transient ischemic attack
- Venous thromboembolism
- Vessel dissection, perforation, or rupture
- Vessel spasm or recoil
- Worsening claudication or rest pain

For Biliary Use

- Abscess
- Additional intervention due to but not limited to:
 - Stent migration
 - Unintentional placement of stent
 - Partial stent deployment
 - Damaged stents
- Bile duct occlusion / obstruction
- Bile duct perforation potentially leading to infection or death
- Cholangitis
- Pancreatitis
- Parenchymal hemorrhage
- Peritonitis
- Sepsis
- Tumor overgrowth at the stent ends

8.0 CLINICIAN USE MANUAL

8.1 Stent Inspection Prior To Use

Prior to using the Omniliink Elite Peripheral Stent System, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent is located between the radiopaque balloon markers. Do not use if any defects are noted.

8.2 Materials Required

- Introducer sheath in the appropriate size and configuration for the selected stent delivery system (refer to packaging label or *Table 1*)
- Two to three 10-20 cc syringes
- 1000 u / 500 cc heparinized normal saline (HepNS)
- One 0.035" (0.89 mm) guide wire of appropriate length
- Contrast diluted 1:1 with normal saline
- One inflation device
- One three-way stopcock
- One torque device (if applicable)
- One guide wire introducer

8.3 Lesion or Stricture Preparation

1. Standard percutaneous technique should be used to place the introducer sheath in the vessel or biliary tree. An appropriate sized (0.035") guide wire should be advanced across the lesion or stricture and into the common vessel or bile duct.
2. Pre-dilate the lesion or stricture with an appropriate size balloon dilatation catheter to closely match the lumen diameter proximal and distal to the lesion or stricture.
3. Withdraw the balloon dilatation catheter leaving the guide wire in place.

8.4 Guide Wire Lumen Flush

1. Remove the protective cover from the tip.
2. Attach the syringe with HepNS to the guide wire port.
3. Flush until fluid exits the distal tip.

8.5 Stent Delivery System Preparation

1. Prepare the inflation device / syringe with diluted contrast medium.
2. Attach the inflation device / syringe to the stopcock; attach to the inflation port.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to the delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled.
Note: If air is seen in the shaft, repeat balloon preparation steps 3 through 5 to prevent uneven stent expansion.
7. If a syringe was used, attach a prepared inflation device to the stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

8.6 Stent Delivery Procedure

1. Wipe the exposed guide wire with heparinized saline to remove residual blood or contrast medium.
2. Fully open the hemostatic valve. Maintain neutral pressure on the inflation device.
3. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion or stricture.

4. Advance the delivery system over the guide wire to the target lesion or stricture. Utilize the radiopaque balloon markers to position the stent across the lesion or stricture; perform angiography to confirm the stent position. If applicable tighten the hemostatic valve.

Note: If during the process of moving the delivery system into position you notice the stent has moved on the balloon, do not deploy the stent. The entire system should be **removed as a single unit**. See *Stent / System Removal - Precautions* section for specific delivery system removal instructions.

5. The stent is now ready to be deployed.

8.7 Stent Deployment Procedure

CAUTION: Refer to the product label for expanded stent outer diameter, nominal deployment pressure and rated burst pressure.

1. Slowly inflate the delivery balloon to low pressure; hold until balloon inflation is observed both proximally and distally to the stent. Continue balloon expansion to the specified stent deployment pressure. Confirm complete expansion of the stent / balloon fluoroscopically. If necessary, the delivery balloon can be used to post-dilate the stent to optimize stent apposition.

Do not exceed the RBP. A larger PTA catheter may be used to dilate stent.

- Do not expand the 4 - 5 mm stent beyond 6 mm maximum.
 - Do not expand the 6 - 7 mm stent beyond 8 mm maximum.
 - Do not expand the 8 - 10 mm stent beyond 11 mm maximum.
2. Deflate the balloon by pulling negative pressure on the inflation device. Ensure that the balloon is fully deflated.

8.8 Removal Procedure

1. Maintain **negative pressure** to allow the balloon to remain fully deflated during removal through the sheath.
2. With the inflation device on **negative pressure**, carefully withdraw the delivery catheter with the guide wire remaining across the lesion or stricture.

Note: Should **unusual resistance** be felt **at any time** during either lesion or stricture access or removal of delivery system post-stent implantation, the entire system should be **removed as a single unit**. See *Stent / System Removal - Precautions* section for specific delivery system removal instructions.

3. Confirm optimal stent apposition using standard angiography techniques. If necessary, post-dilate within the stent. Balloon inflations should closely match vessel / bile duct reference diameter.

ASSURE THAT THE STENT IS NOT UNDERDILATED.

The final stent diameter should match the reference vessel / duct.

9.0 REFERENCES

The physician should consult current literature on current medical practice on balloon dilatation.

10.0 PRODUCT INFORMATION DISCLOSURE

Abbott Vascular has exercised reasonable care in the manufacture of this device. **Abbott Vascular** excludes all warranties, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device as well as factors relating to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond **Abbott Vascular's** control directly affect this device and the results obtained from its use. **Abbott Vascular** shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. **Abbott Vascular** neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

11.0 TRADEMARKS

Omnilink Elite is a trademark of the Abbott Group of Companies.

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©2014 Abbott