
Xpert Pro


Peripheral Self-Expanding Stent System

INSTRUCTIONS FOR USE
НАЧИН НА УПОТРЕБА
NÁVOD NA POUŽITÍ
GEBRAUCHSANWEISUNG
BRUGSANVISNING
INSTRUCCIONES DE USO
KASUTUSJUHIS
MODE D'EMPLOI
KÄYTTÖOHJEET
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
HASZNÁLATI UTASÍTÁS
ISTRUZIONI PER L'USO
NAUDOJIMO INSTRUKCIJOS
LIETOŠANAS INSTRUKCIJAS
GEBRUIKSAANWIJZINGEN
BRUKSANVISNING
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUKCJA UŻYTKOWANIA
INSTRUCȚIUNI DE UTILIZARE
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
BRUKSANVISNING
NÁVOD NA POUŽITÍ
NAVODILA ZA UPORABO
KULLANIM TALİMATLARI

Xpert Pro

Peripheral Self-Expanding Stent System

I. DEVICE NAME AND DESCRIPTION

Xpert Pro Peripheral Self-Expanding Stent System is the brand name. It consists of a self-expanding stent integrated into a delivery system, intended for use in peripheral vasculature or the biliary duct.

The stent

The stent is a self-expanding nitinol stent. The specially treated nitinol material gives the stent properties such as non-crushability. When exposed to body temperature, the stent expands to a preprogrammed diameter. The stent has been specifically designed to show a high degree of flexibility both prior to stent deployment and after implantation in the vessel/duct, and to exert radial strength to establish vessel/duct patency.

The delivery system

The stent (7) is integrated at the distal end of the delivery system inside the outer tube (6), and rests onto an inner tube in which the guide wire channel runs. To facilitate positioning, two radiopaque markers indicate the stent's proximal (8) and distal (9) end. The inner tube is connected distally to an atraumatic soft tip (4) where the guide wire entry is located. Between the outer and inner tube is a middle tube (3), which ends distally at the proximal end of the mounted stent.

The inner and middle tube are connected proximally to a metal shaft (2), which runs through the Tuohy Borst adapter and ends proximally with a luer lock (5) where the guide wire exits. The outer tube is connected proximally to a Tuohy Borst adapter (1). The Tuohy Borst valve is shipped locked. It should be kept locked during handling and introduction of the system to avoid premature deployment of the stent. It is opened by rotating the proximal valve end (10) in a counter clockwise direction. During stent deployment, the metal shaft should be held in a fixed position whereas the opened Tuohy Borst valve together with the outer tube is retracted in a proximal direction over the metal shaft to release the stent.

How the system is supplied

Xpert Pro is supplied sterile and is intended for single use only. This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and/or chemical characteristics introduced under conditions of repeated use, cleaning, and/or resterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or spaces and diminished safety and/or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and/or user.

Storage

Xpert Pro should be stored in a cool, dark and dry place. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that the products are used prior to the use by date on the package label.

MRI Compatibility

The Xpert Pro stent has been shown through non-clinical testing to be MRI conditional (poses no known hazards) immediately following implantation when scanned under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient magnetic field of 2500-Gauss/cm or less
- Normal mode of MR system.

Stent heating was estimated by relating the measured in vitro temperature rises in a 1.5T GE Signa coil and a 3.0T GE Signa HDx MR system to the local specific absorption rates in worst-case locations of limb joints in a human body model. The Xpert Pro stents produced an estimated maximum local temperature rise of 5.4°C.

Different SAR values should be utilized depending on placement of the stent within the body. Placement in the arms: WB SAR should be less than 1 W/kg.

Placement in iliac or SFA: WB SAR less than 2 W/kg for landmark above the umbilicus and 0.5 W/kg for landmark below the umbilicus.

Placement in Tibial artery below the knee: WB SAR less than 2 W/kg for landmark above upper thigh and 0.5 W/kg below upper thigh.

Placement in Biliary: WB SAR less than 2 W/kg.

At the SAR levels listed above, the presence of the Xpert Pro should not expose the patient to added risk of thermal harm versus non-implant conditions during MRI. Nonetheless, it is possible that there will be thermal harm to the patient.

The effect of MRI-related heating for overlapping stents or stents with fractured struts is unknown. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Xpert Pro stent.

II. INDICATIONS

Vascular

Insufficient result after PTA, e.g. residual stenosis, dissection, obstruction due to detached atherosclerotic plaque material or reocclusion.

Biliary

Bile duct obstruction caused by malignant tumors (treated via transhepatic approach).

Malignancy

- Patients presenting pancreatic carcinoma which is compressing the biliary tree
- Patients presenting cholangiocarcinoma which has led to biliary obstruction
- Patients presenting an extra-hepatic biliary obstruction due to gallbladder carcinoma, metastases, or ampullary carcinoma.

III. CONTRAINDICATIONS

Vascular contraindications

Contraindications relevant for PTA, including but not limited to:

- History of coagulation disorders
- Calcified lesions resistant to PTA
- Functionally relevant obstruction of the inflow path, poor outflow or no distal run-off
- Presence of fresh unlysed unorganized thrombus or embolic material
- Bifurcation lesions where the stent placement would obstruct access to a major branch which may in the future benefit from PTA

Biliary contraindications

- All of the customary contraindications to percutaneous transhepatic manipulation of catheters of the same caliber as Xpert Pro, including but not limited to:
 - Ascites
 - Extensive metastatic disease or tumor spread through the liver
 - Intra-arterial chemotherapy

Also including

- stenting a perforated duct, whereby leakage from the duct could be exacerbated by the self-expandable endoprosthesis and leakage could occur across the mesh of the stent
- Failure to pass the stricture with a guide wire

Relative contraindication:

- Stenting across a major bifurcation could lead to compromised future diagnostic or therapeutic procedures

IV. WARNINGS

- Single use only – do not resterilize. Do not use if inner package is damaged or opened.
- Use prior to the "use by" date. Store in a cool, dry and dark place – excessive heat could have major negative influence on the performance of the system.
- Once the expansion of the stent has begun, **the stent cannot be retracted into the delivery system or repositioned.**
- The system must always be introduced, moved or withdrawn over a guide wire (0.014" up to max. 0.018").

V. PRECAUTIONS

- This device should only be used by physicians who have a thorough understanding and experience in the clinical and technical aspects of PTA/stenting.
- When the system is introduced into the vascular system/bile duct, it should only be manipulated under high quality fluoroscopy.
- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass Xpert Pro delivery catheter through a smaller sized sheath introducer than indicated on the label.
- The delivery system is not designed for the use of power injection.
- Do not use with Lipiodol or Ethiodol contrast media.
- Do not advance the Xpert Pro device against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- The stent should deploy easily. Do not release it if unusual force is required, since this indicates a failed system. Use a new system instead.

VI. POTENTIAL COMPLICATIONS

The following complications may occur as a result of PTA/stenting, but may not be limited to:

- Abrupt closure
- Access site hematoma
- Aneurysm
- Angina
- Arrhythmias
- Arteriovenous fistula
- Bleeding complications which may require transfusion
- Cerebral ischemia/transient ischemic attack (TIA)
- Death
- Embolism (air, tissue, thrombotic, systemic or device component)
- Fever/pyrogenic reaction
- Hypersensitivity or allergic reaction to contrast agents and drug reactions
- Hypertension/hypotension
- Infection
- Ischemia, including tissue ischemia, steal syndrome and necrosis
- Leg edema
- Myocardial ischemia or infarction
- Nausea and vomiting
- Neuropathies or nerve injury
- Occlusion
- Organ failure (single, multiple)
- Pain
- Palpitations
- Pseudoaneurysm
- Renal failure/insufficiency
- Restenosis
- Stroke/cerebrovascular accident (CVA)
- Vascular complications, including entry site, which may require vessel repair
- Vascular thrombosis
- Vessel injury, e.g. dissection, perforation
- Vessel spasm
- Complications relevant to transhepatic stent implantations, e.g.
 - Abscess
 - Bile duct occlusion/obstruction
 - Bile duct perforation
 - Bleeding
 - Cholangitis
 - Death
 - External biliary fistula
- Hypersensitivity or allergic reaction (to contrast agents, metal alloys of stent, and drug reactions)
 - Infection
 - Peritonitis

- Pancreatitis
- Sepsis
- Tumor overgrowth at stent ends

VII. PRE-PROCEDURAL PREPARATIONS

Preparation of patient and concomitant medications

Patient preparation and sterile precautions should be the same as for any angioplasty procedure/transhepatic catheter procedure.

The medication is decided by the physician, including precautions to reduce clotting during and after the procedure according to latest scientific guidelines and with respect to the individual patient.

The intervention should be performed in an angiography procedure room and carried out under X-ray control.

Angiography/Cholangiography

Angiography/Cholangiography is performed to map out the lesions and the collateral flow. In cases of thrombus formation, thrombolysis should precede stent deployment using standard methods. Ensure that access vessels are sufficiently patent, or sufficiently recanalized, before proceeding with the intervention.

Pre-dilatation of lesion

Perform a PTA balloon dilatation/bile duct dilatation leaving the guide wire in place across the stenosis.

Check the vessel and stenosis dimensions carefully after dilatation in order to choose the correct stent size.

VIII. SELECTION AND PREPARATION OF DEVICE AND COMPATIBILITY WITH ACCESSORIES

Selection of stent size

The stent diameter should be slightly larger than the diameter of the reference vessel (i.e. healthy vessel adjacent to the lesion)/ dilated bile duct. Refer to labeling for reference vessel diameter information.

The stent length should cover the lesion completely.

The safety and effectiveness of multiple overlapping stents have not been established. However, when multiple stents are required, stent materials should be of similar composition.

Preparation of stent delivery system

Remove the system from the package. Ensure that the sterile barrier (steri-pouch) is not broken or damaged.

Warning: If broken or damaged – do not use the system.

Inspect the system visually for damages that could influence the performance. Inspect that the stent is fully contained within the outer tube. Do not use if the distal marker is visible or the stent is partially deployed. Make sure that the Tuohy Borst valve is locked by rotating in a clockwise direction.

The complete system should be flushed at the Tuohy Borst adapter and the proximal luer lock with heparinised saline solution before use.

- Attach a syringe filled with the heparinised saline solution to the lateral branch of the y-connector and apply positive pressure. Ensure that the Tuohy Borst valve is locked and continue to flush until the heparinised saline solution drips from the distal part of the outer tube.

- Flush the central guide wire lumen at the proximal luer lock.

IX. STENT IMPLANTATION PROCEDURE

Introduction of the system

After the initial balloon dilatation procedure, introduce the system through the introducer over a guide wire (0.014" up to max. 0.018") which has been left positioned through the stenosis/stricture site.

Warning: The system must always be introduced, moved or withdrawn over a guide wire. **Note:** Always use an introducer sheath with a hemostasis valve. At contralateral deliveries, it is recommended to use a long introducer sheath specially designed for contralateral deliveries which facilitates the crossing over the bifurcation.

Advance the system under fluoroscopic guidance and position the stent with the help of the radiopaque markers at the stenosis/stricture.

Caution: Do not advance the Xpert Pro device against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken. Withdraw the system and use a new one.

Deployment of the stent

Ensure correct stent position and prior to stent deployment, unlock the Tuohy Borst valve (10).

Hold the metal shaft in a fixed position and retract the outer tube by slowly pulling the Tuohy Borst valve over the metal shaft towards the proximal end.

Note: During deployment, the proximal metal shaft should be kept in a fixed position as the stent is released by retracting the outer sheath. **Also observe and keep the proximal stent marker in a fixed position for optimal stent placement.**

Caution: Do not push the metal shaft, as this will cause misalignment of the stent and possible vascular/ductal damage.

Warning: Once expansion has begun, the stent **cannot be retracted into the delivery system or repositioned.**

Caution: The stent should deploy easily. Do not release if unusual force is required, since this indicates a failed system. Use a new system instead. Deployment is complete when the Tuohy Borst adapter meets the proximal luer lock.

X. POST STENT IMPLANTATION

Withdrawal and disassembly procedure

After full deployment of the stent, withdraw the delivery system carefully, leaving the guide wire in place. If the tip of the delivery device catches the distal stent filaments, perform gentle movements in order to free it.

Perform routine post-procedural angiography / cholangiography.

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