STERILE EO

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

**Pressure Products, Inc.**  
Customer Service  
1861 N. Gaffey Street, Suite B  
San Pedro, CA 90731 - USA  
Tel +1-310-547-4973  
Fax +1-310-547-4760  
www.pressure-products.com

**Oscor Inc.**  
3816 DeSoto Boulevard  
Palm Harbor, FL 34683 - USA  
Tel +1-727-937-2511  
Fax +1-727-934-9835  
www.oscor.com

**Oscor Europe GmbH**  
Fritz-Vomfelde-Strasse 6  
40547 Düsseldorf, Germany

## Technical Manual

## en-Instructions for Use

This device is intended for one time use only. Read instructions prior to use.

### Indications

For the introduction of various types of pacing leads and catheters.

### Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton. Store in a cool, dark, and dry place.
- Infusion through the side port can be done only after all air is removed from the unit.

### Precautions

- Do not alter this device in any way.
- When assembling the sheath and dilator, care must be taken to insert the dilator tip straight through the center of the valve membrane in order to prevent inadvertent puncturing of the membrane.
- Aspiration and saline flushing of the sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, electrode, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve members resulting in blood flow through the valve. Never advance or withdraw guide wire or sheath when resistance is met. Determine cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the side port only.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

### Use Sterile Technique *A suggested procedure:*

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture as desired.
3. The angle of the needle should be adjusted depending on the Patient's build: shallow in a thin person, deeper in a heavyset person.
4. Insert needle into vessel. The needle position should be verified by blood return.
5. Aspirate the puncture needle using the syringe.
6. Remove the syringe and insert soft tip of guide wire through the introducer needle into the vessel. Advance the guide wire to the required depth. Leave an appropriate amount of guide wire exposed.

### At no time should the guide wire be advanced or withdrawn when resistance is met.



Determine the cause of resistance before proceeding. Fluoroscopic verification of the guide wire's entrance into the vessel is suggested.

7. Hold guide wire in place and remove introducer needle. Do not withdraw the guide wire back into the cannula as this may result in separation of the guide wire. The cannula should be removed first.
8. Insert vessel dilator into sheath and rotate the dilator cap over valve housing to secure the dilator onto sheath assembly.
9. Thread the dilator/sheath assembly over the guide wire.
10. Advance the dilator and sheath together with a twisting motion over the guide wire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guide wire will prevent inadvertently advancing the guide wire entirely into the patient.
11. Once assembly is fully introduced into the vascular system, separate the dilator cap from the sheath

valve housing by rotating the dilator cap off the hub. (see figure 1)

12. Slowly retract the guide wire and dilator, leaving the sheath in position. The hemostasis valve will reduce the loss of blood and the inadvertent aspiration of air through the sheath.
13. Aspirate all air from the sheath valve assembly by using a syringe connected to the side port. Flush the introducer through the side port. If the introducer is to remain in place during catheter positioning and testing, flushing the introducer via the side port periodically with saline is advised.
14. Introduce the catheter through the hemostasis valve/sheath and advance it into position.
15. Flush the sheath with 5cc of saline immediately before peeling the sheath away in order to minimize backbleeding.
16. First, withdraw the sheath from the vessel. Then remove the sheath by sharply snapping the tabs of the valve housing in a plane perpendicular to the long axis of the sheath to split the valve and peel the sheath apart. (see figure 2)

en	Directions for Use	1
cs	Návod k použití	2
da	Brugsanvisning	2
de	Gebrauchsanleitung	2
el	Οδηγίες χρήσης	3
es	Instrucciones para el uso	3
et	Kasutusjuhend	4
fi	Käyttöohjeet	4
fr	Mode d'emploi	4
hu	Használati utasítás	5
it	Istruzioni per l'uso	5
ja	使用説明書	6
ko	사용 방법	6
lt	Naudojimo instrukcijos	7
nl	Gebruiksaanwijzing	7
no	Bruksanvisning	7
pl	Instrukcja użytkowania	8
pt	Instruções de Utilização	8
ru	Инструкции по использованию	9
sk	Pokyny na používanie	9
sv	Bruksanvisning	9
tr	Kullanım Talimatları	10
zh	使用说明	10

The H. Lee SafeSheath™, hemostatic, splitable sheath is protected by Patent Insurance through 2011, USA Patents 5,125,904 and 5,312,355. Canadian Patent 2,112,394. European Patent EP0939685. SafeSheath is a registered trademark of Pressure Products, Inc.

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