Instructions for use

Pusher, pusher with radiopaque ring marker

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GANL Nr 31 Pusher EN/FR_00/24.08.2012

Please note: This medical product may be obtained and used only by trained medical personnel.

1. Description

Pusher and pusher with radiopaque ring marker

2. Contents and packaging

- Polyamide pusher or
- Platinum pusher (PTFE with platinum ring marker):
 - CH 5 = grey colour
 - CH 6 = black colour
 - CH 7 = yellow colour
 - CH 8 = orange colour

Double packaging (length 70 cm also in dispenser with clamp)
Sterilised with ethylene oxide.

3. Intended use

Pushers are used to position ureteral stents precisely in the ureter via a guide wire. The pusher allows the ureteral stent to be advanced forwards. It can also be pulled backwards and rotated if the pusher and ureteral stent are connected to each other (controllable ureteral stents). Instead of ureteral stents, other objects can also be positioned at other places in the urinary tract.

The medical indication corresponds to that of the object being positioned.

4. Contraindications

Contraindications arise from restrictions on the corresponding basic applications.

5. Side effects

There are no known side effects from the use of pushers.

6. Instructions for use

Pushers sold singly are threaded over a guide wire before the objects being positioned (ureteral stents) so that they are behind the object when viewed from the perspective of the wire tip. On pushers with radio markers, the ring marker (silver colour) must be at the front.

During the procedure, the guide wire is advanced to the destination point to which the object is intended to be moved (in the case of ureteral stents, this is the renal pelvis). Under fluoroscopic monitoring, the object is then first threaded onto the guide wire and then, with the aid of the pusher, advanced to its destination point. When the object is positioned correctly, the guide wire is withdrawn while the pusher holds the object in place. The pusher is then also withdrawn.

7. Possible complications and/or risks

- When a hydrophilic wire is used, there is a risk that the pusher will remain stuck to the wire if the wire is no longer adequately hydrated.
- Do not guide the product over or through sharp-edged tools since they can cause it to tear.

Do not treat the products with disinfectant agents, since UROTECH has no information regarding possible material changes that this may cause.

8. Warnings / Precautions

- This product is only sterile if the packaging is undamaged or unopened. For single use only!
- All components must be carefully checked before use. Products that may be damaged must not be used.
- This product should be used only by physicians experienced in the use of such medical devices.
- Check the compatibility of disinfectant agents with the catheter material before using them.
- Ensure that the pusher is compatible with the guide wire.

9. Interaction with other drugs

It is possible that the stent may discolour due to the concomitant use of medications or even as a result of contact with urine. UROTECH has not carried out any tests in this regard, however our many years of experience have documented no indications that this impairs the product's function in any way.

10. Disposal

After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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Medizinische Technologie GmbH