NeoClean Tube

FR

Manuel d'instructions

EN

Instruction Manual

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
The Neoclean Tube is a cost effective, completely disposable alternative to the reusable irrigation tube set in GI Endoscopy. Case after case, our design allows for 24 hour use without the need for reprocessing. Each irrigation tube is equipped with a back-flow valve as a safeguard against potential infection control risks and work in conjunction with connectors for seamless integration with a variety of irrigation pumps and endoscopes.

【Intended Use】

The Neoclean Tube is intended to provide irrigation via sterile water during GI endoscopic procedures when used in conjunction with an irrigation pump (or cautery unit).

【Types and specification】

Product sketch map as following:

1. Tube 2. Cap 3. Adapter 4. Connector

Item Number	Description	Suitable	With CO ₂
41205	Irrigation Tube	Erbe or Medivators Pump	No

【Warnings and Precautions】

- Guidelines for Preventing Cross Contamination:

Asept InMed
Essentiellement simple... et tout simplement Essentiel
 Distributed by: Asept InMed SAS
 9, Avenue Mercure 31130 Quint Fonsegrives.
 Tel 05.62.57.69.00 - Fax 05.62.57.69.01
www.aseptinmed.fr

FARTLEY
 Yangzhou Fartley Medical
 Instrument Technology Co., Ltd
 Beizhou Road, Lidian Town,
 Guangling District,
 Yangzhou, P.R.C., 225106
 China

EC REP Prolinx GmbH
 Brehmstr. 56, 40239
 Düsseldorf Germany
 Tel: 0049 2131 4051968-0
 DE/0000045300

- ◆ ALWAYS prime the auxiliary water channel prior to inserting the GI endoscope into patient.
- ◆ NEVER attach the Neoclean Tube directly to the biopsy channel without also using a secondary backflow valve.
- ◆ ALWAYS close the pinch clip before disconnecting the NeoClean Cap from the endoscope.
- ◆ DO NOT REMOVE THE BACKFLOW VALVE. Backflow valves are permanently attached to Neoclean Tube.
- ◆ The backflow valve attached to the Neoclean Tube must be in place. If the backflow valve is missing, then discard the entire irrigation tube, accessories and water bottle.
- ◆ DO NOT RE-STERILIZE OR REPROCESS

- If any of the above prevention guidelines were not followed, or there is any reason to believe that the Neoclean Tube system was contaminated, immediately discard the entire irrigation tube, accessories and water bottle. Replace with a new Neoclean Tube, accessories, and unopened sterile water bottle.
- Sterility is not guaranteed if package has been opened or damaged. Do not use if packaging or product is damaged.

- Neoclean Tube This disposable medical device must be used within 24 hours of opening the package. It must be disposed of at the latest at the end of this period, irrespective of the number of times it has been used.

【Instructions for Use】

Open the sterile peel pack and remove the Neoclean Tube. Open a bottle of sterile water. Drop the weighted tube into the sterile water bottle and firmly tighten the bottle cap to ensure a secure seal. Remove the label from the pouch and apply it to the bottle of sterile water. Write the date and time on the label.

- 1) Open pump head, insert Neoclean Tube and close pump head. (Note: Ensure that tubing is placed in between tube markers.)
- 2) Connect the backflow valve (distal tip of the Neoclean Tube) to either the auxiliary water connector or biopsy irrigation accessory.
- 3) Prime the tubing prior to use by activating the foot pedal; flush water through the entire GI endoscope.
- 4) **WARNING:** To prevent cross-contamination, always prime the auxiliary water channel prior to inserting the GI endoscope into patient.
- 5) If the water bottle needs to be refilled, use proper aseptic technique.
- 6) When the procedure is finished, turn off the pump unit before disconnecting the Neoclean Tube and sterile bottle.
- 7) Discard the Neoclean Tube and sterile bottle daily..

【Caution】

- ✓ Sterility is guaranteed unless the sterile pouch is damaged or seal is broken.
- ✓ Disposable device for single use.
- ✓ Sterilized by ETO Sterilization.
- ✓ Do NOT re-use or re-sterile.
- ✓ Check the production and expiry date printed on the sterile pouch.
- ✓ Do not store in direct sunlight, at extreme temperature, or in high humidity.
- ✓ Re-use may cause cross-infection.

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Symbol	Explanation	Symbol	Explanation
	Manufacturer		Keep away from sunlight
	Authorized representative in the European		Keep dry
	Date of manufacture		Temperature limit
	Use-by date		Do not re-use
	Batch code		Consult instructions for use
	Sterilized using ethylene oxide		Caution
	Do not use if package is damaged		
	Notified Body, and product meet the basic requirements of MDD93/42/EEC		

File No.: FL/CE-07-02-IFU-D

Rev.: A/0

Issued date: 2017-03-17