IFU-CAAB-INTER-1.1

Cutter Stapler and Reloading Unit for Cutter Stapler

Date of issue: May 04, 2017

Cutter Stapler and Reloading Unit for Cutter Stapler

EN

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY

IMPORTANT

This booklet is NOT a reference to surgical techniques but a supporting material for the use of this product, which was designed, tested and produced merely for single patient use. The product is not reusable so it is shall NOT be reprocessed or resterilized due to the risks of failure, patient injury, contamination as well as patient infection.

DESCRIPTION

The Cutter Stapler CAAB is available in 30, 40, 45, 50, 55, 60, 90 mm staple line lengths for use in different applications. Stapler is featured with perfect osculating capacity. Product and components shall be sterile; the residual amount of ethylene oxide after the sterilization shall lower than 10ug/g.

INDICATIONS

It is applicable to the reconstruction of the digestive tract with the difficulties of exposing the operative eye, the reconstruction of anastomotic stoma after the visceral resections and the closure of an incision or stump.

CONTRAINDICATIONS


1. This medical instrument shall not be applied to the position where the surgical stapler nails are disabled;
2. It shall not be applied to the tissues where the instrument can be closed with great difficulties;
3. It shall not be applied to the ischemic or necrotic tissues;
4. It shall not be applied to the aorta, coronary artery, carotid artery, pulmonary artery, pulmonary vein or any other great vessels, nor superior vena cava, inferior vena cava, common iliac artery, internal iliac artery and vein, external iliac artery and vein and brachiocephalic artery trunk;
5. It shall not be applied to the vaginal vault.

WARNINGS AND PRECAUTIONS

1. The application of this product shall subject to the general principles of stapler;
2. Cutter stapler shall be handled with care, free from any collision or knock. During the disassembly, assembly of cutter staplers, no metal products shall be used for wiping, scratching, knock and other hazard operations.

3. The cutter stapler shall be sealed with aseptic packaging, and once the aseptic packaging is opened, the product shall not be re-used even though after re-sterilization.
4. If case of any missing or unpinning of cutter stapler and its components in the transportation and installation, please recheck the components of stapler cartridge before use.
5. After finishing the anastomosis, it is necessary to check the cutting-off stump or residual tissues; if the incision is complete, it means that the effect of anastomosis is satisfying; if there is any defect or missing, more staplers shall be made on the corresponding incision to ensure the quality of anastomosis.
6. The ex-factory stapler shall be sealed with aseptic packaging, and the residual amount of ethylene oxide after the sterilization shall lower than 10ug/g.
7. This manual is NOT a reference to surgical operation, but a general introduction to the use of this stapler, and the operator shall refer to another related medical literature for surgery.

SCHEMATIC VIEW

- | | | | | | | |
|--------------------------|-------------------------------|--------------------------------|----------------------------------|----------------------------|--------------------------|----------------|
| 1.Cutting stapler device | 2. Cover of stapler cartridge | 3.Stapler cartridge components | 4.Clampi ng plate of compon ents | 5.Push rod of locating pin | 6.Button of locating pin | 7.Reset button |
| 8.Handle | 9.Trigger | 10.Jaw clamping plate | 11.Locati ng pin | 12.Cutter | 13.Stapler cartridge | |

STAPLER SPECIFICATIONS

- T) INSTRUMENT REORDER CODES
- U) RSU REORDER CODES
- V) STAPE LINE LENGTH
- W) COLOR
- X) OPEN STAPLER SIZE
- Y) CLOSED STAPLER SIZE

T	U	V	X	Y	W
CAAB-30	CAZB-30T	30mm	2.5mm	1.0mm	white
	CAZB-30N	30mm	3.5mm	1.5mm	Blue
	CAZB-30D	30mm	4.8mm	2.0mm	Green
CAAB-40	CAZB-40T	40mm	2.5mm	1.0mm	white
	CAZB-40N	40mm	3.5mm	1.5mm	Blue
	CAZB-40D	40mm	4.8mm	2.0mm	Green
CAAB-45	CAZB-45T	45mm	2.5mm	1.0mm	white
	CAZB-45N	45mm	3.5mm	1.5mm	Blue

	CAZB-45D	45mm	4.8mm	2.0mm	Green
CAAB-50	CAZB-50T	50mm	2.5mm	1.0mm	white
	CAZB-50N	50mm	3.5mm	1.5mm	Blue
	CAZB-50D	50mm	4.8mm	2.0mm	Green
CAAB-55	CAZB-55T	55mm	2.5mm	1.0mm	white
	CAZB-55N	55mm	3.5mm	1.5mm	Blue
	CAZB-55D	55mm	4.8mm	2.0mm	Green
CAAB-60	CAZB-60T	60mm	2.5mm	1.0mm	white
	CAZB-60N	60mm	3.5mm	1.5mm	Blue
	CAZB-60D	60mm	4.8mm	2.0mm	Green
CAAB-90	CAZB-90T	90mm	2.5mm	1.0mm	white
	CAZB-90N	90mm	3.5mm	1.5mm	Blue
	CAZB-90D	90mm	4.8mm	2.0mm	Green

INSTRUCTIONS FOR USE

1. Place to-be cut or removed tissue in the jaw of the stapler.

2. Before closing the stapler, it is necessary to set the locating pin in the proper position by pushing the button; and the locating pin can be moved forward by pushing the button of locating pin manually. If the locating pin has not been set in the proper position before closing the stapler, the locating pin will be automatically

in place when the jaw of stapler closes.

3. Hold and squeeze the trigger to be in the semi-closure status, and now the tissue in the jaw can be adjusted.

4. After adjusting the tissue, squeeze the trigger tightly until the jaw of stapler closes. If relax the grip, the trigger will fully rebound to the original position.

5. Fully squeeze the trigger to stitch and cut the tissue, and then continuously squeeze the trigger until the trigger firing does not move, which means that the triggering process finishes.

6. Press the blue reset button at the end of the stapler to open the jaw of the stapler, and then pull the button of locating pin back to the original position.

Note: After moving away the stapler, check the bleeding conditions of suture lines, and if there is mild bleeding, manual suture or electric coagulation hemostasis can be applied.

[Component replacement method]

1. Before removing the components of the used stapler cartridge, make sure that the jaw of the stapler is in the fully open status, and the locating pin has been fully returned back.

2. Pinch the thumb pad beyond the component of the stapler cartridge, then draw it upward to draw out the jaw. Before replacing the components of stapler cartridge, any tissue and/or shaping but still unused suturing nail shall be removed away from the surface of the cutting stapler device, which can ensure the correct shaping of suturing nails when triggered.

3. Replacement of components: Pinch the thumb pad of the components of the unused stapler cartridge, and make the pin hole align with the cutting stapler device. Insert the components of the stapler cartridge into the clamping plate of metal components, and then push down until the components of the stapler cartridge are in the proper position.

Note: Before installing the components of the stapler cartridge, check whether the button of locating pin is in the initial position, and then reset the button of locating pin to the initial position, if not, for normal installation of stapler cartridge components; if the components are not installed in the proper position, or fail to be fully inserted, or fully/partially triggered, the jaw of the stapler may not close normally. Please select the proper components of stapler cartridge based on the thickness of the tissue.

4. Repeat the above-mentioned steps based on the needs.

Note: Each stapler can be replaced with 7 components of the stapler cartridge for eight times of triggering.

5. Suture and cutting operations can be done again according to the above methods from Step 1-4.

STORE AT ROOM TEMPERATURE.

AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

DO NOT EXPOSE TO TEMPERATURES ABOVE 1300F (54 °C) .

From the EO sterilization date, the product is valid for 3 years.

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B.J.ZH.F.PANTHER MEDICAL EQUIPMENT CO., LTD