Tubulure CO2

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
Revision: 0

【Intended Use】

The CO2 Tube with Luer Input is intended to be used with a CO2 insufflator along with clean connecting tube or hybrid clean connecting tube to supply CO2 to a GI endoscope during GI endoscopic procedures.

【Types and specification】

Product sketch map as following:

Example CO2 Tube with Filter for Fujinon®

1. Active Male Luer
2. Tube
3. Air Filter
4. Female connector
5. Fujinon CO₂ Adapter

【Warnings and Precautions】

- Guidelines for Preventing Cross Contamination:
 - **ALWAYS** prime the air/ water channel prior to inserting the endoscope into the patient.
 - **NEVER** turn off the air processor/ CO2 while the clean connecting tube is attached to an endoscope.
 - If any of the above prevention guidelines were not followed, or there is any reason to believe that the clean connecting tube system was contaminated, immediately discard the entire tubing and water bottle. Replace with a new tubing and unopened sterile water bottle.

- **ALWAYS** turn off the CO2 Insufflator between procedures.
- **ALWAYS** inspect the endoscope's o-rings and gaskets (air/ water button, endoscope connection, etc.) for damage. Damaged components can result in less than optimal performance.
- The CO2 tube is to be used only for GI endoscopic procedures. Do not use the tube set for CT Colonography (Virtual Colonoscopy).
- **DO NOT RE-STERILIZE OR REPROCESS.**

- Sterility is not guaranteed if package has been opened or damaged. Do not use if packaging or product is damaged.

- CO2 tube is for 24 hours use (discard daily).

PRECAUTION: If sterile water bottle needs replacing during a procedure, refrain from depressing the air/water valve. There is a potential for backflow to occur if the air/ water valve is depressed while no sterile water bottle is connected to the tubing. If

this occurs, replace the tube immediately.

WARNING:

After use, this product may be a potential biohazard. handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

【Instructions for Use】

- 1) Open the pouch and remove the CO2 source tubing.
- 2) Connect the Active male luer of CO2 tube to the luer connection of clean connecting tube or CO2 tube.
- 3) Connect the adapter of the CO2 tube to the CO2 insufflator.
- 4) Turn on the CO2 insufflator. Verify that the air supply from the GI endoscopic system is turned off.
- 5) Prime the channel and test for gas and water flow prior to insertion of GI endoscope. If the pressure is low, ensure that the water bottle is tightly closed.

WARNING:

To prevent cross-contamination, always prime the auxiliary water channel prior to inserting the GI endoscope into patient.


- ✓ When procedure is completed, close the pinch clip on the clean connecting tube or CO2 tube and disconnect from the GI endoscope.
- ✓ Turn off CO2 insufflator.
- ✓ Discard the CO2 tube daily.

【Caution】

- ✓ Sterility is guaranteed unless the sterile pouch is damaged or seal is broken.
- ✓ Disposable device for single use.
- ✓ Sterilized by ETO Sterilization.
- ✓ Do NOT re-use or re-sterile.
- ✓ Check the production and expiry date printed on the sterile pouch.

- ✓ Do not store in direct sunlight, at extreme temperature, or in high humidity.
- ✓ Re-use may cause cross-infection.

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Symbol	Explanation	Symbol	Explanation
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community		Keep dry
	Date of manufacture		Temperature limit
	Use-by date		Do not re-use
	Batch code		Consult instructions for use
	Sterilized using ethylene oxide		Caution
	Do not use if package is damaged		Notified Body, and product meet the basic requirements of MDD93/42/EEC