ntplast ncranial

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Parte testuale versione 01IFU-02_GH_rev.00

User instructions

Text version 01IFU-02_GH_rev.00

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Texteil Version 01IFU-02_GH_rev.00

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Parte textual versión 01IFU-02_GH_rev.00

Mode d'emploi

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ENGLISH
<p>Intended use ntcranial is a device made of plastic material designed to fix the cranial flap in place in a stable manner, while closing the holes made by the perforator; it consists of a tensor assembled on a base available in two versions (Gap and Hole), and locked by means of a plate descending from above.</p>

Introduction

- The **ntcranial** device should only be implanted by neurosurgeons who are fully aware of the intended use, surgical technique, risks arising from improper use, and the instruments dedicated to the device.
- The **ntcranial** device must be implanted in an operating room under aseptic conditions.
- When implanting **ntcranial** devices, the utmost attention must be paid to the movements.
- The use of **ntcranial** combined with other cranial fixing device is not recommended. If this is the case, it is the surgeon’s responsibility to verify the stability of the fixed flap.

Product information
ntcranial is a long-term implantable device.
The **ntcranial** device consists of:

- a tensor (1a),
- a GAP base (1b) or a HOLE base (1c)
- a locking plate (1d)

The **ntcranial** device is available in two versions:

- 1) **GAP version** (Fig. 2) to be placed along the craniotomy cut.
- 2) **HOLE version** (Fig. 3) to be placed on the holes made by the perforator.

The surgeon chooses the version and number of devices after evaluating the shape and size of the flap.
ATTENTION: **ntcranial** device is indicated for bone flaps ranging between 2.5 mm and 20 mm in thickness. The GAP version is suitable for positioning in craniotomy cuts with a width of between 1.7 mm and that resulting from the use of standard craniotome cutters.
The HOLE version is suitable for use in holes with diameters between 9 mm and 14 mm.

If ncranial devices are used in conditions other than those indicated above, fixing of the flap may prove to be inadequate, as it may be unstable and misaligned with respect to the cranial vault.
The GAP version is suitable for positioning in craniotomy cuts with a width of between 1.7 mm and that resulting from the use of standard craniotome cutters.
The HOLE version is suitable for use in holes with diameters between 9 mm and 14 mm.
If ncranial devices are used in conditions other than those indicated above, fixing of the flap may prove to be inadequate, as it may be unstable and misaligned with respect to the cranial vault.
The **ntcranial** device is made entirely of biocompatible PEEK OPTIMA® plastic, which does not interfere with CT or MRI exams.
WARNING: Do not use the **ntcranial** device on patients who have shown signs of allergic reaction to PEEK.

Sterilisation
The **ntcranial** device (Fig.4) is packaged in double blister packs and supplied sterile.
WARNING: The ncranial device is for single-use only. Do not sterilise again and/or reuse. Reuse of the device may entail the risk of malfunction, breakage, or infection.

Package Contents
The **ntcranial** device is available in two versions: **GAP** and **HOLE**.
The package consists of a cardboard box containing 1 device and the user instructions.
In the GAP version, the box contains the following components already assembled and packaged in double sterile blister packs:

- Tensor, Gap base, locking plate
- 3 traceability labels.

In the Hole version, the box contains the following components already assembled and packaged in double sterile blister packs:

- Tensor, Hole base, locking plate
- 3 traceability labels.

Contraindications
Contraindications include, but are not limited to, the following:

- hypersensitivity or allergies to PEEK.
- latent or active intracranial infections
- fever and leukocytosis
- bone conditions that are not suitable for the application of **ntcranial** device
- cranioplasty performed manually
- bone infections
- bone thickness of less than 2.5 mm or more than 20 mm
- subtoriental craniotomy
- post-traumatic multi-fragmented bone flaps
- patients under 10 years.

Complications
Improper use may cause the breakage of the **ntcranial** device or instability during fixing, with possible dural or brain damage.
The general perioperative risks include:

- infections
- brain injury, neurological deficit
- haematomas at the implantation site
- alterations in wound healing.

Surgical technique
The **ntcranial** device must be implanted following only the instructions in these instructions.
Choice of the number of ncranial devices and configurations, choice of the positioning site and correct application of the devices fall under the sole responsibility of the surgeon.
The selection of eligible patients is the responsibility of the surgeon.

Positioning the devices
1) Prepare, at the surgeon’s request, the number of packaged devices needed to fix the flap.
2) **Make sure the packages are intact:** In the case of damaged packaging or packaging with labels that are no longer legible, do not use the contents. The device must be implanted before the expiration date indicated on the packaging. Check also that the **ntcranial Cutter** is intact
3) Remove the device from the packaging without compromising its sterility. **If any component appears to be defective or damaged, do not use the device.**
4) The **ntcranial** device is already assembled in two configurations: **GAP** (Fig. 5) and **HOLE** (Fig. 6).
WARNING: Do not alter the assembly of the device during use.
WARNING: Do not try to alter the curve of the elastic arms of the tensor by acting directly on it with your hands or tools. Altering the elastic arm curve may hinder proper assembly with the bases.
WARNING: The corners and edges of the **ntcranial** device may pierce any protective barriers (e.g., surgeon’s gloves) or cause patient injury if used improperly. Caution is recommended at each step of the application of the **ntcranial** device.

5) Place the GAP version device along the craniotomy cut in the extradural subtheal site (Fig. 7) and the HOLE version device in the perforator hole in the extradural subtheal site (Fig. 7).
WARNING: Make sure to correctly position the GAP version device. (Fig. 8)
WARNING: Do not place the HOLE version device along the craniotomy cut. Do not place the device in GAP configuration at the perforator hole. If the device is incorrectly positioned for the given configuration, locking the device may be impossible or unstable.
6) It is recommended that **ntcranial** devices should be positioned as evenly as possible along the perimeter of the craniotomy; **do not neglect any side of the craniotomy** (Fig. 9). The number and configurations must be compatible with the shape and area of the flap. **Do not use less than 3 devices to achieve stable fixing.** For standard craniotomies, 3 or 4 devices are generally sufficient. More devices may be required for large craniotomies.
7) After inserting all the **ntcranial** devices, place the bone flap in the craniotomy site.
8) While holding the tensor in place with two fingers, bring the locking plate into contact with the bone by pushing it down with the fingers of your other hand (Fig. 9). Repeat this step for all the devices (Fig. 10).
WARNING: When lowering the plates, it is important to make sure that the bone flap is kept in the correct position.
9) Take the tensor again with two fingers, pulling it upwards to activate the underlying elastic component. At the same time, push the plate down as far as possible using two fingers of your other hand. Compared to the previous step, the locking plate will descend further, consolidating the fixing. Repeat this step for all the **ntcranial** devices.
ATTENTION: Perform the steps above by hand only. It is absolutely forbidden to use tools (such as pliers) that could damage the ncranial device.
10) Once all the **ntcranial** devices are positioned, check the fit and stability of the flap.
11) Remove excess tensor lengths by means of the supplied **ntcranial cutter** (Fig. 11). The cutting of the tensor must be flush with the upper surface of the locking plate.
ATTENTION: During cutting, hold the handle of the tensor with your fingers so that it does not come off accidentally. (Fig. 12).
ATTENTION: To remove the tensor parts, it is recommended to use only the **ntcranial Cutter**.
WARNING: Check that the resulting cut profile is smooth. If necessary, trim. Be careful when doing this; small burrs may pierce the surgeon’s gloves.

Removing the devices
Using the **ntcranial Cutter**, cut the locking plate (Fig. 13). Repeat this step for each **ntcranial** device.

Post-operative warnings

- During the twelve weeks following surgery, the patient should avoid blows, vibrations, pressure and stress on the surgical site, and should not engage in any activity that could expose him/her to such risks.
- Once the bone has completely healed (6 months after surgery), the **ntcranial** device ceases its function, though remaining implanted in the cranial vault.
- If any of the devices should break, it is recommended that the damaged devices be replaced if bone consolidation is not complete yet. Damaged components may migrate from their original location and/or injure the dura mater and brain. In any case, the final decision on whether to replace damaged devices lies with the surgeon.

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