ENGLISH

INTENDED USE

NT CRANIAL CUTTER is a reusable surgical instrument consisting of a surgical steel cutter intended for cutting and removal of the NtCranial medical device and comes with an accessory protective plastic container.

INTRODUCTION

- The NT CRANIAL CUTTER device may exclusively be used by neurosurgeons who perfectly know the intended use, the surgical technique, and the risks arising from improper use of the NtCranial device and the NT CRANIAL CUTTER device.
- The NT CRANIAL CUTTER device must be used in an operating theatre in aseptic conditions.
- Be extremely careful with all movements during use of the NT CRANIAL CUTTER device.

PRODUCT INFORMATION

The NT CRANIAL CUTTER device (*Fig. 1*) is a reusable surgical instrument. The NT CRANIAL CUTTER device consists of a surgical steel cutter whose cutting edges are made in such a way that the NtCranial device can be cut flush with the locking plate. The NT CRANIAL CUTTER device is supplied in an accessory protective plastic container (*Fig. 2*) complete with silicone mat in which the device is to be stored when not used. Store the device in a dry, clean, protected place and away from direct sunlight. Read, keep and follow these instructions for use.

CLEANING, DISINFECTION AND STERILIZATION

The NT CRANIAL CUTTER device is supplied non-sterile and must be subjected to an appropriate cleaning, disinfection and sterilization cycle before each use, especially the first time. The cleaning, disinfection and sterilization cycle must be carried out by trained and qualified personnel. Cleaning and disinfection are essential prerequisites for the effectiveness of sterilization.

WARNING: During cleaning, disinfection and sterilization, comply with the national regulations in force as well as those of the hospital facility. If the equipment and substances described in this guide are not available, the user is responsible for validating his own procedures before using them. Use only validated methods, products and equipment. In particular, the equipment must be subjected to scheduled maintenance and controls. Always comply with the national legislative provisions in force.

For patients with Creutzfeldt-Jakob (CJ) disease, suspected CJ or possible variants of CJ, comply with the national regulations in force as regards sterile preparation of the products. Take all the essential precautions to prevent transmission to other patients, users or third parties. Make sure that you use cleaning agents and disinfectants that do not contain the following constituents: acids (< pH 5 / oxidising acids), alkaline substances (< pH 10), organic solvents, disinfectants containing alcohol, phenols, ammonium, benzene, halogenated hydrocarbons, sodium chloride in high concentrations, oxidants.

Pre-cleaning

Keep the cutter outside the container in the phases described below.

Clean and disinfect the cutter and the container after each use. Remove coarse contamination immediately after use so that it does not dry on the surfaces complicating subsequent cleaning and disinfection. Do not use solutions with a fixing effect or hot water as this would cause the residues to stick to the surfaces with a negative impact on the cleaning result.

Do not immerse the device in saline solution (NaCl) as this would result in corrosion and premature damage of the instrument. Handle the device with care as dropping or carelessly tossing it may damage it. Immerse the cutter and the container in cold water for at least 5 minutes. If possible, clean it with a soft brush in cold water until there are no more visible signs of contamination. Do not use brushes with metal bristles or steel wool. Rinse the cutter joints with a compressed water gun for at least 10 seconds (pulsed method) or an equivalent method. Then immerse the cutter and the container in an ultrasound bath (demineralised water at 40°C, 0.5% alkaline detergent) for at least 15 minutes at a frequency of at least 35 kHz. Remove and rinse with cold water.

Cleaning

Again keeping the cutter outside the container, place both of them in an automatic instrument washer. Start the cycle:

- 1 minute pre-cleaning in cold water
- Draining
- 3 minutes pre-cleaning in cold water
- Draining
- 5 minutes cleaning at 55°C with 0.5% alkaline detergent
- Draining
- 3 minutes neutralization with hot tap water (>40°C) and neutralizer
- Draining
- 2 minutes rinsing with hot tap water (>40°C)
- Draining

At the end of the cycle, the cutter and the container must be free of any kind of residue visible to the naked eye. Carefully visually inspect the cutter and the container before going on to the next phases. Change the cleaning solutions at least once a day.

Disinfection

Again keeping the cutter outside the container, disinfect the cutter and the container using a heat disinfection device and comply with the current national requirements for the A0 value (ISO 15883). Do not overload the tray and make sure that there are no unexposed areas on the cutter and the container.

Dry using the automatic cycle of the device used for disinfection. You can further dry the items by hand using cloths that do not leave residues or sterile compressed air.

Check that there are no marks or residues on the cutter and the container. If you find any trace of residue, repeat the process. Finally, check the integrity of the cutter and the container and the cutter's movements. The device does not require lubrication.

Sterilization

Place the cutter in the accessory container complete with silicon mat. Then wrap the container in a material suitable to maintain the sterility in compliance with ISO 11607 and suitable for sterilization in a steam autoclave.

WARNING: The accessory container is exclusively intended to provide mechanical protection for the cutter during storage and sterilization, but on its own does not maintain the instrument sterile. Before sterilization, it is therefore essential that the container holding the cutter be wrapped in packaging specifically intended to maintain sterility. Remove the device from this packaging only at the time of use, using a procedure that assures that the product will not be contaminated. NT CRANIAL CUTTER must be sterilized by steam using the fractionated vacuum method (in compliance with ISO 13060 / ISO 17665 and the current national requirements).

Make sure that all the external and internal surfaces will be exposed to the sterilizing agent. Do not fill the autoclave compartment over the maximum level. Follow the instructions provided by the

manufacturer. Use a validated steam cycle at a temperature between 132°C and 135°C for at least 5 minutes.

CONTRAINDICATIONS

There are no known contraindications for use of the NT CRANIAL CUTTER device. The surgeon is in any case responsible for evaluating the overall conditions of the patient and deciding whether the instrument can be used.

USING THE DEVICE

Take the NT CRANIAL CUTTER instrument out of the accessory container and check that all its parts are intact. In particular, check that there are no signs of corrosion, blade wear, poor movement, scratches, loosening, bending and/or breakage of one or more parts, which may hamper proper functioning. Do not use the instrument if it appears damaged. Do not attempt to repair the instrument yourself. Return any damaged instruments after cleaning, disinfecting and sterilizing them.

Cutting the excess part of the tensor

Once you have checked fixation stability, cut off the excess parts of the tensor.

- Rest the NT CRANIAL CUTTER blade on the surface of the locking plate;
- Holding the tensor stem with your fingers, make the cut (*Fig. 4*);
- If necessary, trim in the same way;
- Repeat the operation for each tensor.

WARNING: While cutting, hold the tensor stem with your fingers so that it cannot accidentally spring away with the risk of hitting the surgeon or the assistants (*Figs. 3,4*).

WARNING: Do not try and insert the tips of the cutting edges in the locking plate cavities in which the teeth are. Insisting with this action may cause damage to the NtCranial device.

WARNING: It is recommended to exclusively use the NT CRANIAL CUTTER to remove the excess parts of the tensor. Using different cutters may result in an imperfect cut or damage to the NtCranial device.

Removing the NtCranial device (*fig. 5*)

Insert the cutter blade underneath the locking plate and cut it up to the middle. Repeat the operation on the other side. Be careful not to drop any scraps into the operating field. Repeat the operation for each device.

USEFUL LIFE OF THE DEVICE

The device must be replaced in the following cases:

- When there are signs of corrosion, blade wear, scratches, bending and/or breakage of one or more parts
 - When the central pin has loosened or movement is poor
 - When excessive force is required to make the cut
 - When it is difficult to make a precise cut and flush with the locking plate
- The useful life of the device depends on how frequently it is used and the attention paid during its preparation.

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