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
**Package content**  
1 x pRELAX Vasospasm Treatment Device

## Product description

The pRELAX consists of a self-expanding radiopaque Nitinol structure ①, carries one radiopaque marker on its proximal ② and two radiopaque markers on its distal end ③ and is firmly attached to an insertion wire ④.  
The device is stored in compressed form in an introducer sheath (not displayed).  
The pRELAX is introduced into the target vessel through a suitable microcatheter and is deployed in the vessel constricted by vasospasm by withdrawing the microcatheter. The target vessel is dilated by deploying the pRELAX. The pRELAX is removed by pushing the microcatheter in the distal direction over the device left at the respective location; the device is then pulled out of the microcatheter and thereby removed.

## Intended Purpose

The intended use of the pRELAX Vasospasm Treatment Device is to dilate intracranial vessels for acute treatment of vasospasm, e.g. of the internal carotid artery, the middle cerebral artery, the basilar artery and of vessels situated distally of them with a temporary application time < 60 minutes.

## Indications

The pRELAX Vasospasm Treatment Device is designed for mechanical dilation of intracranial arteries for endovascular treatment in the acute phase of vasospasm

- in patients for whom treatment with medicinal vasodilators (e.g. nimodipine, nicardipine, verapamil) or a balloon angioplasty is contraindicated or
- in patients for whom treatment with a triple-H therapy (hypertension, hypervolemia, hemodilution) is contraindicated or
- in patients for whom a treatment with medicinal vasodilators (e.g. nimodipine, nicardipine, verapamil) or a balloon angioplasty, as well as triple-H therapy, has already failed as well as
- as a complement to treatment that has already been initiated with medicinal vasodilators.

## Contraindications

The use of the pRELAX is contraindicated in patients

- for whom the prepared angiography shows that the respective anatomical conditions are not suitable for an endovascular treatment,
- with an established infarction of the entire brain tissue dependent on the vessel to be dilated,
- with a dissection of the target vessel or a potential dissection of the target vessel,
- with atherosclerotic stenosis at the affected segment of the target vessel,
- with an unsupplied aneurysm at the affected segment of the target vessel and
- in whom there may be a pathological increase in the fragility of the constricted vessels (e.g. Ehlers-Danlos syndrome).

## Compatibility

The following pRELAX model and the following microcatheters are compatible with one another and are recommended for the indicated minimum vessel diameters:

pRELAX Device	Insertion wire length	Microcatheter	Min. vessel Ø (relaxed)	Relaxed Ø x usable length
PRELAX-4-20	1.8 m	ID 0.021": - REBAR 18 - TREVO PRO 18 - PROWLER SELECT PLUS - VIA 21 - PHENOM 21	2.5 mm	4 x 20 mm

## Recommended procedure

### Preparing the patient

- Record and document any available information about the time and nature of the onset of the clinical symptoms of vasospasm.
- Record and document any available information on medical history, including current medication.
- As far as possible, inform the patient and document the patient's consent to the intervention, pointing out the possible complications and potential consequences (see Complications). In cases where patients cannot give consent themselves, their relatives should, as far as possible, be asked what they think the patient would want. Otherwise, the rules of emergency care for incapacitated patients apply subject to different institutional or national requirements.

- An unenhanced CT scan should be performed in advance to rule out intracranial hemorrhage or a large infarcted area of the brain. If these changes are present, treatment with the pRELAX is contraindicated.
- The diagnostic angiography and, where applicable, the endovascular treatment should be performed under general anesthesia with neuromuscular relaxation and invasive hemodynamic monitoring. These measures must be configured such that, as far as possible, the start of the endovascular treatment is not delayed. During anesthesia, aim to maintain elevated blood pressure values, which are adapted to the situation, i.e. in the event of constricted vessels. Hypotension must be strictly avoided. Appropriate medicinal measures (anticoagulants, antiplatelets) should be used to reduce the risk of a thrombus forming on the surface of the pRELAX insofar as is permitted by the respective situation.

### Preparing the treatment

- After shaving and disinfecting both groins, a 6F or 8F catheter sheath is inserted, preferably into the right femoral artery.
- For the vessel diagnostics, perform an angiography of the complete anterior and posterior circulation. If possible, avoid magnified views and oblique views. Restrict the contrast agent exposure to the unavoidable minimum.
- The target vessel(s) for the endovascular treatment must be defined. It is important that the device is introduced in sufficiently large target vessels only. Observe the minimum vessel diameters and the catheter compatibility in this regard. If the diameter of the target vessel seems too small for catheterizing with a suitable microcatheter, investigations should be performed to determine the extent to which vascular dilation can be achieved with medicinal measures.

### Inserting the microcatheter

- Never catheterize against any resistance! Avoid catheterizing side branches with a small diameter!**

Insert a suitable microcatheter (cf. compatibility table) into the target vessel using a microguidewire. Using the so-called "road map" technique is recommended. Remove the microguidewire from the microcatheter.

### Inserting the pRELAX

- In the event of major resistance during insertion of the pRELAX, remove the device, and potentially the microcatheter, and recatheterize the vessel.**

After a hemostatic valve has been connected to the microcatheter under continuous pressurized flushing (physiological saline solution), the pRELAX is inserted into the microcatheter as follows:

- Open the hemostatic valve and push the introducer sheath of the pRELAX through the open valve.
  - Carefully close the valve; the introducer sheath will be deaerated by the retrograde entry of the flushing liquid.
  - After the introducer sheath is deaerated in this form, it is then pushed forward as far as it will go into the adapter (hub) of the microcatheter. The introducer sheath is held in place in this position.
  - The pRELAX is then pushed forward from the introducer sheath into the microcatheter using the insertion wire. This process is continued until around 60 cm of the insertion wire is inside the microcatheter. Then the introducer sheath is removed.
  - The device is pushed forward until the marker ("FLUOROSAFE MARKER": distance to distal pRELAX tip approx. 126 cm in compressed condition) on the insertion wire reaches the entry of the hemostatic valve. This procedure does not need to be performed with continuous fluoroscopy because the marker ("FLUOROSAFE MARKER") identifies the position to which the system can be advanced without the pRELAX tip leaving the microcatheter.
- Never push the pRELAX beyond the distal tip of the microcatheter. This may lead to dissection or perforation of the target vessel!**  
The pRELAX is inserted further up to the tip of the microcatheter under continuous fluoroscopy. The distal tip of the device should reach the tip of the microcatheter.

### Deploying the pRELAX

- The deployment process should take place under fluoroscopy to keep the device in its original position.**

Deploy the device completely by carefully and slowly withdrawing the microcatheter until the microcatheter tip is located in a sufficiently proximal position. This is achieved when the distal microcatheter marker is located proximal to the radiopaque proximal pRELAX marker.

The device is ideally deployed in the vessel segment constricted by vasospasm. If longer vessel segments are affected by vasospasm, several deployments of the pRELAX may be performed in sequence.

If the pRELAX is longer than the vessel segment affected by vasospasm, the device should protrude proximally. The picture to the right depicts the radiographic image of a pRELAX.

- Movement of the device must be avoided under all circumstances!**

Wait up to ten (10) minutes until the pRELAX has dilated the constricted vessel segment. Observe the opening behavior of the pRELAX in the target vessel under fluoroscopy and, where possible, without administering contrast agent. The radiopaque coating of the pRELAX can be used as an additional visual check.

### Removing the pRELAX

- The removal process should take place under fluoroscopy to keep the device in its original position.**

Secure the device by carefully and slowly pushing forward the microcatheter until the microcatheter tip is located in a sufficiently distal position. This is achieved when the distal microcatheter marker is located distal to the radiopaque distal pRELAX marker.

- Now slowly pull the pRELAX back through the microcatheter and remove the device once the hemostatic valve of the microcatheter has opened fully.

16. By injecting around ten (10) mL of X-ray contrast agent via the guide catheter, it is now possible to check whether the target vessel has been dilated by using the pRELAX.
- If the target vessel is sufficiently enlarged, check whether other vessels are constricted and decide whether they can be dilated using this method. If this is not the case, remove the guide catheter. After 10-15 minutes, inject the affected supra-aortic vessel to check whether the target vessel is still open.
  - If the target vessel is *not* dilated, a decision must be made regarding how to proceed. This may include: a repeat pRELAX procedure, use of a different mechanical device, medicinal treatment.

*Reuse in the same patient*

17. In the event of a repeat or additional procedure, an intact device can be reused in the same patient up to two (2) times:
- Cleaning: Carefully rinse the pRELAX with sterile water or sterile physiological saline solution; there must be no blood or thrombus residue on the device.
  - Inspection: Examine the device to rule out any damage.
  - Preparation for the repeat application: Insert the proximal end of the insertion wire into the corresponding introducer sheath. Then move the introducer sheath until the device is in front of the tip of the introducer sheath. Now *SLOWLY* and completely pull the device back into the introducer sheath. If the pRELAX cannot be fully retracted, the product must be discarded.
  - The repeat use of the pRELAX takes place as described in steps 10 to 16.

**Precautions**

- Do not use the pRELAX for purposes other than the intended purpose.
- In order to avoid vessel injuries, do not use the device in vessels with a diameter ( $\varnothing$ ) that is less than the value in the compatibility table.
- Use the device only in microcatheters listed in the compatibility table. Microcatheters (or other catheters) with larger internal diameters may lead to mechanical strain on the pRELAX device during its passage and the device or parts of the device may become separated.
- All procedures must be carried out under fluoroscopic monitoring.
- Extensive flushing of the introducer sheath, as described above, is required in order to remove any embedded air bubbles.
- When the device is transferred from the introducer sheath to the microcatheter, make sure that the introducer sheath does not move back proximally. If this happens, the device and/or insertion wire may be bent and/or weakened and therefore become separated from one another.
- If the device is pushed distal to the tip of the microcatheter, the vessel may be dissected or perforated. If a vessel has been perforated, it should be sealed rapidly and completely.
- Forcibly pulling or pushing the insertion wire and twisting the insertion wire may separate the pRELAX from the insertion wire. In such cases, recovery using an appropriate device is recommended. If the recovery is not successful, dual platelet function inhibition must be immediately initiated and maintained.
- In the event of excessive retraction resistance, the microcatheter should be used to splint or recover the device. If the microcatheter cannot be pushed back over the device completely, remove the device and microcatheter together as a unit.
- The pRELAX may be used up to three (3) times in the same patient. The contact with the patient of all applications must be < 60 minutes in total.
- Between applications, the device must be cleaned of all residue, such as blood and thrombus where applicable, and inspected for damage. If the pRELAX or the introducer sheath are damaged or the device cannot be retracted into the introducer sheath, the device must be discarded.
- The pRELAX is a delicate device and requires careful handling. Never push or pull against resistance. Never twist the insertion wire. If necessary, remove the pRELAX via the microcatheter.
- The device must only be used by specialized and trained physicians.
- The pRELAX must not be deployed in vessel implants (e.g. stents or flow diverters).

**General Information**

- Keep away from heat. Store in a cool, dry place.
- Use before the expiry date only, as sterility is not guaranteed otherwise.
- To ensure product integrity and sterility, do not use damaged packaging.
- Check the product for damage before use. Do not use broken or damaged products, as functionality cannot be guaranteed.
- The product is for single use only. The device must not be resterilized or recycled for other patients because it cannot be reliably cleaned by the user.
- The product, blood and thrombus residue and, if applicable, packaging components must be disposed of appropriately in marked containers.

**Complications**

The following and further complications may occur during or despite the use of the pRELAX:

- Air embolism, embolism in distal vessels, vascular occlusion, thrombosis, cerebral ischemia
- Perforation, rupture, dissection and other arterial damage
- Vasospasm, formation of a pseudoaneurysm, intracranial hemorrhage
- Allergic reaction, infection, reperfusion hemorrhage
- Space-occupying infarction, neurologic deficit including consequences of a stroke
- Persistent vegetative state, death

**Symbols and their meanings**

Symbols on the product label:

	Caution
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Batch code
	Sterilized using ethylene oxide
	Do not use if package is damaged
	Non-pyrogenic
	Use-by date
	Catalogue number
	Content
	Keep dry Keep away from sunlight
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Manufacturer
	Material: NiTi (nickel-titanium alloy, nitinol)
	Latex-free
	Phthalate-free
	Radiopaque coating on nitinol structure
	Colored marker (black) on insertion wire

**Limitation of Liability**

femtos GmbH shall not be liable for damage caused by use other than that intended or by reuse of the product.

femtos and pRELAX are registered trademarks of femtos GmbH in the Federal Republic of Germany and other states.

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femtos et pRELAX sont des marques commerciales déposées de femtos GmbH en Allemagne et d'autres États.
🏠: **femtos GmbH** | Universitätsstr. 136 | 44799 Bochum, Germany
Distributor: **phenox GmbH** | Lise-Meitner-Allee 31 | 44801 Bochum, Germany