# INSTRUCTIONS FOR USE

## WALLABY AVENIR™ COIL SYSTEM

**Wallaby Medical, Inc.**

23181 Verdugo Dr., Suite 104A

Laguna Hills, CA 92653

United States of America

WALLABY AVENIR™ COIL SYSTEM INSTRUCTIONS FOR USE	.....	English	.....	1
WALLABY AVENIR™ COIL-SYSTEM GEBRAUCHSANWEISUNG	.....	German/Deutsch	.....	3
SYSTÈME DE COIL WALLABY AVENIR™ MODE D'EMPLOI	.....	French/ Français	.....	5
SISTEMA DI SPIRALI WALLABY AVENIR™ ISTRUZIONI PER L'USO	.....	Italian/Italiano	.....	8
SISTEMA DE ESPIRAL WALLABY AVENIR™ INSTRUCCIONES DE USO	.....	Spanish/Español	.....	10
SISTEMA DE ESPIRAL WALLABY AVENIR™ INSTRUÇÕES DE UTILIZAÇÃO	.....	Portuguese/Português	.....	13
WALLABY AVENIR™-SPIRAALSYSTEEM GEBRUIKSAANWIJZING	.....	Dutch/Nederlands	.....	15
WALLABY AVENIR™ SPOLSYSTEM BRUKSANVISNING	.....	Swedish/Svenska	.....	18
WALLABY AVENIR™ –KOILIJÄRJESTELMÄ KÄYTTÖOHJE	.....	Finnish/Suomi	.....	20
WALLABY AVENIR™ SPOLESYSTEM BRUGSANVISNING	.....	Danish/Dansk	.....	22
WALLABY AVENIR™-SPOLESYSTEM BRUKSANVISNING	.....	Norwegian/Norsk	.....	24
SYSTEM SPIRAL WALLABY AVENIR™ INSTRUKCJA UŻYCIA	.....	Polish/Polski	.....	26
SISTEM TULJAVE WALLABY AVENIR™ NAVODILA ZA UPORABO	.....	Slovenian/Slovenščina	.....	28
SPIRÁLOVÝ SYSTÉM WALLABY AVENIR™ COIL SYSTEM NÁVOD K POUŽITÍ	.....	Czech/Česky	.....	30

SISTEM SPIRALĂ AVENIR™ WALLABY			
Instrucțiuni de utilizare	.....	Romanian/Română	..... 32
WALLABY AVENIR™ SPIRÁLRENDSZER			
HASZNÁLATI ÚTMUTATÓ	.....	Hungarian/Magyar	..... 34
ΣΥΣΤΗΜΑ ΣΠΕΙΡΑΜΑΤΩΝ WALLABY AVENIR™			
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	.....	Greek/Ελληνικά	..... 36
WALLABY AVENIR™ SPIRLES SISTĒMA			
LIETOŠANAS INSTRUKCIJAS	.....	Latvian/Latviski	..... 39
WALLABY СИСТЕМА AVENIR™ COIL			
ИНСТРУКЦИИ ЗА УПОТРЕБА	.....	Bulgarian/Български	..... 41
WALLABY KOILISÜSTEEM AVENIR™			
KASUTUSJUHEND	.....	Estonian/Eesti	..... 44
SUSTAV ZAVOJNICA WALLABY AVENIR™			
UPUTE ZA UPOTREBU	.....	Croatian/Hrvatski	..... 46
CIEVKOVÝ SYSTÉM WALLABY AVENIR™			
NÁVOD NA POUŽÍVANIE	.....	Slovakian/Slovenčina	..... 48
WALLABY AVENIR™ SPIRALĒS SISTEMA			
NAUDOJIMO INSTRUKCIJOS	.....	Lithuanian/lietuvių	..... 50
WALLABY AVENIR™ KOİL SİSTEMİ			
KULLANIM TALİMATLARI	.....	Turkish/Türkçe	..... 52
СПИРАЛЬНАЯ СИСТЕМА WALLABY AVENIR™			
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ	.....	Russian/Русский	..... 54
WALLABY AVENIR™ نظام الملف			
تعليمات الاستخدام	.....	Arabic/عربي	..... 57
SYMBOL GLOSSARY	.....	.....	..... 59

**INSTRUCTIONS FOR USE**

**EN**

**WALLABY AVENIR™ COIL SYSTEM**

**DEVICE DESCRIPTION**

The Wallaby Avenir™ Coil System (ACS) consists of a platinum implantable coil attached to a delivery pusher. The delivery pusher provides a means to deliver the embolization coil to the vascular target treatment site and is designed to mechanically detach the coil from the delivery pusher.


THE ACS is available in the following configurations:

- Complex Framing Coil
- Helical Filling Coil
- Helical Finishing Coil
- Complex Finishing Coil

The ACS frame coils establish the initial framework for embolization of cerebral aneurysms or lesions. Once a frame has been established by one or more Framing coils, additional ACS Framing coil, Filling coils, and Finishing coils can be deployed to provide embolization of the aneurysm or lesion. ACS Filling coils are used to fill an aneurysm or lesion space within the established framing coil. ACS Finishing coils are pliable and are used to complete embolization of the aneurysm neck or lesion. The ACS is offered in various diameters and lengths and the use of the coil size is determined by a physician. The ACS is offered in 10 and 18 compatible systems and are delivered through the following microcatheters with the specified minimum ID:

Coil Type	Microcatheter ID	
	Inches	mm
ACS 10 System	0.017	0.43
ACS 18 System	0.017	0.43

**MODEL NUMBER DESCRIPTION**

**INDICATIONS FOR USE**

The Wallaby Avenir™ Coil System is indicated for the endovascular embolization of intracranial aneurysms and other neurovascular abnormalities, such as arteriovenous malformations and arteriovenous fistulae. The Avenir™ Coil System is also indicated for arterial and venous embolizations in the peripheral vasculature

**POTENTIAL COMPLICATIONS**

Potential complications include, but are not limited to:

- hematoma at the site of entry
- vessel perforation
- aneurysm rupture
- parent artery occlusion
- incomplete aneurysm filling
- emboli
- hemorrhage
- ischemia
- vasospasm
- coil migration or misplacement
- premature or difficult detachment of coil
- clot formation
- revascularization
- post-embolization syndrome
- neurological deficits including stroke and possibly death
- potential complication from X-ray exposure: alopecia, skin reddening, ulceration; ocular cataracts and neoplasia

**REQUIRED ADDITIONAL ITEMS**

- Microcatheter with 2 radiopaque markers
- Guide catheter compatible with microcatheter
- guidewire
- Rotating hemostatic valves (RHV)
- 3-way stopcock
- Sterile saline
- Pressurized sterile saline drip
- Femoral sheath
- IV pole


**CONTRAINDICATION**

There are no known contraindications.

**WARNING AND PRECAUTIONS**

- Rx Only.
- Federal law (USA) restricts this device to sale by or on the order of a physician.
- The ACS is sterile and non-pyrogenic unless the package is opened or damaged.
- The ACS is intended for single use only. Do not re-sterilize and/or reuse the device. Reuse of the device may result in re-infection or cross infection. After use, dispose in accordance with hospital, administrative and/or local government policy.

- Do not use if package is damaged or breached.
- The ACS must be delivered through a reinforced microcatheter with an inner surface coating. Damage to the device may occur and necessitate removal of both the ACS and the microcatheter from the patient.
- High quality digital subtraction fluoroscopic road mapping is mandatory to achieve correct placement of the ACS.
- Advance and retract the ACS device slowly.
- Do not advance the delivery pusher with excessive force. Determine the cause of any unusual resistance and remove the ACS and inspect for damage.
- Remove the entire ACS if excessive friction is noted with a second ACS device. Inspect the microcatheter for damage or kinking.
- If repositioning is necessary, take special care to retract the ACS under fluoroscopy. If the ACS does not move in a coil to delivery pusher one-to-one motion or if repositioning is difficult, the coil may have become stretched and could possibly break. Gently remove and discard the ACS.
- Due to the delicate nature of the ACS, the tortuous vascular anatomy that lead to aneurysms and vessels, and the varying morphologies of intracranial aneurysms, a coil may occasionally stretch while being maneuvered. Stretching is a precursor to potential coil breakage.
- If resistance is encountered while withdrawing a coil that is at an acute angle relative to the microcatheter tip, it is possible to avoid stretching or breaking the coil by carefully repositioning the distal tip of the microcatheter at, or slightly inside, the ostium of the aneurysm. By doing so, the aneurysm and the artery act to funnel the coil back into the microcatheter.
- Delivery of multiple ACS coils is usually required to achieve the desired occlusion of some aneurysms or lesions. The desired procedure endpoint is angiographic occlusion.
- The long-term effect of this product on extravascular tissues has not been established so care should be taken to retain this device in the intravascular space.
- Always ensure at least two ACS coils are available before starting the procedure.
- When initiating detachment by "bend and pull" method of the proximal pusher, ensure that surgical gloves are not punctured.
- Proper positioning of the operator's hands at the center of the bend indicator strips is required to expose the coil release wire.

**Figure 1 Procedure Setup**

1. Flush Line
2. Rotating Hemostatic Valve (RHV)
3. Stopcock
4. Guide Catheter
5. Microcatheter
6. Femoral Sheath
7. 3cm radiopaque marker
8. Distal radiopaque marker
9. Avenir Coil System
10. Aneurysm
11. Avenir Coil System
12. Microcatheter hub

**PREPARATION FOR USE**

1. Refer to the Figure 1 diagram for procedure set-up. Appropriately flush all components of the system.
2. Attach a rotating hemostatic valve (RHV) to the hub of the guiding catheter. Attach a stopcock to the side arm of the RHV and then connect a line for continuous infusion of flush solution.
3. Attach a second RHV to the hub of the microcatheter. Attach a stopcock to the sidearm of the second RHV and connect the flush solution line to the stopcock.
4. Open the stopcock and flush the microcatheter with sterile flush solution and then close the stopcock. To minimize the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate sterile flush solution be maintained into the femoral sheath, the guide catheter and the microcatheter.

**CATHETERIZATION OF THE LESION**

5. Using standard interventional procedure, access the vessel with a femoral sheath and guide catheter. The guide catheter should have an inside diameter (ID) large enough to allow for contrast injection while the microcatheter is in place. This will allow for fluoroscopic road mapping during the procedure.
6. Select a microcatheter with the appropriate inner diameter. After the microcatheter has been positioned inside the lesion, remove the guidewire.

**COIL SIZE SELECTION**

7. Perform fluoroscopic road mapping.
8. Measure and estimate the size of the lesion to be treated.
9. Select the appropriate size coils.
10. Correct coil selection increases ACS effectiveness and patient safety. Occlusive efficiency is, in part, a function of coil compaction and overall coil mass. In order to choose the optimum ACS coil for any given lesion, examine the pre-treatment angiograms. The appropriate ACS coil size should be chosen based upon angiographic assessment of the diameter of the parent vessel, aneurysm dome and aneurysm neck. When assessing aneurysms, the diameter of the first and second coils should never be less than the width of the aneurysm neck or the propensity for the coils to migrate may be increased.

**PREPARATION OF THE ACS FOR DELIVERY**


11. Peel open the pouch by the chevron seal. Retrieve the hoop out of the pouch. Take care not to bend the silicone end fitting in the process.
12. Hold the hoop and pull the end fitting to bring the delivery pusher out. Hold the delivery pusher and pull off the end fitting. Discard the end fitting. Take care not to kink the pusher in the process.
13. Remove the ACS from the packaging hoop by pulling the pusher proximal end until the introducer exits the hoop.
14. At the proximal end of the sheath, push the sheath lock distal to release the pusher from the sheath.
15. Slowly advance the ACS coil out of the introducer and inspect the coil for any irregularities or damage. If any damage to the delivery pusher or coil is observed, DO NOT use the system.
16. Gently immerse the ACS coil and its detachment zone in heparinized saline. Take care not to stretch the coil during this procedure. While remaining in the heparinized saline and holding the introducer sheath pointing vertically downward, gently retract the coil completely back into the introducer sheath.
17. Pull the sheath lock proximal to re-lock the sheath to the pusher.

**INTRODUCTION AND DEPLOYMENT OF THE ACS**

18. Open the RHV on the microcatheter just enough to accept the introducer sheath of the ACS, flush the stopcock until saline flows from the proximal end of the RHV.
19. Insert the introducer sheath of the ACS into the RHV. Stop advancing the introducer sheath just prior to seating the distal tip at the distal end of the microcatheter hub.
20. Lightly tighten the RHV around the introducer sheath.

**WARNING**  
Do not over tighten the RHV around the introducer sheath. Excessive tightening may damage the ACS.

21. Flush the stopcock until saline is seen dripping from the proximal end of the introducer sheath.
22. Loosen the RHV on the microcatheter and seat the introducer sheath into the hub of the microcatheter. Lightly tighten the RHV around the introducer sheath.
23. Push the coil into the lumen of the microcatheter. Use caution to avoid catching the coil on the junction between the introducer sheath and the microcatheter hub.
24. Push the ACS through the microcatheter until the proximal end of the delivery pusher meets the proximal end of the introducer sheath. Loosen the RHV and retract the introducer sheath just out of the RHV. Close the RHV around the delivery pusher. Slide the introducer sheath completely off the delivery pusher. Use care not to kink the delivery pusher.
25. Loosen the RHV, continue to advance the delivery pusher into the microcatheter until the Radio Safe marker approaches the RHV. The Radio Safe marker indicates the implant coil is near the distal tip of the microcatheter. At this time, fluoroscopic guidance must be initiated.
26. Under fluoroscopic guidance, slowly advance the ACS coil out of the tip of the microcatheter. Continue to advance the ACS coil into the lesion until optimal deployment is achieved. Reposition if necessary.
  - a. If the coil size is not correct, remove and replace with another device.
  - b. If undesirable movement of the coil is observed under fluoroscopy following placement and prior to detachment, remove the coil and replace with another more appropriately sized coil. Movement of the coil may indicate that the coil could migrate once it is detached.
  - c. Do not rotate the delivery pusher during or after delivery of the coil into the aneurysm. Rotating the delivery pusher may result in a stretched coil or premature detachment of the coil from the delivery pusher.
  - d. Angiographic assessment should also be performed prior to detachment to ensure that the coil mass is not protruding into the parent vessel.
27. Advance the coil into the desired position until the proximal end of pusher 3cm marker is adjacent to the distal end of the microcatheter proximal marker, forming a "T", see Figure 2. This alignment positions the proximal end of the coil just outside the microcatheter. To minimize the potential risk of aneurysm or vessel rupture, DO NOT advance the delivery system so that the proximal end of the pusher 3-cm marker advances past the distal end of microcatheter proximal marker.
28. Tighten the RHV to prevent movement of the coil.

**Figure 2 Microcatheter and Avenir Coil system pusher distal tip alignment**


1. ACS Pusher 3cm marker
2. Microcatheter proximal marker
3. Microcatheter distal marker
4. ACS Pusher distal marker

29. Verify that the distal shaft of the delivery pusher is not under stress before coil detachment. Axial compression or tension could cause the tip of the microcatheter to move during coil detachment. Catheter tip movement could cause the aneurysm or vessel to rupture.


**WARNING**  
To minimize the potential risk of aneurysm or vessel rupture, DO NOT advance the delivery system so that the proximal end of the 3-cm marker advances past the distal end of microcatheter proximal marker

**DETACHMENT OF THE ACS COIL**


30. When the ACS detachable coil has been placed as desired, proceed with detachment as follows:
  31. Locate the "bend and pull" indicator on the proximal end of the delivery pusher.
  32. Confirm again under fluoroscopy that the proximal end of pusher 3cm marker is just distal to the proximal marker band on the microcatheter.
  33. Place each hand at the center of each bend indicator strips, bend to break the delivery pusher. Once the pusher segments are separated, pull on the proximal end of pusher by at least 1 inch to activate release wire to detach the coil from the delivery pusher. See figure 3-6.
  34. Verify detachment of the coil under fluoroscopy by first loosening the RHV around the delivery pusher, then pulling back slowly on the delivery pusher and verifying there is no movement of the coil. If the coil did not detach, remove and discard the ACS.
  35. Withdraw the delivery pusher from the microcatheter.
  36. Verify the position of the coil angiographically through the guide catheter.
  37. Repeat the above steps if additional coil placement is required.


**Figure 3 Hold center of each bend indicator**

**Figure 4 Bend proximal end of pusher**

**Figure 5 Pusher snaps at center point, and two segments are separated.**

**Figure 6 Pull on the proximal end of pusher to activate release wire and detach implant**

The physician has the discretion to modify the coil deployment technique to accommodate the complexity and variation in the embolization procedure. Any technique modifications must be consistent with the previous described warnings, precautions and patient safety data.

**PACKAGING AND STORAGE**

The ACS is placed inside a protective plastic dispenser hoop and packaged in a pouch and unit carton. The ACS and dispenser hoop will remain sterile unless the package is opened, damaged, or the expiration date is passed. Store in a room with controlled temperature in a dry place.

**SHELF LIFE**

See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

**MRI Conditional**

Non-clinical testing and MRI simulations were performed to evaluate the entire family (i.e., available diameters and lengths) of the ACS. Non-clinical testing demonstrated that the entire family of these embolization coils is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the ACS is expected to produce a maximum temperature rise of 2.2°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the ACS extends approximately 3-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. In addition, the use of a magnetic resonance angiography (MRA) pulse sequence to scan an aneurysm model (i.e., 40% packing density) at 3-T, using an echo time of 3.5-ms, generated a signal void artifact of 9-mm relative to edge of the coils.

**MATERIALS**

The ACS is not made with natural rubber latex or PVC materials

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