Rescue™

Single-Use Grasping Forceps

Instructions For Use

The logo for Anrei, featuring the word "Anrei" in a bold, italicized sans-serif font. A grey swoosh underline is positioned under the "A" and extends under the "n".The logo for rescue, featuring the word "rescue" in a lowercase, rounded sans-serif font. A trademark symbol (TM) is located at the top right of the "e". Below the word "rescue" is the text "Retrieval Devices" in a smaller, plain sans-serif font.

AP30047A Rev A 10-2018

Single-Use Grasping Forceps

Instructions For Use

⚠ ONLY

Warning:

Warning: Use same care as for any other surgical device to avoid punctures or cuts to hands, gloves, etc. Read this document in its entirety prior to use.

Device Description:

The single-use Rescue Grasper consists of an actuating jaw which can be opened and closed. When passed through an endoscope the grasper can be actuated to grasp and/or retrieve.

User Information

The intended user is a physician properly trained in endoscopic procedures.

Contents

Each packaged pouch contains the following:

One Rescue Grasping Forceps

Intended Use/Indications For Use:

Grasping forceps are intended to be used to grasp tissue and/or retrieve foreign bodies, excised tissue, and stents from within the gastrointestinal tract.

Contraindications:

None known within the scope of the indicated use.

Warning:

For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death.

Precautions:

None known.

Adverse Events:

None known.

How Supplied:

Contents supplied sterile using an ethylene oxide (EO) process. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Storage:

Store sterile forceps at room temperature. Do not expose to solvents, radiation, UV light or any other substance that could compromise the sterility of the unit.

Inspection:

- **WARNING:** Inspect the package containing the forceps. If any damage is noticed, do not use the forceps as sterility may be compromised.
- Peel the pouch open and remove the forceps. Visually inspect the forceps for kinks, loose or broken parts. Then form three 6 to 8 inch (15 to 20 cm) loops and open and close the jaws several times to ensure the forceps are operating properly.
- **WARNING:** DO NOT USE if any abnormalities are found or if the forceps are not operating properly.

Operation:

- Close the forceps jaws before inserting the forceps into the working channel of the endoscope. Use short motions when inserting the forceps to avoid kinking.
- If resistance is met when the bending section of the endoscope is angulated, release the bending angle slightly until the forceps pass smoothly. Do not force the forceps into an angulated scope as this could damage the endoscope.
CAUTION: Do not use excessive force to close the jaws as this could damage the forceps.
- To grasp an object, open the jaws carefully while visualizing the procedure through the scope. Advance the forceps to the target item. Close the jaws gently.
- When the grasping operation is complete, close the jaws fully and slowly withdraw the forceps from the scope keeping the jaws closed. If the scope has a forceps elevator, it must be lowered before withdrawing the forceps.
- If an item still remains in the jaws and the jaws will not close completely, do not use force to pull the forceps into the endoscope. Remove the endoscope with the forceps still in place from patient, release the item from the jaws and manually close the jaws before removing the forceps from the scope.
- This is a single-use device and is not designed to be reused or reprocessed. Discard after each procedural use.

Disposal:

- **WARNING:** After use, this product may be a potential biohazard. Handle or dispose of it in accordance with accepted medical practices and applicable local, state and federal laws and regulations.

Warranty:

- Manufacturer warrants that the forceps will be free from defects for use in a single procedure. This warranty does not apply to use beyond the expiration date. This device is not intended to be reprocessed. The user assumes responsibility for any consequence related to the use of reprocessed forceps. Manufacturer's obligations under this warranty are limited to the replacement of the forceps.

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Boston Scientific (Australia) Pty Ltd
Corp. GmbH (Europe)
Anrei Medical (HZ) Co., Ltd.
Boston Scientific Corporation