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Disposable Hemoclip

FR Manuel d'instructions

EN Instruction Manual

【Usage revendiqué】

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Révision : A/0


【Intended Use】

The Disposable Hemoclip is indicated for clip placement within the Gastro-intestinal (GI) tract for the purpose of:

1. Endoscopic marking,
2. Hemostasis for:
 - Mucosal/sub-mucosal defects < 3 cm,
 - Bleeding ulcers,
 - Arteries < 2 mm,
 - Polyps < 1.5 cm in diameter,
 - Diverticula in the colon,
3. Anchoring to affix jejunal feeding tubes to the wall of the small bowel,
4. As a supplementary method, closure of GI tract luminal perforations < 20 mm that can be treated conservatively.

【Types and specification】

Product sketch map as following:

【Warnings】

Contents supplied **STERILE** using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Fartley representative.

For single use only. **DO NOT REUSE, REPROCESS OR RESTERILIZE.** Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

The number of Clips required for hemostasis may vary depending upon anatomical site, histology, lesion type, and patient condition and history. A sufficient quantity of Clips should be prepared in consideration of all of these factors prior to procedure.

【Contraindications】

Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.

- Arteries greater than 2 mm
- Polyps greater than 1.5 cm in diameter
- Mucosal/Submucosal defects greater than 3 cm

Those specific to primary endoscopic procedure to be performed in gaining access to desired site. Those specific to endoscopic hemostasis include, but are not limited to: uncooperative patient, coagulopathy, cricopharyngeal or esophageal narrowing or stricture, and tortuous esophagus.

【Warnings】

Operation of this instrument is based on the assumption that open surgery is possible as an emergency measure if the clip cannot be detached from

the instrument or if any other unexpected circumstance takes place.

- Always have pliers and/or wire cutters ready to cut the delivery system at the handle if the clip cannot be detached

【Potential Complications】

Those associated with gastrointestinal endoscopy and endoscopic hemostasis include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, hematemesis, transient dysphagia, aspiration pneumonia, wound dehiscence, minimal acute inflammatory tissue reaction, transitory local irritation, migration of clip into the bile duct, and anatomy disruption.

Limited studies indicate that the use of clips in the presence of bacterial contamination may increase or prolong infection.

- Re-bleeding may occur if the clips detach within 24 hours.
- Although rates of occurrence are low, recurrent bleeding, ineffective clipping or endoscopic complications could result in the need for surgery.

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The vitro simulative test on single sterile repositionable hemostasis clipping has confirmed that, according to the ASTM F2503 standard, patients with one clip or multiple dispersive clips or even multiple centralized clips which

not touch each other, could get safely scanned under magnetic resonance parameters as follows.

Static Magnetic Field

Static magnetic field of 1.5-Tesla or 3-Tesla. The maximum spatial gradient should be no more than 720 G/cm.

The maximum MRI system report: the averaged specific absorption rate in the whole body level is 2 W/kg.

MRI-Related Heating

After 15 minutes of continuous scan based on the parameters described as above, the increase of temperature caused by the clips is no more than 2°C.

Artifact Information

The simulative test demonstrate that the imaging within about 80 mm from the clips will be distorted when imaging with spin-echo or gradient echo sequence in 3T MRI system.

Caution!

If multiple clips touch each other, the consequences of MRI scanning have not been

adequately studied. We suggest performing MRI scan after clips getting off along with the connective tissues.

【Pre-procedure】

Passage of the Disposable Hemoclip through a retroflexed or tortuous path, may result in the clip separating from the catheter and potentially kinking or damaging the device. If the catheter or over-sheath kinks or becomes damaged during device insertion or passage, do not use it. Call Fartley Customer Service and return the product.

- Applying tangential pressure to an opened or closed clip may result in the clip separating from the catheter and potentially kinking or damaging the device. If the catheter or over-sheath kinks or becomes damaged during device insertion or passage, do not use it. Call Fartley Customer Service and return the product.
- In a difficult scope position, it may be necessary to straighten the endoscope to facilitate the device passage, then reposition scope for treatment. If the catheter or over-sheath kinks or becomes damaged during device insertion or passage, do not use it. Call Fartley Customer Service and return the product.
- The Disposable Hemoclip is supplied sterile. Carefully examine the unit to verify that neither the contents nor the sterile package has been damaged in shipment. DO NOT USE if damaged. Immediately return damaged product to Fartley.
- Refer to package label for minimum channel size required for this device.

- Endoscope must remain as straight as possible when inserting or withdrawing device.
- Clip must remain closed during introduction into, advancement through and removal from endoscope. If clip is open, damage to clip and endoscope may occur.
- Exercising handle while clip is coiled may result in damage to clip.
- Endoscopic hemostasis should be performed under direct endoscopic visualization.
- Clipping hard or severely fibrotic lesions to achieve hemostasis may be more difficult.
- If clip deployment device is used with endoscope in a torqued or retroflexed position, clip deployment difficulties can occur.
- Removing undeployed device through a retroflexed scope can cause detachment of the clip.
- Lesions located in esophagus and lesser curvature of stomach may be difficult to treat with a forward-viewing endoscope.
- Treatment of esophageal varices may require clipping in combination with a sclerosing agent.
- Use of clip in presence of bacterial contamination may potentiate or

prolong infection.

【Instructions for Use】

1.Hemoclip insertion:

A.Remove the Hemoclip Protective shell at the front end of Disposable Hemoclip, and pull backward with the Handle Slip ring to maintain a certain tension to fully close the Hemoclip, as shown in Figures 2 and 3.


Figure 2 Remove the Protective shell

Figure 3 Close the Hemoclip

B. Carefully insert the front end of the Hemoclip into the endoscopic biopsy channel, and advance only 2-3 cm each time.

Note: when this instrument passes through the backward or curved path, it

may cause the Hemoclip to detach from the Outer tube, kink or damage the Hemoclip. Do not use this instrument if the Outer tube is kinked or damaged during the Hemoclip insertion or advance.

2.Clamping and release:

A. When the Hemoclip is in the required position, after finding a suitable position, gently move the distal Slip ring away from the Thumb ring to open the Hemoclip. (as shown in Figure 4) At this time, you can revolve the rotating wheel (knob)) to adjust the position of the Hemoclip (as shown in Figure 5)


Figure 4 Open the Hemoclip


Figure 5 Rotate a
adjust the Hemoclip

B. Close the Hemoclip at the required position and move the Slip ring to the near end. It should be felt that there is resistance on the Handle to gradually close the Hemoclip. (as shown in Figure 6) In addition, the position of the Clip needs to be evaluated before complete closure.

Note: after you feel resistance, do not continue to move the Slip ring to the

near end, otherwise you may not be able to reopen the Hemoclip. If you hear a click, the Hemoclip is released.


Figure 6 Close the Hemoclip

3. At this point, there are two options:

Option 1: you can reopen and position the Hemoclip. (See step 4)

Note: the Hemoclip can be opened and closed for an unlimited number of times before complete closure, and there is no limit of times. It can be used to assist in positioning the Clip at the injury area. The ability to reopen and close may be limited by the clinical environment, the patient's lesion, and other factors.

Option 2: you can permanently close the Hemoclip. (see step 5)

4. To reopen the Hemoclip, and place it in another position. Slowly move the Slip ring to the far end to open the Hemoclip, and then move the Clip to the required position (see step 1)

5. To permanently close the Hemoclip, continue to move the Slip ring

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Yangzhou Fartley Medical
Instrument Technology Co., Ltd
Beizhou Road, Lidian Town,
Guangling District,
Yangzhou, P.R.C., 225106

Prolinx GmbH
Brehmstr. 56, 40239
Düsseldorf Germany
Tel: 0049 2131 4051968-0
DE/0000045300

toward the near end after reaching the resistance point, where you can feel or hear a click, continue to move the Slip ring toward the near end until it reaches the Thumb ring, as shown in Figure 7.


Figure. 7 Permanently close the Hemoclip

6.Gently pull the released Hemoclip outwards to separate the Hemoclip from the conveying device.

7.After moving the Slip ring to the far end, slowly withdraw the Hemoclip from the endoscope.

【Storage】

The single package consists of a sterile paper-plastic bag; one packaging unit consists of paper boxes and the transport package consists of paper cartons. The products should be stored in well-ventilated, dry and non-corrosive rooms.

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Yangzhou Fartley Medical
Instrument Technology Co., Ltd
Beizhou Road, Lidian Town,
Guangling District,
Yangzhou, P.R.C., 225106

Prolinx GmbH
Brehmstr. 56, 40239
Düsseldorf Germany
Tel: 0049 2131 4051968-0
DE/0000045300

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

【Explanation of symbols used on labels and instructions for use】

	Single use only		Do Not Re-Sterilize
	Date of Manufacture		Manufacturer
	Use By		Caution
	Reference		Lot number
	CE Mark		Sterilized by Ethylene Oxide
	Diameter		Working length
	Applied to channel dia.		Authorized Representative in the
	Refer to IFU before operation		Storage in dry environment

	Do Not Use if packaging or product damage is evident.		Pack Quantity
--	--	---	----------------------

【ATTENTION】

Sale, distribution of this product should follow local laws and regulations, used by the order of a physician.

Date: April 2020

Revision: A/0

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Distributed by: Asept InMed SAS