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INSTRUCTIONS FOR USE – PHYSICA ZUK UNICOMPARTMENTAL KNEE SYSTEM**Intended for cemented fixation**

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the available product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

The main goal of a knee joint prosthesis is to reproduce the articular anatomy. The joint prosthesis is intended to reduce pain and give articular mobility to the patient. The degree of pain reduction and mobility depends in part on the pre-operative situation, intra-operative options and on the post-op rehabilitation. Physica ZUK knee system is a unicompartmental knee replacement system: it is designed for the surface reconstruction of individual femoral/tibial condyles. Physica ZUK knee system includes femoral components, tibial baseplates, tibial liners and All-Poly tibial components. These devices are intended for cemented use only.

1.1. MATERIALS

Femoral components are made from Cobalt-Chromium-Molybdenum Alloy, tibial baseplates from Ti-6Al-4V Alloy, articular surfaces and All-Poly components from Ultra-High Molecular Weight polyethylene. All-Poly components have a 316L stainless steel radiographic indicator wire embedded in polyethylene; femoral components and tibial baseplates have a poly methyl methacrylate (PMMA) pre-coating on cement interface surfaces to enhance fixation. This means that, although there is a coating, the use of surgical cement is mandatory to implant the Physica ZUK Implants.

Components	Material
Femoral component	CoCrMo + PMMA
Tibial plate	Ti6Al4V + PMMA
Tibial liner	UHMWPE

All-Poly component

UHMWPE + AISI 316L

Applicable standards

CoCrMo (ISO 5832-4 – ASTM F75) – **UHMWPE** (ISO 5834-2 – ASTM F648) – **AISI 316L** (ISO 5832-1) – **Ti6Al4V** (ASTM F136)

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon. Follow the instructions of the manufacturer of bone cement (PMMA) for preparing the bone and applying the bone cement for prosthesis fixation.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0–50° C / 32–122° F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and sudden changes in temperature. Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or the surface integrity. Careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. Components removed from the package should not be used if they are dropped or suffer other accidental impacts. Devices should not be modified in any way. The device's code and lot number should be recorded in the patient's case history by using the labels included in the components packaging. The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided. Risks associated with reuse of single use devices are:

- infection;
- early or late failure of the device or device fixation;
- lack of appropriate coupling between modular junctions (e.g. taper connections);
- device wear and wear debris associated complications;
- transmission of diseases (e.g. HIV, hepatitis);
- immune system response / rejection.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS

2.1. INDICATIONS

These devices are indicated for patients with:

- painful and/or disabling knee joints due to osteoarthritis or traumatic arthritis;
- previous tibial condyle or plateau fractures with loss of anatomy or function;
- varus or valgus deformities;
- revision of previous arthroplasty procedures.

The Physica ZUK Unicompartmental Knee System is designed for use when load bearing ROM is expected to be less than or equal to 155 degrees. These devices are indicated for cemented use only.

2.2. CONTRAINDICATIONS

Contraindications include: • previous history of infection in the affected joint and/or local/systemic infection that may affect the prosthetic joint; • insufficient bone stock on femoral or tibial surfaces; • skeletal immaturity; • neuropathic arthropathy; • rheumatoid arthritis; • obesity¹; • varus/valgus deformities greater than 15 degrees; • chondrocalcinosis; • osteoporosis; • damage to the articular cartilage of the opposite compartment; • inflammatory synovitis; • eburnated bone in the patellofemoral joint.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight²; • strenuous physical activities (active sports, heavy physical work); • incorrect implant positioning; • medical disabilities which can lead to an unnatural gait and loading of the knee joint; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • overcorrection of the limb.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products should be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques. Physica ZUK articular surfaces are available in six sizes that match exactly the tibial baseplate sizes: they are mated one-to-one (eg., size 1 articular surface is only used with size 1 tibial baseplate); all femoral component can be used with all sizes of articular surface and tibial baseplate components. Physica ZUK components have exactly the same design of Zimmer Unicompartmental Knee System (ZUK) components. Therefore two equivalent versions (ZUK and Physica ZUK) exist for each component. A mixed coupling between ZUK and Physica ZUK components is allowed. Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning. Complications or failures of the total knee replacement are more likely to occur in heavy and highly active patients. The surgeon should perform a careful evaluation of the patient's clinical condition and level of physical activity before performing knee replacement. An inadequate clinical condition and level of physical activities evaluation may adversely affect the outcome of the joint replacement. LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning,

1. According to the definition of the World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 30

2. According to the definition of the World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25

the surgical technique, and product and instrumentation assistance both prior to and during surgery. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can fail or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and may need to be replaced at some time in the future: the longevity and the performance of the implant over time, in fact, can be influenced by the the natural/physiological progression of the patient pathology, the insurgence of co-morbidities and post-operative complications, even if they are not directly implant-related (i.e. infection, pain, stiffness, reduced range of motion). The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors. All-Polyethylene tibial baseplates should be limited to use in low demand (modest weight, activity level) patients with good bone quality. Fat embolism risk is increased with intramedullary instrumentation and/or cement pressurization. Consider venting the femur or tibia. Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person. Do not use this product for other than labeled indications (off-label use); the risk of implant failure is higher with inaccurate component alignment or positioning. Use only instruments and provisionals specifically designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of knee function. Selection of polyethylene components is a matter of physician discretion. Thicker polyethylene components may be needed if the patient is young, heavy, and/or physically active. Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants are available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection as well as the correct seating/placement of the implant is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may affect system performance and survival of the implant. The components must be assembled according to the surgical technique and used only for the labelled indications. Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery. Care must be taken to protect the surfaces involved in the coupling between components; the articular surfaces of the implants should be protected from scratches or any other damage. All component coupling surfaces should be clean and dry before assembly. The stability of component couplings should be verified as described in the surgical technique.

Do not reinsert an articular surface implant that has been inserted previously. There may be visually nondetectable flaws that could reduce the service life of the implant. Do not use any component if damage is found or caused during setup or insertion. With simultaneous bilateral knee surgery, release leg tourniquets ten minutes apart to lessen any lung insult that may occur. Prior to closure, the surgical site should be thoroughly cleansed of bone chips, bone cement.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follow up is recommended to detect possible changes in position or condition of the implant or surrounding tissues. The surgeon should make the patient aware that caution is needed after the unicompartmental knee arthroplasty, especially in the post-operative period. Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear, and/or fracture of the implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

In particular the following precautions should be presented to the patient by the surgeon:

- avoid repeated high weight lifting;
- keep body weight under control, overweight conditions may adversely affect the outcomes of joint replacement;
- avoid sudden peak loads (consequences of activities, such as running and skiing) or movements which can lead to sudden stops or twisting;
- avoid positions that can increase the risk of dislocation.

The following precautions have to be taken:

- avoid notching, scratching, or striking the device;
- the potential for deep sepsis can be minimized by using biocontamination controls;
- continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that can occur in knee arthroplasty include:

- loosening of the prosthetic components;
- prosthesis dislocation and instability;
- damage to the prosthetic implant;
- instability of the system because of inadequate soft tissue balancing;
- dissociation due to incorrect coupling of the devices;
- swelling or infection;
- tissue reaction to the implant material or abrasion;
- local hypersensitivity;
- bone fracture or nerve damage;
- local pain;
- progressive disease in the opposite department;
- periprosthetic fractures;
- temporary or permanent nerve damage;
- flexion contracture;
- limited range of motion;
- inflammation;
- lengthening or

shortening of the leg; • fracture of the devices; • excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles; • wear debris can initiate osteolysis which may result in loosening of the implant; • corrosion of metal components (the significance and long-term implications are uncertain and await further clinical evidence and evaluation); • malalignment of the prosthetic knee components; • additional surgery.

General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, hematomas, systemic allergic reactions and systemic pain. Some adverse effects can lead to death.

4. STERILITY

a. Implants

All implantable components of the Physica ZUK knee system are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . All components are sterilized by radiation. Do not use any component from a package that has been previously opened or appears to be damaged. Do not use implants after the expiration date printed on the label.

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipments.

5. MAGNETIC RESONANCE IMAGING (MRI)

Physica ZUK components have not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near the implant site.

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