Comet™
Pressure Guidewire

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	12
Gebrauchsanweisung	17
Istruzioni per l'uso	22
Gebruiksaanwijzing	27
Instruções de Utilização	32

TABLE OF CONTENTS

WARNING3

DEVICE DESCRIPTION3

 Figure 1. Comet™ Pressure Guidewire3

USER INFORMATION.....3

 Contents3

INTENDED USE.....3

INDICATIONS FOR USE.....3

CONTRAINDICATIONS.....3

WARNINGS.....3

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....4

 Handling and Storage4

OPERATIONAL INSTRUCTIONS.....4

 Preparation5

 Zeroing the Pressure Guidewire5

 Diagnostic Procedure5

 Intervention.....5

 How to Finish the Procedure5

 Technical Specifications6

WARRANTY6

Comet™

Pressure Guidewire

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Comet Pressure Guidewire is a single use, hydrophilic-coated, steerable wire with a custom optical cable. A pressure sensor is mounted approximately 3 cm from the distal end of the radiopaque and shapeable straight tip. For product specifications, including wire diameter, length, and radiopaque tip length, please refer to the product label.

Brachial and femoral markers are located on the proximal segment of the pressure guidewire to aid in estimating the pressure guidewire's position relative to the distal guide catheter tip. The proximal markers are compatible with brachial and femoral guide catheters that are at least 90 cm or 100 cm long, respectively.

USER INFORMATION

Intended users of the pressure guidewire are primarily interventional cardiologists trained in performing endovascular procedures in the sterile field. Catheterization lab technicians and/or nurses trained in assisting endovascular procedures are also intended users for preparation and system utilization.

Contents

Quantity	Material
(1)	Pressure Guidewire
(1)	Torque Device
(1)	Optical Cable

INTENDED USE

The Comet Pressure Guidewire measures blood pressure gradient across coronary lesions during endovascular procedures. The Comet Pressure Guidewire may also be used as a coronary guidewire for interventional treatments.

INDICATIONS FOR USE

The Comet Pressure Guidewire is indicated to direct a catheter through a blood vessel and to measure physiological parameters in the coronary blood vessels.

CONTRAINDICATIONS

The Comet Pressure Guidewire is contraindicated for use in the cerebral vasculature.

WARNINGS

- The Comet Pressure Guidewire should be used only by physicians trained in angiography and percutaneous interventions including transluminal coronary angioplasty (PTCA). Vessel trauma (dissection, perforation, rupture or injury) may result from the improper use of this device.
- Resulting pressure guidewire fractures/separations might require additional percutaneous intervention or surgery.
- Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated.


Figure 1. Comet Pressure Guidewire

- Severe reaction may occur in response to contrast agents that cannot be adequately premedicated.

PRECAUTIONS

1. Remove the pressure guidewire carefully from the carrier tube to reduce the possibility of damage to the distal tip. Refer to the Diagnostic Procedure section.
2. The pressure guidewire should be wiped clean with heparinized saline prior to use to prevent thrombus accumulation.
3. Maintain diligent control of the distal tip at all times during an intervention to avoid vessel dissections and perforations. The pressure guidewire should be manipulated in the body only under fluoroscopy. The pressure guidewire should never be advanced and/or torqued against resistance without first determining the reason for resistance. If ability to torque the pressure guidewire is diminished during procedure and the source of torque loss is not known, remove the guidewire without further torquing.
4. Exercise care in handling the pressure guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, tip or wire separation or other damage to the pressure guidewire.
 - Do not use a pressure guidewire that has been damaged as it may decrease performance characteristics and/or lead to vessel damage.
 - Do not attempt to straighten a pressure guidewire that has been kinked or bent.
 - Do not leave the pressure guidewire in a prolapsed condition, as damage to the pressure guidewire may occur.
 - Wire damage may occur when the pressure guidewire is manipulated in a sharp bend, such as that caused by incomplete coaxial alignment of the delivery catheter and ostium, creating a sharp bend of the pressure guidewire between the delivery catheter and the vessel wall. Sharp bends are more likely to occur when using a less supportive delivery catheter, such as a diagnostic catheter.
5. When crossing a stent, exercise care to avoid entanglement between the pressure guidewire and the stent.
6. Avoid abrasion of the pressure guidewire coating.
 - To avoid damage to the hydrophilic coating, do not withdraw or manipulate the pressure guidewire in a metal cannula or sharp object.
 - Due to the variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, use a different catheter.
 - Excessive tightening of the torque device onto the pressure guidewire may result in abrasion of the coating on the pressure guidewire.
7. Exercise care in handling the optical cable during a procedure to reduce the possibility of accidental breakage, bending, kinking, or other damage.
8. Use only the optical cable provided to connect the pressure guidewire to FFR Link. Use of a different optical cable will produce inaccurate pressure readings.
9. Keep the proximal end of the pressure guidewire clean to ensure proper connection to the optical cable.
10. The accuracy of the diagnostic information is affected by, but not limited to:
 - Failure to achieve maximum coronary and myocardial hyperemia if using FFR (Fractional Flow Reserve) modality.
 - Interventional devices, such as balloon catheters, which are positioned so as to affect the blood flow or guidewires that stretch the vessel.

- Pressure wire positioning relative to the lesion.
- Microvascular resistance.

11. Carefully check and match therapeutic device compatibility to the pressure guidewire prior to use. Refer to instructions supplied with any interventional devices to be used in conjunction with pressure guidewire for their intended uses, contraindications and potential complications.
12. Do not use the pressure guidewire in conjunction with atherectomy catheters. It may damage the pressure guidewire.

ADVERSE EVENTS

Potential adverse events which may result from the use of the device include but are not limited to:

- Abrupt closure
- Allergic reaction
- Embolism
- Exposure to biohazardous material
- Infection
- Prolonged procedure
- Restenosis (reocclusion)
- Spasm
- Stroke/cerebral vascular accident (CVA)/transient ischemic attack (TIA)
- Vascular thrombus
- Vessel trauma (dissection, perforation, rupture or injury)

In addition, when used for interventional procedures:

- Angina or unstable angina
- Arrhythmias
- Cardiac tamponade/pericardial effusion
- Contrast induced renal insufficiency or renal failure
- Death
- Myocardial infarction or ischemia

Some of the above potential adverse events may require additional surgical intervention.

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Do not use if the use by date has been exceeded.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

The following are required for use with the Comet™ Pressure Guidewire:

- iLab™ Polaris Multi-Modality Guidance System or an ANSI/AAMI-BP22 compatible hemodynamic receiving system
- FFR Link
- Hemodynamic Cable Kit

Preparation

1. Ensure the receiving system and FFR Link are connected, powered on, and functioning properly.

Note: Refer to the receiving system User Guide throughout this procedure for instructions on using the system.

2. Remove the sterile pouch from the carton. Verify the device is suitable for use by checking the expiration date and ensuring the packaging and sterile barrier are intact. Do not use if the sterile package is damaged or the product is expired.
3. Use standard sterile technique to open the pouch and remove the tray lid. Leave the pressure guidewire in the carrier tube. Ensure that the pressure guidewire is securely connected to the optical cable.
4. Flush the carrier tube with heparinized saline by injecting into the hub end of the carrier tube to hydrate the hydrophilic segment of the pressure guidewire.

Zeroing the Pressure Guidewire

5. Place the carrier tube on a flat surface.
6. Insert the optical cable connector into the Pd IN port on the FFR Link. Upon recognizing the pressure guidewire, the FFR Link will perform an automatic zeroing.

Note: If the pressure guidewire is not recognized, attempt to disconnect and reconnect the pressure guidewire to the FFR Link. If the system still will not recognize the pressure guidewire, use a new pressure guidewire.

Diagnostic Procedure

7. Carefully remove the pressure guidewire from the carrier tube and inspect it for damage.

Note: Do not use a damaged device. If the pressure guidewire cannot be removed easily, repeat the injection of heparinized saline into the hub end of the carrier tube and attempt to remove the pressure guidewire again.

8. If desired, the tip may be carefully shaped according to standard tip shaping practices. Do not use a shaping instrument with a sharp edge and take care not to damage the sensor during shaping.
9. Open the hemostatic valve and the flush line of the coronary manifold. Insert a guidewire insertion tool through the valve and into the guide catheter.
10. Carefully insert the distal tip of the pressure guidewire through the insertion tool and into the guide catheter.
11. Remove the insertion tool from the hemostatic valve and continue to advance the guidewire. Tighten the hemostatic valve knurled knob so that the valve seals around the guidewire, but does not inhibit intentional wire movement.

Note: The brachial/femoral markers may be used to facilitate advancement.

12. Position the pressure guidewire sensor, which is 3 cm from the tip of the pressure guidewire, just distal to the tip of the guide catheter.
13. Flush the guide catheter with heparinized saline until all contrast is removed.
14. Ensure the aortic pressure transducer is at the height of the heart.
15. Allow the waveforms and trendlines to stabilize.

16. Use the system to equalize the pressure guidewire and aortic pressure transducer to 1.00 +/- 0.01. During equalization, do not disconnect the optical cable from the pressure guidewire or FFR Link.

Note: If the equalization is not successful, repeat the equalization procedure. Refer to receiving system Directions for Use for instructions.

17. Advance the pressure guidewire to the desired location for measurement.

Note: If a different tip shape or pressure guidewire is required, carefully withdraw and remove the pressure guidewire while observing pressure guidewire movement under fluoroscopy.

18. If using FFR modality, use standard technique to induce maximum hyperemic condition.
19. Begin recording pressure using the receiving system. If using FFR modality, record pressure until a steady state of maximum hyperemic condition is reached, or until the hyperemic effect starts to decrease.
20. Use the receiving system to stop recording, and review the data.
21. After every recording, withdraw the pressure guidewire to the equalization position, so that the sensor is just distal to the guide catheter tip, and recheck the pressure equalization. If necessary, repeat steps 16-20.
22. If additional measurements are to be taken from other locations, reposition the pressure guidewire and repeat steps 16-21.
23. If intervention is required, refer to the Intervention instructions below.

Intervention

The Comet™ Pressure Guidewire or a different guidewire may be used for delivery of the interventional device. If using the pressure guidewire for the intervention, perform the following steps:

24. Disconnect the pressure guidewire from the optical cable handle.
25. Remove the torque device.
26. Clean the proximal end of the pressure guidewire with a saline-soaked gauze.
27. Carefully advance the interventional device over the pressure guidewire according to manufacturer's instructions.
28. Perform the intervention using standard practice according to manufacturer's instructions.
29. If additional measurements are needed following intervention, carefully clean the proximal end of the pressure guidewire and reconnect the pressure guidewire to the optical cable handle.

Note: Exercise care to not excessively twist nut on optical cable handle.

30. Verify proper connection on the receiving system.
31. Return to step 16 above.

How to Finish the Procedure

32. When the procedure is complete, remove and discard the pressure guidewire and optical cable.
33. Use the receiving system to review pressure data.
34. Shutdown the FFR Link and receiving system when finished.

Technical Specifications

Parameter	Value
Pressure Accuracy	+/-3% or +/-3 mmHg, whichever is greater
Zero Thermal Effect	0.3 mmHg/°C
Zero Drift	< 3 mmHg/H
Frequency Response	> 25 Hz
Operating Range	-45 mmHg to 300 mmHg

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2018 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd