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MAGNETIC BLACK STAR
Instructions for use

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IFU Nr. 32-A

Symbols see part B

Instructions for use
MAGNETIC BLACK STAR

IFU Nr.32-A Magnetic Black Star EN 02/03.09.2015

Please note: This medical product may only be obtained and used by specialists, trained nursing and medical personnel.

- Description**
Double-J ureteral stent with magnet and accessories for urinary diversion. Retrieval device with magnet (stent removal system for magnetic stent).
- Contents and packaging**
Ureteral stent set
- Ureteral stent with magnet at the bladder end
- Guide wire, teflon-coated (in sets with guide wire)
- Pusher
- Double packaging
- Patient record and instructions for use
Retrieval device (stent removal system for magnetic stent)
- Double packaging
- Disposable catheter with Tiemann tip and magnet at the distal end
- Instructions for use
- Intended use**
Ureteral stents should ensure a urinary diversion from the kidney to the bladder, prevent extravasation of urine, or stent the ureter.
- Indications**
- Stenting of the ureter pre- and postoperatively
- Prophylactic or postoperative urine drainage security after URS
- Counteracting of strictures, also scar silts after ureter silts or silts at the ostium
- Ureteral stenoses of various etiology
- Ureteral strictures, stenoses, kinks, or adhesion
- Postoperative, e.g. after Ureteropelvic surgery or ureterotomy
- Contraindications**
- Serious infections
- Prolonged obstruction
- Contracted bladder
- Urosepsis
- Urge incontinence
- Chronic urinary retention or extrinsic obstruction.
- Cardiac pacemaker
- Possible side effects**
- Foreign body irritation in bladder and renal regions
- Erosions and ulcerations to the ureter
- Dysuria
- Hematuria
- Elevated micturition frequency, nighttime urge to urinate
- Feeling of pressure in the flank, flank pain
- Vesicoureteral reflux
- Incrustation
- Ascending infections
- Instructions**
7.1 Inserting the ureteral stent
a) Pass the guide wire, flexible tip first, into the renal pelvis.
b) When positioning the stent, make sure that the magnet (side of the bladder) is also led over the guide wire.
c) Under fluoroscopic guidance, push the ureteral stent via the guide wire into the renal pelvis using the pusher.
d) When the distal end of the ureteral stent reaches the renal pelvis, do not push the ureteral stent any further. The guide wire is withdrawn while the pusher holds the ureteral stent in position.
Once the guide wire has been removed, the memory effect of the ureteral stent is activated.
- Removing the ureteral stent with retrieval device (stent removal system for magnetic stent). Size "S" for ureteral stents FR/CH 4.8 / Size "L" for ureteral stents FR/CH 6-7**
a) To remove the stent, insert the retrieval device into the bladder retrogradely, b) The connection to the magnet on the ureteral stent is created by rotating the Tiemann tip (in front of the respective ostium).
c) Successful connection between the retrieval device and the stent magnet is indicated with an audible click.
d) After both magnets are connected, the ureteral stent is removed by careful, steady pulling.

- Possible complications and/or risks**
- Dislocation of the ureteral stent.
- Fragmentation of the stent with subsequent surgical or endoscopic removal of the product.
- Encrustation and blockage of the ureteral stent, rendering it non-functional and potentially causing urinary calculi to form in the kidneys.
- Pushing the ureteral stent forward with the guide wire can damage the ureter.
- Retracting the guide wire and pushing it forward again can cause the tip of the guide wire to exit through the holes in the ureteral stent, causing damage to the ureter area.
- Pushing the stent forward too hard with the pusher can damage the bladder wall.
- Due to obstructions in the ureter, it is not possible to advance the guide wire into the kidney.
- Strictures can develop through damage to the ureter.
- Manipulations with the guide wire can trigger ureteral spasms.
- Failure to lead the stent magnet over the wire may result in damage to the urethra.
- Insertion of the stent first with the end for the bladder (with magnet) may result in damage to the urethra.
- If incrustation of the kidney end pigtail is suspected, the stent should not be removed via retrieval device, since the impeded stretching of the pigtail can lead to injury of the ureter. In this case, use a cystoscope.
- Warnings / Precautions**
- This product is only sterile if the packaging is undamaged and unopened.
- For single use only!
- Do not use stents after the date of expiry has passed.
- Do not wipe the product and set components with alcohol/disinfectant agents. This may damage the surface.
- All components must be carefully checked for compatibility and integrity before use.
- Do not damage the stent, pusher, guide wire, and retrieval device with sharp surgical instruments - this increases the risk of breakage. Do not excessively stretch the stent before use. Stretching it too far increases the risk of breakage.
- Correct positioning and complete forming of the pigtail is required to minimize the risk of stent migration. This should be ensured through radiological or contrast agent-assisted fluoroscopy techniques.
- Any air in the lower abdomen of the patient could impair visibility of the ureteral stent in an x-ray.
- Where there are extensive defects in the ureter, a percutaneous nephrostomy should be carried out, since the blind pushing of the stent incurs the risk of enlarging any lesion further.
- In order to prevent loss of function due to encrustation, the max. service life of a magnetic ureteral stent is limited to 30 days.
These are not permanent implants, and it is therefore recommended that the urine status, stent position and patency of the ureteral stent be checked by radiography at regular intervals, and at least every 3-4 weeks, to ensure the early detection of any infection, dislocation, or blockages which might necessitate a stent change.
- Push the guide wire forward carefully and in a controlled manner to avoid damaging the ostium, the ureter, or the parenchyma.
- MRT examinations are not permitted for patients with ureteral stents.**
- Patients with cardiac pacemakers may not have the product implanted due to the magnet.
- Interaction with other drugs**
The ureteral stent could potentially discolor due to the concomitant administration of medication, use of dye-based urine presentation methods, or even as a result of contact with urine.
- Transport and storage conditions**
The products may be transported and stored only in the packaging intended for this. There are no further specific requirements for transport.
Products must be stored dry and protected from direct sunlight, in the temperature range of 5 - 30 degrees Celsius.
- Disposal**
After use, this product may pose a biological hazard. Handling and disposal must be carried out in accordance with recognized medical procedures and be completed pursuant to applicable legal regulations and guidelines.

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