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Yangzhou Farley Medical Instrument Technology Co., Ltd
 Beizhou Road, Lician Town, Guangling District, Yangzhou,
 P.R.C. | 225106
 Tél : +86(0514)87572788 Fax : +86(0514)87268752
www.farley.com

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Sponge protective case for 2 to 6-way infusion manifold and for infusion accessories: Luer connectors and 3-way stopcocks

Please read this leaflet carefully before using this medical device because it contains important information.
If you have more questions, please contact Asept InMed.

1. DESCRIPTION

The range is composed of 3 models of shapes and sizes adapted to the infusion devices to protect:

Reference	Description
200012	Sponge protective case for luer connectors
200006	Sponge protective case for 3-way stopcocks
200000	Sponge protective case for 2 to 6-way infusion manifold

Each case is composed of a red or blue soft shell, containing a foam sealed to the shell by a hot-melt adhesive.

All cases have a closure system with clips or clasp.

Cases for manifolds has an opening on the rear face to associate the manifolds with a clamp to fix on an infusion pole.

2. INDICATIONS

The case NeoProtect'in® is indicated to protect infusion manifolds and connectors.

It acts as a physical barrier that reduces the risk of contact and airborne contamination and prevents accidental mobilisation risks involving connectors.

The use of protective case is particularly recommended in the following cases:

- Paediatric and geriatric units;
- When infusion manifolds or accessories come into contact with patient bedding or any contamination (wound,...);
- When the 3-way stopcocks and manifolds are mobile (patient transfers,...);
- In all cases, to prevent unnecessary handling (by the patient, visitors,...).

3. CONTRAINDICATIONS

Not applicable.

4. WARNINGS

Must be used only by healthcare professionals familiar with good perfusion practices.

Do not use a protective case whose expiry date has passed.

Do not use a protective case whose integrity of the packaging (primary or secondary) is damaged.

Do not leave the case open.

5. USE PRECAUTIONS

Do not to soak the foam; the sponge case physically protects infusion connectors from external contamination and accidental mobilisation risks involving stopcocks and connectors (recommendations SFHH).

Special care is required when installing and then using the protective case (openings / closures) to ensure that the connections of the manifolds or infusion accessories are not unscrewed or mis-connected: screwing and sealing connections as well as the absence of (micro) bubbles in the infusion line before each closure of the protective case.

What to do in case of disconnection: stop the infusion, evaluate the impact on the patient's condition either by the appearance of symptoms of a gas embolism or hypovolemia.

Before handling infusion manifolds and accessories, refer to their leaflets for more information on use instructions, warnings, precautions for use and hygiene precautions.

6. HYGIENE PRECAUTIONS

The level of hygiene must be adapted when using this sterile device as well as sterile infusion manifolds or accessories.

Respect the general rules of hygiene and asepsis of health establishments.

It is recommended to wash your hands before and after use.

It is recommended to wear gloves when administering anticancer drugs.

7. FREQUENCY OF USE

The protective case should be replaced at the same frequency as the infusion line, according to good infusion practices (SFHH and HAS).

The case should not remain in place for more than 4 days (96h).

The case is for single use, it cannot be re-used after replacement of the infusion line.

The case must be replaced whenever its integrity is compromised.

8. STORAGE AND TRANSPORT CONDITIONS

Case should be stored and transported away from heat and moisture.

9. DISPOSABLE CONDITIONS

Used cases must be disposed of in the waste stream of health facilities.

10. MATERIAL

The soft shell is polyethylene, red or blue.

The inner foam is polyurethane polyether.

The glue is hot-melt adhesive.

The retaining clip and the clasp are in polyethylene.

11. PACKAGING AND STERILIZATION

The cases are packaged in single packaging and sterilized with ethylene oxide.

The cases cannot be re-sterilized.

Shelf-life: 5 years (maintaining the sterile state).

12. SYMBOLS

	Manufacturer		Do not re-sterilize
	Warning		Single use
	Special instruction for use		Do not use if packaging is damaged
	Sterilized with ethylene oxide		Reference
	Store away from moisture		Batch number
	Store away from heat		Expiry date

13. ADMINISTRATION INFORMATION

Date of the publication of the leaflet: December 2019 (reference FTL-SP-IFU-01 V/0)

Date of first affixing of the CE marking: December 2019 (CE 0197)

If you have any further questions or need information, please contact us at the following address:

Yangzhou Farley Medical Instrument Technology Co., Ltd
Beizhou Road, Lidian Town, Guangling District, Yangzhou,
P.R.C. | 225106
Tél : +86(0514)87572788 Fax : +86(0514)87268752
www.farley.com

Imported by :
ASEPT INMED
9, Avenue Mercure 31130 Quint Fonsegrives France
Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01
www.aseptinmed.fr

Prolinx GmbH
Address: Brehmstr. 56, 40239, Duesseldorf, Germany
Phone: 0049 2131 4051968-0 Fax: 0049 2131 4051968-9
E-mail: info@eulinx.com DIMDI Code: DE/0000045300