ENGLISH
CAPIO
POLYESTER SURGICAL SUTURES
USP

DESCRIPTION
Braided Polyester Surgical Suture (Green PTFE-Coated size 0) is a sterile nonabsorbable surgical suture composed of Poly(ethylene terephthalate), impregnated with Polytetrafluoroethylene. It is prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component.

Braided Polyester Surgical Sutures meet all requirements established by the United States Pharmacopeia (USP) for Nonabsorbable Surgical Sutures.

INDICATIONS
Braided Polyester Surgical Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, orthopedic and neurological procedures.

CONTRAINDICATIONS
None known.

WARNINGS 
Braided Polyester Surgical Sutures elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyester Surgical Sutures are not absorbed, nor is any significant change in tensile strength retention known to occur [in vivo](#).

As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Braided Polyester Surgical Sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Re-processing of products intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only products may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected, or re-sterilized.

PRECAUTIONS 
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid excessive handling i.e. crushing or crimping resulting from use of surgical instruments such as forceps and needle holders.

Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Strong familiarity with surgical procedures and techniques for nonabsorbable sutures is required before use.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in appropriate container ("sharps"). This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

ADVERSE REACTIONS 
Adverse effects associated with the use of this device include: wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

ENVIRONMENTAL CONDITIONS 
Polyester Suture with needle removed is MR Safe (i.e., an item that poses no known hazards in all MR environments).

HOW SUPPLIED
The Suture is provided sterile, for single use only. The suture is available in a variety of lengths, with and without needles. Finished suture may be packaged in cartons as single packs, multi-packs, or procedure packs.

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