ENGLISH
CAPIO
POLYPROPYLENE
SURGICAL SUTURE
USP EXCEPT FOR DIAMETER

DESCRIPTION

Polypropylene Suture is a nonabsorbable, sterile surgical suture composed of a strand of polypropylene, a synthetic linear polyolefin. Polypropylene Surgical Suture is available dyed (blue, light blue) or undyed (colorless). Polypropylene Suture meets all nonabsorbable surgical suture requirements established by the United States Pharmacopeia (USP) except for suture sizes 7-0 and 8-0 which differ from USP in diameter only.

MAXIMUM SUTURE OVERSIZE DIAMETER (mm) FROM USP	
USP SUTURE SIZE	MAXIMUM OVERSIZE (mm)
Size 7-0	0.015mm
Size 8-0	0.009mm

INDICATIONS

Polypropylene Surgical Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular and neurological procedures, but not for use in ophthalmic procedures.

CONTRAINDICATIONS

None known.

WARNINGS 

Polypropylene Suture elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Polypropylene Suture is not absorbed, nor is any significant change in strength retention known to occur *in vivo*.

As with any foreign body, prolonged contact of this or any other suture with soft solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Polypropylene Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Re-processing of products intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only products may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected, or re-sterilized.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid excessive handling i.e., crushing or crimping resulting from use of surgical instruments such as forceps and needle holders.

Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. The surgeon should employ appropriate suture technique according to their judgment and the standard of care to reduce the risk of suture erosion.

Strong familiarity with surgical procedures and techniques for nonabsorbable sutures is required before use.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in appropriate container ("sharps"). This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, calculus formation in urinary and biliary tracts when prolonged contact with soft solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

ENVIRONMENTAL CONDITIONS

Polypropylene suture with the needle removed is MR Safe (i.e., an item that poses no known hazards in all MR environments.)

HOW SUPPLIED

Polypropylene Surgical Suture is provided sterile for single use only. The suture is available in a variety of lengths, with and without needles. Finished suture may be packaged in cartons as single packs, multipacks, or procedure packs.

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