 CAPIO[®]
Monodek[®] (POLYDIOXANONE) Suture DYED MONOFILAMENT SYNTHETIC ABSORBABLE SUTURES, USP except for oversized diameter.

DESCRIPTION

Monodek[®] Monofilament Synthetic Absorbable Suture is prepared from the polyester, poly(p-dioxanone). The empirical molecular formula of the polymer is (C₄H₆O₃). Monodek[®] suture is provided sterile as follows:

- Undyed (Natural)
- Dyed (violet). The sutures contains less than 0.1% by weight of D & C Violet #2

Monodek[®] suture conforms to all requirements established by the United States Pharmacopeia (USP) for Absorbable Surgical Suture, except for oversized diameter, and the European Pharmacopeia (EP) for diameter for Sterile Synthetic Absorbable Sutures

Suture Size	USP Ave. Diameter Specification (mm) <8>1>	EP Average Diameter Specification (mm) <400>	Maximum Oversized Average Diameter from USP ¹ (mm)
0	0,350–0,399	0,400–0,499	0,1

To calculate the maximum oversized average diameter value subtract the max USP average diameter from the max EP average diameter.

INDICATIONS

Monodek[®] Monofilament Synthetic Absorbable Surgical Sutures are indicated for use in all types of soft tissue approximation, including use in pediatric cardiovascular tissue where growth is expected to occur and ophthalmic surgery. Monodek[®] suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue. These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to six weeks) is desirable.

CONTRAINDICATIONS

These sutures, being absorbable, are not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and are not to be used in conjunction with prosthetic devices, i.e. heart valves or synthetic grafts.

WARNINGS

Store at room temperature.

Avoid prolonged exposure to elevated temperatures.

Two important characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention, and second, the absorption rate (loss of mass).

Monodek[®] Monofilament Synthetic Absorbable Suture has been formulated to minimize the variability of these characteristics and to provide wound support through an extended healing period. The results of in-vitro bench test studies of Monodek[®] Monofilament Synthetic Absorbable Suture in buffered saline solution at 37 degrees Centigrade (body temperature) indicate that approximately 70% of its original straight tensile strength remains after two weeks at body temperature. At four weeks of body temperature exposure in the buffered saline, approximately 50% of its original strength is retained, and at six weeks, approximately 25% of the original strength is retained. This trend in strength loss with time for other polydioxanone sutures is confirmed by several in vivo studies reported in the medical literature. A review of the medical literature indicates that implantation studies in rats show the absorption (mass loss) of these sutures is minimal until about 3 months of implantation and the absorption (mass loss) is essentially complete within six months. Polydioxanone polymer has been found to be non-anitgenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

The use of this suture may be inappropriate in elderly, mal-nourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.

As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching, or distention, or which may require additional support.

As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Monodek[®] Monofilament Synthetic Absorbable Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Re-processing of products intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only products may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected, or re-sterilized.

PRECAUTIONS

The Monodek[®] suture knots must be properly placed to be secure. As with other synthetic sutures, knot security requires the standard surgical technique of flat and square ties with additional throws if indicated by surgical circumstance and the experience of the surgeon.

Avoid excessive handling i.e. crushing or crimping resulting from use of surgical instruments such as forceps and needle holders.

Conjunctival and Vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation and should be removed as indicated.

Subcuticular sutures should be placed as deeply as possible in order to minimize the erythema and induration normally associated with absorption. Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

Discard used needles in appropriate container ("sharps"). This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

ADVERSE REACTIONS

Due to prolonged suture absorption, some irritation and bleed-ing has been observed in the conjunctiva and mild irritation has been observed in the vaginal mucosa.

ENVIRONMENTAL CONDITIONS MR

Monodek[®] suture with the needle removed is MR Safe (i.e., an item that poses no known hazards in all MR environments).

HOW SUPPLIED

The Suture is provided sterile for single use only. The suture is available in a variety of lengths with and without needles. Finished suture may be packaged in cartons as single packs, multi-packs, or procedure packs.

Bondak and Monodek are registered trademarks of Teleflex Incorporated in the US and/or other countries.

Capio is a trademark of Boston Scientific Corporation

Bondak and Monodek are registered trademarks of Teleflex Incorporated in the US and/or other countries.

Capio is a trademark of Boston Scientific Corporation

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Capio ist ein Warenzeichen von Boston Scientific Corporation.