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Transfer Set

REF : APY1100CCG

REF : APY1200CCG

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ASEPT INMED

9, Avenue Mercure 31130 Quint Fonsegrives France

Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01

www.aseptinmed.fr

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Please read this leaflet carefully before using this medical device because it contains important information. If you have more questions, please contact, in France, Asept InMed, the distributor of Adria.med.

1. DESCRIPTION

The transfer set is equipped at one end by a spike (protected by a cap) to perforate the septum of a solution pouch. The other end is mounted with a 3 way stopcock to connect a syringe and the filling port of a portable infuser. A back check valve is positioned on the tube upstream the 3 way stopcock to prevent any back flow of the liquid to transfer. A clamp is placed on the tube to stop the flow at any time of the procedure.

2. INDICATIONS

This device is indicated to the preparation and transfer of drugs and solutions. It allows the transfer of drugs and solutions keeping the system closed. It can also be used connected via Luer Lock connections at all the drug's administration devices.

3. CONTRA-INDICATIONS

No known contraindication.

4. COMPLICATIONS

No known complications.

5. WARNINGS

- Do not use if the product is damaged, present a defect or is not functional ;
- Reuse and cleaning may alter the structural and mechanical properties. It may induce infection or other illness/injury ;
- Adria.med will not be responsible for any direct incidental or consequential damages resulting from products reuse.

6. USE PRECAUTIONS

- The product must be used only by qualified healthcare professionals ;
- The product must be replaced and disposed of as per the facility-approved protocol.

7. HYGIENE PRECAUTIONS

Respect the general hygiene and asepsis rules used by health establishments.

It is recommended to wash your hands before and after use and to disinfect hands with hydro alcoholic solution. The use of sterile gloves is possible strictly respecting the asepsis rules.


8. FREQUENCY OF USE

Single use only.

9. TERMS OF USE

- Prepare the infuser, the 50ml syringe, the transfer set, the solution bag.
- Connect the 50ml syringe at the entry of stopcock of the transfer set.
- Connect the perforator of the transfer set on the solution bag.
- Prime the transfer set aspirating the air towards the syringe
- Remove the stop cock cap and displace the air of the syringe than of the stop cock
- Remove the filling site cap located on the infuser than connect the stopcock on it.
- Realize the filling of the syringe and the successive transfers of the liquid towards the infuser, in order to reach the prescribed volume. It is not necessary to turn the stopcock, thanks to the back check valve (upstream the pouch side). The complete priming of the infuser tube must be checked by temporary and aseptically remove the cap from the tube connexion (patient side) and remplacing it after the dropping has been properly observed.

Example : Assembling of the Transfer Set for the filling of a Neopump infuser

10. STORAGE AND TRANSPORT CONDITIONS

These devices must be conserved, stored and transported avoiding in excessive heat, protect from direct sunlight and moisture.

11. DISPOSAL CONDITIONS

Used devices must be eliminated according to the hygiene and security protocol and to the actual regulation in the waste stream of the health facility or homecare.

12. MATERIAL

The perforator is in Acrylonitrile Butadiene styrene (ABS), the tube is in PVC without DEHP, the clamp is in Polyethylen high density, the stopcock is in Polycarbonate.

13. PACKAGING AND STERILISATION

These devices are conditioned in simple packaging and sterilized with ethylene oxide.
They cannot be re-sterilized; Do not use if packaging is damaged.
Expiring delay: 3 years (maintain performance and sterile state)
Do not use if expired.

14. ADMINISTRATIVE INFORMATION

Date of modification of the leaflet: September 2021
Date of first affixing of the EC marking: 2007

In case of serious incident related to these devices, please inform the manufacturer and ANSM (Competent Authority in France).

If you have any further questions or need information, please contact us at the following address:

ADRIA MED SRL
Via Ignazio Silone 1
65028 Tocco Da Casauria (PE)

Distributed by :

ASEPT INMED
9, Avenue Mercure 31130 Quint Fonsegrives France
Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01
www.aseptinmed.fr