UNI'perf

F Perfuseur

GB Infusion set

MODE D'EMPLOI

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The product should be used only by qualified healthcare professionals.

- Do not re-sterilize. Discard the set after single use.
- Reuse and cleaning of product may alter their structural and mechanical properties. It may lead to infection or other illness/injury.
- The product should be replaced and disposed of as per facility approved protocol or CDC guidelines.
- Store at room temperature, avoid excessive heat, protect from direct sunlight and moisture.
- Poly Medicure Limited will not be responsible for any direct incidental or consequential damages resulting from reuse of product.

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POLY MEDICURE LTD.

Plot No. 104-105, Sector 59, HSIIDC Indl. Area,
Faridabad-121004, (Haryana) INDIA,

INSTRUCTION FOR USE

- Ensure that the airway cap and the flow control clamp are closed.
- Insert the spike to full extent into the fluid container by applying a rotational movement.
- Raise or suspend the fluid container and squeeze the drip chamber to allow the solution to fill around half drip chamber (open the airway cap if the container is rigid).
- Open the flow control clamp allowing solution to displace air entirely in tubing.
- Then close the flow control clamp tightly and connect delicately the end luer lock of the infusion set to the vascular device (lock by carrying out 1/4 of turn).
- Regulate the flow rate by opening the flow control clamp gradually and according to counting of the drops.

CAUTION BEFORE USE

- Ensure the compatibility of the infused active ingredients with the tubing.
- Ensure the integrity of the device and the packaging before use.
- The infusion must be attentively monitored, according to the factors modifying the flow. The counting of the drops is obligatory in order to optimize the adjustment.
- Injection port to be used only for injecting the medicines.

STERILE EO

LATEX FREE

EC REP Obelis s.a
Boulevard Général Wahis 53
1030 Brussels, BELGIUM
Tel: +(32) 2. 732.59.54
E-mail : mail@obelis.net

Imported By :

 **Asept InMed®**
A MEDI-GLOBE COMPANY

9, avenue Mercure, 31130 QUINT FONSEGRIVES FRANCE
tel 33 (0)5 62 57 69 00 fax 33 (0)5 62 57 69 01

Rév. 02. 05/2020

DOC NO. Assept/I-000068	Date :20-1-2017	Dgn. By	Checked By	Approved by
INS SHEET FOR UNIPERF ALL	Size :180x185mm	 (Designer)	(Q.C. Officer)	(Q.A. Incharge)
Colour Break  BLACK	Specifications : 1. 3. 2. 4.			
Rev.	Ammendment		Remarks	Date of Amdt.
(1)	CE No. Changed			22-9-17
(2)	CE No. Changed and remove DEHP logo			28-5-20