**Medtronic**

**ITREL® 3**

**7425**


Neurostimulator  
Neurostimulateur  
Neurostimulator  
Neuroestimulador  
Neurostimulator  
Neurostimolatore  
Nervstimulator  
Νευροδιεγέρτης  
Neurostimulator  
Neuroestimulador

Implant manual • Manuel d'implantation  
• Implantationsanleitung • Manual de implantación  
• Implantatiehandleiding • Manuale per l'impianto  
• Implantationsmanual • Εγχειρίδιο εμφύτευσης  
• Implantationshåndbog • Manual de implante



**USA** Rx only


7425IPG\_fcv.fm 4/27/04 2:22 pm  
UC200xxxxxx EN  
4.625 x 6 inches (117 mm x 152 mm)

**Medtronic Confidential**  
**NeuroImplant\_R01**

221656002 Rev A

### Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Do not reuse

Sterilization: ethylene-oxide gas

Attention: See accompanying documentation.

Use by

Manufacturing date

Storage temperature

Serial number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC and R&TTE Directive 1999/5/EC.

The use of this device might be subject to individual country licensing regimes in Europe.



For USA audiences only

7425IPG\_sym.fm 4/27/04 2:22 pm  
UC200xxxxxx EN  
4.625 x 6 inches (117 mm x 152 mm)

**Medtronic Confidential**  
**NeuroImplant\_R01**

Itrel<sup>®</sup>, Medtronic<sup>®</sup>, and SoftStart<sup>™</sup> are trademarks of Medtronic, Inc.

2

221656002 Rev A

## Table of contents

<b>Device description</b>	<b>5</b>
<b>Package contents</b>	<b>5</b>
<b>Patient identification card and registration</b>	<b>5</b>
<b>Device specifications</b>	<b>5</b>
<b>Declaration of Conformity</b>	<b>9</b>
<b>Instructions for use</b>	<b>9</b>
Verifying neurostimulator operation	9
Connecting the extension to the neurostimulator	10
Implanting the neurostimulator	12
Checking system integrity	12
Completing the implant procedure	13

**Refer to the indications sheet for indications and related information.**

**Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.**



**Refer to System Eligibility, Battery Longevity, Specifications reference manual packaged with the software application card for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.**

**[USA] Refer to the clinical summary booklet packaged with the neurostimulator for information on the clinical study results of the neurostimulation system and individualization of treatment.**

7425IPGTOC.fm 4/27/04 2:22 pm  
UC200xxxxxx EN  
4.625 x 6 inches (117 mm x 152 mm)


**Medtronic Confidential**  
**NeuroImplant\_R01**

4 English 7425 2004-05

221656002 Rev A

## Device description

The Medtronic Model 7425 Itrel 3 Neurostimulator is part of a neurostimulation system for pain therapy.

## Package contents

- Neurostimulator
- Hex wrench
- Product literature
- Warranty card (USA only)
- Registration form
- Patient identification card

## Patient identification card and registration

A patient identification card is packaged with this device. Advise the patient to carry the identification card at all times.

USA The patient identification card packaged with the device is temporary; a permanent card will be mailed to the patient when Medtronic receives the registration form.

The implant registration form registers the device warranties and creates a record of the device in Medtronic's implant data system.

## Device specifications

The neurostimulator is a programmable device that accommodates one extension and delivers one program of stimulation through one or two leads.

A program of stimulation consists of an amplitude, pulse width, and rate delivered to selected electrodes.

**Table 1. Shipping, operating, and power-on-reset (POR) values for the Model 7425 Itrel 3 Neurostimulator<sup>a</sup>.**

Programmable parameter	Shipping	Operating	POR <sup>b</sup>
<b>Amplitude</b>			
<b>Normal resolution</b>	0.0 V	0.1 V steps	0.0 V
Upper limit	0.0 V	10.5 V maximum	0.0 V
Lower limit	0.0 V	0.0 V minimum	0.0 V
<b>Fine resolution</b>	-----	0.05 V steps	-----
Upper limit	-----	6.35 V maximum	-----
Lower limit	-----	0.0 V minimum	-----
<b>Rate</b>	31 Hz	49 values (from 2.1 to 130 Hz)	31 Hz
Upper limit	31 Hz	130 Hz maximum	31 Hz
Lower limit	2.1 Hz	2.1 Hz minimum	2.1 Hz
<b>Pulse width</b>	210 $\mu$ s	increments of 30 $\mu$ s	210 $\mu$ s
Upper limit	210 $\mu$ s	450 $\mu$ s maximum	210 $\mu$ s
Lower limit	60 $\mu$ s	60 $\mu$ s minimum	60 $\mu$ s
<b>Operation mode</b>	Continuous	Continuous or cycling	Continuous
<b>Cycle ON/OFF Time</b>			
W/O SoftStart/Stop	0.1 s	0.1 s to 24 h	0.1 s
With SoftStart/Stop	-----	1 s to 24 h	-----
<b>SoftStart/Stop</b>	Off	1 s, 2 s, 4 s, 8 s, or off	Off
<b>Ramp</b>	Off	15, 20, 25, 30 s, or off at start of stimulation	Off
<b>Dose time</b>	Off	15, 30, 45, 60, 75 m, or off stimulation periods	Off
<b>Dose lockout time</b>	0 h	1.0, 1.5, 2.0, 2.5, 3.0, 3.5, or 4.0 h	0 h
<b>Output On/Off<sup>c</sup></b>	Off	On or off	Off
<b>Magnet control</b>	On	On or off	On

**Table 1. Shipping, operating, and power-on-reset (POR) values for the Model 7425 Itrcl 3 Neurostimulator<sup>a</sup>. (continued)**

Programmable parameter		Shipping	Operating	POR <sup>b</sup>
<b>Polarity</b>	<b>Electrode number</b>			
<b>Program</b>	0	–	OFF, –, or +	–
	1	OFF	OFF, –, or +	OFF
	2	OFF	OFF, –, or +	OFF
	3	+	OFF, –, or +	+
	Case	OFF	OFF, or +	OFF

<sup>a</sup> All measurements are approximate.

<sup>b</sup> Power-on-reset (POR) turns OFF stimulation by resetting the amplitude to 0.0 V and all electrodes to OFF. POR occurs when there is a temporary fluctuation in battery voltage (eg, due to electromagnetic interference during electrocautery or defibrillation) or the battery is depleted. When POR occurs, the serial number is displayed as "?????" on the clinician programmer.

<sup>c</sup> Exposure to a strong source of electromagnetic interference (EMI) can cause the output switch to toggle On or Off. However, since the amplitude setting is 0.0 V, there is no neurostimulator output.

**Table 2. Physical characteristics of the Model 7425 Itriel 3 Neurostimulator.<sup>a</sup>**

Description	Value
Height	55.0 mm (2.2 in)
Length	60.0 mm (2.4 in)
Thickness	
case	10.0 mm (0.4 in)
connector block	10.1 mm (0.4 in)
Weight	42.0 g (1.5 oz)
Volume	22.0 cm <sup>3</sup> (1.3 in <sup>3</sup> )
Battery	2.7 Amp hours, 3.7V lithium-thionyl chloride cell
Storage temperature	-18° to +52°C -0° to + 125°F
Serial Number <sup>b</sup> : Radiopaque Identification (ID)	NAT

<sup>a</sup> All measurements are approximate.

<sup>b</sup> The serial number is the radiopaque ID followed by a number. The clinician programmer displays the entire serial number beginning with the radiopaque ID.

**Table 3. Material of components in the Model 7425 package.**

Component	Material	Material contacts human tissue
Neurostimulator		
Case	Titanium	Yes
Connector block	Polyurethane	Yes
Grommets, seals	Silicone rubber	Yes
Setscrews	Titanium	Yes
Insulation coating	Fluoropolymer	Yes
Adhesive	Silicone adhesive	Yes
Hex wrench		
Handle	Acetal resin	Yes
Shaft	Stainless steel	Yes

## Declaration of Conformity

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment and Directive 90/385/EEC on Active Implantable Medical Devices.

For additional information, contact the appropriate Medtronic representative listed on the inside back cover of this manual.

## Instructions for use

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

### △ Cautions:

- When using sharp instruments near the neurostimulator, be extremely careful to avoid nicking or damaging the case, the insulation, or the connector block. Damaging the neurostimulator may require surgical replacement.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

## Verifying neurostimulator operation

Before opening the sterile neurostimulator package, use the clinician programmer to interrogate the neurostimulator and verify neurostimulator battery status and current settings.

△ **Caution:** Do not implant a neurostimulator if it was dropped onto a hard surface from a height of 30 cm (12 in) or more, because the neurostimulator may be damaged and fail to operate properly.

**Note:** The neurostimulator pocket may be flushed with an antibiotic solution; do not submerge the neurostimulator in fluid.


### Connecting the extension to the neurostimulator

△ **Caution:** Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Wipe the extension connector pins with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution.
2. Make sure the connector block receptacles are dry and clean.
3. Insert the appropriate extension connector pins into the appropriate neurostimulator socket until they are seated fully within the connector block (Figure 1).

**Note:** To retract the setscrews, insert the hex wrench into the self-sealing grommet and rotate the setscrews counterclockwise; however, do not remove the setscrews from the connector block.

△ **Caution:** Do not insert the extension connector into the connector block if the setscrews are not sufficiently retracted. If the setscrews are not retracted, the extension connector pins may damage the setscrews and the extension connector pins will not be seated fully into the connector block.


**Figure 1.** Insert the extension connector pins fully into the neurostimulator.

4. Fully insert the hex wrench (packaged with the neurostimulator) into the self-sealing grommet of the connector block and tighten each setscrew.

**△ Cautions:**

- To avoid undertightening, do not use a torque wrench. Undertightening the neurostimulator setscrews may result in insufficient electrical contact within the connector block, resulting in intermittent or loss of stimulation.
- Be sure the hex wrench is fully inserted into the self-sealing grommet. If the hex wrench is not fully inserted, the setscrew may be damaged, resulting in intermittent or loss of stimulation.
- Before tightening setscrews, ensure that the extension connector pins are inserted into the connector block to prevent damaging the connector block.
- When using the hex wrench, do not overtighten the setscrews. Overtightening may damage the setscrews and prevent setscrew removal.
- Verify that each leaf of the self-sealing grommet is closed after the hex wrench is withdrawn. If fluid leaks through a grommet seal that is not fully closed, the patient may experience shocking, burning, or irritation at the neurostimulator implant location, or intermittent stimulation or loss of stimulation may occur.

**Note:** The hex wrench must be oriented to the same angle as the setscrew (Figure 2).


**Figure 2.** Tightening the setscrews in the self-sealing grommet.

### Implanting the neurostimulator

1. Place the neurostimulator into the subcutaneous pocket with the etched identification side placed outward, away from muscle tissue, and ensure that the extension is not bent sharply.

△ **Cautions:**

- Ensure that the neurostimulator is placed no deeper than 4 cm (1.5 in) below the skin and is parallel to the skin. If the neurostimulator is too deep or is not parallel to the skin, telemetry may be unsuccessful.
- Do not place the etched identification side of the neurostimulator facing inward. Placing the etched side inward could increase the possibility of skeletomuscular stimulation, which the patient may perceive as twitching or burning.
- Do not coil excess extension in front of the etched identification side of the neurostimulator. Wrap excess extension around the perimeter of the neurostimulator (Figure 3) to avoid increasing subcutaneous pocket depth, help minimize potential damage during neurostimulator replacement surgery, help minimize potential kinking of the extension, and minimize interference with telemetry during programming.

**Figure 3.** Wrap excess extension around the perimeter of the neurostimulator.

2. Use the suture holes in the connector block to secure the neurostimulator to the muscle fascia with nonabsorbable silk.

### Checking system integrity

1. Use the clinician programmer to confirm the integrity of the connected system.

△ **Caution:** To use the nonsterile clinician programmer in a sterile field, place a sterile barrier between the patient and the programming head to prevent infection. Do not sterilize any part of the clinician programmer. Sterilization may damage the programmer.


**Note:** The neurostimulator should be in the pocket during system interrogation for integrity to ensure proper readings.

2. Program the basic stimulation parameters and check the battery status; if applicable, check the electrode impedances to rule out a short or open circuit.
3. Complete the stimulation assessment form.



### **Completing the implant procedure**

1. Close and dress all incisions.
2. Ensure that a patient control device is given to the patient.
3. Complete the device tracking and patient registration paperwork and return the documents to Medtronic.


7425IPG\_ch.fm 4/27/04 2:22 pm  
UC200xxxxxx EN  
4.625 x 6 inches (117 mm x 152 mm)

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All Rights Reserved
Medtronic International Ltd.
© Medtronic, Inc. 2004