WALLSTENT™ RX BILIARY

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RÉFÉRENCES

1. Adam, Chetty N, Roddie M, Yeung E, Benjamin IS. Self-expandable stainless steel endoprosthesis for treatment of malignant bile duct obstruction. *AJR* 1991; 156:321-325.
2. Dick R, Gilliams A, Dooley JS, Hobbs KEF. Stainless steel mesh stents for biliary strictures. *Journal of Interventional Radiology* 1989; 4:95-98.
3. Gilliams A, Dick R, Dooley JS, Wallsten H, El-Din A. Self-expandable Stainless steel braided endoprosthesis for biliary strictures. *Radiology* 1990; 174:137-140.
4. LaBerge JM, Doherty M, Gordon RL, Ring EJ. Hilar malignancy: treatment with an expandable metallic transhepatic biliary stent. *Radiology* 1990; 177: 789-792.
5. Lammer J, Klein GE, Kleinert R, Hausegger K, Einspieler R. Obstructive jaundice: use of expandable metal endoprosthesis for biliary drainage. *Radiology* 1990: 789-792.
6. Neuhaus H, Hagemuller F, Griebel M, Classen M. Percutaneous cholangioscopic or transpapillary insertion of self-expanding biliary metal stents. *Gastrointestinal Endoscopy* 1991: 37:31-37.
7. Born P, et al. Initial experience with a new, partially covered wallstent for malignant biliary obstruction. *Endoscopy* 1996; 28:699-702.
8. Prat F, et al. A randomized trial of endoscopic drainage methods for inoperable malignant strictures of the common bile duct. *Gastrointestinal Endoscopy* 1998; 47(1):1-7.
9. Davids P, Groen AK, Rauws E, Tytgat G, Huibregste K. Randomized trial of self-expanding metal stents versus polyethylene stents for distal malignant biliary obstruction. *The Lancet* 1992; 340(Dec 19/26):1488-1492.
10. Schmassman A, et al. WALLSTENTS versus plastic stents in malignant biliary obstruction: effects of stent patency of the first and second stent on patient compliance and survival. *The American Journal of Gastroenterology* 1996; 91(4): 654-659.
11. Yeoh et al., "Comparative costs of metal vs. plastic biliary stent strategies for malignant obstructive jaundice by decision analysis" *Gastrointestinal Endoscopy* 1999; 49(4): Pt.1: 466-477.

Catalog Number
 Número de catálogo
 Numéro de catalogue
 Bestell-Nr.
 Numero di catalogo
 Catalogusnummer
 Referência

Consult instructions for use.
 Consultar las instrucciones de uso.
 Consulter le mode d'emploi.
 Gebrauchsanweisung beachten.
 Consultare le istruzioni per l'uso.
 Raadpleeg instructies voor gebruik.
 Consulte as Instruções de Utilização

Contents
 Contenido
 Contenu
 Inhalt
 Contenuto
 Inhoud
 Conteúdo

EU Authorized Representative
 Representante autorizado en la UE
 Représentant agréé UE
 Autorisierter Vertreter in der EU
 Rappresentante autorizzato per l'UE
 Erkend vertegenwoordiger in EU
 Representante Autorizado na U.E.

Legal Manufacturer
 Fabricante legal
 Fabricant légal
 Berechtigter Hersteller
 Fabbricante legale
 Wettelijke fabrikant
 Fabricante Legal

Lot
 Lote
 Lot
 Charge
 Lotto
 Partij
 Lote

Recyclable Package
 Envase reciclable
 Emballage recyclable
 Wiederverwertbare Verpackung
 Confezione riciclabile
 Recyclebare verpakking
 Embalagem Reciclável

Use By
 Fecha de caducidad
 Date limite d'utilisation
 Verwendbar bis
 Usare entro
 Uiterste gebruiksdatum
 Validade

Australian Sponsor Address
 Dirección del patrocinador australiano
 Adresse du promoteur australien
 Adresse des australischen Sponsors
 Indirizzo sponsor australiano
 Adres Australische sponsor
 Endereço do Patrocinador Australiano

Argentina Local Contact
 Contacto local en Argentina
 Contact local en Argentine
 Lokaler Kontakt Argentinien
 Contatto locale per l'Argentina
 Contactpersoon Argentinië
 Contacto local na Argentina

Brazil Local Contact
 Contacto local en Brasil
 Contact local au Brésil
 Lokaler Kontakt Brasilien
 Contatto locale per il Brasile
 Contactpersoon Brazilië
 Contacto local no Brasil

For single use only. Do not reuse.
 Para un solo uso. No reutilizar.
 À usage unique. Ne pas réutiliser.
 Für den einmaligen Gebrauch. Nicht wieder verwenden.
 Exclusivamente monouso. Non riutilizzare.
 Uitsluitend bestemd voor eenmalig gebruik. Niet opnieuw gebruiken.
 Apenas para uma única utilização. Não reutilize.

Do Not Resterilize
 No reesterilizar
 Ne pas restériliser
 Nicht erneut sterilisieren
 Non risterilizzare
 Niet opnieuw steriliseren
 Não reesterilize

Do not use if package is damaged.
 No usar si el envase está dañado.
 Ne pas utiliser si l'emballage est endommagé.
 Bei beschädigter Verpackung nicht verwenden.
 Non usare il prodotto se la confezione è danneggiata.
 Niet gebruiken als de verpakking is beschadigd.
 Não utilize se a embalagem estiver danificada.

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EC **REP** **EU Authorized
Representative**

Boston Scientific Limited
Ballybrit Business Park
Galway
IRELAND

AUS **Australian
Sponsor Address**

Boston Scientific (Australia) Pty Ltd
PO Box 332
BOTANY
NSW 1455
Australia
Free Phone 1800 676 133
Free Fax 1800 836 666

 **Legal
Manufacturer**

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
USA
USA Customer Service 888-272-1001

 **Do not use if package
is damaged.**

 **Recyclable
Package**

C € 0344

© 2017 Boston Scientific Corporation or its affiliates.
All rights reserved.