# RX Accunet

## Embolic Protection System

CE  
0086

EN	English	HU	Magyar
DE	Deutsch	FR	Français
SV	Svenska	ES	Español
DA	Dansk	SK	Slovensky
PL	Polski	PT	Português
CS	Česky	NL	Nederlands
TR	Türkçe	EL	Ελληνικά
IT	Italiano		

Graphical Symbols for Medical Device Labeling  
 Grafische Symbole zur Beschriftung medizinischer Geräte  
 Symboles graphiques pour l'étiquetage de matériel médical  
 Símbolos gráficos para el etiquetado de productos sanitarios  
 Simboli grafici per l'etichettatura di dispositivi medicali  
 Símbolos Gráficos para Rotulagem de Dispositivos Médicos  
 Grafiska symboler för märkning av medicinsk utrustning  
 Grafische symbolen voor de etikettering van medische hulpmiddelen  
 Grafiske symboler beregnet til mærkning af medicinske produkter  
 Γραφικά σύμβολα για τη σήμανση ιατρικών συσκευών  
 Symbole graficzne do oznaczania sprzętu medycznego  
 Grafikus jelzések orvosi eszközök címkézéséhez  
 Grafické symboly k označení lékařského zařízení  
 Tibbi Cihaz Etiketleri için Grafik Semboller  
 Vysvetlivky ku grafickým symbolom na štítku balenia

**RX AccUNET  
Embolic Protection System**

**CAUTION**

**CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.**

Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects and hazards commonly associated with carotid interventional procedures should use this device.

Refer to the instructions supplied with all interventional devices to be used in conjunction with the RX AccUNET EPS for their intended uses, contraindications and potential complications.

**DEVICE DESCRIPTION**

The RX AccUNET Embolic Protection System (EPS) is comprised of a delivery system and two recovery catheters with different features to address varying patient anatomies and clinical challenges.

**RX AccUNET Delivery System:**

The RX AccUNET EPS guide wire with filter basket is a single use, steerable 0.014" (0.36 mm) guide wire system that comes in 190 cm and 300 cm lengths. Attached to the distal end is a filter basket made of a thin porous membrane supported by a metal basket (Figure 1). The purpose of the filter basket is to capture and recover embolic material produced during interventional procedures. The RX AccUNET EPS guide wire with filter basket is delivered to the target site while contained within the delivery sheath. Once the filter basket is positioned at the target site, the delivery sheath is retracted by peeling it away from the guide wire to allow for full expansion of the filter basket. Radiopaque markers are located on the guide wire, filter basket and delivery sheath to aid in device positioning. Once the filter basket has been deployed, the wire may be used as a standard guide wire to support the delivery and placement of endovascular devices such as balloon dilatation catheters and endovascular stents.


FIGURE 1 - Guide Wire and Filter Basket

**Filter Sizing Recommendation**

Filter size selection is based on reference vessel diameter in the area where the filter basket will be expanded. Measurements of reference vessel diameters may be acquired through:

- Angiogram (Fluoroscopy)
- Quantitative Coronary Angiography (QCA)

Filter Size Fully Expanded	Reference Vessel Diameter Minimum to Maximum Range
4.5 mm	3.25 to 4.0 mm
5.5 mm	4.0 to 5.0 mm
6.5 mm	5.0 to 6.0 mm
7.5 mm	6.0 to 7.0 mm

**RX AccUNET Recovery Catheter - "Shapeable Tip Design"**

The RX AccUNET Recovery Catheter is a rapid exchange (RX) catheter used to remove the RX AccUNET EPS guide wire with filter basket from the vasculature and guiding catheter / sheath while retaining any emboli or particulates collected in the filter basket during the procedure. The RX AccUNET Recovery Catheter has a shapeable tip that is designed to provide torqueability or steerability to the tip of the device. The RX AccUNET Recovery Catheter is available in four sizes and is compatible with the matching four sizes of RX AccUNET Filter Baskets.

**RX AccUNET 2 Recovery Catheter - "Low-Profile, Flexible Design"**

The RX AccUNET 2 Recovery Catheter is a rapid exchange (RX) catheter used to remove the RX AccUNET Embolic Protection System (EPS) guide wire with filter basket from the vasculature and guiding catheter / sheath while retaining any emboli or particulates collected in the filter basket during the procedure. The RX AccUNET 2 Recovery Catheter has a lower profile and is more flexible than the RX AccUNET Recovery Catheter. In addition, the RX AccUNET 2 Recovery Catheter has an extended soft tip and was designed to deflect away from stent struts without the need to shape the distal tip. One size of RX AccUNET 2 Recovery Catheter is used to capture all sizes of RX AccUNET Filter Baskets.

**HOW SUPPLIED**

**Sterile:** Sterilized with electron beam radiation. Non-pyrogenic.

This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

**Contents:** Three pouches containing:

**RX AccUNET Delivery System:** One (1) 0.014" (0.36 mm) guide wire with filter basket, One (1) delivery sheath, One (1) introducer tool, One (1) torque device with peel away adapter, One (1) flushing tool.

**RX AccUNET Recovery Catheter - "Shapeable Tip Design":** One (1) recovery catheter, One (1) catheter clip.

**RX AccUNET 2 Recovery Catheter - "Low-Profile, Flexible Design":** One (1) recovery catheter, One (1) catheter clip.

**Storage:** Store in a dry, dark, cool place.

**INDICATIONS**

The RX AccUNET EPS is indicated for use in patients undergoing percutaneous interventional procedures in carotid arteries to:

- facilitate the placement of diagnostic and therapeutic devices during percutaneous interventional procedures
- capture and remove embolic material that may be released during the procedure

The RX AccUNET EPS is contraindicated for use in:

- Patients in whom anticoagulant and / or antiplatelet therapy is contraindicated.
- Patients with severe vascular tortuosity or anatomy that would preclude the safe introduction of a guiding catheter, sheath, embolic protection system, or stent system.
- Patients with known hypersensitivity to nickel-titanium.
- Patients with uncorrected bleeding disorders.

**WARNINGS**

- The safety and effectiveness of this device as an embolic protection system have not been established in vasculatures outside the carotid arteries.
- The safety and efficacy of the RX AccUNET EPS have not been demonstrated with carotid stent systems other than the over-the-wire or RX AccuLink Carotid Stent System.
- Appropriate antiplatelet and anticoagulation therapy should be administered pre- and post-procedure. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.
- This device is intended for single-use only. Do not reuse. Do not sterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.
- Do not use the product after the "Use By" date specified on the package.
- Maintain the patient's Activated Clotting Time (ACT) > 250 seconds throughout RX AccUNET EPS usage to prevent thrombus formation on the device.
- During positioning, advance the RX AccUNET EPS guide wire with filter basket and the delivery sheath as a unit. Advancing the guide wire independent of the delivery sheath may result in premature filter deployment.
- Always advance or withdraw the RX AccUNET EPS guide wire slowly using fluoroscopy to observe corresponding wire movement and tip movement in particular. Never push, tug, retract or torque a guide wire that meets resistance. If the wire tip becomes entrapped within a lesion or a device, such as a deployed stent, do NOT torque the guide wire. Determine the cause of resistance and take necessary remedial action. Torquing or retracting the guide wire against resistance may damage the wire, cause wire tip separation, or cause vessel trauma. Resistance may be felt and / or observed using fluoroscopy by noting any buckling of the guide wire tip.
- Maintain continuous flush while removing and reinserting devices on the guide wire. Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- When introducing the delivery system, confirm that the wire tip is free within the vessel lumen and is not directed into the vessel wall. Failure to do so may result in vessel trauma. Use the radiopaque marker on the interventional device to confirm position.
- Avoid excessive movement of the RX AccUNET EPS guide wire and filter basket during catheter device exchanges. Excessive movement of the deployed basket may cause vessel trauma or spasm.
- The filter basket must be kept distal to the area to be covered by the stent and proximal to the petrous portion of the carotid artery to avoid excessive forces on the filter basket.
- Allow for and maintain adequate distance between the radiopaque proximal bushing marker on the guide wire with filter basket and the stent delivery system or other compatible interventional devices to avoid potential entanglement.
- Always keep the open filter basket distal to the deployed stent. Do not attempt to pull an open filter basket through the stent. Do not attempt to capture the filter basket by pulling it into the recovery catheter if the recovery catheter tip is in the stent area. Pulling the filter basket into the stent area may lead to stent-filter basket entanglement and / or basket detachment. If filter basket entanglement or detachment occurs, surgical conversion or collapsing the basket with a second stent should be considered.
- Maintain proper guiding catheter / sheath support in the common carotid artery throughout the procedure. If guiding catheter / sheath access cannot be maintained, the case should be discontinued. Failure to maintain proper support of the guiding catheter / sheath can lead to prolapse of the catheter into the aortic arch, resulting in any of the following:
  - Movement of an open filter through an undilated lesion; or
  - Filter-stent entanglement, filter basket detachment and / or proximal movement of the stent; or
  - Filter guide wire breakage.
- Do not rotate either recovery catheter more than 90 degrees in either direction since this can result in the guide wire wrapping around the catheter.
- Use with fixed (passive) hemostatic valves is not recommended.
- If excessive debris is collected in the filter basket such that distal perfusion of dye is significantly reduced or no dye is perfusing past the filter, the RX AccUNET EPS may have reached its maximum capacity to contain emboli. Remove and replace the RX AccUNET EPS. Otherwise it may be difficult to completely recover all embolic debris and the potential for thrombus formation may increase.
- Discard unused RX AccUNET and RX AccUNET 2 Recovery Catheters after completing procedure. Failure to discard unused recovery catheters can result in any or all of the following: the use of a device past its "Use By" date; or the use of an incorrectly sized RX AccUNET Recovery Catheter, which can cause loss of particulates from the filter basket during recovery, filter-stent entanglement, filter basket detachment and / or proximal movement of stent.

**PRECAUTIONS**

- Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and / or inaccurate torque response.
- Do not expose the delivery or recovery systems to organic solvents (e.g. alcohol) as structural integrity and / or function of the device may be impaired.
- Confirm the compatibility of the RX AccUNET EPS with the interventional devices before actual use.
- To avoid air entrapment, do not leave the introducer tool in the valve while advancing the delivery system. Advance the delivery system slowly through the guiding catheter / sheath.
- Ensure that the torque device is secured tightly to the guide wire to avoid difficulties with deployment. If the torque device is not tightened securely on the guide wire, proper deployment may not occur.
- For proper positioning of the filter basket, the vessel distal to the lesion should have an absence of excessive tortuosity and be of adequate length (approximately 4 cm distal to the lesion and proximal to the petrous portion of the vessel).
- Reconfirm filter basket position prior to deployment of the stent to ensure that there is adequate distance between the radiopaque proximal bushing marker on the guide wire with filter basket and the distal end of the desired stent position. Reposition the filter basket as necessary by gently advancing the guide wire. Under fluoroscopy, observe movement of the filter basket to the desired position.
- Observe all guide wire movement in the vessels using fluoroscopic guidance.
- Venous access should be available during carotid stenting in order to manage bradycardia and / or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.
- The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.
- Do NOT deliver the RX AccUNET EPS guide wire with filter basket through any interventional devices other than the RX AccUNET EPS delivery sheath.
- Removal of the RX AccUNET EPS guide wire with filter basket through any interventional devices other than the RX AccUNET and the RX AccUNET 2 Recovery Catheters has not been tested.
- Care must be used when removing the filter basket through a newly deployed stent to maintain filter basket integrity and to avoid disrupting the stent geometry.
- If the RX AccUNET EPS is desired for intervention in subsequent vessels, use a new device.

**POTENTIAL ADVERSE EVENTS**

Based on the literature, and on clinical and commercial experience with the use of embolic protection systems with a carotid stent, the following alphabetical list includes possible adverse events associated with the use of this device:

• Allergic reactions to antiplatelet agents / contrast medium	• Groin hematoma, with or without surgical repair
• Aneurysm	• Hemorrhage, with or without transfusion
• Angina / coronary ischemia	• Hypotension / hypertension
• Arrhythmia	• Infection and pain at insertion site
• Arterial occlusion / thrombosis at puncture site or remote site	• Ischemia / infarction of tissue / organ
• Arteriovenous fistula	• Myocardial infarction (MI)
• Bacteremia or septicemia	• Pain (head, neck)
• Bleeding from anticoagulant or antiplatelet medications	• Pseudoaneurysm (femoral)
• Cerebral edema	• Renal failure / insufficiency
• Cerebral hemorrhage	• Restenosis of stented / dilated vessel
• Cerebral ischemia / transient ischemic attack (TIA)	• Seizure
• Congestive heart failure (CHF)	• Severe unilateral headache
• Death	• Stent / filter entanglement / damage
• Detachment and / or implantation of a component of the system	• Stroke / cerebrovascular accident (CVA)
• Emboli, distal (air, tissue, or thrombotic emboli)	• Total occlusion of carotid artery
• Emergent or urgent endarterectomy surgery (CEA)	• Vessel dissection, perforation, or rupture
• Fever	• Vessel Spasm or recoil
• Filter thrombosis / occlusion	

Any adverse event occurring involving the RX AccUNET EPS should be reported immediately to Abbott Laboratories, Customer Service, at (951) 914-4669.

CLINICIAN USE INFORMATION

Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects and hazards commonly associated with carotid interventional procedures should use this device.

**WARNING: Do not use after the "Use By" date specified on the package.**

**WARNING: The device is intended for single-use only. Do not reuse. Do not resterilize as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.**

**WARNING: Use with fixed (passive) hemostatic valves is not recommended.**

**Materials Required**

- 6F Guiding Sheath or 8F Guiding Catheter (Min. Guiding Catheter / Sheath Size (I.D.) 0.085" / 2.2 mm).
- 10 - 20 cc luer-lock syringe
- Heparinized normal saline (sterile)
- Balloon dilatation catheter (optional)
- ≥ 0.096" (2.44 mm) Rotating Hemostatic Valve (RHV) (optional)
- Acculink (over-the-wire) or RX Acculink Carotid Stent System (optional)
- Doc extension wire (optional)

**CAUTION: Confirm the compatibility of the RX Accunet EPS with the interventional devices before actual use.**

**Filter Basket Size Determination**

Filter size selection is based on reference vessel diameter in the area where the filter basket will be expanded. See Table below for filter basket diameters and recommended reference vessel diameters.

Filter Size Fully Expanded	Reference Vessel Diameter Minimum to Maximum Range
4.5 mm	3.25 to 4.0 mm
5.5 mm	4.0 to 5.0 mm
6.5 mm	5.0 to 6.0 mm
7.5 mm	6.0 to 7.0 mm

**WARNING: Overstretching of the artery may result in rupture and life-threatening bleeding.**

**Inspection Prior to Use**

Inspect the product prior to use. Do not use if the package is open or damaged.

**CAUTION: Guide wires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guide wire carefully for bends, kinks, or other damage. Do not use damaged guide wires. Using a damaged guide wire may result in vessel damage and / or inaccurate torque response.**

**Delivery System Preparation**

**CAUTION: Do not expose the delivery or recovery systems to organic solvents (e.g. alcohol) as structural integrity and / or function of the device may be impaired.**

1. Hold the dispenser hoop and flushing tool securely in one hand.
2. Loosen the RHV on the flush tool.
3. Fill a 10 cc syringe with heparinized normal saline. Holding the flush tool in an upright position (see Figure 2), attach the syringe as shown, and flush gently to remove air from the filter. Observe fluid exiting the RHV.
4. Confirm that no part of the light blue, distal section of the delivery sheath is within the RHV, then tighten the RHV onto the delivery sheath.
5. Release the flushing tool and delivery sheath from the dispenser hooks. Release the proximal end of the wire from its dispenser hook. Grasp the delivery sheath and the flush tool securely together, and remove from the dispenser. Inspect the delivery system for kinks, or other damage.
6. Flush to remove air from the sheath. Observe fluid exiting the proximal end of the delivery sheath.
7. Place the torque device onto the proximal end of the RX Accunet EPS guide wire and tighten.
8. While observing through the flushing tool, pull the filter basket back into the delivery sheath (by pulling proximally on the torque device) until the distal tip of the delivery sheath aligns with the proximal end of the blue filter obturator (see Figures 3 and 4). DO NOT pull the torque device such that the proximal end of the blue filter obturator goes into the distal tip of the delivery sheath. Doing so will cause a reduction in device flexibility and deliverability, and has the potential to damage (crumple) the delivery sheath and adversely affect filter basket delivery.
 

**Note:** Step 8 may be repeated a maximum of two additional times.
9. Loosen the RHV of the flushing tool, and slide the flushing tool off the distal end of the delivery system.
10. Loosen the torque device and position it so that the light blue, proximal end of the delivery sheath is within the central collet tube of the peel away adapter of the torque device. Tighten the torque device to the wire.



Figure 4 Align distal tip of delivery sheath with proximal end of obturator (blue)

**Recovery Catheter Preparation (for both the RX Accunet Recovery Catheter and the RX Accunet 2 Recovery Catheter)**

1. Choose one of the two recovery catheters provided, based on the relevant features of the device and how those features may best address the conditions of the case (e.g., CCA / ICA tortuosity, protruding stent struts, difficulty of aortic access). The RX Accunet Recovery Catheter – "Shapeable Tip Design" has a shapeable tip that is designed to provide torqueability or steerability at the tip of the device. The RX Accunet 2 Recovery Catheter – "Low-Profile, Flexible Design" is more flexible than the RX Accunet Recovery Catheter and has an extended soft tip and lower profile designed to deflect away from stent struts.
2. Remove the selected recovery catheter from the dispenser hoop.
3. Fill a 10 cc syringe with heparinized normal saline.
4. Flush the selected recovery catheter to remove air using one of two methods:
  - Align the tip of the syringe with the distal tip of the recovery catheter and flush to remove air, or
  - Use the flushing tool from the RX Accunet EPS delivery system: Attach a luer lock syringe to the distal end of the flushing tool, and insert the distal 3 cm of the recovery catheter into the proximal end of the flushing tool. Lightly close the RHV of the flushing tool, and flush gently to remove air from the recovery catheter. Do not close the RHV on the clear tip area of the recovery catheter. Observe fluid exiting the exit notch.
5. If using the RX Accunet Recovery Catheter – "Shapeable-Tip Design", the distal 3 cm of the catheter may be shaped as needed to allow for tip deflection.
 

**General guidelines for shaping the RX Accunet Recovery Catheter:**

  - Keep the shaping ribbon on the outside of the curve
  - Place the bend at the proximal part of the shaping ribbon. This section is the area where the orange catheter body transitions to the clear tip area.

**Guiding Catheter / Sheath Support**

Introduce guiding catheter / sheath into vasculature and confirm that it is well positioned and firmly seated prior to introduction of the RX Accunet EPS. Always maintain view of the tip of the guiding catheter / sheath to minimize the likelihood of the catheter backing out or prolapsing into the aortic arch.

**WARNING: Maintain proper guiding catheter / sheath support in the common carotid artery throughout the procedure. If guiding catheter / sheath access cannot be maintained, the case should be discontinued. Failure to maintain proper support of the guiding catheter / sheath can lead to prolapse of the catheter into the aortic arch, resulting in any of the following:**

- Movement of an open filter through an undilated lesion;
- Filter-stent entanglement, filter basket detachment and / or proximal movement of the stent; or
- Filter guide wire breakage

**Lesion Preparation**

**WARNING: Maintain the patient's ACT at > 250 seconds throughout RX Accunet EPS usage to prevent thrombus formation on the device.**

**CAUTION: Venous access should be available during carotid stenting in order to manage bradycardia and / or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.**

If desired, a 2 mm balloon may be used to predilate the lesion.

**Delivery Procedure**

**CAUTION: Observe all guide wire movement in the vessels using fluoroscopic guidance.**

**CAUTION: The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.**

**CAUTION: Do not deliver the RX Accunet EPS guide wire with filter basket through any interventional devices other than the RX Accunet EPS delivery system.**

**WARNING: Use with fixed (passive) hemostatic valves is not recommended.**

1. Use standard percutaneous techniques. Shape the tip of the RX Accunet EPS guide wire, if desired.
2. Carefully insert the guide wire tip of the delivery system into the funnel of the introducer tool and advance the unit completely into the RHV until the delivery sheath has entered the guiding catheter / sheath. Remove the introducer tool.

**CAUTION: To avoid air entrapment, do not leave the introducer tool in the valve while advancing the delivery system. Advance the delivery system slowly through the guiding catheter / sheath.**

**Note:** If using a long sheath with a hemostasis valve, the introducer tool should not be peeled off. Rather, it should be pulled back onto the sheath so that it can be reinserted later to facilitate deployment.

Under fluoroscopy, advance the tip of the delivery system into the selected vessel and across the lesion using the torque device to steer the guide wire. Do not torque the delivery sheath.

**WARNING: During positioning, advance the RX AccUNET EPS guide wire with filter basket and the delivery sheath as a unit. Advancing the guide wire independently from the delivery sheath may result in premature filter deployment.**

**WARNING: When introducing the delivery system, confirm that the wire tip is free within the vessel lumen and is not directed against the vessel wall. Failure to do so may result in vessel trauma. Use the radiopaque marker on the interventional device to confirm position.**

**WARNING: Always advance or withdraw the RX AccUNET EPS guide wire slowly using fluoroscopy to observe corresponding wire movement and tip movement in particular. Never push, auger, retract or torque a guide wire that meets resistance. If the wire tip becomes entrapped within a lesion or a device, such as a deployed stent, do NOT torque the guide wire. Determine the cause of resistance and take necessary remedial action. Torquing or retracting the guide wire against resistance may damage the wire, cause wire tip separation, or may cause vessel trauma. Resistance may be felt and / or observed using fluoroscopy by noting any buckling of the guide wire tip.**

**Note:** A variety of techniques can be used to assist passage of the delivery system if it has difficulty advancing through the lesion. Some options to aid passage of the delivery system are:


- Have the patient rotate her / his neck from side-to-side. This motion may re-orient the carotid artery.
- If the RX AccUNET EPS delivery system cannot cross the lesion, a 2 mm balloon may be used to predilate the lesion.
- If desired, insert a stiff 0.014" guide wire (buddy wire) to straighten the carotid vasculature to aid advancement of the delivery system. Leave the buddy wire in place while advancing the RX AccUNET EPS delivery system. Remove the buddy wire prior to deployment of the RX AccUNET EPS filter basket.

**Filter Basket Positioning and Deployment**

1. Position the delivery system in the artery such that the radiopaque markers are distal to the lesion to be treated.

**CAUTION: For proper positioning of the filter basket, the vessel distal to the lesion should have an absence of excessive tortuosity and be of adequate length (approximately 4 cm distal to the lesion and proximal to the petrous portion of the vessel).**

2. Loosen the torque device from the wire.
3. Feed the proximal, light blue end of the delivery sheath into the adjacent port of the peel away adapter on the torque device (see Figure 5). Advance the torque device until it abuts against the darker blue part of the sheath where the slit in the sheath begins. Tighten down the torque device to the wire.

**Figure 5 –**  
Proximal End of Delivery Sheath Inserted in Peel Away Adapter

4. Remove all slack from the system.

**Note:** If using a long sheath with a hemostasis valve or if using the Abbott Vascular CoPilot Bleedback Control Valve, advance and reposition the introducer tool within the valve to hold it in an open position. This will allow the delivery sheath to move freely within the valve.

5. Stabilize the torque device with one hand (preferably the left) and keep that hand resting over the delivery sheath, distal to the torque device. With the other hand (the right), grab the proximal, light blue end of the delivery sheath that has been fed through the peel away adapter.
6. Under fluoroscopy pull the end of the sheath approximately 3 cm in a horizontal direction away from the patient. The sheath will peel away as it is pulled through the adapter. Observe filter basket expansion.

**CAUTION: Ensure that the torque device is secured tightly to the guide wire to avoid difficulties with deployment. If the torque device is not tightened securely on the guide wire, proper deployment may not occur.**

**Note:** Avoid excessive movement of the filter basket during unsheathing by removing all slack between the patient, guiding catheter / sheath and delivery system.

**Note:** Slight advancement of the guide wire during the unsheathing of the filter basket may be required to avoid excess movement of the filter basket during unsheathing.

7. The filter basket should now be fully expanded within the vessel. While maintaining a stable guide wire position, continue to peel back and retract the delivery sheath from the wire using the peel away feature. When the light blue, distal section exits the guiding catheter, close the RHV. The light blue, distal section will not peel. Once it exits the RHV, remove the torque device and the remaining portion of the delivery sheath.

**Note:** If using a long sheath with a hemostasis valve, remove the introducer tool after the light blue, distal section of the delivery sheath exits the hemostasis valve.

8. Confirm filter basket expansion by checking in two projections for visibility of the radiopaque markers on the frame of the basket and their apposition against the vessel wall. Under fluoroscopy, perform a contrast injection to determine that flow is adequate distal to the filter basket and that the filter basket is in the desired position.

**WARNING: The filter basket must be kept distal to the area to be covered by the stent and proximal to the petrous portion of the carotid artery to avoid excessive forces on the filter basket.**

9. Compatible 0.014" guide wire-based interventional devices may now be backloaded onto the guide wire for treatment of the lesion.

**CAUTION: Reconfirm the filter basket position prior to deployment of the stent to ensure that there is adequate distance between the radiopaque proximal bushing marker on the guide wire with filter basket and the distal end of the desired stent position to avoid stent deployment or entanglement on any part of the filter basket. If it is necessary to reposition the filter basket further distal to the desired stent position, gently advance the guide wire with filter basket. Under fluoroscopy, observe movement of the filter basket to the desired position.**

**Monitoring the Filter Status**

**WARNING: Maintain continuous flush while removing and reinserting devices on the guide wire. Perform all exchanges slowly to prevent air embolism or trauma to the artery.**

**WARNING: Avoid excessive movement of the RX AccUNET EPS guide wire and filter basket during catheter device exchanges. Excessive movement of the deployed basket may cause vessel trauma or spasm.**

1. Keep the position of the guide wire with filter basket steady during the intervention.
2. Check the status of the expanded RX AccUNET EPS filter basket at regular intervals during the interventional procedure.
3. During the case, keep the following three areas in the field of view at all times:
  - Filter basket and / or the proximal bushing marker,
  - Lesion area, and
  - Radiopaque tip of the guiding catheter / sheath. Maintaining a view of the tip of the guiding catheter / sheath will minimize the likelihood of the catheter backing out or prolapsing into the aortic arch.

**WARNING: Prolapse of the guiding catheter / sheath into the aortic arch can result in any of the following:**

- Movement of an open filter through an undilated lesion;
- Filter-stent entanglement, filter basket detachment and / or proximal movement of the stent; or
- Filter guide wire breakage

Inject contrast through the guiding catheter / sheath and observe flow distal to the filter basket.

**WARNING: Allow for and maintain adequate distance between the radiopaque proximal bushing marker on the guide wire with filter basket and the stent delivery system or other compatible interventional devices to avoid potential entanglement.**

**WARNING: If excessive debris is collected in the filter basket such that distal perfusion of dye is significantly reduced or no dye is perfusing past the filter, the RX AccUNET EPS may have reached its maximum capacity to contain emboli. Remove and replace the RX AccUNET EPS. Otherwise, it may be difficult to completely recover all embolic debris and the potential for thrombus release may increase.**

**Filter Basket Recovery (using the RX AccUNET Recovery Catheter or the RX AccUNET 2 Recovery Catheter)**

1. Remove all interventional devices from the guide wire.
2. Backload the selected, prepared recovery catheter onto the proximal end of the wire and advance the system through the open RHV on the guiding catheter / sheath.

**WARNING: Do not rotate the recovery catheter more than 90 degrees in either direction since this can result in the guide wire wrapping around the catheter.**

**WARNING: Use with fixed (passive) hemostatic valves is not recommended.**

**Note:** If using the RX AccUNET Recovery Catheter – "Shapeable Tip Design", the distal 3 cm of the RX AccUNET Recovery Catheter tip may be shaped to allow for tip deflection to maneuver around stent struts or in bends. Tip deflection is accomplished by rotating the proximal luer of the recovery catheter up to 90 degrees in either direction.

3. Under fluoroscopy, carefully advance the selected recovery catheter through the deployed stent.

**Note:** A variety of techniques can be used to assist passage of the recovery catheter if it has difficulty advancing through the deployed stent. These techniques are intended to adjust the bias of the guide wire. Some options to aid passage of the recovery catheter are:

- Have the patient rotate her / his neck from side-to-side. This motion may re-orient the carotid artery.
- Change the position of the guiding catheter or guiding sheath. The new position may either re-orient the entry of, or give better support to, the recovery catheter.
- If using the RX AccUNET Recovery Catheter – "Shapeable Tip Design", the tip shape of the recovery catheter can be modified. If one shape does not pass through the stent, shape the catheter tip in another direction or modify the degree of angulation.
- If stent struts are impeding the advancement of the recovery catheter, post-dilate the stent.
- Insert a guide wire ("buddy wire") to straighten the stented area.

4. If the above techniques are unsuccessful in advancing the recovery catheter through the deployed stent, the alternate recovery catheter should be prepared and used, after withdrawal of the first recovery catheter; i.e., if the RX AccUNET Recovery Catheter will not advance, carefully withdraw it and advance the RX AccUNET 2 Recovery Catheter, or vice versa.


5. Gently advance the selected recovery catheter over the filter basket until the radiopaque tip of the recovery catheter just covers the 4 (four) radiopaque wall apposition markers of the filter basket and causes the basket to collapse into the recovery catheter. The filter basket should not be completely contained within the catheter (see Figure 6).

**WARNING: Always keep the open filter basket distal to the deployed stent. Do not attempt to pull an open filter basket through the stent. Do not attempt to capture the filter basket by pulling it into the recovery catheter if the recovery catheter tip is in the stent area. Pulling the filter basket into the stent area may lead to stent-filter basket entanglement and / or basket detachment. If filter basket entanglement or detachment occurs, surgical conversion or collapsing the basket with a second stent should be considered.**

**Note:** If using the RX AccUNET Recovery Catheter, there is a stop (you will feel resistance) when the tip of the recovery catheter is advanced far enough over the radiopaque wall apposition markers. However, the RX AccUNET 2 Recovery Catheter does not have a stop within the tip of the device. This means that no resistance will be felt if the RX AccUNET 2 Recovery Catheter is advanced too far over the RX AccUNET Filter Basket. See Figure 6 for the proper positioning of both the RX AccUNET and the RX AccUNET 2 Recovery Catheters over the filter basket.

**Note:** If necessary to facilitate recovery of the filter basket, a simultaneous push-pull maneuver may be used. This may be done when the tip of the recovery catheter is completely through the stent and at, or beyond, the radiopaque proximal bushing.

**Note:** Clinical investigators have used other interventional devices, such as compatible guide or balloon catheters, to recover the filter basket in the instance when difficulty was encountered using the RX AccUNET or RX AccUNET 2 Recovery Catheters.

**Figure 6.**  
Recovering the Filter Basket

6. Hold tension on the guide wire and grasp the recovery catheter. Retract the devices together as a single unit with no movement relative to the catheter and guide wire, ensuring radiographically that the filter basket does not re-deploy.

**CAUTION: Care must be used when removing the filter basket through a newly deployed stent to maintain filter basket integrity and to avoid disrupting the stent geometry.**

7. Remove the devices as a unit from the RHV. Ensure that the RHV is fully open when the filter basket is being removed to maintain filter basket integrity.

**Note:** If unable to retract the filter basket through the guiding catheter, stabilize the filter basket and recovery catheter at the guiding catheter tip by tightening down the RHV. Remove the guide wire with filter basket, recovery catheter and guiding catheter as a unit.

**Note:** If using an Abbott Vascular CoPilot Bleedback Control Valve ensure that the CoPilot is held open while the filter basket is being removed to maintain filter basket integrity.

**WARNING: Discard unused RX AccUNET and RX AccUNET 2 Recovery Catheters after completing procedure. Failure to discard unused recovery catheters can result in any or all of the following: the use of a device past its "Use By" date; or the use of an incorrectly sized RX AccUNET Recovery Catheter, which can cause loss of particulates from the filter basket during recovery, filter-stent entanglement, filter basket detachment and / or proximal movement of stent.**

**CAUTION: If the RX AccUNET EPS is desired for intervention in additional vessels, use a new device.**

This product and / or its use are covered by one or more of the following United States Patents: 5,827,321; 6,106,548; 6,177,522; 6,245,090; 6,336,934; 6,402,771; 6,432,122; 6,475,236; 6,514,280; 6,540,722; 6,551,341; 6,565,591; 6,592,607; 6,645,224; 6,669,716; 6,695,813; 6,726,701; 6,752,819; 6,764,506; 6,887,256; 6,913,612; 6,918,921; 6,964,672; 7,014,647; 7,018,372; 7,037,320; 7,094,243; 7,144,408; 7,172,620.  
Additional patents pending.

AccUNET is a trademark of the Abbott Group of Companies.

**RX AccUNET**  
**Embolieschutzsystem (Embolic Protection System – EPS)**

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©2011 Abbott