Emboshield NAV⁶

Embolic Protection System

EN	English	IT	Italiano	HU	Magyar
DE	Deutsch	FI	Suomi	BG	Български
SV	Svenska	FR	Français	RO	Română
DA	Dansk	ES	Español	SK	Slovenčina
PL	Polski	RU	Русский		
CS	Česky	PT	Português		
TR	Türkçe	NL	Nederlands		
NO	Norsk	EL	Ελληνικά		

Graphical Symbols for Medical Device Labeling; Grafische Symbole zur Beschriftung medizinischer Geräte; Symboles graphiques pour l'étiquetage de matériel médical; Símbolos gráficos para el etiquetado de productos sanitarios; Simboli grafici per l'etichettatura di dispositivi medicali; Símbolos Gráficos para Rotulagem de Dispositivos Médicos; Grafiska symboler för märkning av medicinsk utrustning; Grafische symbolen voor de etikettering van medische hulpmiddelen; Grafiske symboler beregnet til mærkning af medicinske produkter; Γραφικά σύμβολα για τη σήμανση ιατρικών συσκευών; Symbole graficzne do oznaczania sprzętu medycznego; Grafikus jelzések orvosi eszközök címkézéséhez; Grafické symboly k označení lékařského zařízení; Grafické symboly na označenie lekárskeho zariadenia

Tıbbi Cihaz Etiketleri için Grafik Semboller; Grafiske symboler for merking av medisinsk utstyr; Графични символи за етикетиране на медицински уреди; Simboluri grafice pentru etichetarea dispozitivelor medicale; Графические символы на этикетках медицинских изделий; Lääketieteellisten laitteiden tuotetarroissa esiintyvät symbolit

Emboshield NAV⁶ Embolic Protection System

1.0 DEVICE DESCRIPTION

The Emboshield NAV⁶ Rapid Exchange (RX) Embolic Protection System is a temporary percutaneous transluminal filtration system designed to be used as a guide wire and to capture embolic material released during an angioplasty and stent procedure within a saphenous vein bypass graft or a carotid artery. The system consists of the following components:

1.1 BareWire Filter Delivery Wire

The BareWire Filter Delivery Wire (Figure 1) is a 0.014" PTFE coated stainless steel wire with a 3 cm platinum / nickel radiopaque distal tip section. The Filter Delivery Wires in the BareWire range have zebra stripes on the proximal end for ease of identification.


Figure 1

The Emboshield NAV⁶ Embolic Protection System is supplied pre-loaded with the BareWire Workhorse Filter Delivery Wire (190 cm). An alternative range of BareWire Filter Delivery Wires are also available for use with the Emboshield NAV⁶ Embolic Protection System (see Table 1). Only the BareWire Filter Delivery Wires shown in Table 1 are intended for use with the Emboshield NAV⁶ System.

Table 1

Description	Length	Product Code
BareWire Distal Access	190 cm	22444-19
BareWire Workhorse	190 cm	22445-19
BareWire Workhorse	315 cm	22445-31
BareWire Support	190 cm	22446-19

1.2 Delivery Catheter

The RX Delivery Catheter (Figure 2) usable length is 135 cm. The crossing profile is either 0.0365" or 0.0415", depending on Filtration Element size. A pull handle is used to deploy the loaded Filtration Element from the pod. Two pairs of silver indicator bands are provided along the catheter shaft. A proximal pair of indicator bands indicates the catheter tip position during advancement through the guide catheter / introducer sheath, and a distal pair of indicator bands indicates the proximity of the RX exit port during catheter retraction. A radiopaque marker band is positioned proximal to the pod.


Figure 2

1.3 Filtration Element

The Filtration Element (Figure 3 and Table 2) consists of a Nylon membrane with an internal Nitinol frame. There are two proximal entry ports and multiple distal perfusion pores. There are proximal and distal marker bands. The Filtration Element is available in two sizes that fit a range of vessel sizes: Small, for use in 2.5 mm - 4.8 mm diameter vessels, and Large, for use in 4.0 mm - 7.0 mm diameter vessels. To aid in visualization, the Filtration Element frame is also radiopaque.


Figure 3

Table 2

Filtration Element	Dim (A) cm	Dim (B) cm	Dim (C) cm	Vessel Size (mm)	Unconstrained Filtration Element (FE) Diameter
Small	1.9	0.2	0.15	2.5 - 4.8	5.0
Large	2.3	0.2	0.15	4.0 - 7.0	7.2

1.4 RX Retrieval Catheter

The RX Retrieval Catheter (Figure 4) has a 139 cm usable length, and a distal radiopaque expansible tip. The maximum outer diameter of the distal tip is 0.067". A handle grip is situated at the proximal end. Two pairs of silver indicator bands are provided along the catheter shaft. A proximal pair indicates the catheter tip position during advancement through the guide catheter / introducer sheath, and a distal pair indicates the proximity of the RX exit port during catheter retraction.


Figure 4

2.0 GENERAL PREPARATION AND USE (ASEPTIC TECHNIQUE)

2.1 Sizing Information

The size of the Filtration Element for placement in the carotid artery or saphenous vein bypass graft is determined by the diameter of the vessel at the intended deployment site. Quantitative Coronary / Carotid Angiography (QCA) is recommended.

2.2 Filter Loading

- Do not move the Delivery Catheter during loading. Movement may cause misalignment of system components.
- Fill the recessed well area with heparinised normal saline (Figure 5) until the Filtration Element is completely immersed (Figure 6).
- Fill the provided syringe with heparinised normal saline in order to flush the Delivery Catheter. Remove air from syringe. Push the syringe tip into the core of the Loading Funnel until it engages the internal diameter of the Loading Funnel (Figure 6). Flush until heparinised normal saline can be seen exiting at the RX exit port (Figure 7).
- Use the syringe with the flushing tip to remove any air bubbles from the Filtration Element. Place the white tip gently inside the Filtration Element and flush (Figure 8).


Figure 5


Figure 6


Figure 7


Figure 8

Load the Filtration Element into the Delivery Catheter as follows:

- Ensure that the Tray is on a level surface. Hold the tray securely with one hand ensuring the well and Filtration Element are clearly visible (Figure 9).
- Unclip the torque device (Figure 10). Ensure that the torque device is securely tightened on the coating free section of the BareWire Filter Delivery Wire (refer to Figure 1). The torque device can be re-tightened if required. Carefully pull the torque device in a controlled manner to retract the BareWire Filter Delivery Wire and observe the Filtration Element being retracted through the Loading Funnel into the Delivery Catheter Pod. The device is fully loaded when the distal marker band on the Filtration Element is at or inside the end of the delivery catheter pod and a definite stop is felt (refer to Figure 11).
- Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.

- Inspect the Delivery Catheter Pod prior to removal from the tray to ensure the Filtration Element is fully loaded (refer to Figure 11).
- Remove the loaded Delivery Catheter and the BareWire Filter Delivery Wire from the tray by gripping the black handle and the Torque Device together and gently pulling them out of the tray (Figure 12).
- **Do not attempt to remove either the Delivery Catheter or the Filter Delivery Wire individually.**
- Inspect the Delivery Catheter to ensure the Filtration Element remains loaded properly into the pod at the end of the Delivery Catheter and that the BareWire Filter Delivery Wire is not damaged.
- Wipe the length of exposed Filter Delivery Wire from the RX port to the torque device with a sterile gauze pad soaked with heparinised normal saline.
- The Delivery System is now ready for use.


Figure 9

Figure 10


Figure 11

Figure 12

2.3 Using Contrast Media

During the insertion of rapid exchange catheters through guide catheters or sheaths, careful handling is required to ensure that air is not drawn into the access device.

It is therefore recommended that flushing of contrast media (or other fluids) is performed before or after insertion of the catheter, but not while the catheter is within the access device.

In the case where a contrast media injection must be performed with the catheter in place, it is essential to ensure that no air is present within the access device prior to injection.

If aspiration is to be performed prior to contrast media injection, it should be performed slowly and steadily at a rate of not more than 0.5 ml (0.5 cc) per second until the user can visually confirm that no further air is entering the aspiration syringe (see to Figure 13). Following all instructions, slow catheter insertion, and good RHV (Rotating Haemostasis Valve) control will help to minimize the risk of air embolization.


Figure 13

2.4 Preparation for Deployment

WARNINGS:

Only use the BareWire Filter Delivery Wires listed in Table 1. Use of other guide wires will lead to loss of the Filtration Element or an inability to retrieve the Filtration Element.

The Emboshield NAV® Embolic Protection System is NOT to be deployed with an access device that uses an integrated leaflet-type valve.

The Emboshield NAV® Embolic Protection System should be used with an RHV with minimum internal diameter of 0.096".

The Emboshield NAV® Embolic Protection System should be used with a guide catheter or sheath access device with an internal diameter equal to or greater than 0.070".

A device or access device that is kinked or damaged in any way should not be used.

Assess all equipment before use.

Maintain a snug seal between the device and the haemostasis valve during catheter insertion.

Catheter insertion should be performed slowly to minimize the risk of air entrapment. Failure to observe this may result in air being drawn into the access device through the RHV.

Perform all device manipulation using fluoroscopic visualization.

Do not torque or pull back on the RX Delivery Catheter during advancement.

Do not pull back on the BareWire Filter Delivery Wire during advancement.

Torquing the BareWire against resistance may cause Barewire damage.

- Confirm the compatibility of the Emboshield NAV® Embolic Protection System with the other interventional devices intended for use in the procedure.
- The BareWire Filter Delivery Wire should be kept moistened throughout the procedure.
- A range of alternative BareWire Filter Delivery Wires is available (see Table 1).
- The Filtration Element can only be positioned proximal to the radiopaque distal section of the Filter Delivery Wire.

Device Introduction

- Use the introducer tool provided to pass the loaded Delivery Catheter and Filter Delivery Wire through the haemostasis valve of the guide catheter or sheath. Indicator bands are positioned 90 cm and 100 cm from the tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the guide catheter / introducer sheath.
- Advance the device through the guide catheter / introducer sheath to a location proximal to the lesion. Advance the Filter Delivery Wire across the lesion.

Use of an Alternative BareWire Filter Delivery Wire

- Load the device with the pre-packed BareWire Filter Delivery Wire as previously described.
- To use an alternative BareWire Filter Delivery Wire, carefully remove the torque device and the Filter Delivery Wire and dispose of the Filter Delivery Wire used to load the device (See Figure 14).
- Remove the replacement BareWire Filter Delivery Wire from its hoop and wipe with a sterile gauze pad soaked with heparinised normal saline.
- Use an introducer tool to pass the wire through the haemostasis valve of the guide catheter or introducer sheath toward the deployment site. **(The BareWire Filter Delivery Wire should be kept moistened throughout the procedure).**
- Advance the BareWire Filter Delivery Wire across the lesion.
- Carefully advance the Delivery Catheter over the BareWire Filter Delivery Wire (Figure 15). The BareWire Filter Delivery Wire will exit the Delivery Catheter at the catheter wire exit port.
- Holding the BareWire Filter Delivery Wire, advance the Delivery Catheter through the haemostasis valve toward the deployment site. Indicator bands are positioned 90 cm and 100 cm from the catheter tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the Guide Catheter. Do not torque or pull back on the RX Delivery Catheter during advancement. Advance the Delivery Catheter to a location proximal to the lesion. Then continue with the Device Placement Instructions below.

Figure 14

Figure 15

Device Placement

- Position the 3 cm radiopaque portion of the Filter Delivery Wire in a suitable section of the vessel distal to the intended site of Filtration Element deployment. **(The Filtration Element can only be positioned proximal to the radiopaque distal section).**
- Carefully advance the Delivery Catheter over the BareWire Filter Delivery Wire to the intended filter deployment site. Do not torque or pull back on the RX Delivery Catheter during advancement.
- Once the Delivery Catheter is distal to the lesion, establish desired stent location and note where the distal edge of the stent will be placed. Using the marker bands on the Filtration Element and Delivery Catheter as guides, position the Delivery Catheter accordingly. Ensure that there is enough distance between the proximal tip of the Filtration Element and the most distal tip of any interventional device to be introduced over the Filter Delivery Wire to avoid contact during the procedure. Ensure that there is enough distance between the proximal end of the 0.014" Radiopaque Distal Tip of the Filter Delivery Wire and the distal tip of the Filtration Element to avoid unnecessary contact during the procedure. **Stent Delivery System / Balloon Catheter distal tips may vary in size and radiopacity. The tip of a balloon catheter or a stent delivery system should not contact the Filtration Element (see Figure 16).**
- Once the Delivery Catheter is positioned, the Filtration Element may be deployed.


Figure 16

2.5 Device Deployment

- Remove the red Locking Clip from the black handle of the Delivery Catheter.
- Hold the black handle steady and pull back on the white pull handle (see Figure 2) to deploy the Filtration Element. Deployment is occurring when the distal marker band on the Catheter moves proximally relative to the Filtration Element marker bands.
- WARNING: A high-pressure contrast media injection may cause Filtration Element movement in the vessel.**
- Remove the Delivery Catheter while maintaining the BareWire Filter Delivery Wire position. Silver indicator bands on the Delivery Catheter shaft help indicate the position of the BareWire Filter Delivery Wire exit port. A single silver band is located 15 cm from the exit port and twin silver bands are located 5 cm from the exit port. Open the haemostasis valve after the twin indicator bands on the hypotube shaft become visible.
- Discard the Delivery Catheter after withdrawal. Perform the required interventions over the BareWire Filter Delivery Wire.

2.6 Device Retrieval

- Remove the Retrieval Catheter tray from the pouch. Remove the Catheter from the tray.
- Flush the catheter lumen by injecting heparinised normal saline through the catheter's distal tip using the flushing tip provided in the tray and a standard syringe. Check that saline is exiting at the RX Port.
- Under fluoroscopy, advance the Retrieval Catheter over the BareWire Filter Delivery Wire. Indicator bands are positioned 90 cm and 100 cm from the tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the guide catheter / introducer sheath.
- Advance through the stent / treated lesion until the radiopaque catheter tip contacts the proximal end of the Filtration Element. Do not torque the Retrieval Catheter during advancement.
- Grip the BareWire Filter Delivery Wire using the torque device provided in the tray, and retract the BareWire Filter Delivery Wire until the proximal radiopaque section of the wire meets the distal end of the Filtration Element.

The standard method of retrieval is as follows:

- Hold the Retrieval Catheter position steady by gripping the shaft close to the RHV.
- Pull on the tightened torque device to retract the BareWire Filter Delivery Wire until the Filtration Element is fully enclosed in the radiopaque expandable tip.
- Complete retrieval of the Filtration Element has been achieved when the distal Filtration Element marker band has been withdrawn to within 2 mm of the distal end of the radiopaque catheter tip (see Figure 17.3).
- Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.**
- Retract the Retrieval Catheter over the BareWire Filter Delivery Wire (or remove Retrieval Catheter and BareWire Filter Delivery Wire as one unit) until the tip of the Retrieval Catheter is adjacent to the tip of the guide catheter or sheath.
- If significant resistance is felt, retract the guide catheter or sheath and the Retrieval Catheter together.
- Remove the Retrieval Catheter from the patient.
- If a rotating haemostasis valve is used, ensure the valve is fully open during removal of the RX section of the Retrieval Catheter. See Figure 17.1, Figure 17.2 and Figure 17.3.
- If an exceptionally large embolic load prevents complete withdrawal of the Filtration Element into the Retrieval Catheter tip, the following method of Retrieval must be employed:**
- Hold the Retrieval Catheter position steady by gripping the shaft close to the RHV.
- Pull on the tightened torque device to retract the BareWire Filter Delivery Wire and Filtration Element into the Retrieval Catheter tip.
- Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.
- WARNING: If the Filtration Element cannot be fully retracted as described above, the BareWire Filter Delivery Wire must be removed with the Retrieval Catheter.** Failure to do so may result in redeployment of the Filtration Element from the Retrieval Catheter tip. The proximal openings of the Filtration Element have been fully contained within the catheter tip once the radiopaque frame has been withdrawn into the radiopaque catheter tip (see Figure 17.2).
- Retract the Retrieval Catheter, Filtration Element and BareWire Filter Delivery Wire as one unit until the tip of the Retrieval Catheter is adjacent to the tip of the guide catheter or sheath.
- Pull the tip of the Retrieval Catheter into and through the guide catheter / introducer sheath. If significant resistance is felt, retract the guide catheter or sheath and the Retrieval Catheter, Filtration Element and BareWire Filter Delivery Wire together.
- Remove the Retrieval Catheter from the patient. If a rotating haemostasis valve is used, ensure the valve is fully open during removal of the RX section and tip of the Retrieval Catheter and retrieved Filtration Element.


Figure 17.1

Figure 17.2

Figure 17.3

3.0 INTENDED USE

3.1 Indications

The system is designed to reduce distal embolization while maintaining blood flow during angioplasty and stent procedures. The system is indicated for use during percutaneous transluminal angioplasty and stenting procedures in saphenous vein grafts and carotid arteries.

3.2 Contraindications

- Inability to cross the lesion with the BareWire Filter Delivery Wire.
- Diffusely diseased vessels where there is no disease-free section in which to deploy the Filtration Element
- Insufficient straight section of vessel distal to the lesion to permit Filtration Element deployment.
- Patients in whom anticoagulant and / or antiplatelet therapy is contraindicated.
- Patients with severe vascular tortuosity or anatomy that would preclude the safe introduction of the Guide Catheter / Introducer Sheath, Barewire, Delivery Catheter, Filtration Element, and / or Retrieval Catheter.
- Patients with a known hypersensitivity to nickel-titanium.
- Patients with uncorrected bleeding disorders.
- Lesions in the ostium of the common carotid artery.

3.3 General Precautions

- Refer to instructions supplied with all interventional devices to be used with the Emboshield NAV® Embolic Protection System for their intended uses, contraindications, and potential complications.
- The Emboshield NAV® System is supplied sterile. Do not use if the package has been opened or is damaged. Carefully inspect system components prior to use to verify that they have not been damaged and that the size, shape and condition are suitable for the procedure for which they are to be used.
- This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.
- Store in a cool, dark, dry area. Do not expose to organic solvents or ionising radiation. Rotate inventory so that devices are used prior to the "Use Before Date" shown on the package label.
- Precautions to prevent or reduce clotting and thrombus should be taken when any interventional device is used. Flush or rinse all devices entering the vascular system with heparinised normal saline prior to use.
- The clinician should be familiar with and experienced in standard techniques of Rapid Exchange percutaneous transluminal angioplasty and stenting and be knowledgeable of the current medical literature concerning the complications of such procedures.
- Torquing the BareWire Filter Delivery Wire against resistance may cause BareWire Filter Delivery Wire damage.

3.4 Specific Precautions / Warnings

The Emboshield NAV® System can only be used with a BareWire Filter Delivery Wire described in Table 1. Use of the system with guide wires other than the BareWire Filter Delivery Wires (see Table 1) will lead to loss of the Filtration Element during the procedure.

- To optimize system performance, it is imperative that careful attention is paid to the preparation of the system using the preparation techniques specified in the Instructions For Use.

- The Emboshield NAV[®] is intended for one time use only. DO NOT resterilize and / or reuse it, as this can potentially result in compromised system performance and risk of cross contamination.
- Do not use the product after the Use by Date specified on the label.
- The size of the Filtration Element should be selected to match the diameter of the vessel at the intended site of deployment. The clinician should take into consideration the slight increase in vessel diameter that may occur after treatment of the lesion.
- The minimum expanded stent Internal Diameter required for retrieval of a large embolic load is ϕ 2.5 mm.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- The minimum guide catheter or sheath internal diameter is printed on the package label. Do not use a smaller guide catheter or sheath than indicated.
- Appropriate anti-platelet, anticoagulant and, if necessary, vasodilator therapy must be used during the procedure. Anticoagulant therapy sufficient to maintain an Activated Clotting Time of at least 250 seconds for the duration of the procedure is recommended.
- Adequate guide catheter or sheath support is required for the introduction of the Delivery Catheter, Retrieval Catheter and all interventional devices to be used during the procedure.
- Maintain a snug seal between the device and the hemostasis valve during catheter insertion. Failure to observe this may result in air being drawn into the access device through the hemostasis valve. Catheter insertion should be performed slowly to minimize the risk of air entrapment.
- During the insertion of Rapid Exchange catheters through guide catheters or sheaths careful handling is required to ensure that air is not drawn into the access device. It is therefore recommended that flushing of contrast media (or other fluids) is performed before or after insertion of the catheter, but not while the catheter is within the access device.
- Do not advance any component of the Emboshield NAV[®] Embolic Protection System against significant resistance. The cause of any resistance should be determined via fluoroscopy and remedial action taken.
- Deployment of the Filtration Element, the subsequent deployment of complementary interventional devices over the BareWire Filter Delivery Wire, and the retrieval of the Filtration Element should only be performed under fluoroscopic observation.
- Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- The Delivery Catheter should not be pulled back quickly in the event that resistance is felt during catheter advancement. Rapid withdrawal of the Delivery Catheter in the presence of resistance may result in premature deployment of the Filtration Element.
- High-pressure contrast injection may cause Filtration Element movement in the vessel.
- Do not attempt to move the Filtration Element after deployment and prior to the start of retrieval.
- Do not torque the RX Delivery and Retrieval Catheters.
- The maintenance of blood flow through the system should be observed throughout the procedure by the use of contrast injection.
- If excessive debris is collected in the Filtration Element such that distal perfusion of contrast media is significantly reduced or no contrast media is perfusing past the filter, the Emboshield NAV[®] System may have reached its maximum capacity to contain emboli. Remove and replace the Emboshield NAV[®] System. Otherwise it may be difficult to completely recover all embolic debris and the potential for thrombus formation may increase.
- Carefully inspect system components prior to use to verify that they have not been damaged and that the size, shape and condition are suitable for the procedure for which they are to be used.
- A device or access device that is kinked or damaged in any way should not be used.
- Confirm the compatibility of the Emboshield NAV[®] Embolic Protection System with the interventional devices before actual use.
- Precautions to prevent or reduce clotting should be taken when any interventional device is used. Flush or rinse all devices entering the vascular system with heparinised normal saline or alternative anticoagulant, prior to use.
- The Emboshield NAV[®] Embolic Protection System must be used with a guide catheter or introducer sheath to maintain adequate support of the BareWire throughout the procedure.
- To reduce the potential for the liberation of emboli during lesion crossing, the system should be carefully manipulated and not advanced against resistance.
- If the Filtration Element moves into the stented vessel segment prior to retrieval, **DO NOT RETRIEVE**. Use the Retrieval Catheter to gently manoeuvre the Filtration Element distally until it is situated in an unstented portion of vessel. Retrieval should then proceed.
- Maintain proper guiding catheter / sheath support in the common carotid artery throughout the procedure. Ensure that there is enough distance between the proximal tip of the Filtration Element and the most distal tip of any interventional device to be introduced over the Filter Delivery Wire to avoid engagement. The tip of a balloon catheter or a stent delivery system should not contact the Filtration Element. Failure to maintain adequate distance could result in inadvertent Filtration Element movement and Stent Delivery System tip / Filtration Element entanglement and / or Filtration Element / Stent entanglement if guide catheter or sheath prolapse occurs.
- Venous access should be available during carotid stenting in order to manage bradycardia and / or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.
- Removal of the Emboshield NAV[®] Filtration Element with any interventional devices other than the Emboshield NAV[®] Retrieval Catheter has not been tested.

- Aneurysm
- Angina / Coronary ischemia
- Arteriovenous Fistula
- Bacteraemia or septicemia
- Bleeding from anticoagulant or antiplatelet medications
- Bradycardia / arrhythmia
- Cerebral edema
- Cerebral hemorrhage
- Congestive Heart Failure
- Death
- Drug reactions
- Embolism (including air and device)
- Emergent or urgent Endarterectomy
- Fever
- Filter thrombosis / occlusion
- Fluid overload
- Groin hematoma, with or without surgical repair
- Hemorrhage or hematoma
- Hemorrhagic stroke
- Headache
- Hypotension
- Hyperperfusion syndrome
- Hypertension
- Infection / sepsis
- Ischemia / infarction of tissue / organ

- Myocardial Infarction
- Other conduction disturbances
- Pain and tenderness
- Pain, infection, or discomfort at the access site
- Pseudoaneurysm
- Renal failure / insufficiency
- Restenosis of the stented artery
- Seizure
- Stent deformation, collapse, fracture, movement of stent, possibly requiring emergency surgery
- Stent / filter entanglement / damage
- Stroke or other neurological complications
- Thromboembolic episodes
- Thrombophlebitis
- Total occlusion of the artery
- Transient ischemic attacks (TIAs)
- Vascular access complications (e.g. loss of pulse, femoral artery pseudoaneurysm and infection)
- Vasospasm
- Ventricular fibrillation
- Vessel dissection, rupture, or perforation
- Vessel thrombosis (partial blockage)
- Unstable angina pectoris

Any adverse event occurring involving the Emboshield NAV[®] Embolic Protection System should be reported immediately to Abbott Vascular, Customer Service.

3.6 Disclaimer of Warranties

There is no express or implied warranty, including any implied warranty of satisfactory quality or fitness for a particular purpose, on the Emboshield NAV[®] System described in this document. Under no circumstances shall Abbott Vascular be liable for any direct, incidental or consequential damages other than as expressly provided by specific law. Description or specifications in Abbott Vascular's printed matter, including this document, are meant solely to generally describe the product at any time of manufacture and do not constitute any express warranties. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Abbott Vascular has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after the device leaves Abbott Vascular's possession. No representative of Abbott Vascular may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

4.0 TRADEMARKS

Emboshield is a trademark of the Abbott Group of Companies.
BareWire is a trademark of the Abbott Group of Companies.

3.5 Potential Adverse Effects

As reported in the literature, the following adverse events are potentially associated with carotid stents and embolic protection systems:

- Abrupt closure
- Allergic reactions
- Loss of all of or part of the Filtration Element

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© 2010, Abbott Laboratories