RX Acculink

Carotid Stent System

EN	English	IT	Italiano	HU	Magyar
DE	Deutsch	FI	Suomi	BG	Български
SV	Svenska	FR	Français	RO	Română
DA	Dansk	ES	Español	SK	Slovensky
PL	Polski	RU	Русский	ET	Eesti keel
CS	Česky	PT	Português	LV	Latviešu
TR	Türkçe	NL	Nederlands	LT	Lietuvių
NO	Norsk	EL	Ελληνικά	SL	Slovenščina

Graphical Symbols for Medical Device Labeling; Grafische Symbole zur Beschriftung von Medizinprodukten; Symboles graphiques pour l'étiquetage des dispositifs médicaux; Símbolos gráficos utilizados en el etiquetado de productos sanitarios; Simboli grafici per l'etichettatura di dispositivi medicali; Símbolos Gráficos para Rotulagem de Dispositivos Médicos; Grafiska symboler för märkning av medicinsk utrustning; Grafische symbolen voor de etikettering van medische hulpmiddelen; Grafiske symboler beregnet til mærkning af medicinske produkter; Γραφικά σύμβολα για τη σήμανση ιατρικών συσκευών

Symbole graficzne do oznaczenia sprzętu medycznego; Grafikus jelzések orvosi eszközök címkézéséhez; Grafické symboly k označení zdravotnického prostředku; Tibbi Cihaz Etiketleri için Grafik Semboller; Grafiske symboler for merking av medisinsk utstyr

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RX Acculink Carotid Stent System

DEVICE DESCRIPTION

The RX Acculink Carotid Stent System includes:

- A pre-mounted self-expanding nickel titanium stent.
- A rapid exchange stent delivery catheter.
- Radiopaque markers on the shaft that mark the stent location, as illustrated in Figure 1.


Figure 1: Tip Detail Showing Stent Location Markers

The catheter is comprised of a retractable sheath covering the stent during delivery, a radiopaque tip, an internal guide wire lumen, a handle assembly with a safety lock (Figure 2), and a pullback handle. The entire system is shown in Figure 3. With the handle in the unlocked position, retracting the pullback handle deploys the stent.


Figure 2: Handle Detail


Figure 3: Delivery System Schematic

The RX Acculink Carotid Stent System is available in several lengths and diameters as listed in Table 1. Stents should always be sized to the internal carotid artery (ICA) reference vessel and should provide stent-to-artery ratios between 1.1:1 and 1.4:1.

For tapered stents, ensure the stent sizing provides both proximal and distal stent apposition. The proximal stent end should be sized to the common carotid artery (CCA) and the distal end sized to the internal carotid artery (ICA). Both stent ends should be sized between the 1.1:1 and 1.4:1 stent-to-artery ratio. See Table 2 for recommended reference vessel diameters.

Table 1 – RX Acculink Carotid Stent System – Stent Diameters

Unconstrained Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility	Guide Compatibility	ICA Reference Vessel Diameter (mm)
5.0	20, 30, 40	6F	8F	3.6 – 4.5
6.0	20, 30, 40	6F	8F	4.3 – 5.4
7.0	20, 30, 40	6F	8F	5.0 – 6.4
8.0	20, 30, 40	6F	8F	5.7 – 7.3
9.0	20, 30, 40	6F	8F	6.4 – 8.2
10.0	20, 30, 40	6F	8F	7.1 – 9.1

Table 2 – RX Acculink Carotid Stent System – Tapered Stent Diameters

Unconstrained Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility	Guide Compatibility	ICA Reference Vessel Diameter (mm)	CCA Reference Vessel Diameter (mm)
6 – 8 Taper	30, 40	6F	8F	4.3 – 5.4	5.7 – 7.3
7 – 10 Taper	30, 40	6F	8F	5.0 – 6.4	7.1 – 9.1

HOW SUPPLIED

Sterile – Sterilized with electron beam radiation. Non-pyrogenic. Do not use if package is open or damaged.

This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents – One (1) RX Acculink Carotid Stent System, one (1) Funnel Introducer.

Storage – Store at room temperature only.

INDICATIONS

The RX Acculink Carotid Stent System is indicated for use in patients with atherosclerotic disease of the carotid arteries.

CONTRAINDICATIONS

The RX Acculink Carotid Stent System is contraindicated for use in:

- Patients in whom anticoagulant and / or antiplatelet therapy is contraindicated.

• Severe vascular tortuosity or anatomy that would preclude the safe introduction of a guide catheter, sheath, embolic protection system or stent.

- The presence of extensive atherosclerotic disease involving the aortic arch and proximal common carotid artery that would preclude the safe introduction of a guide catheter or sheath.
- Patients with known hypersensitivity to nickel titanium.
- Patients with uncorrected bleeding disorders.

WARNINGS

The safety and efficacy of the RX Acculink Carotid Stent System have not been demonstrated with embolic protection systems other than Abbott Vascular's Accunet or Emboshield family of Embolic Protection Systems (EPS). Refer to the Instructions for Use document for the Embolic Protection System that will be used for specific device instructions.

Clinical study results suggest lower event rates when the RX Acculink Carotid Stent System is used in conjunction with an embolic protection device.

Do not use if the temperature indicator is black.

Once deployment is initiated, the stent cannot be repositioned or recaptured.

Stenting across a major bifurcation may hinder or prevent future diagnostic or therapeutic procedures.

When multiple stents are required, stent materials should be of similar composition.

Persons allergic to nickel titanium may suffer an allergic reaction to this implant.

Only physicians familiar with the complications, side effects and hazards commonly associated with carotid stent placement should use this device.

Refer to the instructions for use supplied with any interventional devices to be used in conjunction with the RX Acculink Carotid Stent System for their intended uses, contraindications, and potential complications.

As with any type of vascular implant, infection secondary to contamination of the stent may lead to thrombosis, pseudoaneurysm, or rupture.

In patients requiring the use of antacids and / or H2-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g., aspirin) may be adversely affected. The appropriate antiplatelet and anticoagulation therapy should be administered pre- and post-procedure as suggested in these instructions. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Anticoagulant therapy sufficient to maintain an Activated Clotting Time of at least 250 seconds for the duration of the procedure is recommended.

The stent may cause a thrombus, distal embolization or may migrate from the site of implant down the arterial lumen. Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration. In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.

Overstretching of the artery may result in rupture and life-threatening bleeding.

PRECAUTIONS

Carefully inspect the RX Acculink Carotid Stent System to verify that the device has not been damaged in shipment. Take care to avoid unnecessary handling, which may kink or damage the delivery system.

Use with bleedback control hemostatic valves is not recommended.

For best device performance, the guide wire exit notch should remain within the guiding catheter or sheath.

Ensure the stent system is fully flushed. Do not use the delivery system if flush is not observed exiting at the distal end of the sheath.

Stent Handling – Precautions

- The delivery system has an internal hypotube. Care must be taken to keep the device straight. Do not use if device is kinked.
- This device is intended for single-use only; do not reuse. Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.
- Note the product "Use by" date specified on the package.
- Do not remove stent from its delivery system, as removal may damage the stent and / or lead to stent embolization. Stent system is intended to perform as a system.
- Delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is most important during catheter removal from packaging, mandrel removal, placement over guide wire and advancement through rotating hemostatic valve (RHV) adapter and guiding catheter hub.
- Do not expose the delivery system to organic solvents (e.g., alcohol) as structural integrity and / or function of the device may be impaired.

Stent Placement – Precautions

- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent and may cause acute closure of the vessel requiring additional intervention (carotid endarterectomy, further dilatation, or placement of additional stents).
- Do not attempt to pull a partially expanded stent back through the guiding catheter or sheath; dislodgment of the stent from the delivery system may occur.
- Stent retrieval methods (use of additional wires, snares and / or forceps) may result in additional trauma to the carotid vasculature and / or the vascular access site. Complications may include death, stroke, bleeding, hematoma or pseudoaneurysm.
- If a filter-based embolic protection system (EPS) is used, allow for and maintain adequate distance between the RX Acculink and the EPS to avoid potential filter engagement with the RX Acculink tip and / or filter entanglement with the deployed stent. If filter engagement and / or entanglement or filter detachment occurs, surgical conversion or additional catheter-based intervention may be required.
- If resistance is met during delivery system introduction, the system should be withdrawn and another system used.
- Venous access should be available during carotid stenting in order to manage potential hemodynamic instability (i.e., bradycardia and / or hypotension), by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.
- When catheters are in the body, they should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed.
- The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.
- Prior to stent deployment, remove all slack from the delivery system.
- When more than one stent is required to cover the lesion, or if there are multiple lesions, the distal lesion should be stented first, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent for placement of the distal stent and reduces the chance of dislodging stents that have already been placed.

If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum (approximately 5 mm). In no instance, should more than 2 stents ever overlap.

- The RX Acculink Carotid Stent System must be used with a guiding catheter or introducer sheath to maintain adequate support of the 0.014" (0.36 mm) guide wire throughout the procedure.

Delivery System Preparation

LEAVE THE SAFETY LOCK CLOSED UNTIL THE STENT IS READY TO DEPLOY.

Post Implant – Precautions

- Great care must be exercised when crossing a newly deployed stent with a guide wire, balloon or delivery system to avoid disrupting the stent geometry.
- Through non-clinical testing, the Acculink Carotid Stent has been shown to be MR conditional. Patients may undergo MRI at field strengths of 3.0 Tesla or less, a maximum spatial gradient of 3.3 Tesla/meter and a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of MRI. The Acculink Carotid Stent should not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of stent migration at field strengths higher than 3.0 Tesla or a maximum spatial gradient higher than 3.3 Tesla/meter.
- In this testing, the stent produced a temperature rise of less than or equal to 0.5°C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of MRI. The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.
- MRI quality may be compromised if the area of interest is in the exact same area as or relatively close to the position of the stent.

ADVERSE EVENTS

Based on the literature, and on clinical and commercial experience with carotid stents and embolic protection systems, the following alphabetical list includes possible adverse events associated with use of these devices:

- Allergic reactions to antiplatelet agents / contrast medium
- Aneurysm
- Angina / coronary ischemia
- Arrhythmia
- Arterial occlusion / thrombosis at puncture site or remote site
- Arteriovenous fistula
- Bacteremia or septicemia
- Bleeding from anticoagulant or antiplatelet medications
- Cerebral edema
- Cerebral hemorrhage
- Cerebral ischemia / transient ischemic attack (TIA)
- Congestive heart failure (CHF)
- Death
- Emboli, distal (air, tissue or thrombotic emboli)
- Fever
- Filter thrombosis / occlusion
- Groin hematoma, with or without surgical repair
- Hemorrhage, with or without transfusion
- Hyperperfusion syndrome
- Hypotension / hypertension
- Infection and pain at insertion site
- Ischemia / infarction of tissue / organ
- Myocardial infarction (MI)
- Pain (head, neck)
- Pseudoaneurysm, femoral
- Renal failure / insufficiency
- Restenosis of stented segment
- Seizure
- Severe unilateral headache
- Stent thrombosis / occlusion
- Stroke / cerebrovascular accident (CVA)
- Total occlusion of carotid artery
- Vessel dissection, perforation, or rupture
- Vessel spasm or recoil

CLINICIAN USE INFORMATION

Materials Required

- 8F guiding catheter or 6F introducer sheath compatible with the vascular anatomy. Guiding catheter or sheath should not exceed 100 cm length.
- Rotating hemostatic valve (not recommended for use with bleedback control hemostatic valves)
- Balloon dilatation catheter
- 0.014" (0.36 mm) x 190 cm guide wire
- This device is not compatible with any guide wire larger than 0.014" (0.36 mm).**
- RX Acculink Carotid Stent System
- 1,000 u/500 cc heparinized normal saline (HepNS)
- Torque device
- Guide wire introducer
- Two to three 10 – 20 cc syringes

Inspection Prior to Use

- Inspect the temperature indicator on the inner pouch. **Do not use if BLACK.**
- Remove the stent system from its protective packaging.
Note: Remove the handle from the package prior to removing the shaft from the hoop. Lay the device flat and minimize unnecessary handling. The shaft may kink if not handled carefully. **Do not use if kinked.**
- Ensure that the distal mandrel remains within the inner lumen. Inspect the stent through the delivery system sheath to verify that it has not been damaged during shipment and that the stent does not overlap the proximal marker. Ensure that the stent is fully covered by the sheath.
- Read the specifications on the housing assembly and verify that the stent is the correct diameter and length. Ensure that the dial on the housing assembly is in the locked position. Do not use if any defects are noted.

THIS STENT SHOULD NOT BE REMOVED FROM THE DELIVERY SYSTEM PRIOR TO DEPLOYMENT.

Preparation for Use

Lesion Preparation

- If desired, pre-dilate the lesion with an appropriate size balloon dilatation catheter to a minimum opening of 2.5 mm.
Note: If no pre-dilatation is performed, there must be a minimum luminal opening of 2.5 mm to enable passage of the stent delivery system.
- Maintain guide wire position and withdraw the balloon dilatation catheter.

- Keep the distal mandrel in the guide wire lumen.
- Fill a 10 cc syringe with heparinized normal saline, and inject the saline into the system through the flush port at the proximal end of the housing assembly. Flush will exit the guide wire exit notch first.
- Tightly pinch the guide wire exit notch. Vigorously flush until fluid is observed exiting at both the distal end of the sheath and at the mandrel. Then gently remove the distal mandrel by twisting and pulling as in Figure 4. **DO NOT HOLD THE SHEATH during mandrel removal.** Continue flushing after mandrel removal. If the distal mandrel does not remove easily or flush solution is not observed at the stent sheath junction, do not use the device.
- Keep the device lying flat to avoid kinking in the shaft.


Figure 4: Distal Mandrel Removal
(Hold the tip to remove the distal mandrel.)

Delivery Procedure

- After the pre-dilatation catheter has been removed, **BACKLOAD** the delivery system onto the 0.014" (0.36 mm) guide wire. The guide wire will exit approximately 22 cm from the distal tip.
Note: If using a sheath with a hemostatic valve, the funnel introducer should be placed onto the system prior to backloading of the guide wire.
- Secure the guide wire and sheath position using one hand. Use the other hand to advance the delivery system to the lesion site. Use the radiopaque markers to locate the stent position.
- If applicable, ensure that the RHV remains OPEN and that bleedback is observed.

Stent Deployment

- Place the housing on a stable surface or on the patient's leg. If using a sheath with a hemostatic valve, slide the funnel introducer forward along the shaft of the system and insert it into the valve opening.
- Confirm the stent position angiographically. Adjust the position if necessary. The device is designed to be deployed using one hand. Position the thumb in the textured proximal groove and place two fingers on the pullback handle as shown in Figure 5.
- Turn the safety lock counter-clockwise to the deployment position, symbolized by an open padlock icon . The arrow on the lock will point in the direction the handle will move. ENSURE THAT THE RHV IS OPEN. Remove any slack from the delivery system and reconfirm stent position.
- Ensure that the guide wire and sheath do not move during deployment. Immobilize the guide wire and RHV or sheath by holding them in place with your other hand.
- While pressing down with the thumb to avoid any forward motion, retract the handle.
Note: If significant resistance is encountered during handle pullback before the stent is deployed, relock the handle and remove the device.
- Re-advance the sheath and relock the delivery system before removal into the guiding catheter / sheath. If applicable, remove the funnel introducer from the hemostatic valve, and then remove the delivery system from the patient.
- If necessary, the stent can be post-dilated with a dilatation catheter to ensure good stent apposition. Do not expand the stent past its labeled unconstrained maximum diameter.


Figure 5: Deployment Demonstration
(With the guide position fixed, deploy with one hand. PUSH DOWN on thumb groove and retract pullback handle.)

PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 5,514,154; 5,735,893; 5,759,192; 5,780,807; 6,131,266; 6,325,824; 6,375,676; 6,468,302; 6,537,311; 6,569,193; 6,582,460; 6,599,296; 6,695,862; 6,709,454. Other U.S. patents pending. Foreign patents issued and pending.

Acculink is a trademark of the Abbott Group of Companies.

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©2013 Abbott