StarClose SE

Vascular Closure System

PPL2101550 (9/16/14)

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3. Prior to use, ensure that the StarClose SE Vascular Closure System and sterile packaging have not been damaged during shipment. Examine all components prior to use to verify proper function. Exercise care during device handling to reduce the possibility of accidental device breakage.
4. As with all catheter-based procedures, infection is a possibility. Observe sterile technique at all times when using the StarClose SE Vascular Closure System. Employ appropriate groin management, as per hospital protocol, post-procedure and post-hospital discharge to prevent infection.
5. Use a single wall puncture technique. Do not puncture the posterior wall of the artery.
6. Do not use the StarClose SE Vascular Closure System to close vessels with diameters less than 5 mm.
7. Do not deploy the Clip in areas of calcified plaque.
8. The StarClose SE Vascular Closure System can **ONLY** be used with the StarClose Exchange System (included in the StarClose SE Vascular Closure System packaging).
9. **Do not advance or withdraw the StarClose SE Vascular Closure Device against resistance until the cause of that resistance has been determined.** Excessive force used to advance or torque the StarClose SE device should be avoided, as this may lead to significant vessel damage and / or breakage of the device, which may necessitate interventional and / or surgical removal of the device and vessel repair.

ADDITIONAL CONSIDERATIONS

The StarClose Clip has been shown to be MR Conditional immediately following implantation. A patient with this implant can be scanned safely immediately after clip placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.

In non-clinical testing, the StarClose Clip produced a temperature rise of 0.5°C at maximum MR system-reported whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 Tesla MR system using a transmit / receive body coil.


The MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the StarClose Clip. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary.

CLOSURE PROCEDURE

1. It is necessary to create a 5 – 7 mm skin incision at the sheath site to accommodate the insertion of the Clip Delivery Tube into the tissue tract. This should be done at the beginning of the procedure prior to the administration of anticoagulants and antiplatelet agents, if possible. Blunt dissection by single spread with surgical instrument in the skin incision is recommended, especially in those patients with scar tissue from previous catheterization procedures.
2. Perform a femoral angiogram through the side port of the procedural sheath to determine the location of the arteriotomy site, the vessel size, and the presence of disease (calcified plaque, stenosis, chronic or acute occlusion), tortuosity or presence of arterial wall dissection.
3. Re-prep the access site, placing clean towels around the access site, and wear new sterile gloves prior to handling the Clip Applier and proceed with the closure procedure.
4. Prepare the access site for closure:
 - a. Insert a guide wire and exchange for a StarClose Exchange Sheath.

5. **CLICK 1:** Connect the StarClose Exchange Sheath to the Clip Applier (Figure 3).

Figure 3


- a. Approximate the Clip Applier to the Hub of the StarClose Exchange Sheath. Hold the bottom of the Sheath Hub with the left hand and insert the distal tip of the Flex-guide into the Hub using the right hand.
- b. Slowly advance the Flex-guide down through the Sheath. Make sure not to severely bend the Flex-guide during insertion and advancement.

Note: When connecting the Sheath Hub to the Clip Applier, hold the Hub up at the same angle as the tissue tract above the skin surface. This action allows a firm connection to be made without pushing against the skin.

- c. An audible "click" should be heard. Check to make sure the Hub-to-Clip Applier engagement is secure by gently pulling on the Sheath Hub.

6. **CLICK 2:** Deploy Locator Wings and initiate splitting of Sheath (Figure 4).


Figure 4

- a. Place the left hand on the Stabilizer to stabilize the device at the angle of the tissue tract. Orient the Clip Applier so the front side faces upward (the side with the product logo).
- b. Assume a syringe grip on the proximal end of the device with the right hand by placing the index finger and the middle finger on the proximal posts, and the thumb on the Plunger.
- c. Retract the device out of the tissue 3 to 4 cm with the right hand. Firmly depress the Plunger with the right thumb until the number "2" appears in the number window. This opens and locks the Locator Wings in the expanded position, and advances the Thumb Advancer and the Clip Delivery Tube approximately 3 cm to initiate splitting of the Sheath. The initial splitting of the Sheath can be seen above the skin level.

- d. An audible "click" should be heard. Check to make sure the number "2" is completely visible in the number window. This indicates that the Locator Wings are opened and in the locked position.
7. **CLICK 3:** Advance the Thumb Advancer and deliver the Clip to the access site (**Figure 5**).

Figure 5


- a. Maintain the left hand on the Stabilizer to stabilize the device at the angle of the tissue tract, gently retract the device with the right hand until slight resistance is felt. The goal is to place the Locator Wings against the anterior surface of the arterial wall to locate the access site without applying tension on the artery.
- b. Assume a syringe grip on the distal end of the device with the right hand by placing the index finger on the distal post, the middle two fingers in the loop of the Stabilizer, and the thumb on the Thumb Advancer.
- c. While maintaining apposition of the Locator Wings on the anterior surface of the arterial wall and keeping the Flex-guide straight, advance the Thumb Advancer using the pad of the right thumb until the number "3" appears in the number window. This advances the Clip Delivery Tube over the Flex-guide while splitting the Sheath from the Hub to the distal tip.
- d. An audible "click" should be heard. Check to make sure the number "3" is completely visible in the number window. This indicates that the Thumb Advancer is fully advanced and in the locked position.

Note: Do not advance the Thumb Advancer against excessive resistance.

If excessive resistance is experienced during advancement of the Thumb Advancer, manually collapse the Locator Wings and remove the device following the steps below (**Figure 6**):


Figure 6

- Retract the Thumb Advancer as far proximal as possible until resistance is felt. This is to decrease any interaction between the Clip Delivery Tube and the Flex-guide.
- Slide the Safety Release in the direction of the arrow as seen in **Figure 6** and on the device to collapse the Locator Wings. The Locator Wings are fully collapsed when the number "2" on the Plunger exits the number window completely.
- Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers.
- Provide counter-traction with the left hand on the patient's body, and assertively pull the device out with the right hand.
- It is important not to rock or twist the device as this may cause arterial damage.

8. **CLICK 4:** Deploy Clip to close the arteriotomy (Figure 7).


Figure 7

- While stabilizing the body of the device with the right hand and supporting the Clip Delivery Tube with the left hand, raise the body of the device to a 60 to 75 degree angle.
- Support the Clip Delivery Tube with the left hand and gently push the device down on top of the artery to seat the Clip Delivery Tube on top of the access site. Pulsation of the artery may be felt.
- While maintaining downward pressure and device stability, depress the Deployment Button with the right thumb. An audible "click" should be heard. This releases the Clip to close the arteriotomy and collapses the Locator Wings simultaneously.
- Maintain the downward pressure for 2 – 3 seconds. Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers. Provide counter-traction with the left hand on the patient's body, and remove the device with the right hand.

Note: If resistance is met upon removal of the device, follow the steps below to remove the device (Figure 8):

Figure 8

- Locate the two Access Ports on the lateral and medial sides of the device approximate to the location of the product logo.
- Insert the distal end of the Dilator (or an 18-Gauge Thin-wall Needle) into the Access Ports, one at a time, to release the locking mechanism on the Thumb Advancer.
- An audible "click" should be heard. Check that the number "3" exits the number window completely and the Thumb Advancer is freed from the locked position. Repeat the above steps if necessary.
- Retract the Thumb Advancer as far proximal as possible until resistance is felt. This is to decrease any interaction between the Clip Delivery Tube and the Locator Wings.
- Slide the Safety Release to collapse the Locator Wings and reset the Plunger. The Locator Wings are fully collapsed when the number "2" on the Plunger exits the number window completely.
- Place the left hand on the puncture site in the palm-down position with the Clip Delivery Tube extending up between the index and middle fingers.
- Provide counter-traction with the left hand on the patient's body, and assertively pull the device out with the right hand.
- It is important not to rock or twist the device as this may cause arterial damage.

9. Check for hemostasis.

- Use sterile gauze to firmly depress the area around the puncture site to express any blood in the tissue tract.
- While supporting the patient's leg, have him / her bend his / her knee toward the chest, and have the patient lift his / her head off the table and cough.
- Adjunctive compression can be applied for up to 5 minutes, as needed.
- If hemostasis is not achieved, apply manual compression until hemostasis is achieved.

POST-PROCEDURE PATIENT MANAGEMENT

- Apply an appropriate dressing to the puncture site.
- Assess the insertion site as per hospital protocol.

RECOMMENDATION FOR PATIENT AMBULATION AND DISCHARGE

In determining whether to ambulate or discharge an individual patient, it is important to consider all clinical factors including, but not limited to, anticoagulation regimen, antiplatelet and thrombolytic agents administered, oozing or bleeding from the access site, venous access site hemostasis, the general cardiovascular condition of the patient, anesthetic and medication levels, and the overall clinical condition of the patient.

PRODUCT INFORMATION DISCLOSURE

Abbott Vascular Inc. has exercised reasonable care in the manufacture of this device. Abbott Vascular Inc. excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott Vascular Inc., directly affect this device and the results obtained from its use. Abbott Vascular Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott Vascular Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

HOW SUPPLIED

The StarClose SE Vascular Closure System and accessories are provided sterile and non-pyrogenic in unopened, undamaged packaging. Products are sterilized with ethylene oxide and intended for single use only. This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user. Do not resterilize. The device and the primary packaging do not contain latex. Store in a cool, dry place.

StarClose SE Vascular Closure System contains:

- One (1) StarClose SE Clip Applier
- One (1) 6F Exchange Sheath (11 cm)
- One (1) 6F Dilator
- One (1) 0.038" (0.97 mm) J Tip Guide Wire

PATENTS AND TRADEMARKS

This product and / or its use may be covered by one or more of the following United States Patents: 5,674,231; 5,810,846; 6,197,042; 6,277,140; 6,391,048; 6,461,364; 6,623,510; 6,626,918; 6,632,238; 6,695,867; 6,719,777; 6,749,621; 6,780,197; 6,942,674; 7,001,398; 7,144,411. Other patents pending.

StarClose SE is a trademark of the Abbott Group of Companies.

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©2006, 2013, 2014 Abbott